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LITERATURVERZEICHNIS

- 1 KBV Qualitätsberichte 2011-2019 (Berichtszeiträume 2011-2015).
- 2 Zanetti M et al. HELP LDL apheresis reduces plasma pentraxin 3 in familial hypercholesterolemia. PLoS ONE 2014; 9(7): e101290.
- 3 Mellwig, KP et al. Improved Coronary Vasodilatory Capacity by H.E.L.P Apheresis: Comparing Initial and Chronic Treatment. Therapeutic Apheresis and Dialysis 2006; 10 (6): 510-517.
- 4 Mellwig, KP et al. Heparin-induced Extracorporeal Low density Lipoprotein Precipitation. Therapeutic Apheresis and Dialysis 2003; 7 (3): 365-369.
- 5 Mellwig, KP et al. Myocardial perfusion under H.E.L.P. apheresis measured by positron emission tomography. Z Kardiol 2003 (92); Suppl 3, III/30-III/37.
- 6 Lane DM et al. Selective removal of plasma low density lipoprotein with the H.E.L.P. system: biweekly versus weekly therapy, Atherosclerosis 114 (2), 1995; 203-211
- 7 Heigl F et al. Efficacy, safety, and tolerability of long-term lipoprotein apheresis in patients with LDL- or Lp(a) hyperlipoproteinemia: Findings gathered from more than 36,000 treatments at one center in Germany. Atherosclerosis Supplements Volume 18, May 2015, 154-162.
- 8 Dittich-Riediger J et al. Adverse events of lipoprotein apheresis and immunoadsorption at the Apheresis Center at the University Hospital Dresden. Atherosclerosis Supplements 18, (2015) 45-52.
- 9 Buuren, van F et al. 2012, H.E.L.P. apheresis in hypercholesterolemia and cardiovascular disease: efficacy and adverse events after 8,500 procedures. Clin Res Cardiol Suppl 7:24-30.
- 10 König, W et al., Fibrinogen und kardiovaskuläres Risiko -Pathophysiologie, Klinik und Epidemiologie, Dtsch.Ärztblatt 1994; 91:1563-1567.
- 11 Krobot K et al. Determinants of plasma fibrinogen: relation to body weight, waist-to-hip ratio, smoking, alcohol, age, and sex. Results from the second MONICA Augsburg survey 1989-1990, Survey 1989-1990. Arteriosclerosis and Thrombosis 12:780, 1992.
- 12 Cremer, P et al. Lipoprotein Lp(a) as predictor of myocardial infarction in comparison to fibrinogen, LDL cholesterol and other risk factors: results from the prospective Göttingen Risk Incidence and Prevalence Study (GRIPS) Eur.J.Clin.Invest. 24:444-453, 1994.
- 13 Mellwig KP et al. Acute and Long-Term Effects of LDL Apheresis on Coronary Vasodilation Capacity. 16th Annual Congress of EAMN, Amsterdam 2003.
- 14 Mellwig KP et al. Improvement Of Coronary Vasodilatation Capacity Through Single LDL Apheresis. Atherosclerosis 1998, 139: 173-178.
- 15 Walzl M et al. Effects of ameliorated haemorheology on clinical symptoms in cerebrovascular disease. Atherosclerosis 139 (2), 1998, 385-389.
- 16 Schuff-Werner, P. Clinical long-term results of H.E.L.P.-apheresis. Z Kardiol 2003 (92): Suppl 3, III/28-III/29. 240
- 17 Schuff-Werner, P. Untersuchungen zur hämorheologischen Wirksamkeit der LDL-Apherese. Bibliomed, Medizinische Verlagsgesellschaft mbH, Melsungen 1993.
- 18 Kleophas W et al. Acute effects of extracorporeal LDL cholesterol and fibrinogen elimination on blood rheology and microcirculation. Dtsch Med Wochenschr. 1990 Jan 5; 115(1):3-7.
- 19 Suckfuell M & The Hearing Loss Study Group, Fibrinogen and LDL apheresis in treatment of sudden hearing loss: a randomised multicentre trial. Lancet 2002; 360 (9348): 1811-7
- 20 Stefanutti, C et al. Italian Multicenter Study on Low-Density Lipoprotein Apheresis Working Group 2009 Survey. Therapeutic Apheresis and Dialysis 2013; 17(2):169-178.
- 21 Hagemeyer B, Fibrinogen- und LDL-Apherese zur Behandlung des akuten Hörsturzes im Vergleich zur Standardtherapie nach STENNERT. Dissertation 2013 <http://nbn-resolving.de/urn:nbn:de:bsz:15-qucosa-132656>
- 22 Canis, M et al. Fibrinogen/LDL apheresis is a promising rescue therapy for sudden sensorineural hearing loss. Clin Res Cardiol Suppl (2012) 7:36-40.
- 23 Bianchin, G et al. Treatment with H.E.L.P.-Apheresis in patients suffering from sudden sensorineural hearing loss: A prospective, randomized, controlled study. Laryngoscope. 2010; 120(4):800-7.
- 24 Bianchin, G et al. Role of H.E.L.P.-apheresis in the treatment of sudden sensorineural hearing loss in a group of 230 patients. Acta Otorhinolaryngol Ital 2011; 31:395-398.
- 25 Ali FA et Armogan N Heparin-induced Extracorporeal Lipoprotein Precipitation (H.E.L.P.) Therapy for Dry AMD. Retina Today 2008, 9/10: 72-75
- 26 Moriarty, PM et al. Effect of Low-Density Lipoprotein Cholesterol Apheresis on Blood Viscosity. The American Journal of Cardiology 2004, 93:1044-46.

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B. Braun Deutschland GmbH & Co. KG | Tel. (0 56 61) 9147-70 00 | E-Mail: info.de@bbraun.com | www.bbraun.de