
# Aesculap® Acculan 4

**en Instructions for use/Technical description**

Oscillating saw straight GA336 / Oscillating Sternum saw GA337

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



## Oscillating saw straight GA336 / Oscillating Sternum saw GA337

### Legend

- 1 Product (Oscillating saw straight or Oscillating Sternum saw)
- 2 Trigger (for oscillation frequency control)
- 3 Safety catch
- 4 Battery bay
- 5 Locking pin
- 6 Sterile funnel
- 7 Battery
- 8 Cover
- 9 Cover release element
- 10 Battery removal device
- 11 Saw blade
- 12 Saw head
- 13 Button for saw blade release
- 14 Saw blade coupling

The depictions are schematic only.

### Symbols on product and packages

|                                                                                                                                                                          |                                                                                                                                                                                                                           |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                        | Caution<br>Observe important safety information such as warnings and precautions in the instructions for use.                                                                                                             |
| <br>YYYY-MM                                                                            | Maintenance label<br>Indication of the next maintenance appointment (Date: Year-Month)                                                                                                                                    |
|                                                                                        | Machine-readable, two-dimensional code<br>The code contains a unique serial number which can be used for electronic tracking of the individual instrument. The serial number is based on the global standard sGTIN (GS1). |
|                                                                                        | Date of manufacture                                                                                                                                                                                                       |
|                                                                                        | Follow the instructions for use                                                                                                                                                                                           |
| <br> | Labeling of electrical and electronic devices pursuant to directive 2002/96/EG (WEEE)                                                                                                                                     |
|                                                                                        | Classification Type BF                                                                                                                                                                                                    |
|                                                                                        | Oscillation frequency control                                                                                                                                                                                             |

|                                                                                   |                                                          |
|-----------------------------------------------------------------------------------|----------------------------------------------------------|
|  | Manufacturer's batch designation                         |
|  | Manufacturer's serial number                             |
|  | Manufacturer's article number                            |
|  | Temperature limits during transport and storage          |
|  | Air humidity limits during transport and storage         |
|  | Atmospheric pressure limits during transport and storage |
|  | Operating mode                                           |

## Contents

|       |                                                                                      |    |
|-------|--------------------------------------------------------------------------------------|----|
| 1.    | Applicable to .....                                                                  | 3  |
| 2.    | General information .....                                                            | 3  |
| 2.1   | Intended use .....                                                                   | 3  |
| 2.2   | Main functions and design characteristics .....                                      | 3  |
| 2.3   | Indications .....                                                                    | 3  |
| 2.4   | Absolute contraindications .....                                                     | 3  |
| 2.5   | Relative contraindications .....                                                     | 3  |
| 3.    | Safe handling .....                                                                  | 4  |
| 4.    | Product description .....                                                            | 4  |
| 4.1   | Scope of supply .....                                                                | 4  |
| 4.2   | Components required for operation .....                                              | 4  |
| 4.3   | Operating principle .....                                                            | 4  |
| 5.    | Preparation .....                                                                    | 4  |
| 6.    | Working with the device .....                                                        | 4  |
| 6.1   | System set-up .....                                                                  | 5  |
| 6.1.1 | Connecting the accessories .....                                                     | 5  |
| 6.1.2 | Inserting the rechargeable battery .....                                             | 5  |
| 6.1.3 | Intraoperative battery change .....                                                  | 5  |
| 6.1.4 | Removing the rechargeable battery .....                                              | 5  |
| 6.1.5 | Protection against inadvertent activation .....                                      | 5  |
| 6.1.6 | Attaching and detaching the saw blade .....                                          | 6  |
| 6.2   | Function checks .....                                                                | 6  |
| 6.3   | Safe operation .....                                                                 | 6  |
| 6.3.1 | Operating the product .....                                                          | 6  |
| 7.    | Validated reprocessing procedure .....                                               | 7  |
| 7.1   | General safety notes .....                                                           | 7  |
| 7.2   | General information .....                                                            | 7  |
| 7.3   | Preparations at the place of use .....                                               | 7  |
| 7.4   | Preparation before cleaning .....                                                    | 7  |
| 7.5   | Cleaning/disinfection .....                                                          | 7  |
| 7.5.1 | Product-specific safety instructions for the reprocessing procedure .....            | 7  |
| 7.6   | Manual cleaning with wipe disinfection .....                                         | 8  |
| 7.7   | Automatic cleaning/disinfection with manual pre-cleaning ..                          | 9  |
| 7.7.1 | Manual pre-cleaning with a brush .....                                               | 9  |
| 7.7.2 | Mechanical alkaline cleaning and thermal disinfection .....                          | 9  |
| 7.8   | Inspection, maintenance and checks .....                                             | 10 |
| 7.9   | Packaging .....                                                                      | 10 |
| 7.10  | Steam sterilization .....                                                            | 10 |
| 7.11  | Sterilization for the US market .....                                                | 10 |
| 7.12  | Storage .....                                                                        | 10 |
| 8.    | Maintenance .....                                                                    | 10 |
| 9.    | Troubleshooting list .....                                                           | 11 |
| 10.   | Technical Service .....                                                              | 12 |
| 11.   | Accessories/Spare parts .....                                                        | 12 |
| 12.   | Technical data .....                                                                 | 13 |
| 12.1  | Classification acc. to Directive 93/42/EEC .....                                     | 13 |
| 12.2  | Performance data, information about standards .....                                  | 13 |
| 12.3  | Operating mode .....                                                                 | 13 |
| 12.4  | Environmental conditions .....                                                       | 13 |
| 13.   | Disposal .....                                                                       | 13 |
| 14.   | Distributor in the US/Contact in Canada for product information and complaints ..... | 13 |

## 1. Applicable to

► For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

## 2. General information

### 2.1 Intended use

#### Task/Function

The oscillating saw straight GA336 / Oscillating Sternum saw GA337, combined with the relevant tool, is used for processing hard tissue, cartilage and the like, as well as for bone replacement materials.

#### Application Environment

The product fulfills the requirements for type BF pursuant to IEC/DIN EN 60601-1 and is used in operating rooms in sterile environments of explosion risk areas (such as areas with pure oxygen or anaesthesia gases).

### 2.2 Main functions and design characteristics

|                       |                                                                                         |
|-----------------------|-----------------------------------------------------------------------------------------|
| Oscillation frequency | min. 0 min <sup>-1</sup> to max. 13 000 min <sup>-1</sup>                               |
| Operating mode        | Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1) |
|                       | ■ 30 second application, 30 second pause                                                |
|                       | ■ 7 repetitions                                                                         |
|                       | ■ 30 min cooling time                                                                   |
|                       | ■ Max. Temperature 48 °C                                                                |

Electrical systems generally heat up during continual operation. It is advised to give the system a break after use to cool down, as listed in the table on operating mode.

Heating depends on the tool used and the load. After a certain number of repetitions, the system should cool down. This procedure prevents the system overheating as well as possible injury to the patient or user.

The user is responsible for the use and adherence to the pause sequence described.

### 2.3 Indications

The type and area of application depend on the tool selected.

### 2.4 Absolute contraindications

The product is not licensed for use on the central nervous system or central circulatory system.

### 2.5 Relative contraindications

The safe and effective use of the product greatly depends on influences which can only be controlled by the user. Therefore the specifications provided represent framework conditions only.

Clinically successful use of the product is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

### 3. Safe handling

#### CAUTION

Federal law restricts this device to sale by, or on order of a physician!

**WARNING**

**Risk of injury and material damage if this product is not used as intended!**

► Use the product only for its intended purpose.

**WARNING**

**Risk of injury and damage to property due to improper handling of the product!**

► Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Prior to use, check that the product is in good working order.
- Follow TA022450 "Notes on electromagnetic compatibility (EMC) for Acculan components", see Aesculap Extranet at <https://extranet.bbraun.com>
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
  - Use the product only according to these instructions for use.
  - Follow the safety and maintenance instructions.
  - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.
- Ensure that the electrical installation of the room is consistent with the requirements of IEC/DIN EN.
- Do not operate the product in explosion-hazard areas.
- Sterilize product before use.
- When using the ECCOS fixation system, adhere to the appropriate instructions for use for TA009721, see Aesculap Extranet at <https://extranet.bbraun.com>.

### 4. Product description

#### 4.1 Scope of supply

| Art. no. | Designation                                                                                    |
|----------|------------------------------------------------------------------------------------------------|
| GA336    | Oscillating saw straight                                                                       |
| GA337    | – or –<br>Oscillating Sternum saw                                                              |
| GA678    | Sterile funnel                                                                                 |
| TA014547 | Instructions for use of oscillating saw straight GA336 / Oscillating Sternum saw GA337 (flyer) |

#### 4.2 Components required for operation

- Battery GA676 (charged)
- Sterile funnel GA678
- Cap GA675
- Saw blade (according to indication)

#### 4.3 Operating principle

The product **1** contains an electric motor, which is powered by a replaceable battery **7**.

The charged non-sterile battery **7** is inserted into the product **1** via the sterilized funnel **6** and sealed with the cap **8** so that it is sterile.

The oscillation frequency is controlled electronically and can be continuously regulated with the trigger **2**.

The product **1** has a coupling, which allows for the attachment of various saw blades **11**. These saw blades **11** lock automatically when attached to the product.

The saw blade **11** can be removed by activating a push button **13**.

### 5. Preparation

Aesculap assumes no liability if the following rules are not followed:

- Do not use products from open or damaged sterile packaging.
- Prior to use, inspect the product and its accessories for any visible damage.
- Use the products and their accessories only if they are in perfect technical condition.

### 6. Working with the device

**WARNING**

**Risk of infection and contamination!**

**Product is delivered unsterilized!**

► Sterilize the product before use pursuant to the operating instructions.

**WARNING**

**Risk of injury and material damage due to accidental activation of the product!**

► Products which are not being actively used must be secured against accidental activation (position OFF).

**WARNING**

**Risk of injury and material damage due to inappropriate use of tools!**

► Always follow the safety advice and information given in the instructions for use.  
► When coupling/uncoupling, handle tools with cutting edges with care.

WARNING

**Damage to the product if dropped!**

- Use the products only if they are in perfect technical condition, see Function check

WARNING

**Risk of burns to skin and tissue through blunt tools or if product has not been maintained properly!**

- Use tools only if they are in perfect condition.
- Replace blunt tools.
- Maintain the product properly, see maintenance guide.

## 6.1 System set-up

### 6.1.1 Connecting the accessories

DANGER

**Risk of injury due to unapproved configuration using additional components!**

- Ensure that the classification of all components used is consistent with the classification of the product (such as type BF or type CF).

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

All configurations must fulfill the fundamental standard IEC/DIN EN 60601-1. The person connecting the devices with each other is responsible for the configuration and must ensure that the fundamental standard IEC/DIN EN 60601-1 or relevant national standards.

- Follow the instructions for use of individual accessories.
- Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

### 6.1.2 Inserting the rechargeable battery

CAUTION

**The product does not function when using the wrong battery!**

- Do not use battery GA666 (with red base) or GA346.

- Product 1 with the battery shaft 4 must be turned upwards and the sterile funnel 6 (sterile) attached, see Fig. A.
- Battery 7 (non-sterile) is inserted into the battery shaft 4 by a second (non-sterile) person, see Fig. A.

#### Note

After the battery has been inserted, multiple signals will sound, which indicate that the product is ready to use.

- After insertion of the battery, have the Sterile funnel 6 (non-sterile) removed by a second person.
- Attach the cap 8 (sterile) in such a manner that it clicks into place with both release catches 9.

#### Note

The sterility of the product is only guaranteed with a correctly placed cap.

### 6.1.3 Intraoperative battery change

WARNING

**Risk of injury from hot battery!**

Battery may be hot after use in the machine.

- Remove the battery with the battery removal device and let it cool.

The battery removal device is used for intraoperative battery changes ensuring that sterile conditions are safely maintained.

- Turn the product 1 with the battery shaft 4 upwards.
- Press both safety catches 9 on the cap 8 at the same time and remove the cap 8.
- Attach the sterile battery removal device 10, see Fig. B.
- Shake the product 1 with the attached battery removal device 10 with battery shaft 4 slightly downwards.

The battery 7 will slide easily into the battery removal device 10.

- Give the battery removal device 10 including dead battery 7 to the non-sterile person.
- Insert a charged battery 7, see Inserting the rechargeable battery.

### 6.1.4 Removing the rechargeable battery

WARNING

**Risk of injury from hot battery!**

Battery may be hot after use in the machine.

- Let the battery cool in the machine and then remove.
- or –
- Remove the battery with the battery removal device and let it cool.

CAUTION

**Damage to the battery caused by knocking it against hard objects!**

- Only remove battery by tapping the product on the palm of the hand.

CAUTION

**Damage to, or destruction of the batteries due to processing!**

- Do not sterilize batteries.

Rechargeable battery must be removed after each surgical procedure and prior to processing of the device.

#### Note

The battery removal device can be used for easier battery removal, see Fig. B.

- Turn the product 1 with the battery shaft 4 upwards.
- Press both safety catches 9 on the cap 8 at the same time and remove the cap 8.
- Tightly grasp the battery shaft 4 at the lower end.
- Tap the battery shaft 4 on the palm of the hand until the battery 7 slides out of the battery shaft 4 and can be removed.

### 6.1.5 Protection against inadvertent activation

To prevent inadvertent activation of the product during tool change/alignment of the saw head, the trigger can be locked.

To lock trigger:

- Twist the safety catch 3 to position OFF.

The trigger 2 is blocked and the product 1 cannot be operated.

To unlock trigger:

- Twist the safety catch 3 to position ON.

The trigger 2 is unlocked and the product 1 can be operated.

### 6.1.6 Attaching and detaching the saw blade

WARNING

**Risk of injury when attaching/removing sawblades in the ON position through inadvertent activation of the product!**

- ▶ Only attach/remove saw blades in the OFF position.

#### Coupling the saw blade

- ▶ Secure product 1 against inadvertent activation with the safety catch 3.
- ▶ Insert the saw blade 11 with its coupling end in the slot of saw blade coupling as far as it will go. Ensure that the locking pin engages in the saw blade 11 fenestration, and that the lateral stops of the saw blade 11 rest on the saw blade coupling 14, see Fig. C.
- ▶ If necessary press button for saw blade release/fixation 13.

#### Uncoupling the saw blade

- ▶ Secure product 1 against inadvertent activation with the safety catch 3.
- ▶ Fully press button for saw blade release/fixation 13.
- ▶ Pull out saw blade 11 from the saw blade coupling 14, see Fig. C.

## 6.2 Function checks

The function checks must be carried out prior to each use and after each intraoperative battery change.

- ▶ Make certain there is a battery in the battery compartment.
- ▶ Ensure that the cap has completely clicked into place.
- ▶ Pull on the saw blade to check that it is securely coupled.
- ▶ Check that the cutting edges of the saw blade do not exhibit any mechanical damage.
- ▶ Release product for use (ON position).
- ▶ Briefly operate product at maximum oscillation frequency.
- ▶ Do not use the product if it is damaged or defective.
- ▶ Set aside the product if it is damaged.

## 6.3 Safe operation

WARNING

**Coagulation of patient tissue or risk of burns for patients and user through hot product!**

- ▶ Cool the tool during operation.
- ▶ Keep product/tool out of the reach of patients.
- ▶ Let the product/tool cool down.
- ▶ Use a cloth to protect against burns when changing the tool.

WARNING

**Risk of infection or injury through aerosol formation or particles ejected from the tool!**

- ▶ Use suitable protection (such as waterproof protective clothing, face mask, safety glasses, suction).

WARNING

**Risk of injury and/or malfunction!**

- ▶ Always carry out a function check prior to using the product.

WARNING

**Risk of injury when using the product beyond the field of view!**

- ▶ Apply the product only under visual control.

WARNING

**Risk of injury and damage to the tool/system!**

- The tool may get caught in drapes (such as textiles).
- ▶ Do not let the tool come into contact with drapes (such as textiles) during operation.

WARNING

**Risk of burns to skin and tissue due to blunt saw blades/lack of product maintenance!**

- ▶ Only use perfectly sharp saw blades.
- ▶ Replace blunt saw blades.
- ▶ Maintain the product properly, see maintenance guide.

#### Note

The product drive motor is operated using a magnetic sensor system. In order to prevent inadvertent activation of the motor, the product may not be exposed to any magnetic fields (such as magnetic instrument pads).

### 6.3.1 Operating the product

#### Note

The slight whistling noise when the product is starting up is due to the construction.

- ▶ Activate trigger 2.

The oscillation frequency of the product 4 is continuously variable.

## 7. Validated reprocessing procedure

### 7.1 General safety notes

#### Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

#### Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

#### Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

#### Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemistry was used for validation.

#### Note

If there is no final sterilization, then a virucidal disinfectant must be used.

#### Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bb Braun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

### 7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

On stainless steel, residues containing chlorine or chloride (such as surgical residues, drugs, saline solutions in water for cleaning, disinfection and sterilization) may lead to corrosion (pitting corrosion, tensile corrosion) and thus to the destruction of the product. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Visual material changes (such as fading or color changes in titanium or aluminum). For aluminum, the application/process solution only needs to be pH >8 to cause visible surface changes.
- Material damage (such as corrosion, cracks, breaks, premature aging or swelling).
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see [www.a-k-i.org](http://www.a-k-i.org), link to Publications, Red Brochure – Proper maintenance of instruments.

### 7.3 Preparations at the place of use

- ▶ Remove all attached components from the product (tool and accessories).
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

### 7.4 Preparation before cleaning

- ▶ Prior to the first mechanical cleaning/disinfection process, mount the ECCOS fixations in a suitable tray (e.g. GB243800).
- ▶ Insert the products in the correct position into the ECCOS fixations, see Fig. D.

### 7.5 Cleaning/disinfection

#### 7.5.1 Product-specific safety instructions for the reprocessing procedure

CAUTION

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
  - are approved for plastic material and high-grade steel,
  - do not attack softeners (e.g. in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C during chemical cleaning and/or disinfection.
- ▶ Do not exceed the maximum temperature of 96 °C during thermal disinfection with FD water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

CAUTION

**Damage or destruction of the batteries due to processing!**

- ▶ Protect battery from moisture.

#### Note

The indicated drying time is a guide time only. It must be checked taking into account the specific conditions (e.g. load) and if applicable adjusted.

## 7.6 Manual cleaning with wipe disinfection

| Phase | Step                          | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water quality | Chemicals                             |
|-------|-------------------------------|--------------|------------|--------------|---------------|---------------------------------------|
| I     | Pre-cleaning                  | RT (cold)    | ≥2         | -            | D-W           | Until visually clean                  |
| II    | Cleaning with enzyme solution | RT (cold)    | ≥2         | 0.8          | D-W           | pH-neutral*                           |
| III   | Intermediate rinse            | RT           | ≥5         | -            | D-W           | -                                     |
| IV    | Drying                        | RT           | -          | -            | -             | -                                     |
| V     | Wipe disinfection             | -            | >1         | -            | -             | Meliseptol HBV wipes 50 % Propan-1-ol |
| VI    | Final rinse                   | RT (cold)    | 0.5        | -            | FD-W          | -                                     |
| VII   | Drying                        | RT           | -          | -            | -             | -                                     |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

\* Validated with enzyme cleaner "Cidezyme Johnson & Johnson"

- Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

### Phase I

- Mobilize non-rigid components during cleaning.
- Clean the product under running water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.

### Note

For details on difficult to access surfaces, see *Acculan pre-cleaning and care information TA016000* (available in Aesculap Extranet at <https://extranet.bbraun.com>).

### Phase II

- Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- Spray products with a pH neutral enzyme solution, let soak in for at least 2 minutes and then wipe off.

### Phase III

- Mobilize non-rigid components during cleaning.
- Rinse product under running tap water for at least 5 minutes.
- Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- Contamination should be removed with a lint-free cloth or soft brush moistened with enzyme cleaner.
- Flexible components should be rinsed for 20 seconds with the water pistol (cold water, at least 2.5 bar).
- After manual cleaning, check visible surfaces and areas of flexible components for residues.
- If necessary, repeat the cleaning process (phase 1 to 3).

### Phase IV

- Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

### Phase V

- Wipe all surfaces of the product with a single-use disinfecting wipe.

### Phase VI

- Rinse disinfected surfaces after the prescribed reaction time for at least 1 minute under running demineralized water.
- Drain any remaining water fully.

### Stage VII

- Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

## 7.7 Automatic cleaning/disinfection with manual pre-cleaning

### Note

The cleaning and disinfection device must fundamentally have a tested efficacy (such as FDA approval or CE label pursuant to DIN EN ISO 15883).

### Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

### 7.7.1 Manual pre-cleaning with a brush

| Phase | Step           | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water quality | Chemistry/Note       |
|-------|----------------|--------------|------------|--------------|---------------|----------------------|
| I     | <b>Rinsing</b> | RT (cold)    | –          | –            | D–W           | Until visually clean |
| II    | <b>Brushes</b> | RT (cold)    | –          | –            | D–W           | Until visually clean |

D–W: Drinking water

RT: Room temperature

- Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

#### Phase I

- Mobilize non-rigid components during cleaning.
- Thoroughly clean the product under running water.

#### Phase II

- Mobilize non-rigid components during cleaning.
- Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.

- After manual preparation, check visible surfaces for residue and repeat the pre-cleaning process as needed.

### Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbBraun.com>).

### 7.7.2 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                        | T<br>[°C/°F] | t<br>[min] | Water quality | Chemistry/Note                                                                                                                                                                                                                                                                |
|-------|-----------------------------|--------------|------------|---------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | <b>Pre-rinse</b>            | <25/77       | 3          | D–W           | –                                                                                                                                                                                                                                                                             |
| II    | <b>Cleaning</b>             | 55/131       | 10         | FD–W          | <ul style="list-style-type: none"> <li>■ Concentrate, alkaline: <ul style="list-style-type: none"> <li>– pH ~ 13</li> <li>– &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution <ul style="list-style-type: none"> <li>– pH ~ 11*</li> </ul> </li> </ul> |
| III   | <b>Intermediate rinse</b>   | >10/50       | 1          | FD–W          | –                                                                                                                                                                                                                                                                             |
| IV    | <b>Thermal disinfection</b> | 90/194       | 5          | FD–W          | –                                                                                                                                                                                                                                                                             |
| V     | <b>Drying</b>               | –            | –          | –             | mind. 10 min at max. 120 °C                                                                                                                                                                                                                                                   |

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Insert the product in its proper position in the ECCOS fixation.
- After automatic cleaning/disinfection, check visible surfaces for residues and repeat the cleaning/disinfection process as needed.

## 7.8 Inspection, maintenance and checks

### Note

Aesculap recommends occasional spraying of moving parts (such as knobs, couplings, cover caps) with the Aesculap STERILIT-Power-Systems oil spray.

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Check the product for any damage, abnormal running noise, overheating or excessive vibration.
- ▶ Inspect the saw blades for broken, damaged or blunt cutting edges.
- ▶ Set aside the product if it is damaged.

## 7.9 Packaging

- ▶ Always follow the instructions for use of the respective packaging and storage devices (e.g. instructions for use TA009721 for Aesculap-ECCOS holder system).
- ▶ Insert the products in the correct position into the ECCOS fixations, see Fig. D.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

## 7.10 Steam sterilization

CAUTION

**Damage or destruction of the batteries due to processing!**

- ▶ Do not sterilize the rechargeable battery.

### Note

Remove all attached components from the product (tools, accessories) before sterilization.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Use a validated sterilization method:
  - Steam sterilization using fractional vacuum process
  - Steam sterilizer DIN EN 285 and validated pursuant to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min

When sterilizing multiple products in one steam sterilizer:

- ▶ Ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

## 7.11 Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of  $10^{-6}$ , Aesculap recommends the following parameters:

### Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters\*

| Sterilization method | Temp.         | Time  | Minimum drying time |
|----------------------|---------------|-------|---------------------|
| Prevacuum            | 270 °F/275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

## 7.12 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

# 8. Maintenance

To ensure reliable operation, the product must be maintained in accordance with the maintenance labeling or at least once a year.

YYYY-MM

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

## 9. Troubleshooting list

► Have defective products repaired by Aesculap Technical Service, see Technical Service.

| Issue                                                   | Cause                                                       | Detection                                                                    | Remedy                                                                                                                                                         |
|---------------------------------------------------------|-------------------------------------------------------------|------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product not running                                     | No battery                                                  | No battery in the battery shaft                                              | Insert battery.                                                                                                                                                |
|                                                         | Battery not charged                                         | No signal tone when inserting the battery                                    | Charge battery in charger.                                                                                                                                     |
|                                                         | Battery defective                                           | No signal tone when inserting the battery                                    | Have the manufacturer repair the battery.                                                                                                                      |
|                                                         | Product in locked OFF position                              | Safety catch is in OFF position                                              | Turn the safety catch into the ON position.                                                                                                                    |
|                                                         | Product defective                                           | Product not running                                                          | Have the manufacturer repair the product.                                                                                                                      |
| Product becomes too hot                                 | Excessive use                                               | Product heating                                                              | Follow operating instructions (operating mode).                                                                                                                |
|                                                         | Reprocessing/maintenance carried out incorrectly            | Product heating                                                              | Follow operating instructions (preparation, care).<br>Preventative: lubricate product before every sterilization.                                              |
|                                                         | Fall damage, product defective                              | Product heating                                                              | Have the manufacturer repair the product.                                                                                                                      |
|                                                         | Blunt tool                                                  | Heating of tool and product                                                  | Change tool.                                                                                                                                                   |
| Insufficient power                                      | Product defective                                           | Insufficient product performance                                             | Follow operating instructions (preparation, care).<br>Preventative: lubricate product before every sterilization.<br>Have the manufacturer repair the product. |
|                                                         |                                                             | Large amount of heat generated after a short time                            | Follow operating instructions (operating mode).<br>Have the manufacturer repair the product                                                                    |
|                                                         | Blunt tool                                                  | Tool blades worn out                                                         | Change tool.                                                                                                                                                   |
| Loud running noise                                      | Defective product drive/ball bearings                       | Loud, perceptible noise during operation                                     | Follow operating instructions (preparation, care).<br>Preventative: lubricate product before every sterilization.<br>Have the manufacturer repair the product  |
| Cap cannot be attached/removed                          | Cap not compatible                                          | Cap does not click into place                                                | Use the compatible cap for GA336/GA337.                                                                                                                        |
|                                                         | Cap deformed/defective                                      | Cap is difficult or cannot be mounted/removed                                | Have the manufacturer repair the cap.                                                                                                                          |
|                                                         | Connection to product deformed/defective                    | Cap is difficult or cannot be mounted/removed                                | Have the manufacturer repair the product.                                                                                                                      |
| Battery cannot be inserted/removed                      | Battery not compatible                                      | Battery cannot be completely inserted into the battery shaft                 | Use the compatible battery for GA336/GA337.                                                                                                                    |
|                                                         | Battery deformed/defective                                  | Battery is difficult or cannot be mounted/removed                            | Have the manufacturer repair the battery.                                                                                                                      |
|                                                         | Battery shaft to product deformed/defective                 | Battery is difficult or cannot be mounted/removed                            | Have the manufacturer repair the product                                                                                                                       |
| Sterile funnel or battery removal device not attachable | Sterile funnel or battery removal device not compatible     | Sterile funnel or battery removal device not attachable to the battery shaft | Use compatible sterile funnel or battery removal device for GA336/GA337.                                                                                       |
|                                                         | Sterile funnel or battery removal device deformed/defective | Sterile funnel or battery removal device difficult or not attachable         | Change sterile funnel or battery device.                                                                                                                       |
|                                                         | Battery shaft to product deformed/defective                 | Sterile funnel or battery removal device difficult or not attachable         | Have the manufacturer repair the product.                                                                                                                      |
| Speed control knob not working                          | Product in locked OFF position                              | Safety catch is in OFF position                                              | Turn the safety catch into the ON position.                                                                                                                    |
|                                                         | Speed control knob sticks/is defective                      | Speed control knob not working                                               | Have the manufacturer repair the product.                                                                                                                      |

| Issue                                             | Cause                                         | Detection                                           | Remedy                                                                                      |
|---------------------------------------------------|-----------------------------------------------|-----------------------------------------------------|---------------------------------------------------------------------------------------------|
| Saw blade not moving                              | Loud running noise from motor                 | Defective gearings                                  | Have repaired by manufacturer                                                               |
| Insufficient cutting performance of the saw blade | Saw blade teeth worn                          | Blunt saw blade                                     | Replace saw blade                                                                           |
|                                                   | Battery too weak                              | Insufficient motor power/speed                      | Charge battery                                                                              |
|                                                   | Battery worn/defective                        |                                                     | Replace battery                                                                             |
|                                                   | Saw blade overheated                          | Inadequate chip removal                             | Move the saw back and forth during sawing and ensure sufficient irrigation of the saw blade |
| The saw blade cannot be coupled                   | Coupling components or saw blades deformed    | Deformation                                         | Have repaired by manufacturer                                                               |
|                                                   | -                                             | Incompatible saw blade                              | Insert a compatible Aesculap saw blade, see Accessories/Spare parts                         |
| Saw blade cannot be uncoupled                     | Locking pin not fully releasing the saw blade | Button for saw blade release not pressed completely | Fully press the saw blade release button, see Attaching and detaching the saw blade         |

## 10. Technical Service

**WARNING**

**Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!**

- ▶ **Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.**
- ▶ **Do not modify the product.**

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: [ats@aesculap.de](mailto:ats@aesculap.de)

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042 USA

Aesculap Repair Hotline

Phone: +1 (800) 214 -3392

Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

## 11. Accessories/Spare parts

| Art. no. | Designation                                                                                                 |
|----------|-------------------------------------------------------------------------------------------------------------|
| GA675    | Cover                                                                                                       |
| GA676    | NiMH battery long                                                                                           |
| GA678    | Sterile funnel                                                                                              |
| GA679    | Battery removal device                                                                                      |
| GB600    | STERILIT Power Systems oil spray                                                                            |
| GB498R   | ECCOS fixation for straight machine                                                                         |
| GB487R   | ECCOS fixation for lid                                                                                      |
| GB488R   | ECCOS fixation for sterile funnel                                                                           |
| GB489R   | ECCOS fixation for battery removal device                                                                   |
| TA014546 | Instructions for use of oscillating saw straight GA336 / Oscillating Sternum saw GA337 (A4 for ring binder) |
| TA014547 | Instructions for use of oscillating saw straight GA336 / Oscillating Sternum saw GA337 (flyer)              |

## 12. Technical data

### 12.1 Classification acc. to Directive 93/42/EEC

| Art. no. | Designation              | Class |
|----------|--------------------------|-------|
| GA336    | Oscillating saw straight | Ila   |
| GA337    | Oscillating Sternum saw  | Ila   |

### 12.2 Performance data, information about standards

|                                             |                             |
|---------------------------------------------|-----------------------------|
| Max. power                                  | approx. 250 W               |
| Max. oscillation frequency                  | 13 000 min <sup>-1</sup>    |
| Weight (operational)                        | 1.48 kg ±10 %               |
| Dimensions (L x W x H, ready for operation) | 290 mm x 55 mm x 84 mm ±5 % |
| Applied part                                | Type BF                     |
| EMC                                         | IEC/DIN EN 60601-1-2        |
| Conforming to standard                      | IEC/DIN EN 60601-1          |

The product was subject to 500 preparation cycles as a test by the manufacturer and has passed.

### 12.3 Operating mode

|                |                                                                                                                                                                                      |
|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Operating mode | Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)                                                                                              |
|                | <ul style="list-style-type: none"> <li>■ 30 second application, 30 second pause</li> <li>■ 7 repetitions</li> <li>■ 30 min cooling time</li> <li>■ Max. Temperature 48 °C</li> </ul> |

### 12.4 Environmental conditions

|                      | Operation            | Storage and transport |
|----------------------|----------------------|-----------------------|
| Temperature          | 10 °C to 27 °C       | -10 °C to 50 °C       |
| Relative humidity    | 30 % to 75 %         | 10 % to 90 %          |
| Atmospheric pressure | 700 hPa to 1 060 hPa | 500 hPa to 1 060 hPa  |

## 13. Disposal

### Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

## 14. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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