- GB** Instructions for use/Technical description
- USA** activ C intervertebral disk prosthesis
- D** Gebrauchsanweisung/Technische Beschreibung
- activ C-Bandscheibenprothese
- F** Mode d'emploi/Description technique
- Prothèse de disque intervertébral activ C
- E** Instrucciones de manejo/Descripción técnica
- Prótesis de disco intervertebral activ C
- I** Istruzioni per l'uso/Descrizione tecnica
- Protesi discale activ C
- P** Instruções de utilização/Descrição técnica
- Prótese de disco intervertebral activ C
- NL** Gebruiksaanwijzing/Technische beschrijving
- activ C-tussenwervelschijfprothese
- S** Bruksanvisning/Teknisk beskrivning
- activ C-diskprotes
- RUS** Инструкция по применению/Техническое описание
- Протез межпозвоночного диска актив C
- CZ** Návod k použití/Technický popis
- Protéza meziobratlové ploténky aktiv C
- PL** Instrukcja użytkowania/Opis techniczny
- Proteza krążka międzykręgowego aktiv C
- SK** Návod na použitie/Technický opis
- aktiv C náhrada medzistavcovej platničky
- TR** Kullanım Kılavuzu/Teknik açıklama
- aktiv C disk protezi
- KR** 사용 설명서 / 기술 설명
- 추간체고정보형체 (activ C, 형명개별기재)

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 011995 2017-12 V6 Änd.-Nr. 57994

Aesculap® activ C

activ C intervertebral disk prosthesis

Intended use

The activ C intervertebral disk prosthesis is used for replacing intervertebral disks in the cervical spine (between C3 and C7). It serves to restore the disk height and segmental mobility.

The activ C intervertebral disk prosthesis consists of three components:

- Superior prosthesis plate with spikes for anchoring at the vertebral body
- Inferior prosthesis plate with central anchoring fin for anchoring at the vertebral body
- Polyethylene inlay

The prosthesis plates and the polyethylene inlay together form a ball joint. The polyethylene inlay is anchored in form-fit in the inferior prosthesis plate.

activ C intervertebral disk prostheses are available in six different sizes (XS, S, M, L, XL and XXL) and in up to two different heights (5 mm and 6 mm). The activ C intervertebral disk prostheses are supplied completely pre-assembled.

Materials

The materials used in the implant are listed on the packaging:

- ISODUR®: cobalt-chromium forged alloy CoCrMo according to ISO 5832-12
- PLASMAPORE® surface coating pure titanium acc. to ISO 5832-2
- Ultra-high molecular low-pressure polyethylene acc. to ISO 5834-2

ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for:

- Symptomatic cervical discopathy
- Fresh disk prolapse

Contraindications

Do not use in the presence of:

- Fever
- Acute or chronic vertebral infections of a local or systemic nature
- Systemic diseases and metabolic disorders
- Pregnancy
- Osteoporosis, osteochondrosis, severe osteopenia
- Medical conditions that could preclude a successful implantation (e.g. pre-existing kyphosis)
- Spondylarthrosis
- Chronic degenerative spinal stenosis
- Segmental instability
- Spinal deformity
- Chronic spondylotic myelopathy
- Ossification of the anterior longitudinal ligament (ALL) and/or the posterior longitudinal ligament (PLL)
- Distinct osteophytic and sclerotic changes of the vertebral bodies
- Preoperative immobile segment (mobility less than 2° in flexion/extension)
- Preoperative disk height of the segment to be operated on less than 3 mm
- Narrow, hard disk symptoms
- Poor quality of the articular complex
- Facet joint degeneration
- Vertebral body fracture of the cervical spine
- Tumors of the cervical spine
- Pharmaceutical or other drug abuse; alcoholism
- Foreign body sensitivity to the implant materials
- Inadequate patient compliance
- Cases not listed under indications

Side effects and interactions

- Dislocation, loose, worn or broken implant components
- Loss of fixation, dislocation and migration
- Infections.
- Allergic reactions to the implant materials
- Tissue reaction to implant materials
- Venous thrombosis, lung embolism, cardiac arrest
- Hematomas and wound healing disorders
- Periarticular calcifications and fusions
- Injuries to
 - Nerve roots
 - Spinal cord
 - Blood vessels
 - Organs

Preoperatively existing segmental kyphosis and/or so-called "straight neck" can lead to postoperative kyphotic position of the prosthesis and/or immobility of the segment.

Acute myelopathies can have various reasons. Not all of them can be treated successfully with a prosthesis.

Safety notes

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate aseptics.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Only combine Aesculap modular implant components with each other.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.

WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- Delayed healing can cause implant breakage due to metal fatigue.
- The attending physician shall make any decision with regard to the removal of implant components that have been used, In this case he must be conscious of the potential risks of a further operation as well as of the difficulties of an implant removal.
- Damage to the load-bearing structures of the implant can lead to loosening of components, dislocation, migration, and other severe complications.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.
- Only instruments for activ C intervertebral disk prostheses from Aesculap may be used.
- Follow the instructions for use of trial implants, handle, reamer guide and chisel guide TA011996 and insertion instrument and TA011997.
- Follow the operating manual (no. 031302).

Sterility

- The implant components are supplied pre-assembled in labeled protective packaging.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged

WARNING

Damage to implants caused by processing and resterilization!

- **Do not reprocess or resterilize the implants**

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.
- Select implants according to the indication, preoperative planning and the bone situation found intraoperatively.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The patient is aware of the risks involved in neurosurgery, general surgery, orthopedic surgery and general anesthesia.
- The patient has been informed of the advantages and drawbacks of an activ C intervertebral disk prosthesis, and has been made aware of possible alternative therapeutic methods.
- Excessive load, wear or infection can lead to malfunctions of the activ C intervertebral disk prosthesis.
- The life-span of the implant depends on the patient's body weight.
- The implant components must not be overloaded by extreme strains, hard physical labor or sports.
- Corrective surgery may be necessitated by implant loosening, fracture or loss of correction.
- If revision surgery is carried out, it may not be possible to preserve the segmental movement.
- The patient must undergo medical check-ups of the implant components at regular intervals.

The height indicated on the implant packaging is the total posterior implant height.

The implantation of the activ C intervertebral disk prostheses requires the following application steps:

- Surgical access
- Midline marking
- Disectomy and preparation of the vertebral body end plates
- Application of the distractor and distraction
- Verification of implant size
- Preparation of the fin bed
- Insertion of the activ C intervertebral disk prosthesis

Surgical access

WARNING

Risk of damage to blood vessels, nerves and tissue caused by sharp instruments!

- **Work with extreme care throughout the operation.**

- Establish the surgical access to the affected segment.

Midline marking

WARNING

Mispositioning of the activ C intervertebral disk prosthesis caused by inaccurate marking of the midline!

- **Always mark the midline under X-ray control.**

- Determine the center of the intervertebral disk using the midline marker under X-ray control.
- Apply markings on the superior/ inferior vertebral body.

Dissectomy and preparation of the vertebral body end plates

WARNING

Insufficient anchoring stability of the prosthesis plates due to incomplete removal of cartilage!

- ▶ Remove all cartilage from the vertebral body end plates.

WARNING

Increased risk of migration due to over-preparation of the vertebral body end-plates!

- ▶ Work carefully when preparing the vertebral body end plates.
- ▶ Do not remove the vertebral body end plates.

- ▶ Completely remove the intervertebral disc, using standard instruments.
- ▶ Remove all cartilage from the vertebral body end plates.

Application of the distractor and distraction

WARNING

Compression of the spinal canal and of other posterior elements caused by the distraction forceps inserted too deeply!

- ▶ Always use X-ray control when inserting the distraction forceps.

WARNING

Irritation of the facet joints/strain on spinal cord and nerve roots due to excessive distraction!

- ▶ Extend the distraction forceps slowly and carefully under X-ray control.
- ▶ Observe the tactile response.
- ▶ Use the heights of adjacent (healthy) disks as a guide for the required degree of distraction.

WARNING

Risk of loosening/tear-out of the distraction screws when distracting with the distractor and the distraction screws!

- ▶ Use the distraction forceps for the distraction process.
- ▶ Do not use the distractor and the distraction screws for the actual distraction process.
- ▶ Lock the distractor in extended condition.

- ▶ Place the distractor.
- ▶ Under X-ray control, insert the closed distraction forceps into the intervertebral disk space and extend the distraction forceps.
- ▶ Lock the distractor.
- ▶ Close the distraction forceps and remove it from the intervertebral disk space.

Verification of implant size

WARNING

Compression of the spinal canal and of other posterior elements caused by trial implants inserted too deeply!

- ▶ Use the adjusting wheel to set the depth stop of the trial implant to minimum insertion depth.
- ▶ Always use X-ray control when inserting the trial implants in the intervertebral disk space.

WARNING

Increased risk of migration if the selected implant size is too small!

- ▶ Select the size of the trial implant/implant to achieve maximum coverage of the vertebral body end plates.
- ▶ On the X-ray image, in AP and lateral view, check to ensure that the size, position and alignment of the trial implant are correct.

WARNING

Irritation of the facet joints/strain on spinal cord and nerve roots due to excessive/insufficient distraction!

- ▶ Select the correct trial implant/implant size.
- ▶ Inspect the height of the trial implant in the lateral-view X-ray image.

- ▶ Mount the trial implant of the selected size on the handle, observing the "CRANIAL" and "CAUDAL" markings on the trial implant/handle.
- ▶ Use the adjusting wheel to set the depth stop of the trial implant to minimum insertion depth.
- ▶ Under X-ray control, insert the trial implant with the handle in the intervertebral disk space.
- ▶ Align the markings of the trial implant center to the midline marking.
- ▶ If necessary, correct the depth stop position using the adjusting wheel until the intended position is reached. Ensure that depth stop abuts the anterior edge of the vertebral body.
- ▶ On the intraoperative X-ray image, in AP and lateral view, check the size/height and position of the trial implant.
- ▶ Dismount the handle from the trial implant.

Note

The fin bed can either be prepared by reaming or chiseling.

Fin bed reaming

WARNING

Malpositioning of the activ C intervertebral disk prosthesis due to failure to observe the midline!

- ▶ Check the position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- ▶ Ensure secure positioning of the trial implant during the reaming procedure.

WARNING

Risk of damage to blood vessels, nerves and tissue caused by sharp reamer!

- ▶ Carefully insert the handpiece with the coupled reamer.
- ▶ Do not activate the motor before the handpiece with the coupled reamer is properly positioned in the appropriate guide hole and the respective guide slot of the reamer guide.

WARNING

Risk of damage to blood vessels, nerves and spinal cord caused by reamer applied too deeply!

- ▶ For reaming the holes, apply the reamer only in combination with a correctly selected and positioned trial implant and with the reamer guide coupled.
- ▶ Ensure that depth stop of the trial implant abuts the anterior edge of the vertebral body.
- ▶ Use activ C instruments only.
- ▶ Use instruments according to their intended use.
- ▶ Ream the holes under lateral X-ray control.

- ▶ Mount the reamer guide on the trial implant, observing the "CRANIAL" and "CAUDAL" markings on the reamer guide.
- ▶ Clamp in the activ C reamer in the handpiece.
- ▶ Thread the handpiece with the coupled reamer into the appropriate guide hole and the first guide slot of the reamer guide.
- ▶ Ream the first hole under lateral X-ray control. The reaming depth will be limited by a stop at the reamer guide.
- ▶ Thread the handpiece with the coupled reamer into the appropriate guide hole and the second guide slot of the reamer guide.
- ▶ Ream the second hole under lateral X-ray control. The reaming depth will be limited by a stop at the reamer guide.
- ▶ Remove the handpiece with the coupled reamer, the reamer guide and the trial implant from the intervertebral disk space.
- ▶ Clean the fin bed.

Chiseling the fin bed

WARNING

Malpositioning of the activ C intervertebral disk prosthesis due to failure to observe the midline!

- ▶ Check the position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- ▶ Ensure secure positioning of the trial implant during the chiseling procedure.

WARNING

Risk of damage to blood vessels, nerves and tissue caused by sharp chisel!

- ▶ Carefully insert the chisel holder with the coupled chisel.

WARNING

Risk of damage to blood vessels, nerves and spinal cord caused by chisel applied too deeply!

- ▶ For chiseling the holes, apply the chisel only in combination with a correctly selected and positioned trial implant and with the chisel guide coupled.
- ▶ Ensure that depth stop of the trial implant abuts the anterior edge of the vertebral body.
- ▶ Use activ C instruments only.
- ▶ Use instruments according to their intended use.
- ▶ Chiseling under lateral X-ray control.

- ▶ Mount the chisel guide on the trial implant, observing the "CRANIAL" and "CAUDAL" markings on the chisel guide.
- ▶ Mount the activ C chisel onto the chisel holder.
- ▶ Chisel under lateral X-ray control. The chiseling depth will be limited by a stop at the chisel guide.
- ▶ Remove the handpiece with the coupled chisel, the chisel guide and the trial implant from the intervertebral disk space.
- ▶ Clean the fin bed.

Insertion of the activ C intervertebral disk prosthesis

WARNING

Damage to vertebral body end plates when the activ C intervertebral disk prosthesis is tapped in!

- ▶ Apply proper care when tapping in the activ C intervertebral disk prosthesis.
- ▶ Do not remove the vertebral body end plates.

WARNING

Compression of the spinal canal and of other posterior elements caused by the activ C intervertebral disk prosthesis inserted too deeply!

- ▶ Prior to inserting the activ C intervertebral disk prosthesis, adjust the depth guide to minimum insertion depth.
- ▶ Always introduce the activ C intervertebral disk prosthesis into the intervertebral space under X-ray control.
- ▶ On the X-ray image, in AP and lateral view, check for correct position/alignment of the activ C intervertebral disk prosthesis.

CAUTION

Damage to the activ C intervertebral disk prosthesis due to improper handling!

- ▶ Mount the implant on the insertion instrument immediately after removing it from the implant packaging.
- ▶ Handle with caution the friction surfaces of the prosthesis plates/polyethylene inlay.

- ▶ Put out the implant ready for use.
- ▶ Connect the spacer corresponding to the selected implant height to the insertion instrument, observing the "CRANIAL" and "CAUDAL" markings on the spacer/insertion instrument.
- ▶ Use the adjusting wheel to set the depth stop of the insertion instrument to minimum insertion depth.
- ▶ Mount the intervertebral disk implant on the insertion instrument immediately after removing it from the implant packaging, observing the "CRANIAL" and "CAUDAL" markings. The superior prosthesis plate must abut at the part of the insertion instrument that is marked "CRANIAL", while the inferior prosthesis plate must be at the "CAUDAL" part.
- ▶ Turn the clamping sleeve to clamp in the intervertebral disk prosthesis on the instrument.
- ▶ Remove the disk prosthesis from its packaging, holding it with the insertion instrument.
- ▶ Carefully, under X-ray control, insert the disk prosthesis into the intervertebral disk space. When doing this, the superior prosthesis plate must be oriented towards cranial, the inferior prosthesis plate towards caudal.
- ▶ If necessary, correct the depth stop position using the adjusting wheel until the intended position is reached.
- ▶ Release the clamping sleeve on the insertion instrument and remove the insertion instrument.
- ▶ On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the intervertebral disk implant.

Additional information on Aesculap implant systems can be obtained at any time from B. Braun/Aesculap or from designated B. Braun/Aesculap subsidiaries.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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