
Aesculap Spine

-
- | | |
|------------|----------------------------------|
| GB | Instructions for use |
| USA | activ C insertion instrument |
| D | Gebrauchsanweisung |
| | activ C-Einsetzinstrument |
| F | Mode d'emploi |
| | Instrument d'insertion activ C |
| E | Instrucciones de manejo |
| | Instrumento de inserción activ C |
| I | Istruzioni per l'uso |
| | Strumento inseritore activ C |
| P | Instruções de utilização |
| | Instrumento introdutor activ C |
| NL | Gebruiksaanwijzing |
| | activ C-inzetinstrument |
| CZ | Návod k použití |
| | activ C-nasazovací instrument |

activ C insertion instrument

Legend

- 1 Insertion instrument
- 2 Spacer
- 3 Key (for insertion instrument)
- 4 Repositioning instrument
- 5 Irrigation connector

Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The insertion instrument is used for picking up the activ C intervertebral disk prosthesis from its sterile packaging, inserting it in the prepared intervertebral space, and carrying out position corrections of the prosthesis.

Available sizes

Product no.	Description
FW857R	Insertion instrument with irrigation connector H6
FW863R	Spacer for implant height 5 mm
FW864R	Spacer for implant height 6 mm
FW866R	Insertion instrument with irrigation connector H5
FW867R	Repositioning instrument
FW945R	Key for insertion instrument

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only for its intended purpose, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or using a machine, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Exercise caution when passing products through the working channel (e.g. trocar).

Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

Note

Follow OP-Manual no. O31301 and instructions for use "activ C intervertebral disk prosthesis" TA011995.

Assembling the instrument

Note

Account for the size of the implant when selecting the insertion instrument!

- Remove irrigation connector 5 from the handle end of insertion instrument 1.
- Verify that clamping sleeve 7 on insertion instrument 1 has been loosened.
- Turn clamping sleeve 7 counterclockwise until it is no longer in connection with connecting piece 6 of insertion instrument 1.
- Connect spacer 2 for the intraoperatively selected implant height with insertion instrument 1.
- Insert spacer 2 in the hole of the connecting piece, observing the "CRANIAL" and "CAUDAL" markings on spacer 2 and insertion instrument 1. At the same time keep button 9, "push to clean", at the handle of insertion instrument 1 pressed down.

- Insert spacer 2 down to the stop in insertion instrument 1.


Fig. 1

- Release "push to clean" button 9. By axially pulling at spacer 2, check that the spacer is properly locked in insertion instrument 1.
- Adjust the depth stop to the minimum insertion depth by turning the wheel 8 counterclockwise.

activ C insertion instrument

Coupling the implant

WARNING

Damage to the implant due to improper handling!

- Observe the instructions for use of activ C intervertebral disk prosthesis TAO11995.
- Mount the implant on the insertion instrument immediately after removing it from the implant packaging.
- Handle the friction surfaces of the prosthesis plates and of the polyethylene inlay with care.

- Open the implant packaging.

Note

Do not remove the implant from the packaging!

- Mount insertion instrument **1** on the disk prosthesis with the latter still in its packaging. Observe the "CRANIAL" and "CAUDAL" markings. The superior prosthesis plate (the one with the spikes) must abut the part of insertion instrument **1** marked "CRANIAL", while the inferior prosthesis plate (with the fin) has to abut the part of insertion instrument **1** marked "CAUDAL".
- Turn clamping sleeve **7** clockwise to fixate the intervertebral disk prosthesis on insertion instrument **1**.
- Remove the disk prosthesis from its packaging, holding it with insertion instrument **1**.

Inserting the implant

WARNING

Damage to vertebral body end plates when the intervertebral disk prosthesis is tapped in!

- Apply proper care when tapping in the intervertebral disk prosthesis.
- Carefully, under X-ray control, insert the disk prosthesis into the intervertebral disk space. The superior prosthesis plate (plate with the spikes) must be oriented cranially, the inferior prosthesis plate (the plate with the fin) caudally.
- If necessary impel the implant with light hammer taps on the proximal end of insertion instrument **1**.

WARNING

Compression of the spinal canal and of other posterior elements caused by the intervertebral disk prosthesis inserted too deeply!

- Prior to inserting the intervertebral disk prosthesis, adjust the depth guide to minimum insertion depth.
- Introduce the disk prosthesis into the intervertebral space under X-ray control.
- If necessary, correct the depth stop position with adjusting wheel **8** until the intended position is reached.
- Should repositioning become necessary, screw repositioning instrument **4** onto the proximal end of insertion instrument **1** and carefully, with a slotted hammer, move the implant to its intended position.

Decoupling the implant

Note

To avoid damage to the implant, do not couple an insertion instrument that already has been decoupled once.

- Loosen clamping sleeve **7** on insertion instrument **1** by turning it counterclockwise, and remove insertion instrument **1**.
- If clamping sleeve **7** is difficult to loosen, insert key **3** in one of the lateral drill holes of clamping sleeve **7** and loosen clamping sleeve **7** in this way.
- On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the intervertebral disk implant.

Disassembly

- Uncouple spacer **2** from insertion instrument **1** by pressing "push to clean" button **9** at the handle of insertion instrument **1** and simultaneously pulling it out, axially, towards the working end.
- If applicable, unscrew repositioning instrument **4** from insertion instrument **1** by turning it counterclockwise.
- Prior to processing make sure that clamping sleeve **7** sits loosely on the instrument shaft.

Validated reprocessing procedure

Note

National laws, national and international standards and directives, and product-specific hygiene regulations for reprocessing must be observed.

Note

Machine reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing procedure. The operator/sterile processing technician is responsible for this.

Due to process tolerances, manufacturer's specifications can only serve as a guideline for the processing procedures applied by the individual user.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

activ C insertion instrument

General information

To prevent increased contamination of loaded instrument trays during use, please ensure that contaminated instruments are collected separately and not returned to the instrument tray.

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion of stainless steel. You should therefore leave no more than 6 hours between use and cleaning, not use pre-clean fixing temperatures >45 °C, and ensure that non-fixing disinfectants (aldehyde/alcohol) are used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine in the case of stainless steel.

Residues containing chlorine or chloride – e.g. in surgical residues, drugs, saline solutions and water for cleaning, disinfection and sterilization – may cause corrosion damage to stainless steel (pinholing, stress corrosion), thus rendering the products unusable. These must be removed by rinsing thoroughly with demineralized water and then drying.

Only use process chemicals which have been tested and approved (e.g. VAH/DGHM or FDA-certified or CE marking) and recommended by the chemical manufacturer as being compatible with materials. All the chemical manufacturers' application specifications regarding temperature, concentration and contact time should be strictly observed. Failure to do so may result in the following problems:

- Visible changes to materials e.g. fading or changes in the color of titanium or aluminum. As regards aluminum, visible changes to the surface may already occur at a pH level of >8 in the application/working solution.
- Material damage such as corrosion, cracks, fractures, premature deterioration or swelling.
- Do not use process chemicals which cause stress cracks or brittleness in plastics.
- Clean the product immediately after use.

Please see www.a-k-i.org for more detailed information on hygienically safe reprocessing which is protective of materials and retains instruments' value.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and diminished effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 30 minutes.

Preparation before cleaning

- Disassemble the product prior to cleaning, see Disassembly.

Note

Clamping sleeve 7 must sit loosely on the instrument shaft during cleaning.

Note

To ensure proper cleaning of spacers 2, the PEEK forks must sit loosely (so that they can be turned in both directions) on the thread.

Cleaning/disinfection

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for use e.g. on aluminum, plastic materials and stainless steel,
 - do not attack softeners (e.g. silicone).
 - Observe specifications regarding concentration, temperature and exposure time.
 - Do not exceed the maximum allowable cleaning temperature of 60 °C.
-
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfection.

Manual cleaning/disinfection

- When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.
- After manual cleaning/disinfection, check accessible surfaces visually for residues.
- Where necessary, repeat the cleaning process.

activ C insertion instrument

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
IV	Final rinse	RT (cold)	1	-	DM-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

DM-W: Demineralized water

RT: Room temperature

Phase I

- Completely immerse the instrument in the cleaner/ disinfectant solution. Ensure that all accessible surfaces are moistened.
- Clean the instrument under running tap water with a suitable cleaning brush where necessary for as long as it takes to remove all discernable residues.
- Use a suitable cleaning brush to clean all surfaces which are not visible, e.g. in instruments with concealed crevices, lumens or complex geometries, for at least 1 min or until no further residues can be removed.
- Maneuver non-rigid components, such as adjustable screws, hinges etc during cleaning.
- Then flush these areas thoroughly at least five times with the cleaning solution using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Phase II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Drain any remaining water fully.

Phase III

- Fully immerse the instrument in the disinfectant solution.
- Rinse the lumens at least five times at the beginning of the contact period using a disposable syringe (20 ml) and a suitable irrigation adapter. Ensure that all accessible surfaces are moistened.
- Drain any remaining water fully.

Phase IV

- Rinse/flush the instrument thoroughly (all accessible surfaces).
- Rinse the lumens at least five times using a disposable syringe (20 ml) and a suitable irrigation adapter.
- Drain any remaining water fully.

Phase V

- Dry the instrument with a lint-free cloth or medical compressed air.

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note

For thermal disinfection, always use demineralized water. Ensure that A_0 is $>3\,000$ for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

Note

For mechanical cleaning: Install the irrigation connector on the handle end of the insertion instrument.

activ C insertion instrument

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

Phase I

- Completely immerse the instrument in the cleaner/ disinfectant solution. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all discernable residues have been removed.
- Use a suitable cleaning brush to clean all surfaces which are not visible, e.g. in instruments with concealed crevices, lumens or complex geometries, for at least 1 min or until no further residues can be removed.
- Maneuver non-rigid components, such as adjustable screws, hinges etc during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Phase II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

- Place the instrument in a tray suitable for washing (to ensure that all parts are cleaned).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Place instruments in the tray with their hinges open.

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	VE-W	BBRAUN HELIMATIC CLEANER alkaline with tensides; application solution 0.5%
III	Intermediate rinse	>10/50	1	DM-W	-
IV	Thermal disinfection	90/194	5	DM-W	-
V	Drying	-	-	-	According to disinfectant program

D-W: Drinking water

DM-W: Demineralized water

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. hinges, sliding components and threaded rods) with an oil which is suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I Oil Spray JG600 or STERILIT® I Drip Lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is: dry, clean, operational, and free of damage (e.g. insulation, corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.

activ C insertion instrument

- Assemble dismountable instruments, see Assembling the instrument.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization

Note

The product may only be sterilized in disassembled condition.

Note

The irrigation connector must be removed from the handle end of the insertion instrument prior to sterilization.

- Ensure that the sterilant can reach all external and internal surfaces (e.g. by opening valves and stopcocks).
- Validated sterilization process
 - Disassemble the instrument
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer specified in DIN EN 285 and validated in accordance with DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer: ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F– 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- **Do not modify the product.**

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems, LLC.
3773 Corporate Parkway
Center Valley, PA 18034
USA

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technical Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems, LLC.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

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Aesculap AG