
Aesculap Spine

GB	Instructions for use
USA	activ C handle, reamer guide, trial implants
D	Gebrauchsanweisung activ C-Handgriff, Fräserführung, Probeimplantate
F	Mode d'emploi Poignée, guidage de fraise, implants d'essai activ C
E	Instrucciones de manejo Mango, guía de fresado, implantes de prueba activ C
I	Istruzioni per l'uso Impugnatura, guida per la fresa, impianti di prova activ C
P	Instruções de utilização Cabo, guia de fresa, implantes de ensaio activ C
NL	Gebruiksaanwijzing Activ C-handgreep, freesgeleider, proefimplantaten
CZ	Návod k použití activ C-rukojet', vedení frézy, zkušební implantáty

activ C handle, reamer guide, trial implants

Legend

- 1 Trial implant
- 2 Handle
- 3 Reamer guide
- 4 Mounting block (for reamer handpiece)
- 5 Reamer
- 6 Reamer handpiece
- 7 HiLAN pneumatic motor

Symbols on product

	Attention, see instructions for use
	Position "closed"
	Position "open"
	Ready

Intended use

The trial implants are used in conjunction with the handpiece to verify the size, height and shape of the activ C intervertebral disk prosthesis to be implanted, and in conjunction with the reamer guide to prepare the fin bed for the inferior implant plate of the activ C intervertebral disk prosthesis.

Available sizes

Art. no.	Designation
FW789R	Trial implant flat size XS, H5
FW790R	Trial implant flat size XS, H6
FW791R	Trial implant flat size S, H5
FW792R	Trial implant flat size S, H6
FW793R	Trial implant flat size M, H5
FW794R	Trial implant flat size M, H6
FW795R	Trial implant flat size L, H5
FW796R	Trial implant flat size L, H6
FW797R	Trial implant flat size XL, H5
FW798R	Trial implant flat size XL, H6
FW799R	Trial implant flat size XXL, H5
FW800R	Trial implant flat size XXL, H6
FW870	Handle
FW871R	Reamer guide (with mounting block)
FW874R	Trial implant size XS, H5
FW875R	Trial implant size XS, H6
FW876R	Trial implant size S, H5
FW877R	Trial implant size S, H6
FW878R	Trial implant size M, H5
FW879R	Trial implant size M, H6
FW880R	Trial implant size L, H5
FW881R	Trial implant size L, H6
FW882R	Trial implant size L, H7
FW883R	Trial implant size XL, H5

Art. no.	Designation
FW884R	Trial implant size XL, H6
FW885R	Trial implant size XL, H7
FW886R	Trial implant size XXL, H5
FW887R	Trial implant size XXL, H6
FW888R	Trial implant size XXL, H7
GB771R	Hi-Line XS handpiece angled XL-I
GA529	HiLAN pneumatic motor short with hand control
GA505R	Compressed-air tube Aesculap-Dräger/Aesculap small
GA506R	Compressed-air tube Schrader + diffuser/Aesculap small
GA507R	Compressed-air tube DIN-union/Aesculap small
GE700SU	Hi-Line XS activ C drill XL-I D1, 5

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working tip: Carefully insert the product through the working channel (e.g. trocar).

activ C handle, reamer guide, trial implants

Safe operation

Note

Follow OP manual no. 031302 and instructions for use "activ C intervertebral disk prosthesis" TA011995.

Mounting the trial implant of the selected size on the handle

Note

Height 5 and height 6 trial implants are also available as flat implants marked with the inscription "FLAT".


Fig. 1 Handle for trial implants in open (left) and closed (right) position

- Turn locking pipes **8** of handle **2** to position "open" .

- Mount the trial implant of the selected size, height and shape on handle **2**.
Note markings "CRANIAL" and "CAUDAL" on the handle and on the trial implant.
- Turn locking pipes **8** of handle **2** to position "closed" .
- Check the trial implant for firm seating.

WARNING

Compression of the spinal canal and of other posterior elements caused by trial implants inserted too deeply!

- Prior to insertion of trial implants, adjust the depth stop to minimum insertion depth.
- Insert trial implants in the intervertebral space under X-ray control.
- Use adjusting wheel **10** to set the depth stop of the trial implant to minimum insertion depth.
- Insert the trial implant in the intervertebral space under X-ray control. Ensure correct cranial and caudal alignment of the trial implant.
- Align the markings of the trial implant center to the midline marking.
- If necessary, correct the position of the depth stop with adjusting wheel **10** until the intended position is reached.
Ensure that depth stop abuts the anterior edge of the vertebral body.

WARNING

Increased risk of migration if the selected implant size is too small!

- Select the trial implant or implant size, respectively, to achieve maximum coverage of the vertebral body endplates.
- On the X-ray, in AP and lateral view, check to ensure that the size, position and alignment of the trial implant are correct.

WARNING

Irritation of the facet joints and strain on spinal cord and nerve roots due to excessive or insufficient distraction!

- Select the correct trial implant/ implant height.
- Check for correct height of the trial implant in a lateral X-ray.
- On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the trial implant.
- With the trial implant properly positioned, turn locking pipes **8** to position "open"  and pull off handle **2** carefully in axial direction from the trial implant.

The trial implant remains in the intervertebral space.

Reaming the fin bed/ Mounting the reamer guide

Note

The handle and the reamer guide are fitted with identical locking mechanisms.

- Turn locking pipes **8** of reamer guide **3** to position "open" .
 - Mount reamer guide **3** on the trial implant sitting in the intervertebral space.
Note markings "CRANIAL" and "CAUDAL" on the reamer guide and on the trial implant.
 - Turn locking pipes **8** of reamer guide **3** to position "closed" .
- Reamer guide **3** is fixed to the trial implant.
- Check again, by X-ray imaging, the position of the trial implant.

activ C handle, reamer guide, trial implants

Reaming the fin bed/Mounting the reamer handpiece

WARNING

Damage to nerves, spinal cord and other tissue because of fin bed reamed too deeply due to using too long a reamer handpiece or reamer.

- Only use reamer handpiece Hi-Line XS XL-1 (GB771R).
- Only use activ C reamer (GE700SU).
- Prior to mounting, verify that the safety stop of the pneumatic motor is in position "0".
- Assemble the compressed-air tube, the HiLAN pneumatic motor and reamer handpiece 6 as described and illustrated in the respective instructions for use supplied with the products.

Fig. 2

- Keep pressed locking clip 9 of mounting block 4 and push the clip over reamer handpiece 6, see Fig. 3.


Fig. 3 Mounting block with pressed locking clip

- Slide on mounting block 4 as far as it will go and release locking clip 9, see Fig. 4.

Fig. 4 Mounting block with locking clip engaged

- Pull at mounting block 4 to check whether it is engaged properly.
The arrow marking on locking clip 9 has to point to the pictogram "closed" .
- Mount reamer 5 as described in the instructions for use of reamer handpiece 6.

WARNING

Risk of injury to blood vessels, nerves, spinal cord or other tissue if wrong reamer is used!

- Only use activ C reamer GE700SU.
- Check for technical faults prior to use.
- Do not use bent reamers.
- Do not reuse reamers after first use.

Note

The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk. The fin bed must only be prepared for the inferior implant plate.

WARNING

Mispositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- Make certain that the position of the trial implant does not change during reaming.

WARNING

Risk of injury caused by sharp chisel!

- To avoid injuries to critical blood vessels, nerves and other tissue, as well as injuries to the operating surgeon, always work with appropriate care when a reamer is clamped in the machine.
- Do not activate the motor before the handpiece with the coupled reamer is properly positioned in the guide hole and the guiding slot of the reamer guide.

- **A:** Position reamer handpiece **6** with the coupled reamer through the frontal guide hole of reamer guide **3**, see Fig. 5.
- **B:** Position the laterally protruding pins of mounting block **4** through the rear guiding slots of reamer guide **3**, see Fig. 5.


Fig. 5 Sliding the reamer handpiece into the reamer guide

- With the hand control not actuated, push the safety stop to  .
The system is ready for operation.

activ C handle, reamer guide, trial implants

Note

The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk.

WARNING

Mispositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- Make certain that the position of the trial implant does not change during reaming.

WARNING

Risk of injury to nerves and spinal cord caused by too deep insertion of the reamer!

- Only use activ C instruments.
- Use the instruments according to their intended use.
- For reaming the holes, apply the reamer only in combination with a correctly selected and positioned trial implant and with the reamer guide coupled.
- Take care that the depth stop of the trial implant abuts the anterior edge of the vertebral body.
- Always maintain lateral X-ray control while reaming the holes.

- To ream the first hole under lateral X-ray control: With the motor running, push reamer handpiece **6** at an even speed and with mild forward pressure into the first guiding slot down to the stop.
- Pull back reamer handpiece **6** and, after moving the pins of mounting block **4** to the second guiding slot, repeat the reaming procedure.
The reaming depth will be limited by the stop at reamer guide **3**.
- With the hand control not actuated, shift the safety stop at the pneumatic motor to position "0".
- Cautiously pull out reamer handpiece **6** first from the rear and then from the frontal reamer guide **3**.
- Remove reamer guide **3** and the trial implant from the vertebral disk space.
- Clean the fin bed with nerve hooks from the discectomy set.

Disassembling

Uncoupling the reamer guide/handle from the trial implant

- Turn locking pipes **8** of reamer guide **3** or handle **2**, respectively, to position "open" .
- Pull off the trial implant in axial direction.
- Decouple reamer **5** from reamer handpiece **6** as described in instructions for use TA011267.

Detaching the mounting block from the reamer handpiece

- Hold locking clip **9** of mounting block **4** pressed down (arrow marking points to position "open" ) and pull off the mounting block from reamer handpiece **6**.

Validated processing procedure

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, manufacturer's specifications can only serve as guide values for the processing procedures applied by the individual user.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Single-use products

CAUTION

Possible damage to/or destruction of products, due to processing!

Intended use: These products have been designed for single use only.

➤ **Do not reprocess these products.**

General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Consequently the time interval between application and processing should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures > 45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) should be used.

Excessive doses of neutralizers or basic detergents can cause chemical deterioration and/or fading of the laser engraving on stainless steel.

Residues containing chlorine or chlorides, as contained e.g. in surgical residues, tinctures, medicines, saline solutions, service water used for cleaning, and in cleaning/disinfecting agents, cause corrosion damage (pitting, stress corrosion) and thereby destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals recommended by the chemicals' manufacturers as effective cleaning/disinfecting agents compatible with the materials used in the respective product. Strictly observe all application instructions, e.g. regarding temperatures, concentrations, exposure times, etc. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes or

activ C handle, reamer guide, trial implants

- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use process chemicals that cause stress cracking of plastics such as PPSU, or attack and cause brittleness e.g. of silicone softeners.
- Clean the product immediately after use. Encrustations can be removed by hand, with a soft brush and/or in an ultrasound bath. After that treatment, proceed through the normal steps of processing.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications, Red Brochure – Proper maintenance of instruments.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove visible residues as completely as possible, using a lint-free, single-use cleaning tissue.

Preparation prior to cleaning

- Have the product dry in a disposal container and ready for immediate cleaning and disinfecting within 30 min after use.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use a cleaning and disinfecting agent according to manufacturer's instructions. The cleaning and disinfecting agent – must be approved for the product materials (e.g. aluminum, plastics, high-grade steel) – must not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 60 °C.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Manual cleaning/disinfecting

- Clean hinged or jointed products in open and closed positions.

Note

Clean reamer guide **3** and handle **2** with locking pipes **8** in closed position () and do not move them during cleaning.

- Check visible surface for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfecting and cleaning with brush

FW789R–FW800R, FW874R–FW888R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Final rinse	RT (cold)	0.5	–	FD–W	–
IV	Drying	RT	–	–	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized)

RT: Room Temperature

activ C handle, reamer guide, trial implants

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives, which would damage the surfaces, for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage IV

- Dry the product completely with compressed air.

Manual cleaning with immersion disinfecting and cleaning with brush and ultrasound

FW870/FW871R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting ultrasound cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	1 st Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
IV	2 nd Intermediate rinse	RT	1	-	D-W	-
V	Intermediate drying	-	-	-	-	-
VI	Final rinse	RT (cold)	0.5	-	FD-W	-
VII	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.

activ C handle, reamer guide, trial implants

- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use metal brushes or other abrasives for cleaning the product.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.

Stage V

- Dry the product with compressed air.

Stage VI

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage VII

- Dry the product completely with compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

Manual pre-cleaning with brush

FW789R–FW800R, FW874R–FW888R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	–	D–W	–

D–W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.

- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives that would damage the surfaces for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

activ C handle, reamer guide, trial implants

Manual pre-cleaning with ultrasound

FW870/FW871R

Note

Clean reamer guide **3** and handle **2** with locking pipes **8** in closed position () and do not move them during cleaning.

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting ultrasound cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

FW789R–FW800R, FW870, FW871R, FW874R–FW888R

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.
- Store any products with hinges or joints on the tray in such a way that the joints are kept open.

Note

Clean reamer guide **3** and handle **2** with locking pipes **8** in closed position () and do not move them during cleaning.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	– Concentrate, alkaline: pH = 10.9 <5 % anionic tensides – 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD–W	– Concentrate, acid: pH = 2.6 Basis: Citric acid – 0.2 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD–W	–
V	Thermal disinfecting	94/201	10	FD–W	–
VI	Drying	90/194	40	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized)

activ C handle, reamer guide, trial implants

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts (e.g. joints and links) with a maintenance oil suitable for the respective sterilization process (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Reassemble separable products.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside the product if it is damaged.

Packaging

- Appropriately protect products with fine working tips.
- Secure product with retractor in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

Note

The product may only be sterilized in disassembled condition.

Note

To prevent breakage due to stress crack corrosion, sterilize the instruments with the locking mechanism open, or locked on the first ratchet tooth.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer acc. to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST 46-1993, validated acc. to EN 554/ISO 13683
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270— 275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store processed products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

Technical Service

For service, maintenance and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG