Aesculap® activ L

Aesculap Spine

GB Instructions for use/Technical description
activ L insertion instrument

USA Instructions for use/Technical description
activ L insertion instrument
Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <https://www.aesculapimplantsystems.com/products/instructions-for-use>. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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
Symbols on product and packages

Caution, general warning symbol

Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

Intended use

The activ L insertion instrument is used for implanting the intervertebral disc prosthesis through anterior (0°) or lateral (approx. 45°) approach. The appropriate color-coded and labeled insertion instrument is used according to the implant height.

Available sizes

Art. no.	Size
FW961R	8.5 mm
FW962R	10 mm
FW963R	12 mm
FW964R	14 mm

Safe handling and preparation

Use of the activ L insertion instrument requires precise knowledge about spine surgery and the stabilization and bio-mechanical situation at the spine.

Operative use of the activ L insertion instrument is described in detail in the relevant O.R. manual.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Coupling the implant

- Select the appropriate insertion instrument according to the height of the intervertebral disc prosthesis.
- Completely assemble the intervertebral disc prosthesis.
- Extend spacer 1 as far as possible. If necessary, change the position of the spacer by turning twist cap 11.
- Put the intervertebral disc prosthesis onto clamp 2 of the insertion instrument, observing the marking on the insertion instrument. The inferior prosthesis plate 9 must abut the part of the insertion instrument clamp 2 marked "CAUDAL"; the superior prosthesis plate 7 must abut the part marked "CRANIAL".
- Guide the prosthesis plates on the spacer 1 onto coupling pins 5 and 6 of the clamp 2.
- Make certain that the coupling pins of the clamp fully engage in the coupling holes of the prosthesis plates.
- Turn clamping sleeve 3 clockwise until the intervertebral disc prosthesis is securely couples.

Decoupling the implant

- Release clamping sleeve 3.
- Carefully detach the insertion instrument from the implant.
- If the insertion instrument cannot be detached without applying force: Turn twist cap 11 anticlockwise to move back spacer 1.

Disassembling

- Turn twist cap 11 counterclockwise until it can be removed.
- Hold unlocking button 10 pressed down.
- Pull out spacer 1 to the front.

Assembling

- Insert clamp 1 from the front through clamp 2 and, with unlocking button 10 pressed down, push in the clamp as far as it will go.
- Release unlocking button 10.
- Attach twist cap 11 and screw it on as far as it will go.
- To position spacer 1 as required: Turn twist cap 11 accordingly.

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>.

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open up instruments with hinges.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for, e.g., aluminum, plastics and high-grade steel, according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum permitted cleaning temperature of 60 °C.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfection.
- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Validated cleaning and disinfection procedure

For mechanical cleaning:

Note
Use appropriate (neutral, enzymatic and mild-alkaline) detergents on this product with its aluminum components.

Validated procedure	Specific requirements	Reference
Manual cleaning with ultrasound and immersion disinfection	- Cleaning brush: 30 mm/Ø: 4.5 mm, e.g. TA no. 011944 20 ml disposable syringe - Mount jaws protection on the product. - Keep working ends open for cleaning. - When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. - Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: Chapter Manual cleaning with ultrasound and immersion disinfection
Manual pre-cleaning with ultrasound and brush and subsequent mechanical neutral or mild alkaline cleaning and thermal disinfecting	- Cleaning brush: 30 mm/Ø: 4.5 mm, e.g. TA no. 011944 20 ml disposable syringe - Mount jaws protection on the product. - Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). - Keep working ends open for cleaning. - Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: - Chapter Manual pre-cleaning with ultrasound and brush - Chapter Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D–W	-
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-
D–W:	Drinking water					
FD–W:	Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)					
RT:	Room temperature					
*Recommended:	BBraun Stabimed					

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I
- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - Drain any remaining water fully.

- Phase III
- Fully immerse the product in the disinfectant solution.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

- Phase IV
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
 - Rinse lumens with an appropriate disposable syringe at least five times.
 - Drain any remaining water fully.

- Phase V
- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D–W	-
D–W:	Drinking water					
RT:	Room temperature					
*Recommended:	BBraun Stabimed					

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I
- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical neutral or mild alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	-
II	Cleaning	55/131	10	FD-W	Neutral: - Concentrate: - pH neutral <5 % anionic surfactant 0.5 %* working solution Mildly alkaline: - Concentrate: - pH = 9.5 <5 % anionic surfactant 0.5 % solution
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner neutral

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning/disinfecting process if necessary.

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note

The product may only be sterilized when dismantled.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the instrument
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

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