GB	Instructions for use
USA	Adtec mini instrument system
D	Gebrauchsanweisung
	AdTec mini-Instrumentensystem
F	Mode d'emploi
	Système de mini instruments Adtec
E	Instrucciones de manejo
	Sistema de mini instrumentos Adtec
I	Istruzioni per l'uso
	Sistema di mini strumenti Adtec
P	Instruções de utilização
	Sistema de mini instrumentos Adtec
NL	Gebruiksaanwijzing
	Adtec mini-instrumentensysteem
DK	Brugsanvisning
	AdTec mini-instrumentssystem
N	Bruksanvisning
	AdTec mini-instrumentssystem
S	Bruksanvisning
	Adtec mini instrumentssystem
FIN	Käyttöohjeet
	AdTec-mini-instrumenttijärjestelmä
GR	Οδηγίες χρήσης
	Adtec μίνι σύστημα εργαλείων
CZ	Návod k použití
	Systém nástrojů AdTec mini
PL	Instrukcja użytkowania
	System mini instrumentów Adtec
TR	Kullanım Kılavuzu
	Adtec mini aletler sistemi

B | BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 013248 12/10 V4

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Simbolo CE, em conformidade com a Directiva 93/42/CEE
CE-certificering conform richtlijn 93/42/EEG
CE-markering iht. retningslinje 93/42/EEC
CE-merking ifølge retningslinje 93/42/EWG
CE-märkning i enlighet med direktiv 93/42/EEG
93/42/EEC-standardin mukainen CE-hyväksyntä
Ενδειξη CE σύμφωνα με την Οδηγία 93/42/EEC
CE-označení podle směrnice 93/42/ES
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
93/42/EEG direktifi uyarınca CE işareti

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Technische wijzigingen voorbehouden
Retten til tekniske ændringer forbeholdes
Vi tar forbehold om eventuelle tekniske endringer
Med reservation för eventuella tekniska ändringar
Oikeusteknisistä syistä johtuviin muutoksiin pidätetään
Επιφύλασσόμεθα για τεχνικές αλλαγές
Technické změny vyhrazeny
Zmiany techniczne zastrzeżone
Teknik değişiklikler yapma hakkı saklıdır

A1


A2

B1

B2

Legend

- 1 Handle (complete)
- 2 Moving handle part
- 3 HF pin
- 4 Star wheel
- 5 Actuating lever
- 6 Switch (lubrication point)
- 7 Jaw insert (lubrication point)
- 8 Pressure bar
- 9 Outer tube

Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

- All professional disciplines that apply endoscopy:
- Cutting, preparation, and grasping of tissues
 - Biopsies
 - Sutures

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Exercise caution when passing products through the working channel (e.g. trocar).

WARNING

Danger of injury from burns and/or explosion from flammable gases!
Application of the HF device according to its intended use can involve spark formation.

- Observe the safety notices in the instructions for use of the HF device.

WARNING

Thermal injuries to patients/users due to insufficient insulation of leads to active accessories!

- Adjust the HF device to a setting appropriate to ensure that the maximum peak output voltage does not exceed the accessory voltage rating specified for the product.
- Use the instrument only with an insulated outer tube.

- Adjust the HF power output for the intended operation. Take into account clinical experience or reference values.
- Select the lowest possible HF power output.
- Keep the product's contact surfaces clean during the operation. Wipe off encrusted tissue residues or body fluids, using a moistened swab.

The plug end of the product is fitted with the following connector: Pin, 4 mm, sprung.
Refer to our brochures to find the compatible cable.
The accessory voltage rating of the product is 2 kVp.
The accessory voltage rating must exceed or match the maximum peak output voltage with which the product is operated in combination with an appropriate HF device at an appropriate operating mode/setting (see IEC/DIN EN 60601-2-2).
To avoid HF burns:

- Keep the product's working tip in the user's field of vision whenever HF power is activated.
- Prior to activating the HF device, make certain that the working end of the product is not touching any electrically conductive accessories.
- Prior to each use, visually check the product for: damage to, or surface changes on, the insulation.
- Never put down the product on or next to the patient.
- When using accessories for endoscopy or laparoscopy applications, deactivate the automatic switch-on mode of the HF device.
- Follow the instructions for use of the HF device.

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

WARNING

Risk of injury when applying the product outside the field of view!

- Apply the product only under visual control.

- Open and close the distal end of the jaws:
Open and close the moving handle part 2.

Handle with lock

- Press actuating lever 5 to release the ratchet mechanism.
 - To activate the ratchet mechanism, release actuating lever 5.
- Switch 6 can be used for permanently deactivating the ratchet mechanism.
- Deactivate the ratchet mechanism on handle 1:
Position switch 6 as shown in I.
The switch 6 engages in the front position.
 - Activate the ratchet mechanism on handle 1:
Position switch 6 as shown in II.
The switch 6 engages in the rear position.

Disassembling

Note

For disassembling the product, also see additional series of illustrations A.

- Deactivate the ratchet mechanism on handle 1 with lock:
Position switch 6 as shown in I.
- Detach the sheath from handle 1:
Push star wheel 4 back, hold it there and remove the sheath.
Extract jaw insert 7 with outer tube 9 from handle 1.
- Disassemble the sheath:
Remove outer tube 9 from jaw insert 7.
Rotate pressure bar 8 counterclockwise by 90° and extract it from jaw insert 7.

Assembling

Note

For assembling the product, also see the additional series of illustrations B.

- Deactivate the ratchet mechanism on handle 1 with lock:
Position switch 6 as shown in I.
- Assemble the sheath:
 - Slide pressure bar 8 into jaw insert 7 and rotate it clockwise by 90°.
 - Slide outer tube 9 onto jaw insert 7 down to the positive stop.
- Hold the sheath assembly at jaw insert 7 with one hand.
- With the other hand, hold handle 1 at star wheel 4. The moving handle part 2 must remain freely moveable.
- Push the sheath into handle 1, making certain that the notch on pressure bar 8 is aligned with the mark on star wheel 4.
The moving handle part 2 moves upward.
The sheath engages automatically at the positive stop.
- Test the instrument for proper functioning by opening and closing the jaw elements.

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.
Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures >45 °C nor any fixating disinfectants (active ingredient: aldehyde, alcohol) be used.
Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.
Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Material damage, e.g. corrosion, cracks, fracturing, premature ageing or swelling.
- Do not use process chemicals that would cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.
- Immersion treatment in a 3 % H₂O₂ solution for approx. 5 minutes is a particularly effective and gentle method to dissolve encrustations from HF instruments. Subsequently, the debris can be removed by hand, with a medium-hard brush and/or in an ultrasonic cleaning bath. After this pre-treatment, continue through the normal steps of the processing procedure.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open jointed products.
- Rinse surfaces that cannot be visually inspected, e.g. on products with hidden crevices or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Preparation prior to cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for e.g. aluminum, plastics and high-grade steel, according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 93 °C.

➤ Carry out ultrasound cleaning:

- as an effective mechanical supplement to manual cleaning/disinfecting,
- as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting,
- as an integrated mechanical support measure for mechanical cleaning/disinfecting,
- for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Manual cleaning/disinfecting

- Keep working ends open during cleaning.
- Clean products with movable links with the links opened or mobilized.
- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

1. Handle

Manual cleaning with immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized)
RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surfaces.
- Brush through surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using an appropriate cleaning brush (brush head length: 50 mm/Ø: 3.7 mm, e.g. GK469200).
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, thoroughly (at least five times) rinse through these components with the cleaning/disinfecting solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Make certain that all accessible surfaces are moistened.

Stage IV

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Make certain that all accessible surfaces are moistened.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting

Note
The disinfectant must be of proven effectiveness (e.g. DGHM or FDA approval or CE mark).

Note
Use appropriate (neutral, enzymatic and mild-alkaline) detergents for this product with aluminum components.

Note
For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that Ao is >3 000 for the process.

Note
The disinfectant used for processing must be serviced and checked at regular intervals.

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Keep working ends open during cleaning.
- Ensure that product links and joints are opened when placing the product on the tray.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfector without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER with tensides, application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfector program

D-W: Drinking water
FD-W: Fully desalinated water (demineralized)

2. Sheaths

Manual ultrasonic cleaning and immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	5	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	-	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized)
RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using an appropriate cleaning brush (brush head length: 50 mm/Ø: 3.7 mm, e.g. GK469200).
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, use a 20-ml single-use syringe to rinse thoroughly, for at least 5 times, these difficult to access parts of the product.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Make certain that all accessible surfaces are moistened.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.
- Rinse lumens at least 5 times, using a disposable syringe (20 ml) and an appropriate rinsing adapter.
- Allow water to drip off for a sufficient length of time.

Stage V

- Thoroughly dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

Note
The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note
For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that Ao is >3 000 for the process.

Note
The disinfectant used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Irrigation	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using an appropriate cleaning brush (brush head length: 50 mm/Ø: 3.7 mm, e.g. GK469200).
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

- Machine type: Single-chamber washer/disinfector without ultrasound
- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
 - Connect components with lumens and channels directly to the rinsing port of the injector carriage.
 - Keep working ends open during cleaning.
 - Ensure that product links and joints are opened when placing the product on the tray.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	BBRAUN HELIMATIC CLEANER alkaline with tensides; application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to disinfector program

D–W: Drinking water
FD–W: Fully desalinated water (demineralized)

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Assemble the separable product, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

Note
The product can be sterilized either in disassembled or in assembled condition.

Note
To avoid breakage due to stress crack corrosion, sterilize the instruments with the lock fully open or locked no further than on the first ratchet tooth.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
 - Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.
- To achieve a sterility assurance level of 10^{–6}, Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F–275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Technical Service

WARNING

Risk of injury and/or malfunction!

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-mail: ats@aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 216-3392
Fax: +1 314 895-4420
Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

For accessories and spare parts, see brochure C 28011 or visit <http://www.endoscopy-catalog.com>

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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