
- (GB)** Instructions for use/Technical description
AdTec® needle holder
- (USA)** AdTec® needle holder
- (D)** Gebrauchsanweisung/Technische Beschreibung
AdTec® needle holder
- (F)** Mode d'emploi/Description technique
Porte-aiguilles AdTec®
- (E)** Instrucciones de manejo/Descripción técnica
AdTec® needle holder
- (I)** Istruzioni per l'uso/Descrizione tecnica
AdTec® needle holder
- (P)** Instruções de utilização/Descrição técnica
AdTec® needle holder (porta agulhas)
- (NL)** Gebruiksaanwijzing/Technische beschrijving
AdTec® needle holder
- (S)** Bruksanvisning/Teknisk beskrivning
AdTec® nålförare
- (RUS)** Инструкция по применению/Техническое описание
Иглодержатель AdTec®
- (CZ)** Návod k použití/Technický popis
Držák jehly AdTec®
- (PL)** Instrukcja użytkowania/Opis techniczny
AdTec® needle holder
- (SK)** Návod na použitie/Technický opis
AdTec® needle holder
- (TR)** Kullanım Kılavuzu/Teknik açıklama
AdTec® needle holder

B | BRAUN

SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 011172 2017-07 V6 Änd.-Nr. 57290

CE 0123 – DIR 93/42/EEC

Aesculap® AdTec® needle holder

Legend

- 1 Jaw part
- 2 Lubrication points
- 3 Luer lock connector for cleaning channel
- 4 Screw handle

Symbols on product and packages

	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture

Intended use

The AdTec® needle holders are used for grasping and suturing in laparoscopic operations.
The applier forceps PL406R is used for fixating the counter-pressure plate of a circular stapler.

Available sizes

The AdTec® needle holders available to order are described in brochure C004 and at www.endoscopy-catalog.com

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.
- ▶ To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).

Safe operation

	Risk of injury and/or malfunction! ▶ Always carry out a function check prior to using the product.
	Risk of injury when using the product beyond the field of view! ▶ Apply the product only under visual control.

Closing the jaw piece

- ▶ To close jaw piece 1 and grasp objects without engaging the lock on handle 4, slightly compress the handle.

Engaging the lock

- ▶ To lock handle 4 and jaw piece 1, fully compress the handle until the lock engages automatically with a click.
The instrument is locked now.

Releasing the lock

- ▶ Fully compress handle 4 until the lock is released automatically with a click.

Note

After releasing the lock, the handle must be opened half-way at least before the lock can be engaged again.

Validated reprocessing procedure

General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- ▶ Open up products with hinges.
- ▶ Remove the sealing cap from the Luer lock connector.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

	Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures! ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which – are approved for (e.g. aluminum, plastics, high-grade steel), – do not attack softeners (e.g. in silicone). ▶ Observe specifications regarding concentration, temperature and exposure time. ▶ Do not exceed the maximum allowable temperature of 96 °C.
--	--

Validated cleaning and disinfection procedure

For mechanical cleaning:

Note

Use appropriate (neutral, enzymatic and mild-alkaline) detergents on this product with its aluminum components.

Validated procedure	Specific requirements	Reference
Manual cleaning with ultrasound and immersion disinfection	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Keep working ends open for cleaning. ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection	■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfection

Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D–W	-
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the product in the disinfectant solution.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- ▶ Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- ▶ Rinse lumens with an appropriate disposable syringe at least five times.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D–W	-

D–W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) at the marked lubrication points, using maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

CAUTION

Damage (metal cold welding /friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate both locking parts at the labeled points with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ Lubricate the product at the labeled points 2.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractional vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Validated reprocessing procedure for the US market

Manual cleaning and disinfection

Phase I

- ▶ Prepare washing solution following instructions on solution label using unheated utility water.
- ▶ Open or disassemble instrument/s so solution can reach all instrument surfaces and submerge instrument/s in solution.
- ▶ Brush all surfaces until visibly clean with a soft or medium bristle brush brushing hard to clean areas such as lumens, special rotating areas and hidden surfaces 4 times (back and forth = 1 time).
- ▶ Move moving parts during washing to reach all surfaces and flush lumens 4 times with at least 20cc washing solution using a syringe.

Phase II

- ▶ Rinse by placing instruments under critical water and expose all surfaces of the device under the water flow for 3 times.
- ▶ Move parts during the rinsing to rinse all surfaces and flush lumens 4 times with at least 20cc washing solution using a syringe.

Phase III

- ▶ Use heat or low-linting cloth to dry instruments following final rinse (not required if followed by mechanical cleaning process).

Phase IV

- ▶ Visually inspect instruments for remaining soil or washing solution.
- ▶ Repeat manual washing and/or rinsing steps until instrument is free from soil or washing solution.

Mechanical Cleaning

- Phase I
- Use unheated utility water rinse for at least 1 minute (may be part of automated washer program).
- Phase II
- Using unheated utility water, follow washing instructions for recommended dilution in washer.
 - Open or disassemble instrument/s so washing solution can reach all instrument surfaces and wash tube structure (lumens) using flush attachments in washer if available or use manual method.
 - Use wash cycle for at least 6 minutes.
- Phase III
- Use minimum 180 °F critical water rinse for at least 2 minutes.
- Phase IV
- Dry at 194 °F for 10 minutes
 - or –
 - Use the manual dry method.
- Phase V
- Visually inspect instruments for remaining soil or washing solution.
 - Repeat cleaning steps until instrument is free from soil or washing solution.

Steam sterilization

Steam sterilization parameters			
Sterilization		Validation testing exposure time	Validation testing dry time
Sterilization method	Temperature	Packaged	Length
Terminal sterilization Dynamic air removal (Pre-Vac)	270 °F (132 °C)	4 minutes	15 minutes
Terminal Sterilization Dynamic air removal (Pre-Vac) without dry time			0 minutes

Technical Service

WARNING

Risk of injury and/or malfunction!

- Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Art. no.	Designation
EJ751251	Sealing cap for Luer Lock connector (pack of 20)

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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