
Aesculap® Acculan 4 / ELAN 4 electro

en Instructions for use/Technical description

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USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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











Aesculap® Acculan 4 / ELAN 4 electro

Attachments for small drill GA344 and drill GA844

Legend

- 1 Drill
- 2 Knob (for speed control)
- 3 Safety catch
- 4 Rotating sleeve
- 5 Arrow
- 6 Attachment with collet chuck
- 7 Unlocking sleeve
- 8 Tool coupling
- 9 Attachment with three jaw chuck
- 10 Tightening key
- 11 Tool coupling
- 12 Kirschner wire attachment
- 13 Clamping lever
- 14 Adjustment sleeve
- 15 Kirschner wire chuck
- 16 Drill attachment for radiolucent angle transmission
- 17 Adapter
- 18 Sagittal saw attachment
- 19 Tool attachment with L sagittal tool coupling
- 20 Button for adjusting the tool holder
- 21 L sagittal saw blade
- 22 Piece on saw blade for tool release
- 23 Oil spray adapter
- 24 Reciprocating Saw attachment
- 25 Tool attachment with Reciprocating Saw tool coupling
- 26 Reciprocating Saw blade (e.g. GP550R to GP576R)
- 27 Leaf spring on Reciprocating Saw blade for tool release
- 28 Screw attachment
- 29 Lever

Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
	Machine-readable, two-dimensional code The code contains a unique serial number which can be used for electronic tracking of the individual instrument. The serial number is based on the global standard sGTIN (GS1).
	Manufacturer
	Date of manufacture
	Manufacturer's batch designation
	Manufacturer's serial number
	Manufacturer's article number
	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage

Contents

1.	Applicable to	3
2.	General information	3
2.1	Intended use	3
2.2	Main functions and design characteristics	3
2.2.1	Operating mode	3
2.3	Indications	4
2.4	Absolute contraindications	4
2.5	Relative contraindications	4
3.	Safe handling	4
4.	Product description	5
4.1	Scope of supply	5
4.2	Components required for operation	5
4.3	Operating principle	5
5.	Preparation	5
6.	Working with the device	5
6.1	System set-up	5
6.1.1	Connecting the accessories	5
6.1.2	Protection against inadvertent activation	6
6.1.3	Attaching/removing attachment for drill	6
6.1.4	Attaching and detaching the tool from the attachment	6
6.2	Function checks	7
6.3	Safe operation	8
6.3.1	Operating the product	8
6.3.2	Operating the screw attachment	8
7.	Validated reprocessing procedure	9
7.1	General safety notes	9
7.2	General information	9
7.3	Preparations at the place of use	9
7.4	Preparation before cleaning	9
7.5	Cleaning/disinfection	9
7.5.1	Product-specific safety instructions for the reprocessing procedure	9
7.6	Manual cleaning with wipe disinfection	10
7.7	Mechanical cleaning/disinfection with manual pre-cleaning	11
7.7.1	Manual pre-cleaning with a brush	11
7.7.2	Mechanical alkaline cleaning and thermal disinfection	12
7.8	Inspection, maintenance and checks	12
7.9	Packaging	12
7.10	Steam sterilization	12
7.11	Storage	12
8.	Maintenance	12
9.	Troubleshooting list	13
10.	Technical Service	14
11.	Accessories/Spare parts	14
12.	Technical data	15
12.1	Classification acc. to Directive 93/42/EEC	15
12.2	Performance data, information about standards	15
12.3	Operating mode	17
12.4	Environmental conditions	17
13.	Disposal	17

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

2. General information

2.1 Intended use

Task/Function

The hand-held small drill GA344 and the drill GA844, combined with the appropriate attachment, is used for driving drilling, milling (medullary reaming), sawing and screwing tools.

Application Environment

The product is used in operating rooms in sterile zones outside of the explosion risk zone (such as areas with pure oxygen or anesthesia gases).

2.2 Main functions and design characteristics

Speed	GB870R	0 min ⁻¹ to max. 320 min ⁻¹
	GB871R	
	GB872R	
	GB878R	0 min ⁻¹ to max. 1 250 min ⁻¹
	GB881R	
	GB882R	
	GB884R	
	GB886R	
	GB887R	
	GB891R	0 min ⁻¹ to max. 15 750 min ⁻¹
	GB892R	
	GB896R	0 min ⁻¹ to max. 165 min ⁻¹
	GB897R	
Rotational direction		Right and left rotation, oscillation

2.2.1 Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

Drilling (clockwise/counterclockwise rotation)

- 60 second application, 60 second pause
- 6 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Medullary reaming (clockwise/counterclockwise rotation)

- 30 second application, 30 second pause
- 8 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Drilling (oscillation)

- 15 second application, 15 second pause
- 3 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Saw mode with GB891R

- 30 second application, 60 second pause
- 4 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Saw mode with GB892R

- 30 second application, 60 second pause
- 5 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Screwing with GB896R/GB897R

- 10 second mechanical application, 10 second manual application, 30 second pause
- 30 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Electrical systems generally heat up during continual operation. It is advised to give the system a break after use to cool down, as listed in the table on operating mode.

Heating depends on the tool used and the load. After a certain number of repetitions, the system should cool down. This procedure prevents the system overheating as well as possible injury to the patient or user.

The user is responsible for the use and adherence to the pause sequence described.

2.3 Indications

The small drill GA344 and the drill GA844, combined with the relevant attachment and tool, is used for separating, removing and shaping, hard tissue, cartilage and the like, as well as bone replacement material, to screw bone pins in and out, to drive screw tools and to set transfixion wires.

2.4 Absolute contraindications

The product is not licensed for use on the central nervous system or central circulatory system.

2.5 Relative contraindications

The safe and effective use of electrically driven products greatly depends on influences which can only be controlled by the user. Therefore the specifications provided represent framework conditions only.

Clinically successful use of the product is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

3. Safe handling**⚠ WARNING**

Risk of injury and material damage if this product is not used as intended!

- ▶ Use the product only for its intended purpose.

⚠ WARNING

Risk of injury and damage to property due to improper handling of the product!

- ▶ Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Prior to use, check that the product is in good working order.
- ▶ To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Keep the instructions for use accessible for the user.
- ▶ Always adhere to applicable standards.
- ▶ Ensure that the electrical installation of the room is consistent with the requirements of IEC/DIN EN.
- ▶ Do not operate the product in explosion-hazard areas.
- ▶ Sterilize product before use.
- ▶ When using the ECCOS holder system, adhere to the appropriate instructions for use for TA009721, see Aesculap Extranet at <https://extranet.bbraun.com>

4. Product description

4.1 Scope of supply

Art. no.	Designation
TA014553	Instructions for use, attachments for small drill GA344 and drill GA844

Attachments

One of the following attachments (with corresponding chuck key as applicable) is additionally included with the delivery:

Art. no.	Designation
GB870R	Medullary reaming attachment, large Jacobs chuck – and –
GA031R	Key for large three-jaw chuck
GB871R	Medullary reaming attachment AO large
GB872R	Medullary reaming attachment Hudson/Zimmer
GB878R	Drill attachment for Synthes 511.300
GB881R	Drill attachment large Jacobs chuck – and –
GA031R	Key for large Jacobs chuck
GB882R	Drill attachment small Jacobs chuck – and –
GA062R	Chuck key for small Jacobs chuck
GB884R	Drill attachment AO small
GB886R	Drill attachment Aesculap hexagon
GB887R	Drill attachment dental shaft
GB891R	Sagittal saw attachment
GB892R	Reciprocating Saw attachment
GB894R	Kirschner wire attachment
GB896R	Screw attachment AO small
GB897R	Screw attachment Aesculap hex

4.2 Components required for operation

- Small drill GA344 (ready to use) or drill GA844 (ready to use)
- Tool (depending on indication)

4.3 Operating principle

The attachment can be attached to the drill in three different positions, each 120° difference.

The attachment locks automatically when attached to the drill. The attachment can be removed by operating a rotating sleeve.

A gear in the attachment changes the motor speed.

The attachments have various integrated couplings on the working end, to enable attachment of various tools, with adapter if necessary.

5. Preparation

Aesculap assumes no liability if the following rules are not followed:

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Prior to use, inspect the product and its accessories for any visible damage.
- ▶ Use the products and their accessories only if they are in perfect technical condition.

6. Working with the device

⚠ WARNING

Risk of infection and contamination!

Product is delivered unsterilized!

- ▶ Sterilize the product before use pursuant to the operating instructions.

⚠ WARNING

Risk of injury and material damage due to accidental activation of the product!

- ▶ Products which are not being actively used must be secured against accidental activation (position OFF).

⚠ WARNING

Risk of injury and material damage due to inappropriate use of tools!

- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ When coupling/uncoupling, handle tools with cutting edges with care.

⚠ WARNING

Damage to the product if dropped!

- ▶ Use the products only if they are in perfect technical condition, see Function check

⚠ WARNING

Risk of burns to skin and tissue through blunt tools or if product has not been maintained properly!

- ▶ Use tools only if they are in perfect condition.
- ▶ Replace blunt tools.
- ▶ Maintain the product properly, see maintenance guide.

6.1 System set-up

6.1.1 Connecting the accessories

⚠ DANGER

Risk of injury due to unapproved configuration using additional components!

- ▶ Ensure that the classification of all components used is consistent with the classification of the product (such as type BF or type CF).

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

All configurations must fulfill the fundamental standard IEC/EN/DIN 60601-1. The person connecting the devices with each other is responsible for the configuration and must ensure that the fundamental standard IEC/EN/DIN 60601-1 or relevant national standards.

- Follow the instructions for use of individual accessories.
- Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

6.1.2 Protection against inadvertent activation

To prevent the drill being inadvertently activated when changing a tool or accessory, the knob for speed control can be locked.

Locking the knob for speed control 2:

- Twist the safety catch 3 to position OFF.
The speed control knob 2 is blocked and the product 1 cannot be operated.

Unlocking the knob for speed control 2:

- Twist the safety catch 3 to position ON.
The speed control knob 2 is unlocked and the product 1 can be operated.

Note

For more information on the small drill GA344, see TA014550 and TA014551 (flyer).

For more information on the drill GA844, see TA014436 and TA014437 (flyer).

6.1.3 Attaching/removing attachment for drill

⚠ WARNING

Risk of injury when attaching/removing attachments/tools in the ON position through inadvertent activation of the product!

- **Only attach/remove attachments/tools in the OFF position.**
- Secure drill 1 from inadvertent activation with the safety catch 3, see Protection against inadvertent activation.

Coupling

Note

Aesculap recommends coupling the Reciprocating Saw attachment 24 with the label "TOP" upward, aligned to the drill.

- Push attachment 6/9/12/16/18/24/28 onto the drill 1 until it clicks into place.
- Pull on attachment 6/9/12/16/18/24/28 to check secure connection.

Uncoupling

- Turn sleeve 4 in the direction of the arrow 5 and simultaneously remove the attachment 6/9/12/16/18/24/28 from the drill 1.

6.1.4 Attaching and detaching the tool from the attachment

⚠ WARNING

Risk of injury when attaching/removing attachments/tools in the ON position through inadvertent activation of the product!

- **Only attach/remove attachments/tools in the OFF position.**
- Ensure that the tool connector and attachment type are compatible.

Attachments with quick-mount chuck

GB871R/GB872R/GB884R/GB886R/GB887R/GB896R/GB897R

Coupling:

- Pull back on the unlocking sleeve 7.
- Insert the tool shaft in the correct position up to the stop in the tool holder 8 of the attachment 6.
- Release the unlocking sleeve 7.
Tool is attached.
- Pull at tool to check that it is securely coupled.

Note

Information on operating the screw attachment 28, see Operating the screw attachment.

Uncoupling:

- Pull back the unlocking sleeve 7.
- Remove tool.

Attachments with three-jaw chuck GB870R/GB881R/GB882R

Coupling:

- Open the three jaw chuck with the wrench 10.
- Insert the tool shaft in the correct position up to the stop in the tool holder 11 of the attachment 9.
- Close the three jaw chuck 10 and pull tight.
- Pull at tool to check that it is securely coupled.

Uncoupling:

- Open the three jaw chuck with the wrench 10.
- Remove tool.

Kirschner wire attachment GB894R

⚠ WARNING

Risk of injury when using long guide wires!

- **Use Kirschner wire protection sleeve when using long guide wires.**

Note

The special Kirschner wire chuck is recommended for the placement of guide wires. This quick-action chuck adapter allows quick and easy tensioning of drilling wires.

When using long guide wires, a Kirschner wire protection sleeve must be used with the drill to prevent injury.

For further information on GA344, see TA014550 and TA014551 (leaflet).

For further information on GA844, see TA014436 and TA014437 (leaflet).

The following diameters can be set on the Kirschner wire attachment:

- 0.6 mm to 1.2 mm
- 1.2 mm to 2.2 mm
- 2.2 mm to 3.2 mm

Using a Kirschner wire:

- ▶ Ensure that the clamping lever **13** is in the starting position (deactivated status).
- ▶ Set the adjusting sleeve **14** of the Kirschner wire attachment **12** to the desired diameter range:
 - Push back the adjusting sleeve **14** and turn until the desired diameter range has been reached.
 - Release the adjusting sleeve **14**. Ensure that the adjustment sleeve **14** clicks into place.
- ▶ Insert the Kirschner wire into the Kirschner wire chuck **15** until the desired exposed length is achieved.
The Kirschner wire is held in the intended position in the attachment by an automatic clamping mechanism.

Tension the Kirschner wire:

- ▶ Pull on the clamping lever **13** and hold in the desired position.
The further the clamping lever is pulled back, the higher the tension of the Kirschner wire.

Note

The Kirschner wire will only remain tensioned with the clamping lever pulled back. If the clamping lever is released, it will move back into the starting position and the Kirschner wire is freely mobile.

Drill attachment for radiolucent angle transmission GB878R

The drill attachment allows connecting a radiolucent angle transmission for applying drill holes in the bone under radiographic control.

The drill attachment is only suitable for use of the radiolucent angle transmission 511.300 from the Synthes company.

- ▶ Follow the instructions for the radiolucent angle transmission 511.300 in the operating instructions of the Synthes company.

Attach radiolucent angle transmission on drill attachment GB878R:

- ▶ Push the radiolucent angle transmission up to the stop on attachment **17** of the drill attachment.
- ▶ If necessary, slightly move the radiolucent unit.

Detach radiolucent angle transmission on drill attachment GB878R:

- ▶ Forcefully pull out the radiolucent angle transmission from the drill attachment.

Attach/detach tool in radiolucent angle transmission:

- ▶ Follow the instructions for the radiolucent angle transmission 511.300 in the operating instructions of the Synthes company.

Sagittal saw attachment GB891R

Coupling:

- ▶ Insert saw blade **21** with marking "L" up into the slot of the tool holder **19**, see Fig. A. Ensure that the side stops of the saw blade contact the tool holder.
The tool engages.
- ▶ Pull on the attachment **21** to ensure a secure fit.

Uncoupling:

Note

For easier detachment of the tool, set the tool holder to 45° (see Fig. B) or 90° (see Fig. C), see tool holder adjustment.

- ▶ Move the gold bar **22** on the saw blade **21** slightly downward and hold it to unlock the tool.
- ▶ Remove the saw blade **21** from the tool holder **19**.

Adjusting the tool holder:

- ▶ Press the button for tool holder adjustment **20** and keep it depressed.
- ▶ Turn the tool holder **19** to the desired position (–90°/–45°/0°/45°/90°), see Fig. B and Fig. C.
- ▶ Release the adjustment button of the tool holder **20** and secure the tool holder **19** by clicking it into place. If necessary, slightly move the tool holder back and forth.

Reciprocating Saw attachment GB892R

Coupling:

- ▶ Insert connection side of Reciprocating Saw blade **26** into the slot in the tool holder **25**, see Fig. D.
The Reciprocating Saw blade clicks into place.
- ▶ Pull on the Reciprocating Saw blade **26** to verify a secure fit.

Uncoupling:

- ▶ Gently press the gold piece for tool release on the Reciprocating Saw blade **27** downwards and keep pressed down.
- ▶ Remove the Reciprocating Saw blade **26** from the tool holder **25**.

6.2 Function checks

The function test must be performed before each use and after each intra-operative change of attachment and tool.

- ▶ Check that the attachment is coupled properly. To test this, pull on the attachment.
- ▶ Check that the tool is coupled properly. To do this, pull on the tool.
- ▶ For GB891R: check that the tool attachment has locked into place. To do this, turn the tool attachment.
- ▶ Ensure that the tool blades are not mechanically damaged.
- ▶ Release drill for use (ON position).
- ▶ Run the drill shortly at maximum speed in right and left rotation.
- ▶ Make certain that the rotational direction is correct in each case.
- ▶ For GB896R/GB897R: ensure that there is both mechanical and manual functionality.
- ▶ Check the product for damage, abnormal operation noises, heavy vibration and excessive heat.
- ▶ Do not use the product if it is damaged or defective.
- ▶ Set aside the product if it is damaged.

6.3 Safe operation

WARNING

Coagulation of patient tissue or risk of burns for patients and user through hot product!

- ▶ Do not use drill attachments to operate milling tools.
- ▶ Do not use the medullary milling attachments for acetabulum milling.
- ▶ Cool the tool during operation.
- ▶ Keep product/tool out of the reach of patients.
- ▶ Let the product/tool cool down.
- ▶ Use a cloth to protect against burns when changing the tool.

WARNING

Risk of infection from aerosol formation!

Risk of injuries caused by particles coming loose from the tool!

- ▶ Use suitable protection (such as waterproof protective clothing, face mask, safety gases, suction).

WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!

- ▶ Apply the product only under visual control.

WARNING

Risk of injury and damage to the tool/system!

The rotating tool may get caught in drapes (such as textiles).

- ▶ Do not let the tool come into contact with drapes (such as textiles) during operation.

WARNING

Risk of injury through inadvertent release of the saw blade!

- ▶ Do not press the sagittal saw blade tool release during operation.
- ▶ Do not press the Reciprocating Saw blade tool release leaf spring during operation.
- ▶ Check for secure fixation of the tool after any tool change.

6.3.1 Operating the product

Note

For more information on operation of the small drill GA344, see TA014550 and TA014551 (flyer).

For more information on operation of the drill GA844, see TA014436 and TA014437 (flyer).

- ▶ Start the drill at a moderate speed.
- ▶ Apply moderate pressure to prevent slipping.
- ▶ Do not bend the tool, otherwise there is a danger of breakage.
- ▶ When screwing bone screws and bone pins in and out: ensure correct rotation direction setting.

- ▶ Use the following attachments for drilling:
 - GB881R
 - GB882R
 - GB884R
 - GB886R
 - GB887R
- ▶ Use the following attachments for milling (medullary reaming):
 - GB870R
 - GB871R
 - GB872R
- ▶ Use the following attachment for the placement of Kirschner wires:
 - GB894R
- ▶ Use only the following attachment for using the radiolucent angle transmission 511.300 from the Synthes company:
 - GB878R
- ▶ Use the following attachments for sawing:
 - GB891R
 - GB892R
- ▶ Use the following attachments for driving screw tools:
 - GB896R
 - GB897R

6.3.2 Operating the screw attachment

WARNING

Risk of injury!

- ▶ Do not completely screw bone screws in mechanically.
- ▶ Perform the final rotations and locking of the bone screws manually.

WARNING

Risk of injury and damage to property due to improper handling of the product!

- ▶ Follow the instructions for use of the implant system that is used.

Mechanical operation

- ▶ Activate the knob for speed control 2. Do not activate lever 29 when doing so.

Manual operation with ratcheting function

- ▶ Pull lever 29 to the stop and hold.
- ▶ Move product clockwise.
The rotational movement is transferred to the coupled screw tool.
- ▶ Move product counterclockwise.
The ratcheting function is active. The product can be moved into the starting position without also moving the screw tool.

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bb Braun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

On stainless steel, residues containing chlorine or chloride (such as surgical residues, drugs, saline solutions in water for cleaning, disinfection and sterilization) may lead to corrosion (pitting corrosion, tensile corrosion) and thus to the destruction of the product. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Visual material changes (such as fading or color changes in titanium or aluminum). For aluminum, the application/process solution only needs to be pH >8 to cause visible surface changes.
- Material damage (such as corrosion, cracks, breaks, premature aging or swelling).
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

7.3 Preparations at the place of use

- ▶ Remove all attached components from the product (tool and accessories).
- ▶ Detach the attachment from the drill.
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.4 Preparation before cleaning

- ▶ Before the first automatic cleaning/disinfection: attach ECCOS holders in a suitable tray.
- ▶ Insert the products in the correct position into the ECCOS fixations, see Fig. E.

Kirschner wire attachment GB894R

- ▶ Set adjustment sleeve to the largest Kirschner wire diameter.

7.5 Cleaning/disinfection

7.5.1 Product-specific safety instructions for the reprocessing procedure

⚠ CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastic material and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C during chemical cleaning and/or disinfection.
- ▶ Do not exceed the maximum temperature of 96 °C during thermal disinfection with FD water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

Note

The indicated drying time is a guide time only. It must be checked taking into account the specific conditions (e.g. load) and if applicable adjusted.

7.6 Manual cleaning with wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	≥ 2	-	D-W	Until visually clean
II	Cleaning with enzyme solution	RT (cold)	≥ 2	0.8	D-W	pH-neutral*
III	Intermediate rinse	RT	≥ 5	-	D-W	-
IV	Drying	RT	-	-	-	-
V	Wipe disinfection	-	> 1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol
VI	Final rinse	RT (cold)	0.5	-	FD-W	-
VII	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

* Validated with enzyme cleaner "Cidezyme Johnson & Johnson"

- Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

Phase I

- Flexible components (such as sleeves) will move during cleaning.
- Clean the product under running water, using a suitable plastic cleaning brush until all visible residues have been removed from the surfaces.
- Brush cannulation with cleaning brush TA011944 and difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.

Note

For details on difficult to access surfaces, see *Acculan pre-cleaning and care information TA016000* (available in Aesculap Extranet at <https://extranet.bbraun.com>).

Phase II

- Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- Spray products with a pH neutral enzyme solution, let soak in for at least 2 minutes and then wipe off.

Phase III

- Flexible components (such as sleeves) will move during cleaning.
- Rinse product under running tap water for at least 5 minutes.
- Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- Contamination should be removed with a lint-free cloth or soft brush moistened with enzyme cleaner.
- Flexible components (such as sleeves) and cannulations should be rinsed for 20 seconds with the water pistol (cold water, at least 2,5 bar).
- After manual cleaning, check visible surfaces and areas of flexible components for residues.
- If necessary, repeat the cleaning process (phase 1 to 3).

Phase IV

- Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

Phase V

- Wipe all surfaces of the product with a single-use disinfecting wipe.

Phase VI

- Rinse disinfected surfaces after the prescribed reaction time for at least 1 minute under running demineralized water.
- Drain any remaining water fully.

Phase VII

- Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

7.7 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must fundamentally have a tested efficacy (such as FDA approval or CE label pursuant to DIN EN ISO 15883).

Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

7.7.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry/Note
I	Rinsing	RT (cold)	-	-	D-W	Until visually clean
II	Brushes	RT (cold)	-	-	D-W	Until visually clean

D-W: Drinking water

RT: Room temperature

- Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.
- For Kirschner wire attachment GB894R: Set the adjustment sleeve to the greatest Kirschner wire diameter.

Phase I

- Flexible components (such as sleeves) will move during cleaning.
- Thoroughly clean the product under running water.

Phase II

- Flexible components (such as sleeves) will move during cleaning.
- Brush cannulation with cleaning brush TA011944 and difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.
- After manual preparation, check visible surfaces for residue and repeat the pre-cleaning process as needed.

Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbraun.com>).

7.7.2 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH ~ 13 - <5 % anionic surfactant ■ 0,5 % working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	mind. 10 min at max. 120 °C

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: B Braun Helimatic Cleaner alkaline

- For Kirschner wire attachment GB894R: Set the adjustment sleeve to the greatest Kirschner wire diameter.
- Insert the product in the correct position into the ECCOS fixation, see Fig. E.
- Then connect the interior rinse device to the ECCOS holder and connect to the rinse connector of the cleaning/disinfection device/rinse cart.
- After automatic cleaning/disinfection, check visible surfaces for residues and repeat the cleaning/disinfection process as needed.

7.8 Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- Spray product after every cleaning and disinfection with oil spray adapter 23 GB600880 (green) approx. 2 sec with Aesculap STERILIT Power Systems oil spray GB600, see Fig. F.

Note

Aesculap additionally recommends occasional spraying of moving parts (such as knobs, couplings, cover lids) with the Aesculap STERILIT-Power-Systems oil spray.

- Check the product after each cleaning and disinfection for the following: cleanliness, damage, function, abnormal operation noise, excessive heat or heavy vibration.
- Set aside the product if it is damaged.

7.9 Packaging

- Always follow the instructions for use of the respective packaging and storage devices (e.g. instructions for use TA009721 for Aesculap-ECCOS holder system).
- Insert the products in the correct position into the ECCOS fixations, see Fig. E.
- Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging will prevent a recontamination of the product.

7.10 Steam sterilization

Note

Remove all attached components from the product (tools, accessories) before sterilization.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Use a validated sterilization method:
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer DIN EN 285 and validated pursuant to DIN EN ISO 17665
 - Sterilization using the fractionated vacuum process at 134 °C, holding time 5 min

When sterilizing multiple products in one steam sterilizer:

- Ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

7.11 Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

8. Maintenance

To ensure reliable operation, the product must be maintained at least once a year.

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

9. Troubleshooting list

- Have defective products repaired by Aesculap Technical Service, see Technical Service.

Issue	Cause	Detection	Remedy
Product becomes too hot	Excessive use	Attachment very hot	Follow operating instructions (operating mode).
	Defective attachment drive/ball bearings	Attachment very hot	Follow operating instructions (preparation, care). Preventative: lubricate attachment before every sterilization. Have the manufacturer repair the attachment.
	Fall damage, product defective	Attachment very hot	Have the attachment repaired by the manufacturer.
	Blunt tool	Heating of tool and product	Change tool.
Insufficient power	Product is operating in left rotation mode	Serrated tool run in counterclockwise mode	Use serrated tools only in clockwise mode.
	Attachment defective	Attachment very hot	Follow operating instructions (preparation, care). Preventative: lubricate attachment before every sterilization. Follow operating instructions (operating mode). Have the attachment repaired by the manufacturer.
	Blunt tool	Tool blades worn out	Change tool.
Loud running noise	Defective attachment drive/ball bearings	Loud, perceptible noise during operation	Follow operating instructions (preparation, care). Preventative: lubricate attachment before every sterilization. Have the manufacturer repair the attachment.
Attachment does not attach or detach.	Attachment not compatible	Attachment doesn't engage	Use the compatible attachment for drill.
	Connection on attachment deformed/defective	Attachment is difficult to couple or cannot be coupled	Have the manufacturer repair the attachment.
	Coupling on drill deformed/defective	Attachment is difficult to couple or cannot be coupled	Have the manufacturer repair the drill.
	Sleeve on drill hard to move	Attachment doesn't engage	Turn and hold the drill sleeve, then attach the attachment. Turn the sleeve back if necessary. Preventative: lubricate the sleeve on the drill before each sterilization.
Tool does not attach or detach.	Incompatible tool	Tool does not click into place	Use the tool compatible with the attachment.
	Tool connector deformed/defective	Tool is difficult to couple or cannot be coupled	Use a new tool.
	Coupling on attachment deformed/defective	Tool is difficult to couple or cannot be coupled	Have the manufacturer repair the attachment.
	Tool coupling dirty	Tool is difficult to couple or cannot be coupled	Clean tool or use a new tool. Clean the attachment.
For GB894R: Kirschner wire does not insert into the Kirschner wire chuck	Kirschner wire chuck not set correctly	Diameter of the Kirschner wire is not consistent with the scale setting	Set the adjustment sleeve to the right Kirschner wire diameter.
For GB894R: Kirschner wire not turning	Kirschner wire chuck not set correctly	Diameter of the Kirschner wire is not consistent with the scale setting	Set the adjustment sleeve to the right Kirschner wire diameter.
	Open the clamping lever	Open the clamping lever	Close the clamping lever.
	Clamping lever not pulled back with sufficient force	Kirschner wire not turning	Pull back the clamping lever forcefully and hold it.

Issue	Cause	Detection	Remedy
Tool not moving	Attachment not fully attached to the drill	Attachment can be pulled out of the drill	Correctly attach the attachment and perform a function test.
	Tool not fully attached to the attachment	Tool can be pulled out of the tool coupling	Correctly attach the tool and perform a function test.
	Attachment defective	Drill turning, but attachment is not turning.	Have the manufacturer repair the attachment.
	Drill defective	Drill not turning	Have the manufacturer repair the drill.
	The safety catch of the drill is in the OFF position	Safety catch is in OFF position	Turn the safety catch into the ON position.
	For GB896R/GB897R: lever pulled	For GB896R/GB897R: lever pulled	For GB896R/GB897R: do not pull lever for mechanical operation.
For GB896R/GB897R: changing between mechanical and manual operation is not possible.	Lever mechanism defective	Lever is difficult to move or cannot move. Lever is bent. Lever does not move itself back to the starting position.	Clean and oil attachment. Have the manufacturer repair the attachment. Alternatively, use a screwdriver handle.
	Ratcheting mechanism defective	Lever is difficult to move or cannot move. No ratcheting function in manual operation.	Clean and oil attachment. Have the manufacturer repair the attachment. Alternatively, use a screwdriver handle.
Oil spray adapter not attachable	Oil spray adapter not compatible	Oil spray adapter not attachable	Use the compatible oil spray adapter for GA344.

10. Technical Service

WARNING

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.
- Do not modify the product.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

11. Accessories/Spare parts

Art. no.	Designation
GB257R	ECCOS sieve basket with fixation for GA344
GB262R	ECCOS sieve basket with fixation for GA844
GB600	STERILIT Power Systems oil spray
GB600880	Oil spray adapter for GA344/GA844
TA011944	Cleaning brush
GB715R	ECCOS holder single
GB716R	ECCOS holder triple
GA031R	Key for large Jacobs chuck
GA062R	Chuck key for small three-jaw chuck
TA014552	Instructions for use, attachments for small drill GA344 and drill GA844 (for ring binder)
TA014553	Instructions for use, attachments for small drill GA344 and drill GA844 (flyer)

12. Technical data

12.1 Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GB870R	Medullary reaming attachment, large Jacobs chuck	IIa
GB871R	Medullary reaming attachment AO large	
GB872R	Medullary reaming attachment Hudson/Zimmer	
GB878R	Drill attachment for Synthes 511.300	
GB881R	Drill attachment large Jacobs chuck	
GB882R	Drill attachment small Jacobs chuck	
GB884R	Drill attachment AO small	
GB886R	Drill attachment Aesculap hexagon	
GB887R	Drill attachment dental shaft	
GB891R	Sagittal saw attachment	
GB892R	Reciprocating Saw attachment	
GB894R	Kirschner wire attachment	
GB896R	Screw attachment AO small	
GB897R	Screw attachment Aesculap hex	

12.2 Performance data, information about standards

Rotational direction	Right and left rotation, oscillation
Conforming to standard	IEC/DIN EN 60601-1

The product was subject to 500 preparation cycles as a test by the manufacturer and has passed.

The performance data for the individual attachments can be found in the table below.

Drill attachments

Attachment	Connector	Dimensions Length x $\varnothing$ [mm] $\pm 5\%$	Weight [g] $\pm 10\%$	Max. speed [min ⁻¹], approx.	Max. torque [Nm], approx.	Cannulation [mm]
GB878R	for Synthes 511.300	111 x 27	181	1 250	1.5	3.2
GB881R	Large Jacobs chuck $\varnothing$ 0.5 mm to 7.4 mm	89 x 29	163	1 250	1.5	3.2
GB882R	Small Jacobs chuck $\varnothing$ 0.5 mm to 4 mm	76 x 22	100	1 250	1.5	3.2
GB884R	AO small	60 x 22	66	1 250	1.5	2.8
GB886R	Aesculap hex	57 x 24	67	1 250	1.5	3.2
GB887R	Dental	56 x 22	50	1 250	1.5	1.8

Medullary milling attachments

Attachment	Connector	Dimensions Length x $\varnothing$ [mm] $\pm 5\%$	Weight [g] $\pm 10\%$	Max. speed [min ⁻¹], approx.	Max. torque [Nm], approx.	Cannulation [mm]
GB870R	Large Jacobs chuck $\varnothing$ 0.5 mm to 7.4 mm	105 x 29	230	320	5	3.2
GB871R	AO large	88 x 29	170	320	5	3.2
GB872R	Hudson/Zimmer	86 x 25	150	320	5	3.2

Saw attachments

Attachment	Connector	Dimensions Length x $\varnothing$ [mm] $\pm 5\%$	Weight [g] $\pm 10\%$	Max. oscillation frequency [min ⁻¹], approx.	Cannulation [mm]
GB891R	L sagittal saw blade	74 x 33.5	179	15 750	–
GB892R	Reciprocating Saw blade	135 x 24.5	166	15 750	–

Kirschner wire attachment

Attachment	Connector	Dimensions L x W x H [mm] $\pm 5\%$	Weight [g] $\pm 10\%$	Max. speed [min ⁻¹], approx.	Cannulation [mm]
GB894R	Kirschner wire $\varnothing$ 0.6 mm to 3.2 mm	92 x 27 x 124	179	1 250	3.2

Screw attachments

Attachment	Connector	Dimensions L x $\varnothing$ x H [mm] $\pm 5\%$	Weight [g] $\pm 10\%$	Max. speed [min ⁻¹], approx.	Cannulation [mm]
GB896R	AO small	95 x 28 x 129.5	235	165	1.7
GB897R	Aesculap hex	92 x 28 x 129.5	231	165	1.7

12.3 Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

Drilling (clockwise/counterclockwise rotation)

- 60 second application, 60 second pause
- 6 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Medullary reaming (clockwise/counterclockwise rotation)

- 30 second application, 30 second pause
- 8 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Drilling (oscillation)

- 15 second application, 15 second pause
- 3 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Saw mode with GB891R

- 30 second application, 60 second pause
- 4 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Saw mode with GB892R

- 30 second application, 60 second pause
- 5 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Screwing with GB896R/GB897R

- 10 second mechanical application, 10 second manual application, 30 second pause
- 30 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

12.4 Environmental conditions

	Operation	Storage and transport
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

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