



SPINE SURGERY

### AESCULAP<sup>®</sup> TSPACE<sup>®</sup> PEEK

TRANSFORAMINAL LUMBAR INTERBODY FUSION SYSTEM SURGICAL MANUAL

## AESCULAP® THORACOLUMBAR SPINE

## PROTECTING AND PRESERVING SPINAL STABILITY

Modern lifestyle has resulted in increasing physical inactivity among people all over the world. Of the many medical problems associated with this, spinal disorders are among the most critical. This is even more significant as the spinal column is one of the most important structures in the human body.

It supports and stabilizes the upper body and is the center of our musculoskeletal system, which gives the body movement. Our work in the field of spine surgery is dedicated to protecting the spinal column and preserving its stability. We support spine surgeons with durable, reliable products and partner services for reliable procedures and good clinical outcomes (1-7).

Our philosophy of sharing expertise with healthcare professionals and patients allows us to develop innovative implant and instrument systems that help to preserve stability and stabilize the cervical and thoracolumbar spine.

### RELIABLE PARTNER IN SPINE SURGERY

Discover our comprehensive product portfolio by clicking on this area.



www.bbraun.com/spine-surgery

### AESCULAP® TSPACE® PEEK CONTENT



### А **GENERAL INFORMATION** В SURGICAL MANUAL B.1. POSITIONING OF THE PATIENT AND INCISION MARKING **B.2. FASCIAL INCISION** B.3. EXPOSURE AND BLUNT DISSECTION OF THE PARASPINAL MUSCLES B.4. INTRODUCTION OF THE SPINE CLASSICS **RETRACTOR SYSTEM B.5. INSERTION OF PEDICLE SCREW** B.6. REMOVAL OF FACET JOINT B.7. OPENING OF THE DISC AND REMOVAL OF **DISC MATERIAL** B.8. RESTORATION OF DISC HEIGHT B.9. CLEANING OF THE INTERVERTEBRAL SPACE **B.10. DETERMINATION OF IMPLANT SIZE USING** TRIAL IMPLANTS **B.11. IMPLANT INSERTION B.12. FINAL IMPLANT POSITIONING B.13. APPLICATION OF ROD AND SET SCREW B.14. PEDICLE SCREW POSITIONING ON THE** CONTRA-LATERAL SIDE IMPLANT & INSTRUMENT OVERVIEW С

A | GENERAL INFORMATION

### PHILOSOPHY

TSPACE<sup>®</sup> PEEK implants are used for stabilization of the lumbar and thoracic spine through transforaminal approach, monosegmental and multisegmental.

Always use TSPACE<sup>®</sup> PEEK implants in conjunction with an internal fixator.

### TSPACE® PEEK IS DESIGNED TO DELIVER

- PRIMARY STABILITY (8).
- RESTORATION OF THE NATURAL DISC HEIGHT AND LORDOSIS (9).



A | GENERAL INFORMATION

## IMPLANT MATERIAL





The material used is biocompatible PEEK-OPTIMA<sup>®</sup>. PEEK stands for PolyEtherEtherKetone. PEEK-OPTIMA<sup>®</sup> polymer complies with ISO, ASTM and USP standards for use as a medical implant material.

The use of PEEK-OPTIMA<sup>®</sup> as an orthopaedic device material enjoys increased popularity due to the material's special combination of characteristics (9). Its properties include radiolucency, high mechanical strength, high fatigue resistance, a low wear factor and biocompatibility (10-14).

The intrinsic radioscopic transparency of the material provides permeability on X-rays and CT scans, allowing to visualize bone growth adjacent to the implant. This enables a quick and simple assessment of the bone structure and progress towards bone fusion (11). To verify the position of PEEK implants on radioscopic images, non-radiolucent tantalum markers were integrated serving as location indicators (Fig. 1/2).

Of particular interest is the modulus of elasticity of PEEK-OPTIMA<sup>®</sup>, which is similar to that of cortical bone. This modulus of elasticity may reduce implant subsidence and allow for improved bone growth (11, 15).

In vitro results of PEEK-OPTIMA® test specimens show a high long-term material stability after oxygen aging.\* These results correspond with extensive biocompatibility investigations for PEEK-OPTIMA® proving the material suitable for use as a longterm implant (10, 11).

## IMPLANT FEATURES











#### POSITION VERIFICATION DESPITE X-RAY TRANSPARENCY

- The radiolucency of PEEK-OPTIMA<sup>®</sup> enables assessment of the bone structure and progress towards bone fusion (11).
- **I** X-ray pins facilitating implant positioning and localization.

#### IMPLANT DESIGN

- Anatomical shape and serrated profile aim for an implant fit and high primary stability (8).
- Option of filling with bone or bone substitute to enhance bone bridging.

#### IMPLANT VARIETY

Adequate range of sizes to enable the choice of implant size to fit the patient.

#### INSTRUMENT DESIGN

Specifically designed and clearly arranged instruments.

B | SURGICAL MANUAL



#### B.1. POSITIONING OF THE PATIENT AND INCISION MARKING

- A minimally invasive approach requires the patient to be placed on a radiolucent table which allows for AP views of the various anatomic structures (Fig. 1).
- The appropriate position of the longitudinal incision (4-5 cm in length) is determined by using a C-arm. The intended skin incision is marked paraspinally on the right and respectively on the left side (Fig. 2).



#### **B.2. FASCIAL INCISION**

A slightly curved fascial incision 1.5 cm from the mid-line is performed. This allows a firm hold of the speculum and counter retractor, facilitating the exposure of the individual segment (Fig. 3/4).



B.3. EXPOSURE AND BLUNT DISSECTION OF THE PARASPINAL MUSCLES

After splitting of the thoracolumbar fascia a blunt dissection of the paraspinal muscles is performed with the fingertip. In accordance with the palpatory finding, a correction of the skin incision is still possible, as the muscle retractor should be introduced as vertically as possible and in the direction of the interlaminar space. The length of the retractor is selected by using the index finger (Fig. 5/6).



B.4. INTRODUCTION OF THE SPINE CLASSICS RETRACTOR SYSTEM

- Spine Classics retractor system: see brochure 011402
- I The muscle retractor is introduced with closed blades and with the handle in the longitudinal direction. It is then turned 90° and expanded (Fig. 7/8).

B | SURGICAL MANUAL



#### B.5. INSERTION OF PEDICLE SCREW

Using the standard technique, the pedicle screws are inserted, e. g. AESCULAP<sup>®</sup> Ennovate<sup>®</sup> (Fig. 1).

For more information please visit www.bbraun.com/ennovate

#### INFORMATION

Thoracolumbar pedicle screw system shown in images is the AESCULAP\*  $\mathsf{S}^{\mathtt{4*}}$  Spinal System.



#### B.6. REMOVAL OF FACET JOINT

A complete unilateral facetectomy should be considered on the side targeted for the implant insertion, for which an osteotome may be used. The inferior articular process of the facet joint is resected first, then the subjacent superior articular process is resected (Fig 2/3).



#### B.7. OPENING OF THE DISC AND REMOVAL OF DISC MATERIAL

- To open the disc a small window is cut into the annulus (Fig. 4).
- Rongeurs are used to remove the opened annulus (Fig. 5).
- Posterior osteophytes are removed by using Kerrisons.



#### B.8. RESTORATION OF DISC HEIGHT

- The desired restoration of the natural disc height can be set using the distractors. They are available in heights from 7-17 mm in 2 mm increments and can be attached to the t-handle in the set.
- I The distractor must be inserted horizontally and then rotated. Rotating clockwise the distractors are blunt and allow for a blunt height restoration maneuver. Rotating counterclockwise the distractors have a specially designed sharp rim to allow for removal of disc material (Fig. 6).

B | SURGICAL MANUAL



#### B.9. CLEANING OF THE INTERVERTEBRAL SPACE

- I The disc space is cleared using rongeurs, bone curettes and box curettes (Fig. 1).
- I The bone rasps are used for endplate preparation.
- Alternatively, the box curettes can be used.

#### INFORMATION

Make certain that the endplates of the neighboring vertebral bodies are not weakened, in order to minimize the risk of migration.

Make certain that the implant bed is properly prepared to avoid damage to the implant when it is driven in.

Use the nerve root retractors to protect the dura during insertion.

#### B.10. DETERMINATION OF IMPLANT SIZE USING TRIAL IMPLANTS

- The TSPACE<sup>®</sup> PEEK trial implants are available in heights from 7 17 mm in 2 mm increments.
- Using the slap hammer the desired trial implant is inserted (Fig. 2).

#### INFORMATION

The trials are essential to ensure the correct implant size to be used.





#### **B.11. IMPLANT INSERTION**

■ The TSPACE<sup>®</sup> PEEK implant should be filled with bone or bone substitutes by using the packing block.

#### INFORMATION

- **I** Do not use force during filling to avoid implant damaging.
- Mount the implant on the insertion instrument hand-tight as far as it will go.
- When inserting the implant into the intervertebral space, avoid canting and levering, and take care to maintain an alignment parallel to the endplate.
- Do not use excessive force when mounting or implanting the implant.
- I It is recommended to place bone material harvested from the facet joint around the TSPACE<sup>®</sup> PEEK implant.
- I The insertion guide supports the implantation of the TSPACE<sup>®</sup> PEEK implant.
- The TSPACE<sup>®</sup> PEEK implant is inserted into the disc space using the implant inserter.
- The TSPACE<sup>®</sup> PEEK inserter has a clamp mechanism.

B | SURGICAL MANUAL



#### **B.12. FINAL IMPLANT POSITIONING**

- Using the straight or angled impactor the implant is rotated 90° to achieve the final positioning (Fig. 1/2).
- **I** X-ray control to verify the implant positioning.
- I It is recommended to put bone material harvested from the facet joint around the TSPACE<sup>®</sup> PEEK implant.







#### B.13. APPLICATION OF ROD AND SET SCREW

- Final assembly of the thoracolumbar pedicle screw system (Fig. 3).
- Compression is applied to the pedicle screws to support the contact area between the TSPACE® PEEK implant and the endplates.
- I Final tightening of the pedicle screws and closure.

B.14. PEDICLE SCREW POSITIONING ON THE CONTRA-LATERAL SIDE

I Thoracolumbar pedicle screw system is applied on the contralateral side (Fig. 4).

To view our full portfolio please visit www.bbraun.com

#### INFORMATION

Thoracolumbar pedicle screw system shown in images is the AESCULAP\*  $\mathsf{S}^{\scriptscriptstyle 4^{\scriptscriptstyle 0}}$  Spinal System.

C | IMPLANT & INSTRUMENT OVERVIEW



Article No.	Size (Length x Width x Height )	Angle
FJ687P	26 x 11 x 7 mm	5°
FJ689P	26 x 11 x 9 mm	5°
FJ691P	26 x 11 x 11 mm	5°
FJ693P	26 x 11 x 13 mm	5°
FJ695P	26 x 11 x 15 mm	5°
FJ697P	26 x 11 x 17 mm	5°
FJ317P	30 x 11 x 7 mm	5°
FJ319P	30 x 11 x 9 mm	5°
FJ321P	30 x 11 x 11 mm	5°
FJ323P	30 x 11 x 13 mm	5°
FJ325P	30 x 11 x 15 mm	5°
FJ327P	30 x 11 x 17mm	5°
FJ707P	34 x 11 x 7 mm	5°
FJ709P	34 x 11 x 9 mm	5°
FJ721P	34 x 11 x 11 mm	5°
FJ723P	34 x 11 x 13 mm	5°
FJ725P	34 x 11 x 15 mm	5°
FJ727P	34 x 11 x 17 mm	5°

### FJ015 – TSPACE<sup>®</sup> PEEK INSTRUMENTATION

TSPACE<sup>®</sup> PEEK – Preparation Tray

INSTRUMENTS	Article No.	Description	Quantity
-	FJ646R	T-Handle for distractors	1
	FJ647R	Distractor, 7 mm	1
	FJ649R	Distractor, 9 mm	1
	FJ651R	Distractor, 11 mm	1
	FJ653R	Distractor, 13 mm	1
	FJ655R	Distractor, 15 mm	1*
	FJ657R	Distractor, 17 mm	1*
	FJ658R	Straight osteotome, 8 mm	1
	FJ679R	Left angled bone curette, 45°	1
	FJ680R	Right angled bone curette, 45°	1
	FJ698R	Left angled bone curette, 20°	1*
	FJ699R	Right angled bone curette, 20°	1*
	FJ681R	Curette, straight	1
	FJ682R	Left angled curette, 45°	1
	FJ683R	Right angled curette, 45°	1
	FJ702R	Left angled curette, 20°	1*
	FJ703R	Right angled curette, 20°	1*
	FJ685R	Left angled bone rasp, 45°	1
	FJ686R	Right angled bone rasp, 45°	1
	FJ704R	Left angled bone rasp, 20°	1*
	FJ705R	Right angled bone rasp, 20°	1*
	FJ641P	Tray for prep. instr.	1

\*This article is optional. Recommended container for FJ641P: JK442/JN442 and lid JK489.

C | IMPLANT & INSTRUMENT OVERVIEW

### FJ015 – TSPACE<sup>®</sup> PEEK INSTRUMENTATION

TSPACE<sup>®</sup> PEEK – Implantation Tray

INSTRUMENTS	Article No.	Description	Quantity
	FJ051R	Retractor S	1
	FJ052R	Retractor M	1
ISAN -	FJ053R	Retractor L	1
	FJ054R	Retractor XL	1
	FJ666R	Slap hammer for trials	1
	FJ667R	TSPACE <sup>®</sup> trial, 7 mm	1
	FJ669R	TSPACE <sup>®</sup> trial, 9 mm	1
7 3x85x1x73x838xxx11 5	FJ671R	TSPACE <sup>®</sup> trial, 11 mm	1
	FJ673R	TSPACE <sup>®</sup> trial, 13 mm	1
	FJ675R	TSPACE <sup>®</sup> trial, 15 mm	1*
	FJ677R	TSPACE <sup>®</sup> trial, 17 mm	1*
	FJ660R	TSPACE <sup>®</sup> PEEK insertion instr.	1
	FJ661R	TSPACE <sup>®</sup> insertion guide	1
	FJ662R	TSPACE <sup>®</sup> impactor, straight	1
	FJ663R	TSPACE <sup>®</sup> impactor, angled	1
Table a Mar 2 Ann	FJ664R	TSPACE <sup>®</sup> PEEK packing block	1
	FF913R	Punch	1
	FJ642P	Tray for implantation instruments	1

#### REFERENCES

- Tropiano P, Bronsard JJ, Louis C, Tallet JM, Sauget Y. Three column stabilisation through a posterior appraoch with a titanium PLASMAPORE<sup>®</sup> intervertebral block (PROSPACE<sup>®</sup>). Radiological and clinical study after 4 years. Rivista di Neuroradiologia. 1999;12(Suppl 1):89–94.
- (2) Kroppenstedt S, Gulde M, Schonmayr R. Radiological comparison of instrumented posterior lumbar interbody fusion with one or two closed-box PLASMAPORE<sup>®</sup> coated titanium cages. Follow-up study over more than seven years. Spine. 2008;33(19):2083-8.
- (3) Kreinest M, Schmahl D, Grutzner PA, Matschke S. Radiological Results and Clinical Patient Outcome After Implantation of a Hydraulic Expandable Vertebral Body Replacement following Traumatic Vertebral Fractures in the Thoracic and Lumbar Spine: A 3-Year Follow-Up. Spine (Phila Pa 1976). 2017 Apr 15;42(8):E482-E489.
- (4) Takeuchi M, Yasuda M, Niwa A, Wakao N, Nakura T, Osuka K et al. PLASMAPORE<sup>®</sup>-coated titanium cervical cages induce more rapid and complete bone fusion after anterior cervical discectomy and fusion as compared to noncoated titanium cage. World Neurosurgery. 2014;82(3/4):519-22.
- (5) Vanek P, Bradac O, Konopkova R, de Lacy P, Lacman J, Benes V. Treatment of thoracolumbar trauma by short-segment percutaneous transpedicular screw instrumentation: prospective comparative study with a minimum 2-year follow-up. J Neurosurg Spine. 2014;20:150-6.
- (6) Beisse R. Endoscopic surgery on the thoracolumbar junction of the spine. Eur Spine J. 2006;15:687-704.
- (7) Finger T, Bayerl S, Onken J, Czabanka M, Woitzik J, Vajkoczy P. Sacropelvic fixation versus fusion to the sacrum for spondylodesis in multilevel degenerative spine disease. Eur Spine J. 2014;23:1013-20.
- (8) Mechanical testing of the AESCULAP® TSPACE® PEEK 1st Generation implant evaluating shear resistance. Tuttlingen, 2006. The AESCULAP® TSPACE® PEEK 1st Generation implant was tested in a shear force test following internal standards to evaluate the risk of implant dislocation related to shear forces occurring in situ. For this test, a synthetic bone-equivalent polyurethane test block was used. A quasi-static shear force was applied to the implant. The implant was axially loaded during the whole test with a constant force to simulate a clinically relevant

loading situation. Based on the maximum shear load obtained in the test, the risk of implant dislocation due to shear forces at a physiological load situation is unlikely to occur in the case of the AESCULAP<sup>®</sup> TSPACE<sup>®</sup> PEEK 1st Generation implant.

- (9) Kurtz SM, Devine JN. PEEK biomaterials in trauma, orthopedic, and spinal implants. Biomaterials 2007; 28(32):4845-69.
- (10) Morrison C, Macnair R, MacDonald C, Wykman A, Goldie I, Grant MH. In vitro biocompatibility testing of polymers for orthopaedic implants using cultured fibroblasts and osteoblasts. Biomaterials 1995; 16(13):987-92.
- (11) Invibio<sup>®</sup> Biomaterial Solutions. PEEK-OPTIMA<sup>®</sup> Natural Typical Material Properties; 2013. Available from: URL: https://invibio.com/materials/peek-optima-natural.
- (12) Chen Y, Wang X, Lu X, Yang L, Yang H, Yuan W et al. Comparison of titanium and polyetheretherketone (PEEK) cages in the surgical treatment of multilevel cervical spondylotic myelopathy: a prospective, randomized, control study with over 7-year follow-up. Eur Spine J 2013; 22(7):1539-46.
- (13) Macnair R, Rodgers EH, MacDonald C, Wykman A, Goldie I, Grant MH. The response of primary rat and human osteoblasts and an immortalized rat osteoblast cell line to orthopaedic materials: comparative sensitivity of several toxicity indices. J Mater Sci Mater Med 1997; 8(2):105-11.
- (14) Dr. Stuart Green, Keith Cartwright. An investigation into the effect of accelerated oxygen ageing on the properties of PEEK-OPTIMA<sup>®</sup>. March 2004. (B.DoCS-ID 21870982) Tensile, flex and impact specimens in PEEK-OPTIMA<sup>®</sup> were tested according to appropriate ISO standards (ISO 527, ISO 178 and ISO 180). A comparison was made between control PEEK-OP-TIMA<sup>®</sup> and aged PEEK-OPTIMA<sup>®</sup> to investigate the effect of accelerated oxygen ageing on the material properties. Accelerated ageing was conducted exposing the specimens 40 days to 70°C oxygen at 5 bars pressure. The results show no significant effect on the mechanical properties of the PEEK polymer with the aged and control specimens showing similar values. The retention of good mechanical properties after the intense ageing cycle demonstrates that PEEK-OPTIMA<sup>®</sup> is very resistant to oxygen ageing.
- (15) Kuhn JL, Goldstein SA, Choi K, London M, Feldkamp LA, Matthews LS. Comparison of the trabecular and cortical tissue moduli from human iliac crests. J Orthop Res 1989; 7(6):876-84.

### AESCULAP<sup>®</sup> – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "AESCULAP" and the product trademarks "Ennovate", "S<sup>4</sup>" and "TSPACE" are registered trademarks of Aesculap AG. "PEEK-OPTIMA" is a registered trademark of Invibio Biomaterial Solutions.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.