## BiCONTACT®

|            |                                                                                                                                                                                                                                                                                                                                                                           |
|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>GB</b>  | <b>Instructions for use</b><br>Hip endoprosthesis shaft                                                                                                                                                                                                                                                                                                                   |
| <b>USA</b> | <b>Instructions for use</b><br>Hip endoprosthesis shaft<br><i>Note for U.S. users</i><br><i>U.S. cleared device instructions for use can be found at</i><br><i>www.aesculapimplantsystems.com "Products &amp; Services".</i><br><i>If you do not have internet access please call your local Aesculap</i><br><i>representative or customer service at 1-800-258-1946.</i> |
| <b>D</b>   | <b>Gebruiksaanweisung</b><br>Hüftendoprothesenschaft                                                                                                                                                                                                                                                                                                                      |
| <b>F</b>   | <b>Mode d'emploi</b><br>Tige fémorale prothétique                                                                                                                                                                                                                                                                                                                         |
| <b>E</b>   | <b>Instrucciones de manejo</b><br>Vástago para endoprótesis de cadera                                                                                                                                                                                                                                                                                                     |
| <b>I</b>   | <b>Istruzioni per l'uso</b><br>Stelo di endoprotesi coxo-femorale                                                                                                                                                                                                                                                                                                         |
| <b>P</b>   | <b>Instruções de utilização</b><br>Haste para endoprótese de quadril                                                                                                                                                                                                                                                                                                      |
| <b>NL</b>  | <b>Gebruiksaanwijzing</b><br>Heupprotheseschaft                                                                                                                                                                                                                                                                                                                           |
| <b>DK</b>  | <b>Brugsanvisning</b><br>Hofte-endoproteseskraft                                                                                                                                                                                                                                                                                                                          |
| <b>N</b>   | <b>Bruksanvisning</b><br>Hoftendoproteseskraft                                                                                                                                                                                                                                                                                                                            |
| <b>S</b>   | <b>Bruksanvisning</b><br>Höftendoproteseskraft                                                                                                                                                                                                                                                                                                                            |
| <b>FIN</b> | <b>Käyttöohjeet</b><br>Lonkan sisäproteesin varsi                                                                                                                                                                                                                                                                                                                         |
| <b>GR</b>  | <b>Οδηγίες χρήσης</b><br>Στέλεχος ενδοπρόθεσης ισχίου                                                                                                                                                                                                                                                                                                                     |
| <b>CZ</b>  | <b>Návod k použití</b><br>Dřik kyčelní endoprotézy                                                                                                                                                                                                                                                                                                                        |

# B | BRAUN

SHARING EXPERTISE

**Aesculap AG**  
Am Aesculap-Platz  
78532 Tuttlingen  
Germany  
Phone +49 7461 95-0  
Fax +49 7461 95-2600  
www.aesculap.de

TA-Nr. 010238 01/09 Änd.-Nr. 31521/32607

CE marking according to directive 93/42/EEC  
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG  
Marquage CE conforme à la directive 93/42/CEE  
Identificación CE en conformidad con la directriz 93/42/CEE  
Marchio CE conforme alla direttiva 93/42/CEE  
Símbolo CE, em conformidade com a Directiva 93/42/CEE  
CE-certificering conform richtlijn 93/42/EEG  
CE-markering iht. retningslinie 93/42/EEC  
CE-merking ifølge retningslinje 93/42/EWG  
CE-märkning i enlighet med direktiv 93/42/EEG  
93/42/EEC-standardin mukainen CE-hyväksyntä  
Ενδειξη CE σύμφωνα με την Οδηγία 93/42/EEC  
CE-označení podle směrnice 93/42/EHS

Technical alterations reserved  
Technische Änderungen vorbehalten  
Sous réserve de modifications techniques  
Sujeto a modificaciones técnicas  
Con riserva di modifiche tecniche  
Salvo alterações técnicas  
Technische wijzigingen voorbehouden  
Retten til tekniske ændringer forbeholdes  
Vi tar forbehold om eventuelle tekniske endringer  
Med reservation för eventuella tekniska ändringar  
Oikeusteknisistä syistä johtuviin muutoksiin pidätetään  
Επιφύλασσόμεθα για τεχνικές αλλαγές  
Technické změny vyhrazeny

## Intended use

- Component of a human hip endoprosthesis: Hip endoprosthesis shaft
- Combination with Aesculap hip endoprosthesis components
- Implant range: BiCONTACT®, BiCONTACT® S, H, SD, or N and BiCONTACT® Revision.
- Implantation without bone cement, for PLASMAPORE®-coated shaft
- Implantation with bone cement, for uncoated shaft

## Materials

The materials used in the implant are listed on the package.

- ISOTAN®F titanium forged alloy Ti6Al4V according to ISO 5832-3
- ISODUR®F Cobalt-Chromium forged alloy CoCr29Mo according to ISO 5832-12
- ISOTAN®P pure titanium, ISO 5832-2
- PLASMAPORE® surface coating made of pure titanium acc. to ISO 5832-2 or PLASMAPORE®µ-CaP surface coating made of pure titanium acc. to ISO 5832-2, with additional calcium phosphate coating

BiCONTACT®, ISOTAN®, ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

## Indications

Use for acute disorders of the hip that cannot be treated with other therapies:

- Degenerative osteoarthritis
- Rheumatic arthritis
- Fractures of the joint
- Femoral head necrosis

**Risk of aseptic implant loosening when using bone cement with uncoated BiCONTACT® prosthesis shafts made of ISOTAN®.**

- Do not use uncoated BiCONTACT® prosthesis shafts made of ISOTAN® (titanium forged alloy ISO 5832-3) for implantations with bone cement.

Exceptional indications: Patients with proven foreign-body sensitivity for the alloy elements nickel, cobalt or chromium where a prosthesis shaft anchoring with bone cement is indicated at the same time.

## Indications in cases of interventions to replace hip endoprosthesis shafts, for the shaft type BiCONTACT® Revision

- Cement-free revision operations
- Periprosthetic fractures
- Subprosthetic fractures
- Treatment of large bone defects on the proximal femur

## Contraindications

**Implant failure due to excessive body weight of the patient and mainly distal fixation of the small BiCONTACT® implants!**

- Avoid distal implant fixation by appropriate femoral bone preparation.
- Bore up the distal marrow cavity, using flexible marrow cavity drills.
- Ensure sufficiently deep prosthesis seat through proximal-medial support.
- Carry out intraoperative check of the prosthesis seat, using an imaging device.

In case of distal implant fixation observe maximum body weight of the patient: 50 kg for NK709T and NK710T (BiCONTACT® SD sizes 9 and 10), 60 kg for NJ208T, NJ209T, NJ210T (BiCONTACT® D sizes 8,9 and 10) and 65 kg for NK210T (BiCONTACT® Revision SD size 11).

Do not apply in the presence of:

- Possibility of reconstructive interventions for the treatment of the joint illness, e.g. osteotomic resetting
- Acute or chronic infections near the joint, or systemic infections
- Secondary diseases that could influence joint implant functionality
- Systemic diseases and metabolic disturbances
- Acute osteoporosis or osteomalacia
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumor in the region of implant fixation
- Bone malformations, defective axial positioning, or other bone conditions that rule out implantation of a hip-joint prosthesis
- Anticipated excessive load on the joint implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

## Contraindications for interventions to replace hip endoprosthesis shafts

The replacement of an endoprosthetic implant is a difficult intervention subject to individual conditions.

- To arrive at a decision on whether to provide such replacement, the operating surgeon must assess and, if necessary, consider the contraindications with regard to the success of the intervention.

## Side-effects or adverse interactions

- Changes in position, loosening, wear and tear on, or fracture of implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Venous thrombosis, lung embolism, and cardiac arrest
- Tissue reaction to implant materials
- Injury to vessels and nerves
- Hematoma and impaired wound healing

- Periarticular calcification
  - Decreased joint mobility and flexibility
  - Arthralgia and decreased tolerance for exercise
- When using the prosthesis shaft BiCONTACT® Revision:
- Secondary operation to remove the distal locking screw (if applicable)

## Safety information

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- General risk factors associated with surgical procedures are not described in the present instructions for use.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Under no circumstances should damaged components or surgically excised components be used.
- Implants once used must not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the lot number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keep the patient well informed.
- Damage to the weight-bearing bone cement and/or bone structures can give rise to loosening of the components, fracture of the bone or implant, as well as to other grave complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Only combine modular implant components with the appropriate Aesculap hip endoprostheses.
- Observe material, friction coupling diameter and cone specifications.
- Note further restrictions on combined implants.
- Avoid damaging the implant, especially in the neck or cone region, by applying instruments (e.g. HF surgical devices) close to the implant.

**Risk of implant component breakage due to combination with implant components from other manufacturers!**

- Use Aesculap implant components only.

## Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are sterilized by irradiation (minimum dose 25 kGy).
- Store implant components in their original packaging and only remove them from such packaging immediately before use.
- Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiration dates or whose packaging is damaged.

## Sterilization method and parameters

- Sterilize with steam, taking note of the following:
  - The sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, and validated in conformance with EN ISO 17665 or EN 554/ISO 13683). In case of application of the fractionated vacuum procedure the sterilization has to be carried out for a minimum of 5 minutes at 134 °C and at 2 bar pressure.

Implants made of non-coated metal may be individually resterilized in steam following appropriate preliminary cleaning and disinfecting.

- Make sure the implant mechanism is not damaged during cleaning and resterilization.

Prosthesis cone protection caps that are labeled accordingly can be resterilized with steam.

- If the resterilization of ceramic implants should become necessary, follow the directions in the relevant Aesculap instructions for use.

**The responsibility for resterilization of implant components lies with the end-user. Such resterilization is expressly prohibited by law in certain countries (e.g. France).**

- Never resterilize or reuse surface-coated implants (PLASMAPORE® and PLASMAPORE®µ-CaP).

## Application

The operating surgeon draws up an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bones
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components must be ready to hand
- Operating conditions must be highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, must be complete and in working condition.
- The operating surgeon and operating room team must be thoroughly conversant with the operating technique, as well as the range of implants and instruments to be applied; complete information on these subjects must be readily available at the workplace.

- Those performing operations must be thoroughly conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant articles by medical specialists from the professional literature.
- In case of an ambiguous preoperative situation and if an implant is already present in the joint region to be provided for, information was sought from the manufacturer.

The operative procedure has been explained to the patient, and the latter's understanding of the following information has been documented:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint can become loose owing to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the body weight and the load put on the joint.
- The prosthetic joint must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant becomes loose.
- In the event corrective surgery is performed, under certain circumstances it may not be possible to restore joint mobility and flexibility.
- The revision of a hip endoprosthesis shaft is a complex joint-replacement operation; it is generally inferior to the primary joint replacement.
- The patient must have their physician carry out follow-up examinations of the prosthetic joint at regular intervals.

The implantation and implantation site are prepared in the following way:

- Having performed the head osteotomy and opened the marrow cavity, finish the femur with BiCONTACT®-specific A- and B-osteoprofilers.
- Carry out intraoperative orientation of the implant components. When doing this, pay attention to the resection line and the trochanter wing.
- Select the implant according to B-osteoprofiler last introduced in the correct position.
- Before inserting the implants, carry out a test reposition and examine joint mobility, joint stability and leg length.

**Bone fractures in the implant site reduce the viability of the implant anchorage.**

- Avoid bone fractures through careful operative techniques.
- Treat bone fractures with appropriate intra- and postoperative therapies.
- Handle the implant components properly.

- Do not under any circumstances allow the implant surface to be damaged.
- Ensure that the cone sizes of the prosthesis shafts and the prosthesis heads match (see cone size on the implant packaging, e.g. 12/14 or 8/10).
- Only remove the protective cap from the prosthesis cone immediately before setting the prosthesis head in place.
- Rinse, clean and dry the outside cone of the shaft and, if necessary, the inside cone of the prosthesis head prior to setting the prosthesis head in place.
- Couple the prosthesis head and cone at room temperature only. If necessary, allow the implant to cool to room temperature.
- Prepare and apply bone cement as described in the instructions for use provided by the cement supplier.
- If a distal centralizer is used with the cemented BiCONTACT® prosthesis shaft: Follow the instructions for use for the Aesculap centralizer.
- Before closing the wound, ensure, if necessary by means of an image converter, that the implant components are correctly positioned.
- In order to forestall abnormal wear and tear on the prosthesis: Remove any loose bone cement or bone chips before closing the wound.

When carrying out interventions to replace hip endoprosthesis shafts of the type BiCONTACT® Revision, the following instructions must be observed:

- Prepare the marrow cavity with A and B reamers if you are going to remove the implant through a proximal access (if necessary with an additional bone window). Pay attention to the reamer marking and size according to the planned prosthesis size and length.
- When preparing the marrow cavity, work the reamers manually, taking care to remove as little bone substance as possible. If necessary, finish the marrow cavity preparation with B-osteoprofilers and conclude by implanting a BiCONTACT® revision shaft.
- If the primary stability is not sufficient (axial fit and rotation), carry out distal locking. Use special aiming or/and image amplification device for this purpose.
- Use Aesculap 5-mm locking screws for locking. Follow the instructions for use of the locking screws.

Transfemoral access:

- When operating through a transfemoral access, lock the BiCONTACT® revision shaft.
- Use A-reamers to prepare the distal fitting of the prosthesis. When doing this, keep to a distal anchoring length of at least 100 mm.

**Risk of locking screw breakage due to permanent distal force transmission!**

- Remove the locking screws after 1 to 2 years (reference period), depending on bone remodeling and integration of the revision implant in the proximal bone fraction. The appropriate timing of this intervention depends on the individual patient's situation; it can be later in cases of delayed bone remodeling.

*Note*

*The BiCONTACT® revision shafts are locked for the primary stabilization of the prosthesis shaft in cases of severe bone defects. The distal force transmission via the locking screws is designed to be a device for temporary stabilization.*

Further information on Aesculap implant systems is readily available from Aesculap or the Aesculap office responsible.

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