BiCONTACT® MS

GB	Instructions for use Hip endoprosthesis shaft
USA	Instructions for use Hip endoprosthesis shaft <i>Note for U.S. users</i> <i>U.S. cleared device instructions for use can be found at</i> <i>www.aesculapimplantsystems.com "Products & Services".</i> <i>If you do not have internet access please call your local Aesculap</i> <i>representative or customer service at 1-800-258-1946.</i>
D	Gebrauchsanweisung Hüftendoprothesenschaft
F	Mode d'emploi Tige fémorale prothétique
E	Instrucciones de manejo Vástago para endoprótesis de cadera
I	Istruzioni per l'uso Stelo di endoprotesi coxo-femorale
P	Instruções de utilização Haste para endoprótese de quadril
NL	Gebruiksaanwijzing Heupprotheseschacht
DK	Brugsanvisning Hofte-endoproteseskraft
N	Bruksanvisning Hofteproteseskraft
S	Bruksanvisning Höftendoproteseskraft
FIN	Käyttöohjeet Lonkan sisäproteesin varsi
GR	Οδηγίες χρήσης Στέλεχος ενδοπρόθεσης ισχίου
CZ	Návod k použití Dřik kyčelní endoprotézy

B|BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 012639 06/08

0123

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Símbolo CE, em conformidade com a Directiva 93/42/CEE
CE-certificering conform richtlijn 93/42/EEG
CE-markering iht. retningslinje 93/42/EEC
CE-merking ifølge retningslinje 93/42/EWG
CE-märkning i enlighet med direktiv 93/42/EEG
93/42/EEC-standardin mukainen CE-hyväksyntä
Ενδειξη CE σύμφωνα με την Οδηγία 93/42/EEC
CE-označení podle směrnice 93/42/EHS

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Technische wijzigingen voorbehouden
Retten til tekniske ændringer forbeholdes
Vi tar forbehold om eventuelle tekniske endringer
Med reservation för eventuella tekniska ändringar
Oikeusteknisistä syistä johtuviin muutoksiin pidätetään
Επιφύλασσόμεθα για τεχνικές αλλαγές
Technické změny vyhrazeny

Intended use

- Component of a human hip endoprosthesis: Hip endoprosthesis shaft
- Combination with Aesculap hip endoprosthesis components
- Implant range: BiCONTACT® MS
- Implantation without bone cement

Materials

The materials used in the implant are listed on the package.

- ISOTAN®F titanium forged alloy Ti6Al4V according to ISO 5832-3
 - PLASMAPORE® surface coating, pure titanium according to ISO 5832-2
- BiCONTACT®, ISOTAN® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for acute disorders of the hip that cannot be treated by other therapies:

- Degenerative osteoarthritis
- Femoral head necrosis

Contraindications

Do not use in the presence of:

- Possibility of reconstructive interventions for the treatment of the joint illness, e.g. osteotomic resetting
- Acute or chronic infections near the joint, or systemic infections
- Secondary diseases that could influence joint implant functionality
- Systemic diseases and metabolic disorders
- Acute osteoporosis or osteomalacia
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Bone malformations, defective axial positioning, or other bone conditions that rule out implantation of a hip-joint prosthesis
- Anticipated excessive load on the joint implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side-effects and adverse interactions

- Changes in position, loosening, wear and tear on, or fracture of implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Venous thrombosis, lung embolism, and cardiac arrest
- Tissue reaction to implant materials
- Injury to blood vessels and nerves
- Hematomas and wound healing disorders
- Periarticular calcification
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise

Safety notes

- Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- The implant components are tested and licensed for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not use damaged or surgically excised components under any circumstances.
- Implants that have already been used must not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Only combine modular implant components with the appropriate Aesculap hip endoprostheses.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damaging the implant, especially in the neck or cone region, by applying instruments (e.g. HF surgical devices) close to the implant.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).

Danger to the patient and possible loss of implant functionality due to resterilization!

- **Do not resterilize or reuse surface-coated implants (PLASMAPORE®) under any circumstances.**
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The operating procedure has been explained to the patient, whose understanding of the following information has been documented:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint insert can come loose due to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the body weight and the load put on the joint.
- The prosthetic joint must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may become necessary if the implant loosens.
- In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The revision of a hip endoprosthesis shaft is a complex joint-replacement operation; it is generally inferior to the primary joint replacement.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

The implantation and implantation site are prepared in the following way:

- Following the femoral head osteotomy and after opening the medullary cavity, prepare the femur with the BiCONTACT®-specific A and B osteoprofilers.
- Carry out intraoperative orientation of the implant components. When doing this, pay attention to the resection line and the trochanter wing.
- Select the implant according to B-osteoprofiler last introduced in the correct position.
- Before inserting the implants, carry out a test reposition and examine joint mobility, joint stability and leg length.

Bone fractures in the implant site reduce the viability of the implant anchorage.

- **Avoid bone fractures through careful operative techniques.**
- **Treat bone fractures with appropriate intra- and postoperative therapies.**
- **Handle the implant components properly.**
- Do not under any circumstances allow the implant surfaces to be damaged.
- Verify that the cone sizes of the prosthesis shafts are the same as the cone sizes of the prosthesis heads (see cone size on the implant packaging, e.g. 12/14 or 8/10).
- Only remove the protective cap from the prosthesis cone immediately before setting the prosthesis head in place.
- Rinse, clean and dry the outside cone of the stem and, if necessary, the inside cone of the prosthesis heads prior to setting the prosthesis head in place.
- Couple the prosthesis head and cone at room temperature only. If necessary, allow the implant to cool to room temperature.
- Before closing the wound, ensure, if necessary by means of an image converter, that the implant components are correctly positioned.
- To avoid abnormal wear and tear on the prosthesis: Remove any loose bone cement or bone chips before closing the wound.

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

TA-Nr.: 012639 06/08

Risk of implant component breakage due to combination with implant components from other manufacturers!

- **Use Aesculap implant components only.**

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