- GB** Instructions for use/Technical description
CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm
- CDN** Instructions for use/Technical description
CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm
- USA** Instructions for use/Technical description
CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm
Note for U.S. users
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
- D** Gebrauchsanweisung/Technische Beschreibung
CRANIOFIX2-Titanklammern Ø 11 mm, Ø 16 mm, Ø 20 mm
- F** Mode d'emploi/Description technique
Clamps en titane CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- E** Instrucciones de manejo/Descripción técnica
Pinzas de titanio CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- I** Istruzioni per l'uso/Descrizione tecnica
Clamp in titanio CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- P** Instruções de utilização/Descrição técnica
Grampos de titânio CRANIOFIX2 de Ø 11 mm, Ø 16 mm, Ø 20 mm
- NL** Gebruiksaanwijzing/Technische beschrijving
CRANIOFIX2-titaniem klemmen Ø 11 mm, Ø 16 mm, Ø 20 mm
- S** Bruksanvisning/Teknisk beskrivning
CRANIOFIX2-titanklamrar Ø 11 mm, Ø 16 mm, Ø 20 mm
- RUS** Инструкция по применению/Техническое описание
Титановые клипсы CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- CZ** Návod k použití/Technický popis
Titanové svorky CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- PL** Instrukcja użytkowania/Opis techniczny
Zaciski tytanowe CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- SK** Návod na použitie/Technický opis
Titánové svorky CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- TR** Kullanım Kılavuzu/Teknik açıklama
CRANIOFIX2 titanyum klemensler Ø 11 mm, Ø 16 mm, Ø 20 mm
- KR** 사용 설명서 / 기술 설명
CRANIOFIX2 티타늄 클램프 Ø 11mm, Ø 16mm, Ø 20mm


B | BRAUN

SHARING EXPERTISE


Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company


TA-Nr. 011200 09/16 V6 Änd.-Nr. 49759


A

B

C

D

Aesculap®

CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm

Legend

Fig. A

CRANIOFIX2 titanium clamps FF490T/FF490T-UNI, Ø 11 mm

CRANIOFIX2 titanium clamps FF491T/FF491T-UNI, Ø 16 mm

CRANIOFIX2 titanium clamps FF492T/FF492T-UNI, Ø 20 mm

Fig. B

1 CRANIOFIX2 titanium clamps FF491T/FF491T-UNI, Ø 16 mm

CRANIOFIX2 titanium clamps FF492T/FF492T-UNI, Ø 20 mm

2 CRANIOFIX2 titanium clamps FF490T/FF490T-UNI, Ø 11 mm

Fig. C

1 Correct positioning of the CRANIOFIX2 titanium clamps FF491T/FF491T-UNI/FF492T/FF492T-UNI: Burr hole Ø 12/15 mm

2 Incorrect positioning of the CRANIOFIX2 titanium clamps FF491T/FF491T-UNI/FF492T/FF492T-UNI: Craniotomy gap

Fig. D

1 CRANIOFIX2 titanium clamps

2 New craniotomy line

3 Old craniotomy line

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Product description

The CRANIOFIX2 titanium clamp system from Aesculap consists of sterile CRANIOFIX2 titanium clamps for single use and reusable neurosurgical instruments for the application of titanium clamps.

Intended use

The CRANIOFIX2 titanium clamp system from Aesculap is used for the fixation of craniotomized bone flaps and fractures to the neurocranium.

The CRANIOFIX2 titanium clamp FF490T/FF490T-UNI, Ø 11 mm, is inserted in the craniotomy gap (see Fig. B) or in the fracture gap.

The CRANIOFIX2 titanium clamp FF491T/FF491T-UNI, Ø 16 mm, is inserted in a burr hole with Ø 12 mm (see Fig. C) or is used to bridge defects in case of fractures.

The CRANIOFIX2 titanium clamp FF492T/FF492T-UNI, Ø 20 mm, is inserted in a burr hole with Ø 15 mm (see Fig. C) or is used to bridge defects in case of fractures.

Materials

The materials used in the implant are listed on the packaging:

- Ti6Al4V acc. to ISO 5832-3

Indications

The CRANIOFIX2 titanium clamp system from Aesculap is used for the fixation of craniotomized bone flaps and fractures to the neurocranium.

Contraindications

- Patients with hypersensitivity to metals or allergies to the implant materials, inflammations in the region of the implant site or insufficient compliance and co-operation of the patient
- Bone conditions that rule out the application of CRANIOFIX2 titanium clamps
- Artificial cranial bone flaps (excluding acrylic patient-specific cranial implants from the company Anatomics Pty Ltd, intended for fixation with CRANIOFIX2 titanium clamps)
- Bone tumors in the area supporting the implant
- Degenerative bone diseases
- Missing dura mater
- Application in the facial skull (viscerocranium) and in the orbital or skull-base region
- Relative contraindications, e.g. infections, osseous systemic diseases or metabolic disorders, including medical or surgical conditions, which could impede the success of the implantation.

Side effects and interactions

The following side-effects and/or interactions are referred to as part of the legal obligation to inform. These are procedure-specific, non product-specific and not limited to infections, nerve damage, hematomas and wound healing disturbances. The patient must be informed of these risks and other risks associated with neurosurgery, general surgery and the administration of anesthetics.

The following risks may occur in connection with the use of artificial bone flaps:

- Unstable or sinking artificial bone flap
- Damage to the artificial bone flap
- Surgical complications of a different nature

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or lack of asepsis.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically removed components.
- Implants already used must not be reused.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- Implants produce moderate MRI artifacts.

- Inappropriate fixation of craniotomized bone flaps or any combination with other fixation methods (sutures, mini plates or clamp systems from other manufacturers) can result in the formation of steps, shifts of the calvarium and loosening and breakage of implant components. For maximum stability of a craniotomized bone flap, the use of three CRANIOFIX2 titanium clamps in a triangular formation are recommended. Larger bone flaps or a skull fracture with multiple fragments may require a larger number of CRANIOFIX2 titanium clamps to ensure sufficient stability.
- During application of CRANIOFIX2 titanium clamps, friction between the upper titanium disc and the pin indentation can lead to spark formation. Appropriate safety precautions must be taken to prevent this. Due to the risk of fire, surgical interventions with this product must not be carried out in rooms where flammable anesthetics or other flammable gases, liquids, objects or oxidizing agents are stored.
- After the operation, exceptionally heavy mechanical strains on the fixated bone flap must be avoided.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Application

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information about product training, please contact your national B. Braun/Aesculap agency.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

Preparing and handling the implants

Damage to/destruction of the titanium clamps due to use an incompatible applicator!

- Use CRANIOFIX2 applicator (FF494R) only for CRANIOFIX2 titanium clamps (FF490T/FF490T-UNI–FF492T/FF492T-UNI).
- Use CRANIOFIX applicator (FF102R, FF106R and FF107R) only for CRANIOFIX titanium clamps (FF099T–FF101T).

Prior to the surgical intervention, the surgeon must prepare an operation plan to ensure that there is a sufficient supply of sterile CRANIOFIX2 titanium clamps available, that the operation will be performed in aseptic conditions, that the required CRANIOFIX2 applicator and all other CRANIOFIX instruments are complete and ready for use, and that the surgeon and operating team are familiar with the instruments and the surgical technique. Attention must be paid to an osteoplastic incision for the craniotomy.

The correct handling of the CRANIOFIX2 titanium clamps before and during the operation is an important prerequisite for the success of the bone fixation.

The CRANIOFIX2 titanium clamps are individually packed in protective packaging and sterilized by gamma radiation.

- The CRANIOFIX2 titanium clamps must be kept in their original packaging until immediately before use. Prior to use, check the expiry date and make sure that the sterile packaging is intact.
- Check to make certain that the surfaces of the CRANIOFIX2 titanium clamps are not bent or damaged.
- CRANIOFIX2 titanium clamps are intended for single use only and must not be reused under any circumstances.
- Damaged or surgically removed CRANIOFIX2 titanium clamps must not be reused.
- To avoid spark formation, immerse the CRANIOFIX2 titanium clamps in a sterile liquid prior to application.

DANGER

Risk of injury if the upper disc of the CRANIOFIX2 titanium clamps becomes loose due to incorrect handling!

- ▶ Make certain that at the application site of the CRANIOFIX2 titanium clamps the cranial bone and artificial bone flap are of the same thickness.

DANGER

Risk of injury due to artificial bone flaps fixated with CRANIOFIX2!

- ▶ Fixate artificial bone flap only with CRANIOFIX2 if the artificial bone flap is stable against rotation in the cranial bone defect.
- ▶ Before fixating the artificial bone flap with CRANIOFIX2, make certain that the bone flap is supported on the skull and that it does not sink in without fixation.
- ▶ Do not fixate multipiece artificial bone flap with CRANIOFIX2.
- ▶ Fixate artificial bone flap only with CRANIOFIX2 if the artificial bone flap is made from a fatigue-resistant material.
- ▶ Pay attention to the fixation points in the artificial bone flap for the pins of the CRANIOFIX2 titanium clamps.

Application of the CRANIOFIX2 titanium clamps

The **CRANIOFIX2** titanium clamps Ø 16 mm from Aesculap may only be applied in the burr hole Ø 12 mm (see Fig. C) or used to bridge defects in case of fractures.

The **CRANIOFIX2** titanium clamps Ø 20 mm from Aesculap may only be applied in the burr hole Ø 15 mm (see Fig. C) or used to bridge defects in case of fractures.

The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. For maximum stability of a bone flap approx. 20 cm² in size, it is recommended to apply three **CRANIOFIX2** titanium clamps at equal distance from each other in the craniotomy gap or burr hole. Figure 2 shows a typical application of a craniotomized bone flap.

Brief description of operative procedure:

- ▶ Position the three **CRANIOFIX2** titanium clamps at an equal distance from each other on the bone flap.
- ▶ Slide the lower disc between the dura mater and skull cap.
- ▶ Insert the bone flap in its original position.
- ▶ When fixating the artificial bone flap, position the pins of the **CRANIOFIX2** titanium clamps in the fixation points intended for this purpose in the artificial bone flap.
- ▶ Grasp the first **CRANIOFIX2** titanium clamp with the holding forceps. To ensure a secure fixation and to prevent damage to the pin indentation, the pin of the **CRANIOFIX2** titanium clamp must be held in the recess in the jaw of the holding forceps intended for this purpose.
- ▶ Position the **CRANIOFIX2** applicator over the pin of the **CRANIOFIX2** titanium clamp.
- ▶ Compress the handle of the **CRANIOFIX2** applicator to move the upper disc along the pin indentation, closer to the lower disc. The spring leaves will audibly engage in the pin serrations of the **CRANIOFIX2** titanium clamp.
- ▶ Only remove the holding forceps immediately before closing the **CRANIOFIX2** titanium clamp.
- ▶ The **CRANIOFIX2** applicator has an automatic strain relief mechanism, which ensures that the **CRANIOFIX2** applicator disengages when the maximum allowable application force is reached. After the **CRANIOFIX2** applicator has disengaged, the handle should be pressed down completely several times until no further audible engaging of the spring reeds in the pin indentation of the **CRANIOFIX2** titanium clamp can be heard.

CAUTION

The upper disc can loosen after application if the pin is not cut off properly!

- ▶ After tightening the upper disc, cut the excess pin off at a right angle to the pin's axis using the CRANIOFIX cutting forceps.
- ▶ Be sure that the indentation above the upper disc remains completely intact.

The same procedure is used for the other **CRANIOFIX2** titanium clamps.

- ▶ Remove the remaining excess part of the pin with the cutting forceps. When doing so, the excess pin of the **CRANIOFIX2** titanium clamp must be held onto, to prevent it springing into the operating site in an uncontrolled manner.
- ▶ In the event of intraoperative removal of **CRANIOFIX2** clamps, remove these with **CRANIOFIX2** extraction forceps by loosening the upper disc.

Removing the CRANIOFIX2 titanium clamps in a repeat intervention

In rare cases in which a repeat intervention is necessary, the surgeon can leave the **CRANIOFIX2** titanium clamps in position on the bone flap. Aesculap recommends that the surgeon makes a second incision around the original clamp positions. If this is not possible (e.g. for anatomical reasons), Aesculap recommends a second incision within the original clamp positions. The incision must in no way be made through a **CRANIOFIX2** titanium clamp. Decisions regarding the approach for the repeat intervention fall under the responsibility of the surgeon (see Fig. D).

Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact us directly:

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-0
Fax: +49 7461 95-2600
Website: www.aesculap.de

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Patient ID

The Patient ID enables a written record of important information concerning the implanted **CRANIOFIX2** titanium clamps and the surgical intervention to be kept:

- Patient data
- Hospital data
- Date of provision
- Attending physician and hospital
- Art. no. and batch no. of the **CRANIOFIX2** titanium clamps implanted

To facilitate post-operative radiological examinations, the patient ID should be given to every patient. It can be requested from Aesculap separately. Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact Aesculap directly.

Accessories

Art. no.	Designation
FF490T	CRANIOFIX2 titanium clamps Ø 11 mm
FF490T-UNI	CRANIOFIX2 titanium clamps Ø 11 mm
FF491T	CRANIOFIX2 titanium clamps Ø 16 mm
FF491T-UNI	CRANIOFIX2 titanium clamps Ø 16 mm
FF492T	CRANIOFIX2 titanium clamps Ø 20 mm
FF492T-UNI	CRANIOFIX2 titanium clamps Ø 20 mm
FF494R*	CRANIOFIX2 applicator, cannot be disassembled
FF105R*	CRANIOFIX holding forceps
FF103R*	CRANIOFIX cutting forceps
FF104R*	CRANIOFIX removal forceps
FF094P	CRANIOFIX storage tray

* Regarding the operation, function and reprocessing of this product, follow instructions for use TA011161.

TA-Nr. 011200 09/16 V6 Änd.-Nr. 49759

Aesculap®

CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm

Legend

Fig. A

- CRANIOFIX2 titanium clamp FF490T/, Ø 11 mm
- CRANIOFIX2 titanium clamp FF491T, Ø 16 mm
- CRANIOFIX2 titanium clamp FF492T, Ø 20 mm

Fig. B

- 1 CRANIOFIX2 titanium clamp FF491T, Ø 16 mm
- CRANIOFIX2 titanium clamp FF492T, Ø 20 mm
- 2 CRANIOFIX2 titanium clamp FF490T, Ø 11 mm

Fig. C

- 1 Correct positioning of the CRANIOFIX2 titanium clamps FF491T/FF492T: Burr hole Ø 12/15 mm
- 2 Incorrect positioning of the CRANIOFIX2 titanium clamps FF491T/FF492T: Craniotomy gap

Fig. D

- 1 CRANIOFIX2 titanium clamps
- 2 New craniotomy line
- 3 Old craniotomy line

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Product description

The CRANIOFIX2 titanium clamp system from Aesculap consists of sterile CRANIOFIX2 titanium clamps for single use and reusable neurosurgical instruments for the application of titanium clamps.

Intended use

The CRANIOFIX2 titanium clamp system from Aesculap is used for the fixation of craniotomized bone flaps and fractures to the neurocranium.

The CRANIOFIX2 titanium clamp FF490T, Ø 11 mm, is inserted in the craniotomy gap (see Fig. B) or in the fracture gap.

The CRANIOFIX2 titanium clamp FF491T, Ø 16 mm, is inserted in a burr hole with Ø 12 mm (see Fig. C) or is used to bridge defects in case of fractures.

The CRANIOFIX2 titanium clamp FF492T, Ø 20 mm, is inserted in a burr hole with Ø 15 mm (see Fig. C) or is used to bridge defects in case of fractures.

Materials

The materials used in the implant are listed on the packaging:

- Ti6Al4V acc. to ISO 5832-3

Indications

The CRANIOFIX2 titanium clamp system from Aesculap is used for the fixation of craniotomized bone flaps and fractures to the neurocranium.

Contraindications

The following list shows some, but not all, conceivable contraindications.

- Systemic illnesses and metabolic disorders, hypersensitivity to metals or allergies to the implant materials, inflammations in the region of the implant site, insufficient compliance and co-operation of the patient
- Bone conditions that rule out the application of CRANIOFIX2 titanium clamps
- Artificial cranial bone flaps
- Bone tumors in the area supporting the implant
- Degenerative bone diseases
- Missing dura mater
- Application in the facial skull (viscerocranium) and in the orbita or skull-base region

Side effects and interactions

The following side-effects and/or interactions are referred to as part of the legal obligation to inform. These are procedure-specific, non product-specific and not limited to infections, nerve damage, hematomas and wound healing disturbances. The patient must be informed of these risks and other risks associated with neurosurgery, general surgery and the administration of anesthetics.

The following risks may occur in connection with the use of artificial bone flaps:

- Unstable or sinking artificial bone flap
- Damage to the artificial bone flap
- Surgical complications of a different nature

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or lack of asepsis.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically removed components.
- Implants already used must not be reused.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- Implants produce moderate MRI artifacts.

- Inappropriate fixation of craniotomized bone flaps or any combination with other fixation methods (sutures, mini plates or clamp systems from other manufacturers) can result in the formation of steps, shifts of the calvarium and loosening and breakage of implant components. For maximum stability of a craniotomized bone flap, the use of three CRANIOFIX2 titanium clamps in a triangular formation are recommended. Larger bone flaps or a skull fracture with multiple fragments may require a larger number of CRANIOFIX2 titanium clamps to ensure sufficient stability.
- During application of CRANIOFIX2 titanium clamps, friction between the upper titanium disc and the pin indentation can lead to spark formation. Appropriate safety precautions must be taken to prevent this. Due to the risk of fire, surgical interventions with this product must not be carried out in rooms where flammable anesthetics or other flammable gases, liquids, objects or oxidizing agents are stored.
- After the operation, exceptionally heavy mechanical strains on the fixated bone flap must be avoided.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Application

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information about product training, please contact your national B. Braun/Aesculap agency.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

Preparing and handling the implants

Damage to/destruction of the titanium clamps due to use an incompatible applicator!

- Use CRANIOFIX2 applicator (FF494R) only for CRANIOFIX2 titanium clamps (FF490T–FF492T).
- Use CRANIOFIX applicator (FF102R, FF106R and FF107R) only for CRANIOFIX titanium clamps (FF099T–FF101T).

Prior to the surgical intervention, the surgeon must prepare an operation plan to ensure that there is a sufficient supply of sterile CRANIOFIX2 titanium clamps available, that the operation will be performed in aseptic conditions, that the required CRANIOFIX2 applicator and all other CRANIOFIX instruments are complete and ready for use, and that the surgeon and operating team are familiar with the instruments and the surgical technique. Attention must be paid to an osteoplastic incision for the craniotomy.

The correct handling of the CRANIOFIX2 titanium clamps before and during the operation is an important prerequisite for the success of the bone fixation.

The CRANIOFIX2 titanium clamps are individually packed in protective packaging and sterilized by gamma radiation.

- The CRANIOFIX2 titanium clamps must be kept in their original packaging until immediately before use. Prior to use, check the expiry date and make sure that the sterile packaging is intact.
- Check to make certain that the surfaces of the CRANIOFIX2 titanium clamps are not bent or damaged.
- CRANIOFIX2 titanium clamps are intended for single use only and must not be reused under any circumstances.
- Damaged or surgically removed CRANIOFIX2 titanium clamps must not be reused.
- To avoid spark formation, immerse the CRANIOFIX2 titanium clamps in a sterile liquid prior to application.

DANGER

Risk of injury if the upper disc of the CRANIOFIX2 titanium clamps becomes loose due to incorrect handling!

- Make certain that at the application site of the CRANIOFIX2 titanium clamps the cranial bone and artificial bone flap are of the same thickness.

DANGER

Risk of injury due to artificial bone flaps fixated with CRANIOFIX2!

- Fixate artificial bone flap only with CRANIOFIX2 if the artificial bone flap is stable against rotation in the cranial bone defect.
- Before fixating the artificial bone flap with CRANIOFIX2, make certain that the bone flap is supported on the skull and that it does not sink in without fixation.
- Do not fixate multipiece artificial bone flap with CRANIOFIX2.
- Fixate artificial bone flap only with CRANIOFIX2 if the artificial bone flap is made from a fatigue-resistant material.
- Pay attention to the fixation points in the artificial bone flap for the pins of the CRANIOFIX2 titanium clamps.

Application of the CRANIOFIX2 titanium clamps

The **CRANIOFIX2** titanium clamps Ø 16 mm from Aesculap may only be applied in the burr hole Ø 12 mm (see Fig. C) or used to bridge defects in case of fractures.

The **CRANIOFIX2** titanium clamps Ø 20 mm from Aesculap may only be applied in the burr hole Ø 15 mm (see Fig. C) or used to bridge defects in case of fractures.

The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. For maximum stability of a bone flap approx. 20 cm² in size, it is recommended to apply three **CRANIOFIX2** titanium clamps at equal distance from each other in the craniotomy gap or burr hole. Figure 2 shows a typical application of a craniotomized bone flap.

Brief description of operative procedure:

- Position the three **CRANIOFIX2** titanium clamps at an equal distance from each other on the bone flap.
- Slide the lower disc between the dura mater and skull cap.
- Insert the bone flap in its original position.
- When fixating the artificial bone flap, position the pins of the **CRANIOFIX2** titanium clamps in the fixation points intended for this purpose in the artificial bone flap.
- Grasp the first **CRANIOFIX2** titanium clamp with the holding forceps. To ensure a secure fixation and to prevent damage to the pin indentation, the pin of the **CRANIOFIX2** titanium clamp must be held in the recess in the jaw of the holding forceps intended for this purpose.
- Position the **CRANIOFIX2** applicator over the pin of the **CRANIOFIX2** titanium clamp.
- Compress the handle of the **CRANIOFIX2** applicator to move the upper disc along the pin indentation, closer to the lower disc. The spring leaves will audibly engage in the pin serrations of the **CRANIOFIX2** titanium clamp.
- Only remove the holding forceps immediately before closing the **CRANIOFIX2** titanium clamp.
- The **CRANIOFIX2** applicator has an automatic strain relief mechanism, which ensures that the **CRANIOFIX2** applicator disengages when the maximum allowable application force is reached. After the **CRANIOFIX2** applicator has disengaged, the handle should be pressed down completely several times until no further audible engaging of the spring reeds in the pin indentation of the **CRANIOFIX2** titanium clamp can be heard.

CAUTION

The upper disc can loosen after application if the pin is not cut off properly!

- After tightening the upper disc, cut the excess pin off at a right angle to the pin's axis using the CRANIOFIX cutting forceps.
- Be sure that the indentation above the upper disc remains completely intact.

The same procedure is used for the other **CRANIOFIX2** titanium clamps.

- Remove the remaining excess part of the pin with the cutting forceps. When doing so, the excess pin of the **CRANIOFIX2** titanium clamp must be held onto, to prevent it springing into the operating site in an uncontrolled manner.
- In the event of intraoperative removal of **CRANIOFIX2** clamps, remove these with **CRANIOFIX2** extraction forceps by loosening the upper disc.

Removing the CRANIOFIX2 titanium clamps in a repeat intervention

In rare cases in which a repeat intervention is necessary, the surgeon can leave the **CRANIOFIX2** titanium clamps in position on the bone flap. Aesculap recommends that the surgeon makes a second incision around the original clamp positions. If this is not possible (e.g. for anatomical reasons), Aesculap recommends a second incision within the original clamp positions. The incision must in no way be made through a **CRANIOFIX2** titanium clamp. Decisions regarding the approach for the repeat intervention fall under the responsibility of the surgeon (see Fig. D).

Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact us directly:

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-0
Fax: +49 7461 95-2600
Website: www.aesculap.de

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Patient ID

The Patient ID enables a written record of important information concerning the implanted **CRANIOFIX2** titanium clamps and the surgical intervention to be kept:

- Patient data
- Hospital data
- Date of provision
- Attending physician and hospital
- Art. no. and batch no. of the **CRANIOFIX2** titanium clamps implanted

To facilitate post-operative radiological examinations, the patient ID should be given to every patient. It can be requested from Aesculap separately. Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact Aesculap directly.

Accessories

Art. no.	Designation
FF490T	CRANIOFIX2 titanium clamps Ø 11 mm
FF491T	CRANIOFIX2 titanium clamps Ø 16 mm
FF492T	CRANIOFIX2 titanium clamps Ø 20 mm
FF494R*	CRANIOFIX2 applicator, cannot be disassembled
FF105R*	CRANIOFIX holding forceps
FF103R*	CRANIOFIX cutting forceps
FF104R*	CRANIOFIX removal forceps
FF094P	CRANIOFIX storage tray

* Regarding the operation, function and reprocessing of this product, follow instructions for use TA011161.

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