GB

Instructions for use/Technical description
Ennovate Spinal System Implants

USA

Instructions for use/Technical description
Ennovate Spinal System Implants
Note for U.S. users
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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
TA-Nr. 014887 2017-07 V6 Änd.-Nr. 56458

General advisory note

This document describes general prerequisites for use of the named product. For procedural instructions, please refer to the operating manual. If using additional components, please also refer to their respective instructions for use.

Intended use

The Ennovate Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar, thoracic and sacral spine.

Special instruments must be used for implanting these components, as well as for the distraction, compression and reduction of the thoracic and sacral spine.

Indications

Use for severe disorders of the spine that cannot be treated by conservative therapies:

- Fractures
- Pseudarthrosis or delayed healing
- Dislocation
- Degenerative instability
- Post-discectomy syndrome
- Spondylolisthesis
- Kyphosis
- Lordosis
- Scoliosis
- Stenosis
- Post-trauma instability
- Tumors
- Osteoporosis (with the simultaneous use of bone cement and fenestrated screws)

Absolute contraindications

Do not use in the presence of:

- Severe damage to the bone structures of the spine that could prevent the stable implantation of the implant components; for example, osteopenia, severe osteoporosis, Paget's disease, bone tumors etc.
- Metabolic or degenerative metabolic bone diseases that could compromise the stable anchoring of the implant system
- Suspected allergy or sensitivity to the implant materials
- Acute or chronic vertebral infections of a local or systemic nature
- Poor patient compliance or limited ability to follow medical instructions, particularly in the post-op phase, including with regard to the restrictions on range of movement in terms of physical exercise and occupational activity
- Cases not listed under indications

Relative contraindications

In the following circumstances, use of the implant system could represent an increased clinical risk and therefore requires precise, individual assessment by the surgeon:

- Medical or surgical conditions that could negatively impact the success of the implantation, including wound healing disorders
- Conditions that could subject the spine and implants to excessive pressure; for example, pregnancy, obesity, neuromuscular diseases or disorders
- Generally poor condition of the patient; for example, drug or alcohol addition

Side effects and interactions

The general risks of surgery are assumed known and are therefore not described in these instructions for use.

Potential risks and clinical consequences associated with use of the implant system include, but are not limited to:

- Implant failure resulting from excessive load:
 - Warping or bending
 - Loosening
 - Breakage
- Poor stabilization, requiring subsequent revision surgery or surgery to remove the implant system
- Failed or delayed fusion
- Infection
- Fractured vertebral body or bodies
- Injury of adjacent structures, in particular:
 - Nerve roots
 - Spinal cord
 - Blood vessels
 - Organs
- Changes in the physiological curvature of the spine
- Neurological disorders, including impairment of the gastrointestinal, urological and/or reproductive system
- Pain or indisposition
- Bursitis
- Decreased bone density due to load avoidance
- Bone atrophy/fracture above or below the treated spine section
- Limited physical ability
- Persistence of symptoms that were to be treated by the implantation
- Paresis
- Subsequent surgery to remove the implant system
- Metal/foreign body allergy or hypersensitivity
- Dural tears
- Meningitis

Application

Documentation

The operating surgeon draws up an operation plan and documentation covering:

- The selection and dimensioning of the implant components (including article, lot and, if applicable, serial numbers)
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

Patient information

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The patient must be informed about possible side effects and interactions of the product with the body, and relevant rules of behavior.
- The risks of not observing these rules must be explained to the patient.
- The patient must be informed and educated about the aspects and consequences described in the Side effects and interactions section, and the behaviors associated with these.
- The patient must undergo regular medical follow-up examinations.
- Aesculap is not liable for damage in connection with non-observance of the physician's instructions by the patient.
- The patient must be advised that smoking presents an increased risk of pseudarthrosis and diffuse disc degeneration. The progressive degeneration of adjacent segments as a result of smoking could lead to later clinical failure.
- The patient must be advised that being overweight places such excessive pressure on the implant that the prospect of failure becomes likely.
- The patient must be advised that strenuous physical activity in the form of heavy lifting, excessive muscle strain, exaggerated twisting, repeated bending and bending over or running should be avoided until the bones have healed fully.
- The patient must be instructed to keep pressure off the implant post surgery in order to prevent the onset of clinical problems associated with implant failure.

Implantation

Inserting the Ennovate Spinal System polyaxial/monoaxial screws and rods

Please refer to the operating manual for the required information about implanting the implant system. Please note the following in particular:

- All requisite implant components must be ready to hand.
- All the implantation instruments, including the special Aesculap implant system instruments, must be complete and in working order.

WARNING

Incorrect application can result in trauma to the spinal cord and nerve roots!

Severe or fatal hemorrhage as a result of vascular injury!

- ▶ Align and screw in the instrument and, if applicable, pedicle screw under X-ray control or with the aid of a navigation system.
- ▶ Take care not to breach the anterior and lateral cortices when screwing in.
- ▶ Ensure that the great vessels are not eroded, punctured during implantation or damaged as a result of implant breakage or migration after implantation.
- ▶ Take care to avoid pulsatile erosion of the vessels because of close apposition of the implants.

- ▶ Only allow axial load on the screws during screw insertion. Avoid induction of lateral forces during screw insertion.
- ▶ Do not bend metal implants, except for the rods of the Ennovate Spinal Systems.
- ▶ Do not bend the rods back or bend them excessively.
- ▶ Only use the bending instruments of the Ennovate Spinal System instrument set for bending the rods.

WARNING

Inadequate fixation due to incorrect rod position!

- ▶ Always position rods with hexagon and tip in such a way that the hexagon or tip is outside the implant head.
- ▶ Adjust the curvature of the rod to the anatomical situation.
- ▶ Do not bend the rod in the region of the screw head.
- ▶ Make certain and check that the screw head is at right angles to the rod.
- ▶ Make certain and check that the rod is positioned securely in the groove base.

- ▶ To position the set screw, always use the dedicated set screw driver.
- ▶ For tightening and loosening the set screw, always use the dedicated screw driver and counter torque.
- ▶ Always tighten the set screws with the dedicated torque wrench handle.

WARNING

Risk of injury due to improper mounting of the set screw!

- ▶ Position the set screw correctly.
- ▶ Make certain that the rods are correctly positioned on the floor of the groove.
- ▶ Using the torque wrench handle, tighten the set screw until the torque wrench handle release is audibly triggered.

WARNING

Implant damage due to overtightening the set screw or insufficient fixation due to undertightening the set screw!

- ▶ Check the recalibration/maintenance date on the torque wrench handle maintenance label prior to each use.
- ▶ When tightening or loosening the set screws, make certain that the working end of the torque wrench handle shaft 10 Nm is positioned entirely within the socket of the set screws.
- ▶ Apply the torque slowly and without sudden movement until the torque wrench handle release is audibly triggered.
- ▶ Use a rigid screw driver to loosen the screw connection.

WARNING

Implant head expansion caused by tightening the set screw without application of counter torque!

- ▶ Always use the counter torque when tightening the set screw.

- ▶ If the screw head comes loose during revision of the larger polyaxial screws (7.5 mm or larger), remove the remaining screw component with the set screw driver.

Ennovate Spinal System Inserting the rod-to-rod connectors/cross connectors

Follow the steps below to implant the Ennovate Spinal System rod-to-rod connectors/cross connectors:

- ▶ Using the dedicated cross connector holder or rod-to-rod connector inserter, place the rod-to-rod connectors/cross connectors on the rod at the site calculated by the operating surgeon.

CAUTION

Rod-to-rod connector/cross connector is unusable if set screws are removed or not present!

- ▶ Do not unscrew the set screws from the rod-to-rod connector/cross connector.
- ▶ Do not swap set screws with other rod-to-rod connector/cross connector components.
- ▶ Only use rod-to-rod connectors/cross connectors if all the set screws are present.

- ▶ Before tightening the set screws, ensure that the rod-to-rod connector/cross connector is positioned correctly.

WARNING

Inadequate fixation due to incorrect rod position!

- ▶ Always position the rod-to-rod connector/cross connector in such a way that the rod is completely inserted into the rod-to-rod connector/cross connector.
- ▶ Always position rods in such a way that the end of the hexagon and tip of the rod are located outside the clamping area of the rod-to-rod connector/cross connector.

- ▶ Always use the dedicated torque wrench shaft 5 Nm in combination with the dedicated counter torque to tighten and loosen the set screw.

WARNING

Risk of injury if rod-to-rod connector is positioned incorrectly!

- ▶ Always follow the correct sequence of steps when tightening the set screw.
- ▶ Check correct positioning.

- ▶ Always tighten the set screws with the dedicated torque wrench handle 5 Nm.
- ▶ With the rod-to-rod connectors open, use the dedicated torque wrench handle 5 Nm to first tighten the set screw marked 1 and then the set screw marked 2.

WARNING

Implant damage due to overtightening the set screw or insufficient fixation due to undertightening the set screw!

- ▶ Check the recalibration/maintenance date on the torque wrench handle maintenance label prior to each use.
- ▶ When tightening or loosening the set screws, make certain that the working end of the torque wrench handle shaft 5 Nm is positioned entirely within the socket of the set screws.
- ▶ Apply the torque slowly and without sudden movement until the torque wrench handle release is audibly triggered.
- ▶ Tighten the set screws with the torque wrench handle 5 Nm.
- ▶ Position the rod-to-rod/cross connector correctly.
- ▶ Make certain that the rods are completely inserted into the rod-to-rod/cross connector.
- ▶ Make certain that there is no soft tissue caught in the rod-to-rod/cross connector.
- ▶ Use a rigid screw driver to loosen the screw connection.

- ▶ When tightening or loosening, only subject the set screws to axial pressure. Avoid exerting lateral forces when tightening or loosening the set screws.

Ennovate Spinal System fixing the fenestrated screws with bone cement

To increase anchoring stability, the fenestrated Ennovate screw can be fixed with bone cement. The following must be observed when cementing:

Note

Do not start cementing until all the screws have been aligned and the rods have been adapted; however, cementing should be completed before locking the fenestrated Ennovate screw!

WARNING

Implant damage caused by the application of excessive shear forces with the injection cannula!

- ▶ Do not subject the injection cannula to transverse loads.
- ▶ Do not bend the injection cannula.

Note
Once the fenestrated Ennovate screw has been cemented, it is no longer possible to change the implant position, use the reduction lever, or make corrections!

Note
After the fenestrated Ennovate screw has been cemented, explantation of the screw can prove difficult, especially in cases of osteopenia.

- ▶ Prior to applying the Ennovate bone cement, fully insert and screw the injection cannula into the screw head without exerting appreciable force.
- ▶ Check the bone cement for proper texture and viscosity before positioning the applicator on the injection cannula. Observe the bone cement manufacturer's instructions!
- ▶ Only cement the fenestrated Ennovate screw with a cement suitable for screw augmentation.

WARNING

Risk of injury from cement leakage!

- ▶ Always apply imaging techniques for control during the preparation and positioning of the implants.
- ▶ Make certain that the slots of the fenestrated Ennovate screw are positioned fully inside the bone.
- ▶ Take care not to perforate the pedicle cortex.
- ▶ Do not penetrate the opposite cortex and ensure that the opposite cortex is undamaged.
- ▶ Do not insert or cement the fenestrated Ennovate screw in the fractured or partially resected vertebral body. Carry out preoperative thin-section CT control.
- ▶ Ennovate Insert and screw the injection cannula into the screw head without exerting appreciable force. Make certain that the axes of the pedicle screw and the injection cannula are aligned with each other.
- ▶ Cement application must always be accompanied by high-quality, image-producing X-ray screening. Use no more than 2 ml of cement for each pedicle screw.

WARNING

Risk of injury due to material incompatibility!

- ▶ Only use application systems that have been approved for use with bone cement.

WARNING

Contamination of the locking thread/insufficient fixation when removing the application system due to incorrect viscosity of the cement!

- ▶ Make certain the cement has the correct viscosity.

WARNING

Wound contamination due to insufficient curing of cement! Risk of set screw thread adhesion, resulting in system instability!

- ▶ Observe the correct curing time.
- ▶ For monitoring the cement texture, prepare a reference volume of cement together with the cement to be applied.
- ▶ Make certain the implant system is not exposed to any load before the cement has reached its final strength.

▶ Wait for the cement to cure completely before removing the injection cannula.
Further information about Aesculap implant systems is available at any time from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Safety advice

Users
Appropriate clinical training as well as a theoretical and practical command of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.
Aesculap is not responsible for complications caused by:

- Incorrect indication or implant selection
- Incorrect combination of implant components, particularly combination with components of other manufacturers
- Incorrect surgical technique
- Exceeding the limitations of the treatment method or non-observance of essential medical precautions

The user has a duty to consult with the manufacturer if the preoperative situation regarding use of the implants is ambiguous:

- Implants that have been used before must not be reused.
- When shaping the rod, the operating surgeon must avoid scoring or scratching the product and bending it back and forth. Bending the screws will significantly shorten their service life and could cause premature implant failure.
- The attending physician shall make any decision with regard to the removal of implant components that have been used.
- Postoperatively, individual patient information, as well as moderate mobility and muscle training is of particular importance.

Product
The following warnings, precautions and side effects must be understood by the surgeon and explained to the patient:

- Implants that are damaged or have already been used once must not be used.
- Do not combine stainless steel implant components with implant components made of titanium as this accelerates the corrosion process in stainless steel, which can lead to implant failure as a result of material fatigue.
- The implant components have been tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Delayed healing can cause implant breakage due to metal fatigue.
- Implant components can still become bent, broken or loose despite complete bone hardening.
- Observe the instructions for use of the individual Aesculap implant components and instruments and the valid operating manual.

CAUTION

Risk of implant component breakage due to combination with implant components from other manufacturers!

- ▶ Only use Aesculap implant components from the Ennovate System.

WARNING

Risk of screw breakage caused by pedicle screws used for spondylolisthesis!

- ▶ Support the stabilization measure with an intercorporeal fusion (ALIF, PLIF or TLIF).

WARNING

Risk of more complicated revision!

- ▶ When using bone cement, use only the fenestrated screws intended for this purpose.
- ▶ There is an increased risk of more complicated revision if using fenestrated pedicle screws with a screw diameter greater than 8.5 mm.

- Damage to the load-bearing structures of the implant can lead to loosening of components, dislocation, migration, and other severe complications.
- To ensure successful application, the operating surgeon must undergo appropriate product-specific training prior to initial use.
- In case of questions relating to the indication or product, further information is available at any time from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

Information about durability and reusability

Sterility

Sterile implant components

- Implant components with article numbers ending in **TS** (e.g. SY510**TS**) are individually packed in labeled protective packaging.
- Implant components with article numbers ending in **TS** are gamma sterilized.
- ▶ Store implant components with article numbers ending in **TS** in their original packaging. Do not remove them from their original and protective packaging until immediately before implantation.
- ▶ Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- ▶ Do not use implant components that are past their expiration date or whose packaging is damaged.

WARNING

Damage to implants with article numbers ending in TS caused by reprocessing and resterilization!

- ▶ Do not reprocess or resterilize implants with article numbers ending in **TS**.

Non-sterile implant components

- Implant components with article numbers ending in **T** (e.g. SY510**T**) are supplied in an unsterile condition.
- Implant components with article numbers ending in **T** are packed individually.
- Implants delivered in an unsterile condition are delivered pure.
- Follow the sterilization procedure recommended in ISO 8828 for hospitals for all components.
- ▶ Store implant components with article numbers ending in **T** in their original packaging. Do not remove them from their original and protective packaging until immediately before reprocessing.
- ▶ Use suitable storage equipment for reprocessing, sterilization and sterile supply.
- ▶ Use a suitable tray for cleaning/disinfection, see TA015247.
- ▶ Use the system storage device for sterilization and sterile provision.
- ▶ Ensure that the implant components in their implant system storage devices do not come into contact with each other or with instruments.
- ▶ Ensure that the implant components are not damaged in any way.

Prior to initial sterilization and subsequent resterilization, the implant components must be cleaned using the following validated reprocessing procedure:

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.
The recommended chemical was used for validation.

Note
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Validated reprocessing procedure

WARNING

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- ▶ Handle the implants with new gloves only.
- ▶ Keep the implant system storage devices covered or closed.
- ▶ Process implant system storage devices separately from instrument trays.
- ▶ Clean implants must not be processed together with contaminated implants.
- ▶ Process the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.
- ▶ Do not reuse surgically contaminated implants!

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- ▶ Do not reprocess implants that have been directly or indirectly contaminated with blood.

Manual cleaning/disinfection

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry
I	Ultrasound cleaning	RT (cold)	5	2	D-W	Concentrate, aldehyde-free, phenol-free, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Concentrate, aldehyde-free, phenol-free, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

DW: Drinking water
FD-W: Fully desalinated water (demineralized, microbiological, at least of drinking water quality)
RT: Room temperature
*Recommended: B Braun Stabimed

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for 5 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Do not clean with brushes or other abrasives that would damage the product surface.

Phase II

- ▶ Rinse all accessible product surfaces under running water.
- ▶ Drain any remaining water fully.

Phase III

- ▶ Fully immerse the instrument in the disinfectant solution.

Phase IV

- ▶ Rinse all accessible product surfaces under running water.
- ▶ Drain any remaining water fully.

Phase V

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air).

For implant components with article numbers ending in **T**, which are resterilized, please note the following:

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfecting machine must be tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Ultrasonic cleaning	RT (cold)	5	2	D-W	Concentrate, aldehyde-free, phenol-free, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

DW: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Phase I
- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for 5 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - ▶ Do not clean with brushes or other abrasives that would damage the product surface.

- Phase II
- ▶ Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	D [°C/°F]	t [min]	Water quality	Chemistry
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH ~ 13 - <5 % anionic surfactant working solution 0.5% - pH ~ 11*
III	Intermediate rinse	> 10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

DW: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Immediately set aside damaged or inoperative products.
- ▶ Place the product in its holder or on a suitable tray.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Make sure that the packaging will prevent a recontamination of the product during storage.

Sterilization method and parameters

- ▶ Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Product properties

The materials used in the implant are listed on the packaging:

- ISOTAN®: titanium forged alloy Ti6Al4V acc. to ISO5832-3
- ISODUR®: wrought cobalt-chromium alloy CoCrMo according to ISO 5832-12

The titanium implants are anodized with a colored oxide layer. Slight changes in coloration may occur, but do not affect the implant quality. ISOTAN® and ISODUR® are trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Warnings

Warnings in this document that highlight risks to the patient, user and/or product are displayed as follows:

WARNING

Danger of death or serious injury.

CAUTION

Danger of minor injury or product damage.

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