GB Instructions for use/Technical description

USA Notes on the electromagnetic compatibility (EMC) for Acculan 3Ti charger GA677

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Aesculap®

Notes on the electromagnetic compatibility (EMC) for Acculan 3Ti charger GA677

Contents

1.	Scope.....	0
2.	Guidance and manufacturer's declaration.....	0
2.1	Classification of medical electrical equipment acc. to CISPR 11.....	0
2.1.1	Groups.....	0
2.1.2	Classes.....	0
2.1.3	Set-up environment/set-up location.....	0
2.2	Electromagnetic emission.....	1
2.3	Electromagnetic immunity.....	2
2.4	High frequency conducted disturbance - IEC 61000-4-3.....	5
2.5	Separation distances.....	6
3.	Technical Service.....	6
4.	Distributor in the US/Contact in Canada for product information and complaints.....	6

Group 2

Not relevant for Aesculap products.

2.1.2 Classes

The class to which medical electrical equipment belongs depends on the electromagnetic environment at the site where the equipment is operated.

Class A

Devices are of **Class A** if they are identified in the instructions for use as "CISPR Class A".

Class B

Devices are of **Class B** if no particular classification is specified in the instructions for use, or if they are identified there as "CISPR Class B".

2.1.3 Set-up environment/set-up location

The equipment is intended for use in a specific environment. It may be used in areas in which medical treatments are performed with HF surgical devices.

- ▶ Do not put device directly next to other devices or stack them with other devices, except with or on Aesculap devices.
- ▶ When operation directly next to other devices or must be stacked with other devices: observe the device to ensure proper operation in this configuration.

1. Scope

This information on electromagnetic compatibility (EMC) applies to the following medical electrical device:

Art. no.	Designation
GA677	Acculan 3Ti charger

2. Guidance and manufacturer's declaration

2.1 Classification of medical electrical equipment acc. to CISPR 11

According to CISPR 11, all medical electrical equipment is divided into groups and classes. Different guidelines concerning electromagnetic emissions apply, depending on which group or class the equipment belongs to.

2.1.1 Groups

The group any medical electrical equipment belongs to depends on the way in which it uses RF energy.

Group 1

Group 1 is the group of "equipment and systems which generate or use RF energy exclusively for their internal functions", such as camera systems, light sources and motor systems. Group 1 also comprises "equipment that feed energy into the patient, but not energy in the form of high-frequency electromagnetic energy", e.g. infusion equipment.

RF surgery devices also belong to group 1 since, according to IEC 60601-2-2, they are assessed with the RF output inactive.

2.2 Electromagnetic emission

- ▶ Always observe guidelines according to the group and class the equipment belongs to.
- ▶ For details regarding the meaning of different groups and classes, see Classification of medical electrical equipment acc. to CISPR 11 in the previous section.

Guidance and manufacturer's declaration - Electromagnetic emission

The Acculan 3Ti charger is intended for use in the electromagnetic environment specified below. The user of the equipment should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - Guidance
RF emissions CISPR 11	Complies with group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Harmonic emissions acc. to IEC 61000-3-2	Consistent with Class A (IEC 61000-3-2)	-
Voltage fluctuations/flicker emissions acc. to IEC 61000-3-3	Complies	-
RF emissions acc. to CISPR 11	Class B	The device is suitable for use in all facilities including those in a residential area, which are directly connected to the public power supply that also supplies buildings used for residential purposes.

Notes on the electromagnetic compatibility (EMC) for Acculan 3Ti charger GA677

2.3 Electromagnetic immunity

Guidelines and Manufacturer Explanation – Electromagnetic Compatibility – IEC 60601-1-2, 3rd edition

The Acculan 3Ti charger is intended for use in the electromagnetic environment specified below. The user of the Acculan 3Ti charger should ensure that it is used in such an environment.

Immunity test	Compatibility test level IEC 60601-1-2, 3rd Edition	Compliance level	Electromagnetic environment – Guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV ±8 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Rapid transient electrically conducted disturbance acc. to IEC 61000-4-4	±2 kV for power supply lines ±1 kV for power supply lines	±2 kV ±1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge acc. to IEC 61000-4-5	±1 kV voltage external conductor - external conductor ±2 kV voltage outer conductor- protective earth	±1 kV ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines acc. to IEC 61000-4-11	<5 % U_T^1 (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Complies	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Acculan 3Ti charger requires continued operation during power mains interruptions, it is recommended that the Acculan 3Ti charger be powered from an uninterruptible power supply or a battery.
Magnetic field at power supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note 1 U_T is the AC mains voltage prior to application of the test level..

Guidelines and Manufacturer Explanation – Electromagnetic Compatibility – IEC 60601-1-2, 4th edition

The Acculan 3Ti charger is intended for use in the electromagnetic environment specified below. The user of the Acculan 3Ti charger should ensure that it is used in such an environment.

Immunity test	Compatibility test level IEC 60601-1-2, 4th Edition	Compliance level	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV ±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst acc. to IEC 61000-4-4	±2 kV for power supply lines ±1 kV for power supply lines Modulation 100 kHz	±2 kV ±1 kV Modulation 100 kHz	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	±1 kV voltage outer conductor- outer conductor ±2 kV voltage external conductor - earth	±1 kV ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short term interruptions and fluctuations in power supply acc. to IEC 61000-4-11	0 % U_T^1 at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° for 1/2 periods 0 % U_T ; one period 70 % U_T 25/30 periods Single phase: at 0° 0 % U_T ; 250/300 periods	Complies	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Acculan 3Ti charger requires continued operation during power mains interruptions, it is recommended that the Acculan 3Ti charger be powered from an uninterruptible power supply or a battery.
Magnetic field at power supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note 1 U_T is the AC mains voltage prior to application of the test level.

Aesculap®

Notes on the electromagnetic compatibility (EMC) for Acculan 3Ti charger GA677

Guidance and manufacturer's declaration – Electromagnetic immunity – For equipment that are not life supporting

The Acculan 3Ti charger is intended for use in the electromagnetic environment specified below. The user of the Acculan 3Ti charger should ensure that the Acculan 3Ti charger is used in such an environment.

Immunity test	Test level IEC 60601-1-2	Compliance level	Electromagnetic environment – Guidance ²
Conducted RF acc. to IEC 61000-4-6	3 V _{eff} 150 kHz to 80 MHz	3 V _{eff} 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the Acculan 3Ti charger, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance¹: $d = 1,2 \times \sqrt{P}$
	6 V _{eff} 80 % AM at 1 kHz for ISM bands ^c	6 V _{eff} 80 % AM at 1 kHz for ISM bands ^c	
Transmitted RF acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	$d = 1,2 \times \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \times \sqrt{P}$ for 800 MHz to 2.7 GHz

Where **P** is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and **d** is the recommended separation distance in meters [m].

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range^b.

Interference may occur in the vicinity of equipment marked with the following symbol.

Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Acculan 3Ti charger is used exceeds the applicable RF compliance level above, the Acculan 3Ti charger should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Acculan 3Ti charger.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

c The ISM bands (en: frequency bands used for Industrial, Scientific and Medical purposes) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz, 13.553 MHz to 13.567 MHz, 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz. The amateur radio frequency bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

2.4 High frequency conducted disturbance - IEC 61000-4-3

Test frequency range (MHz)	Band (MHz)	Service	Modulation	Max. output (W)	Distance (m)	Conducted disturbance test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GRMS 460, FRS 460	FM ±5 kHz deviation	2	0.3	28
710	704-787	LTE-Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810						
870	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE-Band 5	Pulse modulation 18 Hz	2	0.3	28
930						
1 720						
1 845	1 700-1 990	GSM 1 800, CDMA 1 900, GSM 1 900, DECT, LTE- Band 1, 3, 4, 25, UMTS	Pulse modulation 217 Hz	2	0.3	28
1 970						
2 450	2 400-2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2 4500, LTE- Band 7	Pulse modulation 217 Hz	2	0.3	28
5 240	5 100-5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5 500						
5 785						

Aesculap®

Notes on the electromagnetic compatibility (EMC) for Acculan 3Ti charger GA677

2.5 Separation distances

Recommended separation distances² between portable and mobile RF communications equipment and the Acculan 3Ti charger.

The Acculan 3Ti charger is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the Acculan 3Ti charger can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Acculan 3Ti charger as recommended below, according to the maximum output power of the communications equipment.

PORTABLE HF communication devices (radio devices), including their ACCESSORIES (such as antenna cables and external antennas) should not be used at a distance of less than 30 cm from the lines of the Acculan 3Ti charger. Non-observance may lead to a reduction of the device performance features.

Rated maximum output of transmitter [W]	Separation distance ¹ according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1,2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \times \sqrt{P}$	800 MHz to 2.7 GHz $d = 2,3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters [m] can be determined using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

Note 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

3. Technical Service

► For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95 -1601
Fax: +49 (7461) 14 -939
E-Mail: ats@aesculap.de

Or in the US:

Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042 USA
Aesculap Repair Hotline
Phone: +1 (800) 214 -3392
Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

4. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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