**GB** Instructions for use/Technical description  
**USA** Notes on electromagnetic compatibility (EMC) for Acculan components

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### Contents

1.	Scope.....	0
2.	Guidance and manufacturer's declaration.....	0
2.1	Classification of medical electrical equipment acc. to CISPR 11.....	0
2.1.1	Groups.....	0
2.1.2	Classes.....	1
2.1.3	Set-up environment/set-up location.....	1
2.2	Electromagnetic emission.....	1
2.3	Electromagnetic immunity.....	2
2.4	High frequency conducted disturbance - IEC 61000-4-3.....	4
2.5	Separation distances.....	5
3.	Technical Service.....	5
4.	Distributor in the US/Contact in Canada for product information and complaints.....	5

## 1. Scope

These instructions on electromagnetic compatibility (EMC) apply to the following medical electrical components of the Acculan systems.

### Acculan 3Ti System

Art. no.	Designation
GA668	Oscillating saw, straight
GA670	Dermatome
GA671	Small drill
GA672	Drill and reamer
GA673	Oscillating saw
GA674	Reciprocating saw
GT668	TPLO saw

### Acculan 4 System

Art. no.	Designation
GA330	Drill and reamer
GA331	Oscillating saw
GA334	Reciprocating Saw attachment
GA336	Oscillating saw, straight
GA337	Oscillating Sternum saw
GA340	Dermatome
GA341	Dermatome 0.1 mm
GA344	Small drill
GT666	TPLO saw

### Batteries

Art. no.	Designation
GA346	Li-Ion battery short
GA666	NiMH battery short
GA676	NiMH battery

#### Note

The use of cables other than those provided by or recommended by Aesculap may lead to an increased emission or reduced interference immunity.

## 2. Guidance and manufacturer's declaration

### 2.1 Classification of medical electrical equipment acc. to CISPR 11

According to CISPR 11, all medical electrical equipment is divided into groups and classes. Different guidelines concerning electromagnetic emissions apply, depending on which group or class the equipment belongs to.

#### 2.1.1 Groups

The group any medical electrical equipment belongs to depends on the way in which it uses RF energy.

### Group 1

Group 1 is the group of "equipment and systems which generate or use RF energy exclusively for their internal functions", such as camera systems, light sources and motor systems. Group 1 also comprises "equipment that feed energy into the patient, but not energy in the form of high-frequency electromagnetic energy", e.g. infusion equipment.

RF surgery devices also belong to group 1 since, according to IEC 60601-2-2, their emission is only considered with the output switch inactive.

### Group 2

Not relevant for Aesculap products.

#### 2.1.2 Classes

The class to which medical electrical equipment belongs depends on the electromagnetic environment at the site where the equipment is operated.

##### Class A

Devices are of **Class A** if they are identified in the instructions for use as "CISPR Class A".

##### Class B

Devices are of **Class B** if no particular classification is specified in the instructions for use, or if they are identified there as "CISPR Class B".

## 2.2 Electromagnetic emission

- ▶ Always observe guidelines according to the group and class the equipment belongs to.
- ▶ For details regarding the meaning of different groups and classes, see Classification of medical electrical equipment acc. to CISPR 11 in the previous section.

### Guidance and manufacturer's declaration – Electromagnetic emission

The Acculan components are intended for use in the electromagnetic environment specified below. The user of the equipment should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – Guidance
RF emissions CISPR 11	Complies with group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions acc. to CISPR 11	Class B	The EQUIPMENT is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

## Notes on electromagnetic compatibility (EMC) for Acculan components

### 2.3 Electromagnetic immunity

#### Guidelines and Manufacturer Explanation – Electromagnetic Compatibility – IEC 60601-1-2, 3rd edition

The Acculan components are intended for use in the electromagnetic environment specified below. The user of the Acculan components should ensure that it is used in such an environment.

Immunity test	Compatibility test level IEC 60601-1-2, 3rd Edition	Compliance level	Electromagnetic environment – Guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV ±8 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Magnetic field at power supply frequency (50/60 Hz) acc. to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note 1  $U_T$  is the AC mains voltage prior to application of the test level.

#### Guidelines and Manufacturer Explanation – Electromagnetic Compatibility – IEC 60601-1-2, 4th edition

The Acculan components are intended for use in the electromagnetic environment specified below. The user of the Acculan components should ensure that it is used in such an environment.

Note: IEC 60601-1-2, 4th edition is only applicable if explicitly indicated in the product-specific instructions for use.

Immunity test	Compatibility test level IEC 60601-1-2, 4th Edition	Compliance level	Electromagnetic environment – Guidance
Electrostatic Discharge (ESD) acc. to IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV ±15 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Magnetic field at power supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note 1  $U_T$  is the AC mains voltage prior to application of the test level..

**Guidance and manufacturer's declaration – Electromagnetic immunity – For equipment that are not life supporting**

The Acculan components are intended for use in the electromagnetic environment specified below. The user of the Acculan components should ensure that it the Acculan components are used in such an environment.

Immunity test	Test level IEC 60601-1-2	Compliance level	Electromagnetic environment – Guidance <sup>2</sup>
Transmitted RF acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Acculan components, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance<sup>1</sup>:</b></p> <p><b>d = 1,2 x √P</b> 80 MHz to 800 MHz</p> <p><b>d = 2,3 x √P</b> for 800 MHz to 2.7 GHz</p> <p>Where <b>P</b> is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and <b>d</b> is the recommended separation distance in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup> should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol.</p> 
Note 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Acculan components are used exceeds the applicable RF compliance level above, the Acculan components should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Acculan components.		
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		
c	The ISM bands (en: frequency bands used for Industrial, Scientific and Medical purposes) between 0.15 MHz and 80 MHz are 6.765 MHz to 6.795 MHz, 13.553 MHz to 13.567 MHz, 26.957 MHz to 27.283 MHz and 40.66 MHz to 40.70 MHz. The amateur radio frequency bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.		

# Aesculap®

## Notes on electromagnetic compatibility (EMC) for Acculan components

### 2.4 High frequency conducted disturbance - IEC 61000-4-3

Test frequency range (MHz)	Band (MHz)	Service	Modulation	Max. output (W)	Distance (m)	Conducted disturbance test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GRMS 460, FRS 460	FM ±5 kHz deviation	2	0.3	28
710	704-787	LTE-Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE-Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1 720	1 700-1 990	GSM 1 800, CDMA 1 900, GSM 1 900, DECT, LTE- Band 1, 3, 4, 25, UMTS	Pulse modulation 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400-2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2 4500, LTE- Band 7	Pulse modulation 217 Hz	2	0.3	28
5 240	5 100-5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5 500						
5 785						

## 2.5 Separation distances

Recommended separation distances<sup>2</sup> between portable and mobile RF communications equipment and the Acculan components.

The Acculan components intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of the Acculan components can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Acculan components as recommended below, according to the maximum output power of the communications equipment.

Portable HF communication devices (radio devices), including their accessories (such as antenna cables and external antennas) should not be used at a distance of less than 30 cm from the lines of the Acculan components. Non-observance may lead to a reduction of the device performance features.

Rated maximum output of transmitter [W]	Separation distance <sup>1</sup> according to frequency of transmitter [m]		
	150 kHz to 80 MHz $d = 1,2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \times \sqrt{P}$	800 MHz to 2.7 GHz $d = 2,3 \times \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters [m] can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

Note 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## 3. Technical Service

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

### Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 (7461) 95 -1601  
Fax: +49 (7461) 14 -939  
E-Mail: ats@aesculap.de

Or in the US:

Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042 USA  
Aesculap Repair Hotline  
Phone: +1 (800) 214 -3392  
Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

## 4. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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