

B. Braun Biosurgicals

Histoacryl®

Tissue adhesive



Histoacryl® – The real thing

Histoacryl®

The real thing – now even simpler and faster

For decades medical professionals around the world have put their trust into Histoacryl®, the first surgical tissue adhesive based on cyanoacrylate. The successful application of Histoacryl® is described in more than 1.100 publications. Now the classic product for the closure of skin wounds is available

- in a new ampoule with a twist-off tip
- with increased content of n-butyl-2-cyanoacrylate
- with the indication for skin closure and for sclerosation therapy of esophageal and fundal varices

New twist-off tip

Histoacryl®

now even simpler and faster



- Simply twist off the ribbed tip of the ampoule
- Then, as usual, Histoacryl® can be applied precisely with the slim cannula

Translucent Histoacryl® L

for facial applications



To be on the safe side, it is recommended to use Histoacryl® L, a version which does not contain dye, for facial applications.

Classic advantages



Well-known advantages of Histoacryl®:

- fast wound closure
- superior strength (ref. bar chart on page 5)
- simple and precise dosage due to slim cannula and blue dye
- applying just one layer is enough
- choice between blue and translucent version
- showering possible

Why tissue adhesive?


- no suturing – no pain
- saves time and costs (no local anaesthetics, no suture removal, no second consultation)
- antibacterial film protects the wound
- no wound dressing necessary

Histoacryl®

The real thing – now even simpler and faster

Tensile strength of Histoacryl®:

Comparative testing in accordance with ASTM F 2458-05, unpublished data




Histoacryl® for closure of endoscopic incisions

The Histoacryl® ampoule is supplied in sterile condition; therefore, Histoacryl® is the ideal tissue adhesive for use in the OR.

For example, it may be used for the closure of endoscopic incisions, as described in:

Rosin D et al. (2001) Closure of laparoscopic trocar site wounds with cyanoacrylate tissue glue: a simple technical solution. J Laparoendosc Adv Surg Tech 11(3):157-9.




Sclerosation of esophageal and fundal varices with Histoacryl®

The formation of esophageal and fundal varices is a common and dangerous consequence of portal hypertension. Effective sclerosation of these varices is possible with Histoacryl®.

Varically injected Histoacryl® polymerizes intravasally to form a plastic cylinder and effects immediate obturation or thrombosing of the vessel. The procedure is described in detail in the relevant literature, for instance in the following articles:

- Binmoeller K. F., Soehendra, N. (1995)
Superglue: the answer to variceal bleeding and fundal varices? Endoscopy 1995; Vol.27: 392-396.
- Gotlib J.-P. (1990)
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive. Can J Gastroenterol 4 (9): 637-638.



Product Range

Description	Color	Article No.
5 x 0.5 ml Histoacryl®	blue	105 0052
10 x 0.5 ml Histoacryl®	blue	105 0044
5 x 0.2 ml Histoacryl®	blue	105 0036
10 x 0.2 ml Histoacryl®	blue	105 0028
5 x 0.5 ml Histoacryl® L	translucent	105 0060

■ Literature is available upon request.


Wound closure in a minute:




Clean the wound and turn open the ampoule



Apply Histoacryl® sparingly, do not allow tissue adhesive to flow into the wound



Adapt wound edges (wounds under tension and wounds longer than 3 cm require additional suturing)



Keep wound edges aligned for about one minute – that's all!



AESCULAP®

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen
Germany

Phone +49 7461 95-0
Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.