- GB** **Instructions for use/Technical description**
BiCONTACT® hip endoprosthesis stem
- USA** **Instructions for use/Technical description**
BiCONTACT® hip endoprosthesis stem
Note for U.S. users
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
- D** **Gebruiksaanweisung/Technische Beschreibung**
BiCONTACT® Hüftendoprothesenschaft
- F** **Mode d'emploi/Description technique**
Tige d'endoprothèse de hanche BiCONTACT®
- E** **Instrucciones de manejo/Descripción técnica**
Vástago de endoprótesis de cadera BiCONTACT®
- I** **Istruzioni per l'uso/Descrizione tecnica**
Stelo per endoprotesi coxo-femorale BiCONTACT®
- P** **Instruções de utilização/Descrição técnica**
Haste de endoprótese coxo-femoral BiCONTACT®
- NL** **Gebruiksaanwijzing/Technische beschrijving**
BiCONTACT®-heupprotheseschacht
- S** **Bruksanvisning/Teknisk beskrivning**
BiCONTACT®-höftendoproteschkäft
- RUS** **Инструкция по применению/Техническое описание**
Ножка эндопротеза тазобедренного сустава BiCONTACT®
- CZ** **Návod k použití/Technický popis**
Dřík kyčelní endoprotézy BiCONTACT®
- PL** **Instrukcja użytkowania/Opis techniczny**
Trzpień endoprotezy stawu biodrowego BiCONTACT®
- SK** **Návod na použitie/Technický opis**
BiCONTACT® telo bedrovej endoprotézy
- TR** **Kullanım Kılavuzu/Teknik açıklama**
BiCONTACT® kalça endroportez şaftları
- KR** **사용 설명서 / 기술 설명**
BiCONTACT® 인공 고관절 스템

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 010238 02/17 V6 Änd.-Nr. 56126

Aesculap® BiCONTACT®

BiCONTACT® hip endoprosthesis stem

Intended use

The implant is used

- as a component of a human hip endoprosthesis: Hip endoprosthesis stem
- to combine with Aesculap hip endoprosthesis components
- for implantation without bone cement with PLASMAPORE® or PLASMAPORE®-μ-CaP-coated stem
- for implantation with bone cement, for uncoated stem

The implant range consists of BiCONTACT®, BiCONTACT® S, H, MS, SD, D, N, E or EH and BiCONTACT® Revision.

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_F titanium forged alloy Ti6Al4V acc. to ISO 5832-3
 - ISODUR®_F cobalt-chromium forged alloy CoCrMo according to ISO 5832-12
 - PLASMAPORE® surface coating pure titanium acc. to ISO 5832-2
 - PLASMAPORE®-μ-CaP surface coating pure titanium acc. to ISO 5832-2 with additional calcium phosphate coating
- BiCONTACT®, ISOTAN®, ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for acute disorders of the hip that cannot be treated by other therapies:

- Degenerative osteoarthritis
- Rheumatic arthritis
- Fractures of the joint
- Femoral head necrosis

WARNING

Risk of aseptic implant loosening when using bone cement with uncoated BiCONTACT® prosthesis shafts made of ISOTAN®!

- Do not use uncoated BiCONTACT® prosthesis shafts made of ISOTAN® (titanium forged alloy ISO 5832-3) for implantations with bone cement.

Exceptional indications: Patients with proven foreign-body sensitivity for the alloy elements nickel, cobalt or chromium, coinciding with indication for prosthesis stem anchoring with bone cement.

Indications for using stem type BiCONTACT® Revision in surgery to replace hip endoprosthesis stems

- Cement-free revision operations
- Periprosthetic fractures
- Subprosthetic fractures
- Treatment of large bone defects on the proximal femur

Contraindications

Absolute contraindications

Do not use in the presence of:

- Joint diseases that can be treated by reconstructive surgery (e.g. displacement osteotomy)
- Acute or chronic infections near the joint, or systemic infections
- Secondary diseases that could influence joint implant functionality
- Acute osteoporosis or osteomalacia
- Poor bone quality as well as osseous deformations, diseases in the area of the implant anchors, which could affect the primary or future status of anchor stability of the joint replacement.
- Known hypersensitivity to the implant material.
- Prosthesis heads with XXL neck length in combination with BiCONTACT® H cement-free, SD, D, N, E or EH and BiCONTACT® Revision.
- All of the areas of application not listed under the indications.

WARNING

Implant failure due to excessive body weight of the patient and mainly distal fixation of the small BiCONTACT® implants!

- Avoid distal implant anchoring by appropriate bone preparation.
- Drill open the distal marrow cavity, using flexible marrow cavity drills.
- Ensure sufficiently deep prosthesis seat through proximal-medial support.
- Carry out intraoperative check of the prosthesis seat, using an imaging device.

Observe the maximum allowable body weight of the patient for distal implant fixation: 50 kg for NK709T and NK710T (BiCONTACT® SD sizes 9 and 10), 60 kg for NJ208T, NJ209T, NJ210T (BiCONTACT® D sizes 8, 9 and 10), NJ628T (BiCONTACT® EH sizes 8) and 65 kg for NK210T (BiCONTACT® Revision SD size qq11).

Relative contraindications

The following conditions, individual or combined, can lead to delayed healing or compromise the success of the operation:

- Anticipated excessive load on the joint implant
- Pronounced osseous deformities
- Disruption of the bone metabolism
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Expected lack of patient cooperation

Use of the products requires individual, critical assessment by the surgeon.

Even the change of an endoprosthetic implant is a complicated procedure with individual preconditions which must be assessed accordingly by the surgeon.

Side effects and interactions

Always consider current known possible risks, including resulting clinical consequences regarding the use of the product:

- Changes in position, loosening, wear and tear of, or fracture of implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Thromboses, embolisms
- Tissue reaction to implant materials
- Injury to the surrounding tissue, including nerve and vascular damage
- Hematomas and wound healing disorders

- Periarticular calcification
- Reduced joint mobility and flexibility
- Limited resilience of the joint and joint pain

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Combine modular implant components only with suitable Aesculap hip endoprostheses.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damage to the implant, especially in the neck or cone region, caused by instruments (e.g. HF surgical devices) applied close to the implant.

CAUTION

Risk of implant component breakage due to combination with implant components from other manufacturers!

- Use Aesculap implant components only.

WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

WARNING

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the patient's body weight and the strain to which the joint is subjected.
- The prosthetic joint must not be overstrained by extreme loads, hard physical labor or sporting activities.
- Corrective surgery may become necessary if the implant loosens.
- In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The revision of a hip endoprosthesis stem is a complex joint replacement operation.
- Generally, the revision joint replacement will be inferior in its performance to the primary joint replacement.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

The implantation and implantation site are prepared in the following way:

- Following the femoral head osteotomy and after opening the medullary cavity, prepare the femur with the BiCONTACT® specific A and B osteoprofilers.
- Carry out intraoperative orientation of the implant components. When doing this, pay attention to the resection line and the trochanter wing.
- Select the implant according to B-osteoprofiler last introduced in the correct position.
- Before inserting the implants, carry out a test reposition and examine joint mobility, joint stability and leg length.

WARNING

Bone fractures in the implant bed will affect the anchoring of the implants!

- ▶ Avoid bone fractures by applying careful surgical techniques.
- ▶ Treat bone fractures by appropriate intraoperative and postoperative therapies.
- ▶ Observe correct application of the implant components.

- ▶ Do not under any circumstances allow the implant surfaces to be damaged.
- ▶ Make certain the prosthesis shafts and heads are of the same cone size (see cone size printed on implant packaging, e.g. 12/14).
- ▶ Only remove the protective cap from the prosthesis cone immediately before setting the prosthesis head in place.
- ▶ Prior to setting the prosthesis head in place, rinse, clean and dry the exterior cone of the stem, and if necessary, the inner cone of the prosthesis head.
- ▶ Position prosthesis head with a slight turn under axial pressure and fixate it with a hammer blow using a plastic (not metal) head punch.
- ▶ Prepare and apply bone cement as described in the user instructions provided by the cement supplier.
- ▶ Follow the instructions for use of Aesculap centralizers.
- ▶ Before closing the wound, ensure, if necessary by means of an image converter, that the implant components are correctly positioned.
- ▶ To prevent abnormal wear and tear of the prosthesis: Remove any loose bone cement or bone chips before closing the wound.

When carrying out interventions to replace hip endoprosthesis stems of the type BiCONTACT® Revision, the following instructions must be observed:

- ▶ Prepare the marrow cavity with A and B reamers if you are going to remove the implant through a proximal access (if necessary with an additional bone window). Pay attention to the reamer marking and size according to the planned prosthesis size and length.
- ▶ When preparing the marrow cavity, work the reamers manually, taking care to remove as little bone substance as possible. If necessary, finish the marrow cavity preparation with B-osteoprofilers and conclude by implanting a BiCONTACT® revision stem.
- ▶ If the primary stability (axial fit and rotation), carry out distal locking. Use special aiming or/and image amplification device for this purpose.
- ▶ Use Aesculap 5-mm locking screws for this step. Observe the instructions for use of the locking screws.

Transfemoral access:

- ▶ When operating through a transfemoral access, lock the BiCONTACT® revision stem.
- ▶ Use A-reamers to prepare the distal fitting of the prosthesis. When doing this, keep to a distal anchoring length of at least 100 mm.

WARNING

Risk of locking screw breakage due to permanent distal force transmission!

- ▶ Remove the locking screws after 1 to 2 years (reference period), depending on bone remodeling and integration of the revision implant in the proximal bone fraction. The appropriate timing of this intervention depends on the individual patient's situation; it can be later in cases of delayed bone remodeling.

Note

The BiCONTACT® revision stems are locked for the primary stabilization of the prosthesis stem in cases of severe bone defects. The distal force transmission via the locking screws is designed to be a device for temporary stabilization.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

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