GB	Instructions for use Hydrolift®
USA	Instructions for use Hydrolift® <i>Note for U.S. users</i> <i>U.S. cleared device instructions for use can be found at</i> <i>www.aesculapimplantsystems.com "Products & Services".</i> <i>If you do not have internet access please call your local Aesculap</i> <i>representative or customer service at 1-800-258-1946.</i>
D	Gebrauchsanweisung Hydrolift®
F	Mode d'emploi Hydrolift®
E	Instrucciones de manejo Hydrolift®
I	Istruzioni per l'uso Hydrolift®
P	Instruções de utilização Hydrolift®
NL	Gebruiksaanwijzing Hydrolift®
RUS	Инструкция по применению Hydrolift®
CZ	Návod k použití Hydrolift®
PL	Instrukcja użytkowania Hydrolift®
SK	Návod na používanie Hydrolift®

B | BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 012603 02/11 Änd.-Nr. 42327

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Simbolo CE, em conformidade com a Directiva 93/42/CEE
CE-certificering conform richtlijn 93/42/EEG
Маркировка CE согласно директиве 93/42/ЕЭС
CE-označení podle směrnice 93/42/ES
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
Označenie CE podľa smernice 93/42/EHS

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Technische wijzigingen voorbehouden
Мы оставляем за собой право вносить технические изменения.
Technické zmeny vyhradeny
Zmiany techniczne zastrzeżone
Technické zmeny vyhradené

Intended use

- Ventral monosegmental and multisegmental stabilization of the lumbar and thoracic spine.
- Complete or partial replacement by implants of one or more vertebral bodies.
- For anterior and anterior-lateral approach from an implantation height in excess of 21 mm.
- For posterior approach from an implantation height in excess of 29 mm.

System component		
pre-assembled implants	15 typical sizes	SV001T-SV018T SV030T-SV032T
additional endplates	S, M, L	SV019T-SV024T

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_F titanium forged alloy Ti6Al4V acc. to ISO 5832-3
- PLASMAPORE® surface coating, pure titanium acc. to ISO5832-2
- MP35N CoNiCrMo forged alloy acc. to ISO 5832-6
- PEEK-OPTIMA®
- Silicone

ISOTAN® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

PEEK-OPTIMA® is a registered trademark of Invibo Ltd, Lancashire FY5 4QD, UK.

Indications

Surgically installed implants serve to support normal healing processes. They are not supposed to replace normal body structures or to support permanently loads occurring in cases where healing did not take place.

Use for:

- Fractures
- Anterior correction of spinal deformities
- Degenerative disc disorder (spondylolisthesis, spondylolysis, spinal stenosis)
- Degenerative instability
- Deformities
- Post-trauma instability
- Spinal tumors

Contraindications

Do not use in the presence of:

- Multisegmental fusions with more than 2 vertebral bodies
- Fever
- Infection
 - acute
 - systemic
 - in the spine
 - local
- Pregnancy
- Acute osteopenia
- Medical or surgical conditions that could negatively affect the success of the implantation
- Foreign body sensitivity to the implant materials
- Inadequate patient compliance
- Cases not listed under indications

Side-effects and adverse interactions

Implantation failure due to:

- Overload
- Warping or bending
- Loosening
- Breakage
- Inadequate fixation
- Failed or delayed fusion
- Infection
- Fractured vertebral body or bodies

Injuries to:

- nerve roots
- spinal cord
- blood vessels
- organs

Safety notes

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Only combine Aesculap modular implant components with each other.
- Do not use damaged or surgically excised components under any circumstances.

- Implants that have already been used must not be reused.
- Delayed healing can cause implant breakage due to metal fatigue.
- The physician in charge decides whether an implant component implanted should be removed.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Damage to the load-bearing structures of the implant can lead to loosening of components, dislocation and migration and other severe complications.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiry date or whose packaging is damaged.

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Application

The operating surgeon draws up an operating plan that specifies and appropriately documents the following steps:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
 - Operating conditions are highly aseptic
 - The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
 - The surgeon and the operating team are thoroughly conversant with the surgical technique and the available range of implants and instruments; appropriate information materials are complete and ready to hand.
 - Those performing operations must be thoroughly conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant articles by medical specialists from the professional literature.
 - The manufacturer has been consulted if the preoperative situation was ambiguous or if implants were found in the operating site.
- The intervention has been explained to the patient, whose consent concerning the following information has been documented:
- The patient is aware of the risks involved in neurosurgery, general surgery, orthopedic surgery and general anesthesia
 - The patient has been informed of the advantages and disadvantages of a Hydrolift® implant, and has been made aware of possible alternative methods of treatment.
 - The special instrument set for installing the implant components and for the distraction, compression and reduction of the lumbar and thoracic spine has been prepared and is ready.
 - In the case of delayed or incomplete fusion, the implants can break and loosen due to high loads.
 - The life-span of the implant depends on the patient's body weight.
 - The implant components must not be overloaded by extreme strains, hard physical labor or sports.
 - Surgical revision may become necessary in cases of implant loosening, implant breakage and correction loss.
 - The patient must undergo medical check-ups of the implant components at regular intervals.

Implantation of the Hydrolift® implant requires the following steps:

- To avoid internal stress on, and weakening of, the implant: avoid scoring or scratching of the implant components.
- Do not deform metal implants.

Discectomy and preparation of the vertebral body endplate

- 

Insufficient anchoring stability due to incomplete removal of cartilage!

 - Completely remove the cartilage layer from the vertebral bodies.
- 

Increased risk of migration due to over-preparation of the vertebral body endplates!

 - Avoid aggressive preparation of the endplates.
 - Do not remove the vertebral body endplates.
- Using standard instruments, remove just enough disc material so that the trial implant can be inserted in the disc compartment.
 - Remove all cartilage from the vertebral body endplates.
 - Carry out partial or complete corporectomy of the damaged vertebral body.

Verification of implant size

- 

Risk of injury due to use of an unsuitable implant size!

 - Only use Hydrolift® implants for posterior approach from an implantation height in excess of 29 mm.

- Select Hydrolift® implants according to indication, preoperative planning and intraoperative findings.
 - When selecting the implant, consider the sizes of the vertebral body endplates, the resected space between the vertebral body endplates and the intended distraction distance.
 - To select the implant, use trial implants for implant endplates and body. The height of the trial implants for the implant endplates corresponds to the height of the implant endplate with spikes. The length of the trial implants for the implant body corresponds to the minimal length of the implant with spikes.
 - Always use the largest possible trial implant for selecting the implant size.
 - Use radiographic control to ensure that the trial implants are positioned centrally. The holes represent the rotation axes of the implant endplates.
 - If necessary, replace implant endplates unilaterally.
- The pre-mounted endplates of the Hydrolift® implant can be replaced by a different size of the additional endplates that are available, see Replacing the pre-mounted implant endplate.

Filling the hydraulic applier

- 

Risk of infection due to using any other medium than sterile 0.9 % NaCl solution!

 - Only use sterile 0.9 % NaCl solution for filling the hydraulic applier.
- Follow the instructions for use of the hydraulic applier.
 - Aspirate the sterile 0.9 % NaCl solution into the hydraulic applier and hold the hydraulic applier with the tube vertically upright to remove any air bubbles.
 - If necessary, connect the extension hose to the hydraulic applier for complete de-aeration.
 - Connect the hydraulic tube to the hydraulic applier or the extension hose, respectively, and de-aerate completely.
 - Top up the applier with hydraulic fluid, if necessary.

Presetting the endplate angle

- 

Risk of dislocation and vascular injury caused by endplate clamping opened to both sides!

 - Preset and lock the angle of the implant endplates.
 - Distract with only one endplate clamp opened.
- 

Inaccessible endplate clamping due to rotation of the endplates!

 - The endplate clamping must be locked prior to implantation through an anterior or anterolateral approach.
 - Use the appropriate torque wrench for locking the clamps.

- Adjust the endplate angles with the adjusting device.
- For the anterior or anterolateral approach, rotate the implant endplates appropriately before locking them in position.
- Lock the endplate clamps with the torque wrench.

Implantation

- 

Insufficient distraction force due to inadequate de-aeration of the system!

 - De-aerate the hydraulic applier, the extension hose and the hydraulic tube completely.
 - Completely compress the implant.
- Couple and lock the Hydrolift® insertion instrument at the selected, pre-adjusted implant.
 - Use the insertion aid to facilitate the intercostal approach for thoracoscopic implantation and to avoid tissue damage.

- 

Incorrect anatomic positioning of the implant endplates. Risk of vascular injury caused by clamping screw positioned ventrally!

 - Always position the Hydrolift® implant in such a way that the clamping screw is in a dorsal position.
 - When replacing the endplates, make certain that the endplate clamping screw is in a dorsal position and accessible from above.

- Note the position of the implant endplates.

- 

Vascular damage caused by protruding implant!

 - Take care that the Hydrolift® implant is positioned centrally.

Distraction

- 

Implant breaking into the vertebral body endplates!

 - Apply proper care when distracting the Hydrolift® implant.
 - If necessary apply radiographic monitoring.
 - Be aware of the tactile feedback of the hydraulic applier.
 - Avoid overdistraction.
- Insert the hydraulic tube through the Hydrolift® insertion instrument and couple it to the implant.
 - Distract the implant by actuating the hydraulic applier.

Locking the axial clamping mechanism and the endplate clamps

- Open the endplate clamps, if applicable, distract and lock the endplate clamps again.

Risk of dislocation and vascular injury caused by endplate clamping opened to both sides!

- Preset and lock the angles of the implant endplates.
- Distract with only one endplate clamp opened.

Correction loss due to insufficient distraction!

- Lock the Hydrolift® implant only at the required distraction.
- Ensure that the hydraulic applier is under pressure when the implant is locked.

Correction loss due to insufficient axial clamping and clamping of the implant endplates!

- Always use the appropriate torque wrench for locking the axial clamping mechanism and the endplate clamps (tightening torque for axial clamping: 12 Nm; tightening torque for endplate clamping: 5 Nm).

Implant breakage caused by excessive tightening torque applied to the axial clamping mechanism or the endplate clamps!

- Always use the appropriate torque wrench for locking the axial clamping mechanism and the endplate clamps.
- Do not exceed the tightening torques of 12 Nm for axial clamping and 5 Nm for endplate clamping.

Implant breaking out or malpositioning due to torque transfer to the spine!

- Always use the Hydrolift® insertion instrument as counter torque when locking the axial clamping mechanism and the endplate clamps.
- Make certain that the Hydrolift® insertion instrument is safely coupled/locked to the implant.

- Lock the axial clamping mechanism.
- Lock the endplate clamps.

Removing the connection port

- Reduce the pressure at the hydraulic applier.
- Remove the hydraulic tube from the insertion instrument.
- Remove the connection port, using the connection uncoupling instrument.

Removing the insertion instrument

- Unlock, uncouple and remove the Hydrolift® insertion instrument.

Inadequate stabilization of the spine due to insufficient instrumentation!

- Use the Hydrolift® implant only in combination with an additional stabilizing system.

Replacing the pre-mounted implant endplate

Intraoperative locking of the endplate clamps impossible due to incorrect screw placement!

- Mount the endplate clamping screws in such a way that they are accessible from the direction of the insertion instrument.

Insufficient clamping stability of the endplate due to incorrect mounting of the endplate inner ring!

- The slot of the inner ring and the slot of the endplates must be in alignment with each other.

Incorrect anatomic positioning of the implant endplates. Risk of vascular injury caused by clamping screw positioned ventrally!

- Always mount the endplate in such a way that the endplate clamping screws are on the same side.
- Always mount the endplates with the endplate clamping screws positioned posterior.

- Adapt and lock the implant in the Hydrolift® insertion instrument.
- Loosen the endplate clamping screw with the appropriate torque wrench.
- Slightly rotate the pre-mounted endplate to remove it from the implant body.
- Align the slot of the inner ring of the endplate with the slot of the endplate.
- Install new endplate on the implant body by slightly rotating the new endplate.
- Fully screw in the clamping screw in the new endplate.
- Pre-adjust the endplate with the adjusting device and lock it with the appropriate torque wrench.

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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