Aesculap Neurosurgery

Zerlegbare Kerrison Knochenstanzen

DEMONTAGE

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Detachable Kerrison Bone Punches

DISASSEMBLY

I Remove spring from outer part of handle.

II Close handle completely and pull down the inner part of the handle.

III Remove sliding part backwards.


Note

- Bone punches with a gold-colored distal working tip or a gold-colored inner handle have a thin footplate. These punches are for the removal of small pieces of bone or soft tissue only.
- Clean the new instrument either manually or mechanically prior to the initial sterilization.
- Inspect the instrument after each cleaning and disinfecting cycle and before use in order to be sure it is clean, functioning properly, not damaged and does not have any loose, bent, broken, cracked, worn or fractured components. Do not use the instrument if it is damaged or defective.
- The bone punches can be sterilized both in the disassembled and in the assembled state.
- Please follow the single instruction steps for disassembly and assembly of the bone punch!

Cleaning, Disinfection and Sterilization

After each use:

- Clean the contaminated instrument as soon as possible. Perform manual or mechanical cleaning. Follow the manufacturer's instructions. When carrying out mechanical cleaning, place instrument in wire baskets suitable for cleaning (avoid watermarks).
- Detachable instrument should be disassembled, see "Disassembly", dry discharge is recommended. If the instrument is reprocessed while wet, use an active cleaning disinfectant. Prior to mechanical cleaning and disinfecting rinse the instrument thoroughly in clear running water.
- Place the instrument in a suitable active cleaning disinfectant solution in such a way that all surfaces, interior surfaces, lumina, and openings are covered. Follow the disinfectant manufacturer's instructions.
- After disinfecting chemically, always rinse the item thoroughly in plenty of clear running water. Follow the disinfectant manufacturer's instructions.
- Remove encrusted material with a soft nylon brush. Do not use harsh cleaning agents or metal brushes. Clean lumina and channels with a soft round nylon brush. Use brushes that are an appropriate diameter for the lumina they are being used to clean.
- Carry out the final rinse in distilled or fully demineralized water.
- Dry the instrument with a soft, absorbent, lint-free cloth. Dry lumina and channels with compressed air.


Lubricate the defined places on the instruments (instruments jaw and red toggle switch) with Aesculap Sterilit Oil e. g. JG598 after each cleaning cycle.

- Moving parts such as hinges and latches should be lubricated with a sterilizable and moisture-permeable light weight machine oil (such as AEscULAP STERILIT maintenance oil JG598).
- The sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, and validated in conformance with EN 554/ISO 13683). In case of application of the fractionated vacuum procedure, the sterilization has to be carried out for a minimum of 5 minutes at 274 °F/134 °C and at 2 bar pressure.

Storage

- Keep instrument in an appropriate storage tray and store in a dry (e.g. storage rack for bone punches JF120R), clean, and safe place.

ASSEMBLY

I Put on sliding part and push it forward.

II Press up inside part of handle and open handle.

III Insert spring into outer part of handle.

For further informations please see our brochure:

Proper Maintenance of Instruments

C21711 Kerrison Detachable bone punches

C08502 Kerrison Noir

C61202 Kerrison bone punches Long version


C61302 Kerrison bone punches Half sizes

C61402 Kerrison bone punches Ejector

C81902 Aesculap Universal Storage System

www.aesculap-neuro.com

Aesculap AG & Co. KG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de