- GB** **Instructions for use/Technical description**
Metha® short stem prosthesis
- USA** **Instructions for use/Technical description**
Metha® short stem prosthesis
Note for U.S. users
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
- D** **Gebrauchsanweisung/Technische Beschreibung**
Metha® Kurzschaftprothese
- F** **Mode d'emploi/Description technique**
Prothèse à tige courte Metha®
- E** **Instrucciones de manejo/Descripción técnica**
Prótesis de vástago corto Metha®
- I** **Istruzioni per l'uso/Descrizione tecnica**
Protesi a stelo corto Metha®
- P** **Instruções de utilização/Descrição técnica**
Haste protésica curta Metha®
- NL** **Gebruiksaanwijzing/Technische beschrijving**
Metha®-korteschaftprothese
- S** **Bruksanvisning/Teknisk beskrivning**
Metha® Kortskafprotes
- RUS** **Инструкция по применению/Техническое описание**
Короткая протезная ножка Metha®
- CZ** **Návod k použití/Technický popis**
Protéza s krátkým dřikem Metha®
- PL** **Instrukcja użytkowania/Opis techniczny**
Proteza z krótką szypką Metha®
- SK** **Návod na použitie/Technický opis**
Protéza s krátkym driekom Metha®
- TR** **Kullanım Kılavuzu/Teknik açıklama**
Metha® kısa şaft protezi
- KR** **사용 설명서 / 기술 설명**
Metha® Short stem 임플란트

Aesculap® Metha®

Metha® short stem prosthesis

Intended use

The implant is used

- as a component part of a human hip endoprosthesis: hip endoprosthesis stem,
- in combination with other Aesculap hip endoprosthesis components,
- for implantation without bone cement.

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_F titanium forged alloy Ti6Al4V acc. to ISO5832-3
 - PLASMAPORE®_{μ-CaP} surface coating pure titanium acc. to ISO 5832-2 with additional calcium phosphate coating
- ISOTAN® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for acute disorders of the hip that cannot be treated by other therapies:

- Degenerative osteoarthritis
- Rheumatic arthritis
- Femoral head necrosis

Contraindications

Absolute contraindications

Do not use in the presence of:

- Joint diseases that can be treated by reconstructive surgery (e.g. displacement osteotomy)
- Existing infections near the joint, or systemic infections
- Secondary diseases that could be expected to influence the joint implant functionality
- Insufficient quality of the bone or osseous deformities, or disorders in the implant fixation area which could impair the stability of the fixation of the joint replacement primarily or over time
- Known hypersensitivity to the implant materials
- Cups whose implantation, based on ISO 21535, result in a limited range of motion (less than 80° in flexion/extension or less than 60° in abduction/adduction), e.g. Bipolar cup
- Prosthesis heads with neck length XL in combination with short stem prosthesis size 0 and 1
- Prosthesis heads with neck length XXL
- All areas of application not mentioned under 'Indications'

Relative contraindications

The following conditions, individual or combined, can compromise the success of the operation:

- Expected overload on the joint implant, especially due to high patient weight or intense physical strain and activity
- Severe osseous deformities
- Prevents axial misalignments
- Disorders of bone metabolism
- Pharmaceutical or other drug abuse; alcoholism
- Expected lack of patient cooperation

Use of the product in these cases requires the individual, critical assessment of the doctor.

Side effects and interactions

- Dislocation, loosening, wear, corrosion and fracture of implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Injury to blood vessels and nerves
- Hematomas and wound healing disorders
- Periarthral calcification
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components have been tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Combine modular implant components only with suitable Aesculap hip endoprostheses.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damage to the implant, especially in the neck or cone region, caused by instruments (e.g. HF surgical devices) applied close to the implant.

The implant has not been tested for safety and compatibility in an MRI environment. It has not been tested for heating, movements or image artifacts in MRI examinations. Performing an MRI scan of a patient with this implant may cause injury to the patient.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Risk of implant component breakage due to combination with implant components from other manufacturers!

- Use Aesculap implant components only.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the body weight and the load put on the joint.
- The prosthetic joint must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may become necessary if the implant loosens.
- In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

The implantation site is prepared in the following way:

- Following a high femoral neck osteotomy, open the femur with the implant-specific awl.
- Introduce the awl until it reaches the lateral corticalis.

Bone fractures/perforations in the implant bed will adversely affect the implant fixation!

- Avoid bone fractures by applying careful surgical techniques.
- Treat bone fractures by appropriate intraoperative and postoperative therapies.
- Handle the implant components properly.

- Gradually prepare the implant bed with the implant-specific form rasps (starting with the smallest size).
- Check and, if necessary, correct the implant position, depth and antetorsion.
- Check the stem size against the form rasp that was introduced last in the correct position.
- Attach the rasp trail adapter and perform trial reduction with the appropriate trial heads; change rasp trial adapter or trial heads as appropriate.
- Check the mobility/range of movement and stability of the joint and the length of the leg.
- Remove the trial implants and the rasp.
- Insert short stem with specific insertion instruments.
- Select a prosthesis ball according to the trial heads.
- Verify that the cone size of the prosthesis stems are the same as the cone size of the prosthesis heads (see cone size on the implant packaging, e.g. 12/14).
- Rinse, clean and dry the outside cone of the stem and, if necessary, the inside cone of the prosthesis spheres prior to setting the prosthesis head in place.
- Do not remove the protective cap until immediately before positioning the prosthesis head.
- Attach the prosthesis ball and prosthesis cone to each other at room temperature only. If necessary, allow the implant to cool down to room temperature.
- Apply the prosthesis head with a slight turn under axial pressure and fix in place by hammering with a plastic (not metal) head impactor.
- To prevent abnormal wear and tear of the prosthesis: Remove all obvious bone cement residues and bone chips before closing the wound.

Further information on Aesculap implant systems is available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

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