- GB

Instructions for use
 Metha® short stem prosthesis
- USA

Instructions for use
 Metha® short stem prosthesis
 Note for U.S. users
 This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapImplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
- D

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Intended use

The implant is used

- as a component of a human hip endoprosthesis: Hip endoprosthesis stem
- to combine with Aesculap hip endoprosthesis components
- for implantation without bone cement

The short stem is available as either modular or non-modular implant.

The modular short stem consists of:

- Prosthesis stem
- Modular cone adapter

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®: Titanium forged alloy Ti6Al4V acc. to ISO 5832-3
 - ISODUR®: Cobalt-chromium forged alloy CoCrMo according to ISO 5832-12
 - PLASMAPORE®: Ti-CaP surface coating pure titanium acc. to ISO 5832-2 with additional calcium phosphate coating
- ISOTAN®, ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for acute disorders of the hip that cannot be treated by other therapies:

- Degenerative osteoarthritis
- Rheumatic arthritis
- Femoral head necrosis

Contraindications

Do not use in the presence of:

- Joint diseases that can be treated by reconstructive surgery (e.g. displacement osteotomy)
- Acute or chronic periarticular or systemic infections
- Secondary diseases that could influence joint implant functionality
- Systemic conditions and metabolic disorders influencing the functioning of the joint implant
- Acute osteoporosis or osteomalacia
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Bone malformations, defective axial positioning, or other bone conditions that rule out implantation of a hip-joint prosthesis preserving the femoral neck
- Expected overload on the joint implant, especially due to high patient weight or major physical strains and activities
- Medicinal or other drug abuse or alcohol addiction affecting the functionality of the joint implant, especially through temporary loss of control or impairment of the implant integration in the bone
- Expected lack of patient cooperation in dealing with the joint implant, especially during rehabilitation, and non-compliance with doctors' recommendations
- Foreign body sensitivity to the implant materials
- Skeletal immaturity
- Neuromuscular diseases impairing the affected extremity
- Cups whose implantation, in analogy to ISO 21535:2002, result in a limited range of movement (less than 80° in flexion/extension or less than 60° in abduction/adduction), e.g. Bipolar cup
- Prosthesis heads with neck length XL in combination with short stem prosthesis size 0 and 1
- Prosthesis heads with neck length XXL

Side effects and interactions

- Dislocation, loosening, wear, corrosion and fracture of implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Injury to blood vessels and nerves
- Hematomas and wound healing disorders
- Periarticular calcification
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Combine modular implant components only with suitable Aesculap hip endoprostheses.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damage to the implant, especially in the neck or cone region, caused by instruments (e.g. HF surgical devices) applied close to the implant.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

CAUTION

Risk of implant component breakage due to combination with implant components from other manufacturers!

- Use Aesculap implant components only.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

WARNING

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the body weight and the load put on the joint.
- The prosthetic joint must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may become necessary if the implant loosens.
- In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

The implantation site is prepared in the following way:

- Following a high femoral neck osteotomy, open the femur with the implant-specific awl.
- Introduce the awl until it reaches the lateral corticis.

WARNING

Bone fractures/perforations in the implant bed will adversely affect the implant fixation!

- Avoid bone fractures by applying careful surgical techniques.
- Treat bone fractures by appropriate intraoperative and postoperative therapies.
- Handle the implant components properly.

- Gradually prepare the implant bed with the implant-specific form rasps (starting with the smallest size).
- Check and, if necessary, correct the implant position, depth and antetorsion.
- Check the stem size against the form rasp that was introduced last in the correct position.
- Attach the rasp trail adapter and perform trial reduction with the appropriate trial heads; change rasp trial adapter or trial heads as appropriate.
- Check the mobility/range of movement and stability of the joint and the length of the leg.
- Remove the trial implants and the rasp.
- When using the modular short stem, couple the prosthesis stem and the modular cone adapter prior to insertion.

WARNING

Implant malfunction due to insufficiently fixated cone adapter!

- Fixate the modular cone adapter in the short stem by applying a hammer blow and using instrument ND401R.
- Securely fixate the modular cone adapter in the short stem.

- Insert the short stem (modular/non-modular) with specific insertion instruments.

WARNING

Loosening of the cone adapter during insertion from applying blows to the Metha® stem!

- Also perform the insertion procedure by applying blows to the cone adapter with the cone protection cap applied.

Optional procedure

Trial reduction with inserted modular short stem

- ▶ Remove the trial cone adapter only after the stem component has been impelled.
- ▶ Insert the stem component and punch it into the implant bed with well-dosed blows.
- ▶ If necessary, change the trial cone adapter, making certain that the arrow points to "MEDIAL".
- ▶ Carry out a trial reduction.
- ▶ Check the mobility/range of movement and stability of the joint and the length of the leg.
- ▶ Remove the trial implants.

WARNING

Fracture and increased wear of the modular implant components, or loosening of the implant connection, due to soiled or damaged joining surfaces!

- ▶ Carefully rinse the inside cone of the stem and the joining surfaces of the cone adapter; clean with cleaning swabs (ND619 or ND622) and dry.
- ▶ Check the joining surfaces for any damage.
- ▶ Do not join modular implant components more than once.
- ▶ Always rinse, clean and dry the cone prior to joining.

- ▶ Select the implant cone adapter according to the fitting trial cone adapter.
- ▶ Carefully insert the implant cone adapter in the short stem already implanted, making certain that the arrow points to "MEDIAL".
- ▶ Fixate the implant cone adapter with hammer taps on the appropriate insertion instrument.

Concluding steps for applying (modular/non-modular) short stems

- ▶ Select a prosthesis ball according to the trial heads.
- ▶ Verify that the cone size of the prosthesis stems are the same as the cone size of the prosthesis heads (see cone size on the implant packaging, e.g. 12/14).
- ▶ Rinse, clean and dry the outside cone of the stem and, if necessary, the inside cone of the prosthesis spheres prior to setting the prosthesis head in place.
- ▶ Do not remove the protective cap until immediately before positioning the prosthesis head.
- ▶ Attach the prosthesis ball and prosthesis cone to each other at room temperature only. If necessary, allow the implant to cool down to room temperature.
- ▶ Put the prosthesis ball in its position.
- ▶ To prevent abnormal wear and tear of the prosthesis: Remove all obvious bone cement residues and bone chips before closing the wound.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

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