- (GB)** Instructions for use/Technical description
Metha rasp handle
- (USA)** Gebrauchsanweisung/Technische Beschreibung
Metha-Raspelhandgriff
- (D)** Mode d'emploi/Description technique
Porte-râpe Metha
- (F)** Instrucciones de manejo/Descripción técnica
Mango de raspador Metha
- (E)** Istruzioni per l'uso/Descrizione tecnica
Impugnatura per raspa Metha
- (I)** Instruções de utilização/Descrição técnica
Cabo de fresa Metha
- (P)** Gebruiksaanwijzing/Technische beschrijving
Metha-rasphandgreep
- (NL)** Bruksanvisning/Teknisk beskrivning
Metha-rasphandtag
- (S)** Инструкция по применению/Техническое описание
Рукоятка для рашпиля Metha
- (RUS)** Návod k použití/Technický popis
Rukojeť rašple Metha
- (CZ)** Instrukcja użytkowania/Opis techniczny
Uchwyt raszpili Metha
- (PL)** Návod na používanie/Technický opis
Rukoväť rašple Metha
- (SK)** Kullanım Kılavuzu/Teknik açıklama
Metha raspa el tutamağı
- (TR)**

B | BRAUN
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 011592 02/13 V6 Änd.-Nr. 46355

CE 0123 – DIR 93/42/EEC

Technical alterations reserved


Fig. 1


Fig. 2


Fig. 3

Aesculap®

Metha rasp handle

Legend

- 1 Rasp handle
- 2 Hammer plate
- 3 Handle part
- 4 Push button
- 5 Spring
- 6 Locking screw
- 7 Rasp

Symbols on product and packages

Caution, general warning symbol

Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

The rasp handle is used in hip-endoprosthetic interventions. It holds prosthesis-specific rasps, which are driven in and out during the preparation of the femoral marrow cavity

Available sizes

Different variants of the rasp handle are available for different operating techniques, e.g.:

Designation		Art. no.
lateral access	- straight	NF180R
	- left angled	NF141R
	- right angled	NF142R
posterior access	- straight	NF144R
	- left angled	NF142R
	- right angled	NF141R
anterior access	- straight	NF140R
	- left angled	NF138R
	- right angled	NF139R

Safe handling and preparation

CAUTION
Federal law restricts this device to sale by, or on order of a physician!

Note
Use of the Metha rasp handle requires precise knowledge of hip endoprosthetics and the biomechanical situation at the hip!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Inserting a rasp into the rasp handle

- Push and hold down push button 4.
- Slide rasp 7 in the handle part 3, see Fig. 1. making certain that the notch in the rasp connector points towards the push button.
- Release pushbutton 4.
- Check that the rasp is firmly seated: verify that
 - there is no gap between the endplates,
 - the rasp cannot be rotated against the handle.

Driving the rasp in and out

WARNING

Risk of injury and/or damage to the product due to incorrect application of the hammer!

- When driving the rasp in and out, using a hammer, only hit the hammer plate of the metha rasp handle.

Disassembling

- Push in push button 4 at least halfway down.
- Rotate push button 4 by approx. 90° in any direction and release the pushbutton.
The locking mechanism becomes detached.
- Remove push button 4 and spring 5.

Note
Never unscrew lateral locking screw 6 on the rasp handle. This screw has been sealed with a special adhesive during the manufacturing process; it must not be removed.

Assembling

- Set spring 5 on the push button 4.
- Slide push button 4 into handle part 3, making sure that the notch on the pushbutton is aligned as shown in Fig. 3.
- Push in push button 4 at least halfway down.
- Rotate push button 4 clockwise by approx. 90°, or turn the notch on pushbutton 4 to the front.

Validated reprocessing procedure

General safety instructions

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary. Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for, e.g., aluminum, plastics and high-grade steel, according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum permitted cleaning temperature of 120 °C.

Validated cleaning and disinfection procedure

Note
Sterilization should always be favored over disinfection under all circumstances. Only if final sterilization is not possible should disinfection with a virucidal agent be preferred.

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection	Cleaning brush: for example TA no. 006874 20 ml disposable syringe	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: Chapter Manual pre-cleaning with ultrasound and brush Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

Risk to patients!
► The product must only be cleaned mechanically!

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-
D-W:		Drinking water				
RT:		Room temperature				
*Recommended:		BBraun Stabimed				

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I
- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: pH = 13 <5 % anionic surfactant 0.5 % working solution pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!
► Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

Risk of injury and/or malfunction!
► Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:

Aesculap Implant Systems LLC
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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