
Aesculap Neurosurgery

-
- GB** Instructions for use/Technical description
MINOP instruments
 - USA** MINOP instruments
 - D** Gebrauchsanweisung/Technische Beschreibung
MINOP-Instrumente
 - F** Mode d'emploi/Description technique
Instruments MINOP
 - E** Instrucciones de manejo/Descripción técnica
Instrumentos MINOP
 - I** Istruzioni per l'uso/Descrizione tecnica
Strumenti MINOP
 - P** Instruções de utilização/Descrição técnica
Instrumentos MINOP
 - NL** Gebruiksaanwijzing/Technische beschrijving
MINOP-instrumenten
 - PL** Instrukcja użytkowania/Opis Techniczny
Instrumenty MINOP

MINOP instruments

Legend

- 1 Star wheel
- 2 Handle
- 3 Button
- 4 Moveable handle part
- 5 Shaft complete (outer tube, inner tube, working insert)
- 6 Jaw piece
- 7 Working insert
- 8 Leaf springs
- 9 Inner tube
- 10 Outer tube
- 11 Spring elements (tactile feedback)

Symbols on product and packages

Symbol	Explanation
	Caution: general warning symbol Caution: see documentation supplied with the product

Intended use

The dismountable Aesculap MINOP instruments (hollow-shaft mini instruments) are used for cutting, dissecting, grasping and removing of biopsies, depending on the type of working end.

They are introduced into the body by means of access instruments (e.g. trocar).

System components

Designation	Art. no.
MINOP micro scissors sharp/sharp	FF385R
MINOP micro scissors blunt/blunt	FF386R
MINOP biopsy spoon forceps	FF387R
MINOP grasping and dissecting forceps	FF388R
MINOP surgical micro forceps	FF389R

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

WARNING

Risk of injury caused by incorrect operation of the product!

- **Attend appropriate product training before using the product.**
- **For information regarding such training, please contact your national B. Braun/Aesculap agency.**
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.

- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Exercise caution when passing products through the working channel (e.g. trocar).

Safe operation

Risk of injury and/or malfunction!

- **Always carry out a function check before using the product.**

Risk of injury when applying the product outside the field of view!

- **Apply the product only under visual control.**

Note

The spring elements give the user some tactile feedback just before the instrument emerges from trocar FF399R into the operating site. The functionality of this tactile feedback mechanism is only ensured when using the product with trocar FF399R.

Disassembling

The MINOP instruments are disassembled into the following components:

- Handle 2
- Outer tube 10
- Inner tube 9
- Working insert 7


Fig. 1

- To separate handle 2 and working insert 7, see Fig. 1:
 - Fully open moveable handle part 4.
 - Push and hold down button 3.
 - At the same time, pull down moveable handle part 4.


Fig. 2

- Separate shaft 5 from handle 2, see Fig. 2:
 - Pull star wheel 1 back towards handle 2 as far as it will go.
 - Hold star wheel 1 at its positive stop and extract shaft 5.

MINOP instruments

Assembling


Fig. 3


Fig. 5


Fig. 4

- To disassemble shaft 5, see Fig. 3:
- Remove outer tube 10 and inner tube 9 from working insert 7.
 - Remove inner tube 9 from working insert 7, see Fig. 4.


Fig. 6

- To assemble shaft 5, see Fig. 5:
- Slide inner tube 9 with leaf springs 8 towards the working end onto working insert 7 until leaf springs 8 engage.
 - Slide outer tube 10 over inner tube 9 onto working insert 7 as far as it will go, see Fig. 6.
 - Ensure that the arrows of the circumferential marking on outer tube 10 point towards handle 2.


Fig. 7


Fig. 8

Fig. 9

- To connect shaft **5** and handle **2**:
 - Ensure that outer tube **10** is pushed over inner tube **9** with working insert **7** to the positive stop, see Fig. 7.
 - Slide back and hold star wheel **1** at its positive stop, see Fig. 8.

- Hold the shaft **5** at its closed working end.
- Slightly turn shaft **5** and handle **2** during insertion until they connect to each other.
- Allow star wheel **1** to slide forward, see Fig. 9. Ensure that shaft **5** is securely locked in handle **2** and cannot be removed even by pulling.


Fig. 10

- To connect handle **2** and shaft complete **5**:
 - Fully open handle **2**.
 - Hold working insert **7** and push up moveable handle part **4** until button **3** engages, see Fig. 10.
 - Check that moveable handle part **4** is securely locked.

MINOP instruments

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, the manufacturer's specifications can only serve as an approximate guide for assessing the processing procedures applied by the individual operator/processors.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures >45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process can already cause visible surface changes.
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.

- Do not use process chemicals that cause stress cracking or brittleness of plastics.
- Clean the product immediately after use. Immersion treatment in a 3 % H_2O_2 solution for approx. 5 minutes is a particularly effective and gentle method to dissolve encrustations from HF instruments. Subsequently, the debris can be removed by hand, with a medium-hard brush and/ or in an ultrasonic cleaning bath. After that treatment, proceed through the normal steps of processing.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open jointed products.
- Rinse surfaces inaccessible to visual inspection, e.g. on products with hidden crevices or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Preparation prior to cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/Disinfecting

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to manufacturers' instructions. The cleaning/disinfecting agents must
 - be approved for e. g. aluminum, plastics and high-grade steel
- Observe specifications regarding concentration, temperature and exposure time.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting,
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting,
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting,
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

MINOP instruments

Manual cleaning/disinfecting

- Keep working ends open during cleaning.
- Clean products with movable links with the links opened or mobilized.
- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	0.5	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make sure that all accessible surfaces are moistened.
- Clean the product under running water, if necessary with a suitable cleaning brush (TA no. 012889), until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometry, for at least 1 min or until no more residues can be removed. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly (at least five times) rinse through these components with the cleaning/disinfecting solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution. Make sure that all accessible surfaces are moistened.

Stage IV

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that A_0 is $>3\,000$ for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

MINOP instruments

Manual pre-cleaning with brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Irrigation	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make sure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e. g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush (TA no. 012889).
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfector without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Keep working ends open during cleaning.
- Ensure that product links and joints are opened when placing the product on the tray.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER ALCALINE with tensides; application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfectant program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Assemble the separable product, see Assembling.
- Check for compatibility with associated products.
- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.

MINOP instruments

Packaging

- Appropriately protect products with fine working tips.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

Note

The product can be sterilized either in disassembled or in assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterilcontainer (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F– 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap Sterilcontainer cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Technical Service

Risk of injury and/or malfunction!

- **Do not modify the product.**

WARNING

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Art. no.	Designation
C60901	Assembly/Disassembly poster (German)
C60902	Assembly/Disassembly poster (English)
FF399R	Trocar
FF432R	Handle (only, for FF385R-FF489R)
FF433R	Outer tube (only, for FF385R-FF489R)
FF435R	MINOP micro scissors sharp/sharp
FF436R	MINOP micro scissors blunt/blunt
FF437R	MINOP biopsy spoon forceps
FF438R	MINOP grasping and dissecting forceps
FF439R	MINOP surgical micro forceps
TA012889	Cleaning brush

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG