GB
USA

Instructions for use

MINOP®/Paediscscope Disposable Introducer for Neuroendoscopy

D

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
CE

0123

CE marking according to directive 93/42/EEC

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MINOP®/Paediscscope Disposable Introducer for Neuroendoscopy

Legend

- 1 Obturator
- 2 Sheath
- 3 Obturator and sheath locked assembly
- 4 Locking flange
- 5 Locking knob
- 6 Sheath handles
- 7 Locking threads
- 8 Handles

Symbols on product

	Sterilization using ethylene oxide
	Not for reuse in terms of the intended application defined by the manufacturer
	Use by
	Attention, see instructions for use

Intended use

The purpose of the Disposable Introducer is to obtain and maintain a temporary pathway to the ventricular system of the patient. The Disposable Introducer is designed to be split lengthwise and peeled down to the skull level of the patient, accommodating different depth requirements. The Disposable Introducer consists of two main parts: Obturator and sheath, which are locked to an assembly.

Available sizes

Disposable Introducer	Size (Ø inside of sheath)
FH603SU	10 French (3.2 mm)
FH604SU	19 French (6.1 mm)

Disposable Introducer	Typical usage with
FH603SU	Paediscscope PF010A
FH604SU	MINOP® 6 mm trocar FF399R

Note
These products DO NOT contain latex.

Safe handling and preparation

- Caution**
Federal law restricts this device to sale by or on order of a physician!
The product has been EO-sterilized and is supplied in sterile packaging. The product must not be reused.
- Read, follow and keep the instructions for use.
 - Use the product only in accordance with professional standards and practices, see Intended use.
 - Do not use products taken from open or damaged sterile packaging.
 - Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Immediately set aside the product if it is found damaged.
 - Do not re-sterilize the product.
 - Do not use products that are past their expiration date.

Safe operation

- Ensure that the obturator 1 is locked into the sheath 2 by twisting the locking flange 4 of the obturator, thereby locking it into the locking threads 7 of the sheath. Compare with picture 3.
The tip of the obturator extends from the sheath.
There is no gap between the tip of the sheath and the outer circumference of the obturator.
- Determine the depth the introducer has been inserted by means of the graduation on the side of the obturator.
- Insert the Disposable Introducer to the desired depth, unlock and remove the obturator, being sure to maintain its sterility.
- Insert the scope carefully, making sure that the sheath does not advance.
The scope allows for visual confirmation that the sheath is in the desired location.
- If the sheath is not in the desired location, carefully remove the scope and sheath.

- Reinsert and lock the obturator into the sheath.
- Repeat ventricular insertion. The tip of the sheath should be maintained inside the ventricular system after removal of the obturator.
This reduces the likelihood of the tip to cause hemorrhaging due to the movement or advancement of the scope.

The circular tip of the sheath could cause haemorrhaging!

- **Never insert or advance the sheath into the brain without having the obturator locked into the place.**
- **Extend the obturator completely and lock it into the sheath.**

- Two sterile persons are required to affix the sheath.
- One person has to hold the depth of the sheath constant.
- The other firmly grips the sheath handles 6.
- Pulling apart the handles 8 will result in the splitting of the sheath lengthwise.
- Split the sheath down to a level just above the skull.
- Firmly affix the split sheath in place by suturing or stapling down the two split halves of the sheath.
This has now established the pathway to the ventricle.
- When you now insert anything through the sheath, always be cognizant that the sheath must not move, so tight fitting objects should be avoided.

Single-use products

Possible damage or destruction of products, due to processing!
Intended use: These products have been designed for single-use only.

- **Do not reprocess these products.**

Technical Service

For service, maintenance or repairs, please contact your national B. Braun/ Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420
Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA
Phone: +1 800 282-9000
www.aesculapusa.com

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