

Aesculap Neurosurgery

-
- | | |
|------------|---|
| GB | Instructions for use/Technical description |
| USA | MINOP trocars |
| D | Gebrauchsanweisung/Technische Beschreibung
MINOP-Trokare |
| F | Mode d'emploi/Description technique
Trocarts MINOP |
| E | Instrucciones de manejo/Descripción técnica
Trocars MINOP |
| I | Istruzioni per l'uso/Descrizione tecnica
Trocac MINOP |
| P | Instruções de utilização/Descrição técnica
Trocates MINOP |
| CZ | Návod k použití/Technický popis
Trokar MINOP |

MINOP trocars

Legend

- 1 MINOP trocar (therapeutic), FF399R
- 2 Mono obturator
- 3 MINOP trocar (diagnostic), FF398R
- 4 MINOP trocar (diagnostic), FF397R
- 5 Knurled nut (obturator)
- 6 Obturator (for irrigation or drainage channel, respectively)
- 7 Obturator (for irrigation/drainage channels)
- 8 Guideway
- 9 Obturator (for optics and working channel)
- 10 Locking pin
- 11 Release button
- 12 Irrigation or drainage channel
- 13 Obturator (for optics channel)
- 14 Endoscope
- 15 Electrode
- 16 Instrument

Symbols on product and packages

Symbol	Explanation
	Caution: General warning symbol Caution, see documentation supplied with the product

Intended use

The MINOP trocars are used for endoscopic examinations and therapies in the central nervous system.

Available sizes

Designation	Art. no.
MINOP trocar (diagnostic)	FF397R
MINOP trocar (diagnostic)	FF398R
MINOP trocar (therapeutic)	FF399R

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Note

For the MINOP trocar system, there are specially developed Aesculap endoscopes available. For reasons of compatibility, the MINOP trocars may only be used with these special Aesculap endoscopes!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- For monopolar contact coagulation, make certain that the output power at the HF device does not exceed 50 W.
- To avoid burns when using the ventriculoscope in combination with HF electrodes, ensure that HF current is activated only under visual control.
- When using the HF electrode in combination with a flexible instrument, ensure that the working tips are kept as far apart as possible.
- Set the HF device to the lowest possible output power required for the intended effect.

- Prior to each use, visually check the product for: damage to, or surface changes on, the insulation.
- Never put down the product on or next to the patient.
- Follow the instructions for use of the HF device.

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

WARNING

Risk of injury when applying the product outside the field of view!

- Apply the product only under visual control.

WARNING

Risk of injury caused by protruding 30°-endoscope in locked position!

- Insert the endoscope only when the trocar is in its final position.

WARNING

Risk of burns due to high temperature of the instrument-side end of the optical cable of a light source!

- Apply proper care when operating the light source.

CAUTION

Damage to the endoscope due to incorrect handling or operation!

- Insert the endoscope in the MINOP trocar only if the trocar is not deformed, bent or kinked.

Note

The endoscopes offered for application with the MINOP system must be used with a light source equipped with a spare bulb!

Mounting the camera

Note

Mount the camera only with the endoscope locked in the working trocar!

- Lock the endoscope in the MINOP trocar.
- Mount the camera.

Demounting the camera

Note

Extract the endoscope from the working trocar only after the camera has been demounted from the endoscope!

- Demount the camera from the endoscope.
- Extract the endoscope from the MINOP trocar.

Disassembling

Removing the obturator (for optics and working channel)

- Press the release button and carefully remove the obturator from the MINOP trocar.

Removing the obturator (for irrigation/drainage channels)

- Turn the knurled nut by 90° and carefully remove the obturator from the MINOP trocar.

Assembling

Inserting the endoscope/obturator

- Carefully insert the endoscope/obturator in the MINOP trocar, making sure that the locking pin of the endoscope/obturator is on the same side as the release button of the MINOP trocar.
- Insert the endoscope/obturator down to the positive stop.

MINOP trocars

Inserting the obturators (for irrigation/drainage channels)

Note

Be careful not to commute the obturators. Mind the markings, "L" (= left-hand obturator) and "R" (= right-hand obturator)!

- Hold the obturators at the distal end and carefully slide them into the respective channels, down to the positive stop. Ensure that the guideway on the obturator is aligned with the pin in the main body of the MINOP trocar.
- Turn the knurled nut by 90°.

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray. Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures >45 °C nor any fixating disinfectants (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process solution can already cause visible surface changes.

- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing chemicals (e.g. H_2O_2) for mechanical cleaning, as these can cause bleaching or layer loss.
- Do not use process chemicals that would cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

- Use appropriate cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove the sealing cap from the Luer lock connector.
- Rinse surfaces that cannot be visually inspected, e.g. on products with hidden crevices or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Preparation prior to cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/Disinfecting

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions. The cleaning and disinfecting agents must
 - be approved for e.g. aluminum, plastics and high-grade steel,
 - not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

MINOP trocars

Manual cleaning/disinfecting

- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual ultrasonic cleaning and immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	5	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush.

- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, use a 20-ml single-use syringe to rinse thoroughly, for at least 5 times, these difficult to access parts of the product.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Make certain that all accessible surfaces are moistened.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.
- Rinse lumens at least 5 times, using a disposable syringe (20 ml) and an appropriate rinsing adapter.
- Allow residual water to drip off sufficiently.

Stage V

- Thoroughly dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note

Use appropriate (neutral, enzymatic and mild-alkaline) detergents for this product with aluminum components.

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that A_0 is $>3\ 000$ for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

MINOP trocars

Manual pre-cleaning with ultrasound and brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Irrigation	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Remove jaws protection.
- Clean the product with the appropriate cleaning brush (TE654202) until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush.
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, use a 20-ml single-use syringe to rinse thoroughly, for at least 5 times, these difficult to access parts of the product.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical neutral or mild-alkaline cleaning and thermal disinfection

Machine type: Single-chamber washer/disinfector without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Keep working ends open during cleaning.
- Ensure that product links and joints are opened when placing the product on the tray.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Neutral: BBRAUN HELIMATIC CLEANER neutral, pH-neutral, process solution 0.5 % Mild-alkaline: - concentrate, pH = 9.5 - < 5 % anionic tensides - 0.5 % solution
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfector program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning process if necessary.

MINOP trocars

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

Packaging

- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).

- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

Note

The product can be sterilized either in disassembled or in assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Disassemble the product.
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F–275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.
- Store single-use products in sterile packaging in dust-protected conditions in a dry, dark and temperature-controlled room.

Technical Service

WARNING

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Art. no.	Designation
EJ751200	MINOP sealing caps for side channels (black)
PE184A	0°-optics; Ø 2.7 mm/0°/180 mm
PE204A	30°-optics; Ø 2.7 mm/30°/180 mm
PF893800	MINOP cleaning brush

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG