- (GB)** Instructions for use/Technical description
MINOP optics and angled optics
- (USA)** Gebrauchsanweisung/Technische Beschreibung
MINOP- und Winkeloptiken
- (D)** Mode d'emploi/Description technique
Optiques MINOP et optiques coudées
- (F)** Instrucciones de manejo/Descripción técnica
Ópticas MINOP y ópticas anguladas
- (E)** Istruzioni per l'uso/Descrizione tecnica
Ottiche MINOP e angolari
- (I)** Instruções de utilização/Descrição técnica
Endoscópios MINOP e angulares
- (P)** Gebruiksaanwijzing/Technische beschrijving
MINOP- en hoekendoscopen
- (NL)** Bruksanvisning/Teknisk beskrivning
MINOP- och vinkelendoskop
- (S)** Инструкция по применению/Техническое описание
MINOP и угловая оптика
- (RUS)** Návod k použití/Technický popis
MINOP a úhlové optiky
- (CZ)** Instrukcja użytkowania/Opis techniczny
Optyki MINOP i optyki kątowe
- (PL)** Návod na používanie/Technický opis
Optika MINOP a uhlové optiky
- (SK)** Kullanım Kılavuzu/Teknik açıklama
MINOP ve açılı optikler
- (TR)**

B | BRAUN
SHARING EXPERTISE


Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 013238 08/12 V6 Änd.-Nr. 44392

CE 0123 – DIR 93/42/EEC

Technical alterations reserved

Aesculap®
MINOP optics and angled optics

Legend

- 1 Distal window
- 2 Sheath
- 3 Optical cable connection
- 4 Illumination surface of the optical cable connection
- 5 Eyepiece housing
- 6 Eyepiece window
- 7 Inscription ring
- 8 ACMI adapter (fixed)
- 9 Wolf adapter
- 10 Storz adapter
- 11 Color ring

Symbols on product and packages

	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture
	Steam sterilization at temperatures up to 134 °C
	Low temperature and plasma sterilization
	High-Definition resolution

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

MINOP optics and angled optics are used for visualizing the ventricles and other structures in the brain, as well as structures at the spine.

The following MINOP and angled optics are approved for direct contact with the central nervous system:

- PE184A
- PE204A
- PE486A
- PE506A
- PE526A

Color coding

The viewing direction of the endoscope is indicated by color ring 11 at optical cable connection 3.

Colour	Viewing direction
Green	0°
Red	30°
Yellow	70°

Available sizes

The rigid, autoclavable neuroendoscopes are available in the following sizes and designs:

- Sheath Ø 2.7 mm and 4.0 mm
- Angled endoscopes

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

This manual contains instructions for the preparation, reprocessing, and disposal of the endoscope. It does not contain information on the actual application of the endoscope.

Risk of injury from defective endoscopes!

- Only use an endoscope if it is in perfect condition.

The products are supplied in an unsterile condition and must be reprocessed and sterilized prior to use.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).

- Make certain a reserve product is available.
- Do not, under any circumstances, put down the product on the patient or on the surgical drape covering the patient.

The following table shows the approved combinations of endoscopes, sheaths, and working trocars:

Insert the endoscope	Sheath and working trocar
PE184A	FF397R, FF398R, FF399R, FH601R
PE204A	FF397R, FF398R, FF399R, FH620R

Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Function checks

- Check the optical function of the endoscope. The image has to be clear and distinct.
- Ensure that the distal window 1, the eyepiece window 6 and the illumination surface 4 of the optical cable connection 3 are not nicht cloudy, dirty or scratched.
- Hold optical cable connection 3 of the endoscope against the light and check that the light fibers at the distal end light up evenly.
- Inspect the sheath for dents, bends, and scratches.

Use of the endoscope

Burns caused by high-frequency current (HF current)!

- When using an endoscope and HF electrodes at the same time, be cautious to activate the HF current under visual control only.
- Ensure that the active electrode is never in the immediate vicinity of electrically conductive components (e.g. trocar, endoscope).

Burns suffered by patient or user from exposure to high light intensity!

- Make certain the distal end of the endoscope or optical cable connection does not touch human tissue or any flammable or heat-sensitive materials while the light source is active.
- Do not put down the endoscope on the patient.
- Do not touch the distal end of the endoscope and the optical cable connection.
- Adjust the light source to the minimum required power for optimal illumination of the endoscopic image.
- Only use light sources of up to 300 W electric power.

Damage to the endoscope by bending the endoscope sheath!

- Do not bend the endoscope sheath.
- Use the endoscope with its appropriate sheaths and working trocars only.
- Always hold the endoscope at the eyepiece housing but never at the sheath.

- Use endoscopes only with halogen light sources equipped with a spare lamp or with xenon light sources.
- Observe the respective manufacturer's instructions when using the endoscope in combination with other equipment, instruments, or optical cables.
- When using the endoscope in combination with any electromedical device, make sure that the respective BF conditions are adhered to (insulated, floating patient application piece).

Disassembling

- Unscrew adapter 9 or 10, respectively, from the endoscope.

Assembling

- Screw adapter 9 or 10, respectively, onto the endoscope.

Validated reprocessing procedure

General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

Only clean medical products can be sterilized safely and effectively. Cleaning is therefore of special importance in the reprocessing procedure.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive dosages of neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking on stainless steel becoming unreadable either visually or by machine.

Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes in materials, e.g., fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org Rubrik Veröffentlichungen Rote Broschüre – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse surfaces that are not accessible to visible inspection (preferably with demineralized water), using a disposable syringe, for example.
- Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for rigid endoscopes,
 - are approved for use e.g. on aluminum, plastic materials and stainless steel,
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum permitted cleaning temperature of 55 °C.

CAUTION

Damage to the optical system caused by loosening of connections during ultrasound cleaning!

- Do not clean the endoscope by ultrasonic treatment under any circumstances.

- Do not use oxidizing chemicals (e.g. H₂O₂), which could cause bleaching/layer loss of the product.
- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.
- Only use bactericidal, fungicidal, and virucidal disinfectants.
- Preferably use thermal disinfecting processes.
- Following the chemical disinfection, thoroughly rinse the product with clear running water. Always adhere to the manufacturer's instructions.

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Cleaning	35-45 95-113	5	0,8	D-W	Cidezyme/Enzol
II	Intermediate rinse	RT (cold)	3 x 1	-	D-W	-
III	Disinfection	20-25 68-77	12	0,55	D-W	Cidex OPA (process solution, 0.55% ortho-phthalaldehyde)
IV	Final rinse	RT (cold)	3 x 2	-	FD-W (sterile)	-
V	Drying	RT	-	-	-	-

- Phase I**
- Completely immerse the instrument in the cleaner/disinfectant solution. Ensure that all accessible surfaces are moistened.
 - Clean the product immersed in the solution, using a soft cloth or, if necessary, a suitable cleaning brush, until all visible residues have been removed from the surfaces.
 - Use a suitable cleaning brush to clean all surfaces that are not visible, e.g.in products with concealed crevices, lumens, or complex geometries, for at least 1 minute or until no further residues can be removed. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop.
 - Thoroughly rinse these components with the cleaning disinfectant solution (at least five times), using a disposable syringe (20 ml).
 - Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion

- Phase II**
- Completely rinse/ rinse through the product (all accessible surfaces) 3 times for at least 1 min. Mobilize non-rigid components, e.g. set screws, joints, slides, etc. 3 times in each direction to the positive stop. Use fresh water for every rinse.
 - Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
 - Drain any remaining water fully.

- Phase III**
- Fully immerse the product in the disinfectant solution. Ensure that all accessible surfaces are moistened.
 - During disinfection, mobilize non-rigid components, such as set screws, joints etc., 3 times in each direction to the positive stop.
 - Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml). Ensure that all accessible surfaces are moistened.

- Phase IV**
- Following disinfection, rinse/ rinse through the product completely (all accessible surfaces) 3 times for at least 2 minutes. Use fresh water for every rinse.
 - During the final rinse, mobilize non-rigid components, such as set screws, joints etc., 3 times in each direction to the positive stop.
 - Thoroughly (at least 5 times) rinse through surfaces inaccessible to visual inspection, e.g. in products with hidden crevices or lumens (e.g. empty/working channel) or complex geometries, using a single-use syringe (20 ml).
 - Drain any remaining water fully.

- Phase V**
- Dry the product with a soft, lint-free cloth.
 - Dry inaccessible areas with medical compressed air (P_{max} = 5 bar).

Mechanical cleaning/disinfecting

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Note
Only use cleaning agents that are suitable for this product class (rigid endoscope).

- Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).
- Products with a sheath diameter of ≤4 mm must be cleaned and sterilized in special Aesculap optics trays only.
- Process the product in a cleaning/disinfecting machine. Follow the instructions of the washing machine's manufacturer.
- Avoid abrupt cooling of the product (e.g. in water).

Mechanical alkaline cleaning and thermal disinfecting

Phase	step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Pre-rinse	<25/77	2	D-W	-
II	Cleaning	55/131	10	D-W	Neodisher MediClean forte 0.5 % (5 ml/L) pH > 10
III	Rinse I	>10/50	1	D-W	-
IV	Rinse II	>10/50	1	FD-W	-
V	Thermal disinfection	90/194	5	FD-W	-
VI	Drying	-	-	-	According to disinfectant program

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) at the marked lubrication points, using maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Inspect surfaces, cavities, lumens and openings for visible debris.
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- For products with locking mechanism (e.g. MINOP): Check the locking mechanism for smooth movement.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Sterilization

Note
Alternating sterilization by different sterilization processes can result in damage to materials and adapters of the product.

- The following processes can be used to sterilize the endoscope:
- Steam sterilization
 - Sterrad® sterilization process: using Sterrad® sterilizers Sterrad® 50, Sterrad® 100S, or Sterrad® 200

CAUTION

Damage to the optical system caused by flash sterilization!

- Do not flash sterilize the endoscope.
- Do not expose the endoscope to temperatures above 134 °C.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Do not autoclave damaged products.
- Protect the products against mechanical impacts.

Steam sterilization

Note
The product can be steam-sterilized either in disassembled or in assembled condition.

Note
Only autoclave products carrying the symbol for steam sterilization on the inscription ring 7.

- ▶ Validated sterilization process
 - Steam sterilization using the fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterrad® sterilization process Sterrad® 50, Sterrad® 100S, Sterrad® 200

Note
The Sterrad® sterilization process can cause cosmetic changes to the product. These changes will, however, not affect its functionality.

Note
The product must be disassembled prior to sterilization through the Sterrad® process.

- ▶ Disassemble the product.
- ▶ Only sterilize endoscopes carrying the SDS symbol for low-temperature and plasma sterilization on the inscription ring 7 by means of the Sterrad® sterilization process 50, 100S, 200.
- ▶ Unscrew dismountable parts of the endoscope, see Disassembling.
- ▶ Sterilize using the Sterrad® sterilization process 50, 100S, 200. When doing so, observe the following: Follow the Sterrad® system manufacturer's instructions. Use of a biological indicator is recommended to confirm effective sterilization.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom)			
Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.
- ▶ Prior to dispatching the product for repairs:
 - Have the product cleaned and disinfected or sterilized, and mark it as "disinfected" or "sterilized", respectively.
 - Pack the endoscope in such a way that it will be protected against transport damage.

Service addresses

AESCULAP FLEXIMED GmbH
Robert-Bosch-Strasse 1
79211 Denzlingen / Germany
Phone: +49 (0) 7666 / 9321-0
Fax: +49 (0) 7666 / 9321-580
E-Mail: techserv@aesculap-fleximed.de
www.aesculap-fleximed.de
Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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