Neuro-Patch®

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Neuro-Patch®

Intended Use

Neuro-Patch® is used in neurosurgery as dura mater replacement.

Indications

Use:

- For covering cerebral and cerebellar dura defects
- For cerebral decompression surgery when there is elevated intracranial pressure
- For covering spinal dura defects
- For spinal decompression surgery

Contraindications

Do not use:

- In infected regions
- In open cerebrocranial traumata
- In open spina bifida
- In case of known hypersensitivity against implant materials; for fixation materials please note the corresponding instructions for use
- In any application area that is not mentioned in "Indications"

Risks, Side Effects, Interactions

Potential complications that the manufacturer is currently aware of:

- Infection
- Cerebrospinal fluid leakage
- Adhesion
- Foreign body reaction

Please note:

The points mentioned above include potential clinical consequences.

No risks, side effects and interactions as a result of comorbidities of the patient have been identified.

Application

- Choose implant size suitable for the closure of the defect.
- Cut Neuro-Patch® according to the application situation, in order to embed with as little stress as possible.
- Fix Neuro-Patch® with nonabsorbable suture material (polyester, polypropylene).
- The use of atraumatic round-bodied needles allows suturing without causing great damage to the implant.
- An additional seal with fibrin glue can be used.

Safety and warning notes

a. Users

- For successful use of this product the following knowledge is necessary:
 - Appropriate clinical training
 - Theoretical and practical knowledge about the required surgical technique
 - In vivo characteristics of Neuro-Patch®
- Aesculap is not responsible for complications correlated to:
 - Incorrect indication, incorrect choice of the implant size or incorrect cutting of the implant
 - Other applications and combinations with products than mentioned in this instructions for use
 - Exceeding the limits of the method of treatment or noncompliance with fundamental medical conditions

b. Product

Please note: During the application of Neuro-Patch® in combination with bone cement, chemical damage of the patch material, depending on the application situation, can not be excluded.

c. Shelf life, Storage and Reuse

- Do not use damaged product packages.
- Discard any pieces of the implant that are not used intraoperatively.
- Do not resterilize Neuro-Patch®.
- Do not use Neuro-Patch® if the expiration date has passed.
- Neuro-Patch® is stored at 25 ± 5 °C.

Product characteristics

Neuro-Patch® is a fine fibred, microporous fleece manufactured from highly purified, aliphatic polyester urethane. The porous surface structure of the fleece supports the rapid immigration of fibroblasts. The patch is anchored into the tissue by the collagen that is secreted. Thin neomembranes which surround the fleece used as a dura substitute are produced. Because of their texture, these neomembranes do not form any expectable material-dependent adhesions with the brain. Furthermore, the risk of leucocytes infiltrating the membrane, which could indicate a chronic inflammation reaction or a hyperallergic reaction, could be reduced.

Neuro-Patch® is well tolerated by tissue and is biostable.

Many years of market experience as well as the conducted long-term trials have shown that even after years of implantation no loss of substance and no change in the physical and chemical characteristics occur. Histological examinations confirm these results.

Product range

Neuro-Patch® is supplied EO sterilized for single use only in double peel packages and different sizes.

Package of 1 piece		Package of 2 pieces	
Size	REF	Size	REF
12 x 14 cm	1064002	4 x 5 cm	1064045
6 x 14 cm	1064010	2 x 10 cm	1064053
4 x 10 cm	1064037	1.5 x 3 cm	1064061
2 x 10 cm	1064122		
6 x 8 cm	1064029		
8 x 9 cm	1064020		
5 x 6 cm	1064040		
4 x 5 cm	1064110		
1.5 x 3 cm	1064123		

Symbols used on labeling

	Do not reuse!
	Use until Year, Month, Day
	Date of manufacture
	Temperature limitation
	Sterile unless package is opened or damaged. Method of Sterilization: Ethylene oxide
	CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.
	Manufacturer
	Lot Number
	Attention. See Instructions for Use
	Catalogue number
	Size

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Neuro-Patch®

Product description

Neuro-Patch® is a microporous fabric, manufactured from specially polyester urethane and used as a dura substitute. The fine-fiber microstructure is characterized by intercommunicating pores and many openings on the surface. This allows immigration of the body's own connective tissue cells. The device is provided sterile, for single use only.

Composition

Neuro-Patch® consists of specially purified aliphatic polyester urethane. It is EO sterilized.

Indications for use

Neuro-Patch® is indicated as a dura mater substitute in neurological procedures for soft tissue reconstruction of damaged, impaired, or missing tissue, such as for:

- Closing cerebral and cerebellar dural defects, such as after excision of tumors, operations on the posterior cranial fossa and after the removal of basal meningiomas;
- Decompressive surgery in the presence of increased intracranial pressure;
- For closure of spinal dural defects after removal of spinal tumors and after spinal trauma; and
- Spinal decompressive surgery.

Contraindications

Do not implant in infected areas, open cranial trauma or open spina bifida.

How supplied

The product is supplied sterile, non-pyrogenic and packaged as single or double pieces. The products are available in various dimensions ranging from 1,5 x 3 cm sections to as large as 6 x 14 cm pieces. Specific sizes are noted by catalog number below.

Precautions

Before using Neuro-Patch®, the surgeon should be familiar with the surgical technique, special applications and the properties of Neuro-Patch® in vivo.

Warnings

Neuro-Patch® should not be used together with bone cement as the patch material could be damaged depending on the application.
Neuro-Patch® should only be used if the package is intact. Visually inspect the packaging; do not use product if the package is found opened, punctured, torn or tampered with as sterility may be compromised.
Do not resterilized, as the sterilization process may alter the structure of Neuro-Patch® and thus its properties in vivo.
U.S.A. Federal Law restricts this device to sale by or on the order of a physician.

Possible adverse reactions

Hypersensitivity or an immune response to Neuro-Patch® may occur. The physician should closely monitor the patient for such a reaction and treat accordingly.
Postoperative adjacent tissue adhesions in situ may occur.

Directions for use

The correct size of Neuro-Patch® should be chosen and cut as closely as possible to fit the defect.
The Neuro-Patch® should be fixed with non-absorbable suture material (polyester, polypropylene). To improve suture security, non-absorbable suture material (polyester, polypropylene) is recommended. The use of atraumatic round-bodied needles will allow suturing without significant damage to the implant.

Storage

Store the Neuro-Patch® box in a clean, dry and protected area with protection from extreme changes in temperature and humidity. Storage conditions should not exceed temperature of 25 ± 5 °C. (68 °F – 86 °F)

Product range

Packs of 1 piece		Packs of 2 pieces	
Dimensions	REF	Dimensions	REF
12 x 14 cm	1064002	4 x 5 cm	1064045
6 x 14 cm	1064010	2 x 10 cm	1064053
4 x 10 cm	1064037	1,5 x 3 cm	1064061
6 x 8 cm	1064029		
8 x 9 cm	1064020		
5 x 6 cm	1064040		

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

Symbols used on labelling

	Do not reuse!
	Use until Year + Month + Day
	Date of manufacture
	Temperature limitation
	Sterile unless package is opened or damaged. Method of Sterilization: Ethylene oxide
	CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.
	Manufacturer
	Batch Number
	See Instructions for Use!
	Cat. No.
	Dimensions

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany