
# Aesculap Spine

(GB)

## Instructions for use

S<sup>4</sup> Element Instruments - S<sup>4</sup> Element Rod  
Persuader FW285R

(D)

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S<sup>4</sup> Element Instrumente - S<sup>4</sup> Element  
Stabandrückzange FW285R

(F)

## Mode d'emploi

Instruments S<sup>4</sup> Element - Pousse-tige  
S<sup>4</sup> Element FW285R

(E)

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### Symbols on product and packages

| Symbol                                                                            | Explanation                                                                             |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
|  | Caution, general warning symbol<br>Caution, see documentation supplied with the product |
|  | Date of manufacture                                                                     |

### Intended use

FW285R is intended for the introduction of S<sup>4</sup> Element screws into proper placement, and introduction of S<sup>4</sup> Spinal system rods into proper location.

### Safe handling and preparation

#### CAUTION

**Federal law restricts this device to sale by or on order of a physician!**

- Make certain that the product and its accessories are operated and used only by persons with the requisite training, expertise or experience.
- Read, follow and keep safe the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product thoroughly, either by hand or mechanically, prior to initial sterilization.
- Store the new or unused product in a dry, clean and protected place.
- Prior to each use, visually inspect the product for: loose, bent, broken, cracked, worn, or fractured parts.
- Do not use the product if it is damaged or defective. Immediately set aside the product if it is damaged.
- Immediately replace damaged components with original spare parts.

### Safe operation

**Risk of injury and/or product malfunction!**

- Always carry out a function check before using the product.

### Using the S<sup>4</sup> Element Rod Persuader

**Damage to the screw head/to the instrument and/or difficulty with instrument function and removal due to non-engagement or not fully seating the instrument onto the screw head!**

- Ensure that the positive stop of the S<sup>4</sup> Element Rod Persuader is fully engaged with the S<sup>4</sup> Element Screw head prior to actuating the instrument.
- Place the S<sup>4</sup> Element Rod Persuader over the head of the S<sup>4</sup> Element screw and press down onto the screw body until the internal stop is fully seated on the screw.
- Do not use the S<sup>4</sup> Element Rod Persuader as a counter torque device.
- Removing the S<sup>4</sup> Element Rod Persuader from S<sup>4</sup> Element Screw heads:  
Release front actuation lever to move freely, depress the release thumb lever on the rear of the instrument. The actuation lever should spring forward and the pusher ring on the working end should slide proximally.

## Validated processing procedure

### Note

*Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.*

### Note

*For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.*

### Note

*Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.*

### Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure.

The user/processor is responsible for the validation.

### Note

*Up-to-date information on processing can be found on the Aesculap Extranet at [www.aesculap-extra.net](http://www.aesculap-extra.net)*

## General notes

To avoid unnecessary, excessive contamination of the entire instrument tray, care must be taken at the stage of application that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated surgical residues can make cleaning more difficult or ineffective, and can cause corrosion of stainless steel products. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehyde, alcohol) should be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions, in terms of visual reading and machine-readability, on stainless steel surfaces.

Residues containing chlorine or chlorides, e.g. surgical residues or residues from medicines, saline solutions and service water used for cleaning, disinfecting and sterilizing will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be adequately rinsed with fully desalinated water, followed by thorough drying.

## S<sup>4</sup> Element Instruments – S<sup>4</sup> Element Rod Persuader FW285R

Only tested and approved process chemicals (e.g. VAH/ DGHM or FDA approval or CE mark) that are, according to the chemical manufacturers' recommendations, compatible with the product's materials may be used for processing the product. All application parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times must be observed under all circumstances. Failure to do so can have the following adverse effects:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing process chemicals (e.g. H<sub>2</sub>O<sub>2</sub>), as these can cause bleaching or layer loss.
- Do not use process chemicals that cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.

Further detailed advice on hygienically safe and material-preserving/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org)

- Use suitable cleaning/disinfecting agents if the product is to be put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, thoroughly rinse the product with running water.

### Preparations at the place of use

- Rinse surfaces that cannot be visually inspected, e.g., on products with hidden gaps or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible residues as completely as possible, using a lint-free wet wipe.
- Have the product transferred, in dry condition and stored in a closed disposal container, for cleaning and disinfecting within 6 h.

### Cleaning/Disinfecting

CAUTION

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- Use cleaning/disinfecting agents according to the manufacturer's instructions that
  - are approved for high-grade steel
  - do not attack softeners (e.g., silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Carry out ultrasound cleaning:
  - as an effective mechanical supplement to manual cleaning/disinfecting.
  - as a pre-cleaning procedure for products with encrusted residues, prior to mechanical cleaning/disinfecting.
  - as an integrated mechanical supplement to mechanical cleaning/disinfecting.
  - as an additional process to remove residues left on products after mechanical cleaning/disinfecting.

## Manual cleaning/disinfecting

- Keep working ends open during cleaning.
- Clean products featuring movable links with the links opened or being moved.
- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- Inspect visible surfaces for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

## Manual cleaning with immersion disinfection

| Stage | Step               | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                         |
|-------|--------------------|--------------|------------|--------------|------------------|--------------------------------------------------|
| I     | Cleaning           | RT<br>(cold) | >15        | 2            | D-W              | BBraun Stabimed; aldehyde-, phenol- and QAV-free |
| II    | Intermediate rinse | RT<br>(cold) | 1          | —            | D-W              | —                                                |
| III   | Disinfecting       | RT<br>(cold) | 15         | 2            | D-W              | BBraun Stabimed; aldehyde-, phenol- and QAV-free |
| IV    | Final rinse        | RT<br>(cold) | 1          | —            | FD-W             | —                                                |
| V     | Drying             | RT           | —          | —            | —                | —                                                |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room Temperature

## **S<sup>4</sup> Element Instruments – S<sup>4</sup> Element Rod Persuader FW285R**

### **Stage I**

- Fully immerse the product in the cleaning/disinfecting solution for at least 15 minutes. Make certain all accessible surfaces are moistened by the solution.
- Clean the product under running tap water with a suitable cleaning brush, until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g., in products with hidden gaps, lumens or complex geometries, for at least 1 minute or until no more residues can be removed, using a suitable cleaning brush (TA011944/TE654202, GK469200).
- Mobilize non-rigid components such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning/disinfecting solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

### **Stage II**

- Rinse (through) the entire product (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Allow residual water to drip off sufficiently.

### **Stage III**

- Fully immerse the product in the disinfecting solution. Make certain all accessible surfaces are moistened by the solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during disinfection.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Make certain that all accessible surfaces are moistened.

### **Stage IV**

- Rinse (through) the entire product (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens at least 5 times, using a disposable syringe (20 ml) and an appropriate rinsing adapter.
- Allow residual water to drip off sufficiently.

### **Stage V**

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

## Mechanical cleaning/disinfecting with manual pre-cleaning

### Note

The disinfectant must be of tested and approved effectiveness (e.g., DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

### Note

The disinfectant used for processing must be serviced and checked at regular intervals.

### Note

For thermal disinfection, always use fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least). Ensure that  $A_0$  is  $>3\,000$  for the process.

## Manual pre-cleaning with ultrasound and brush

| Stage | Step                | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                         |
|-------|---------------------|--------------|------------|--------------|------------------|--------------------------------------------------|
| I     | Ultrasonic cleaning | RT<br>(cold) | >15        | 2            | D-W              | BBraun Stabimed; aldehyde-, phenol- and QAV-free |
| II    | Rinse               | RT<br>(cold) | 1          | —            | D-W              | —                                                |

D-W: Drinking water

RT: Room Temperature

### Stage I

- Clean the product in the ultrasonic cleaning bath (frequency 35 kHz) for at least 15 minutes. Make certain that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g., in products with hidden gaps, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush (TA011944/TE654202, GK469200).
- Mobilize non-rigid components such as set screws, links, etc. during cleaning.

- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

### Stage II

- Rinse (through) the entire product (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

## S<sup>4</sup> Element Instruments – S<sup>4</sup> Element Rod Persuader FW285R

### Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfector without ultrasound

- Place the product on a tray suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the special rinsing attachment of the injector carriage.
- Keep working ends open during cleaning.
- Position the product, with all links and joints open, on the tray.

| Stage | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical/Note                                                                       |
|-------|----------------------|--------------|------------|------------------|-------------------------------------------------------------------------------------|
| I     | Pre-rinse            | <25/77       | 3          | D-W              | —                                                                                   |
| II    | Cleaning             | 55/131       | 10         | FD-W             | BBRAUN HELIMATIC CLEANER<br>alkaline with tensides, applica-<br>tion solution 0.5 % |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | —                                                                                   |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | —                                                                                   |
| V     | Drying               | —            | —          | —                | According to disinfecter pro-<br>gram                                               |

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

## Inspection, maintenance and checks

### CAUTION

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- **Prior to function checks, lubricate moving parts (e.g., joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g., for steam sterilization: Aesculap STERILIT® spray JG600 or STERILIT® I drip lubricator JG598).**
- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g., broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that are still not clean.
- Check the product for proper functioning.
- Immediately set aside damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

## Packaging

- Sort the product into its appropriate storage device or put it on a suitable tray.
- Package trays appropriately for the sterilization process (e.g., in Aesculap sterile containers).
- Make certain that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

## Sterilization

- Make certain that the sterilizing agent will be in contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- Validated sterilization process
  - Steam sterilization through fractionated vacuum process
  - Steam sterilizer compliant with DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

## S<sup>4</sup> Element Instruments – S<sup>4</sup> Element Rod Persuader FW285R

### Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of  $10^{-6}$ , Aesculap recommends the following parameters:

| Aesculap Orga Tray/sterile container (perforated bottom)<br>Minimum cycle parameters* |                   |       |                     |
|---------------------------------------------------------------------------------------|-------------------|-------|---------------------|
| Sterilization method                                                                  | Temp.             | Time  | Minimum drying time |
| Pre-vacuum                                                                            | 270 °F–<br>275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

### WARNING for the US market

**If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.**

### Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

### Technical Service

**WARNING**

**Risk of injury and/or product malfunction!**

- **Do not modify the product.**

For service, maintenance or repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

### Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-1602  
Fax: +49 7461 16-5621  
E-mail: [ats@aesculap.de](mailto:ats@aesculap.de)

Or in the US:

Aesculap Implant Systems, LLC  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood, MO 63042  
Aesculap Repair Hotline  
Phone: +1 800 214-3392  
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

## **Disposal**

- Always adhere to national regulations when disposing of or recycling the product, its components and its packaging!

## **Distributor in the US/Contact in Canada for product information and complaints**

Aesculap Implant Systems, LLC  
3773 Corporate Parkway  
Center Valley, PA 18034  
USA

## S<sup>4</sup> Element Instrumente – S<sup>4</sup> Element Stabandrückzange FW285R

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#### Aesculap AG