Aesculap Spine

-
- GB** Instructions for use/Technical description
 - USA** S⁴ Element Instruments
 - D** Gebrauchsanweisung/Technische Beschreibung
S⁴ Element-Instrumente
 - F** Mode d'emploi/Description technique
Instruments S⁴ Element
 - E** Instrucciones de manejo/Descripción técnica
Instrumental S⁴ Element
 - I** Istruzioni per l'uso/Descrizione tecnica
Strumenti S⁴ Element


Fig. 1


Fig. 2


Fig. 3

S⁴ Element Instruments

Legend

- 1 S⁴ Element Monoaxial Screwdriver (FW276R)
- 2 S⁴ Element Polyaxial Screwdriver (FW277R)
- 3 S⁴ Element Countertorque Instrument (FW283R)

Symbols on product and packages

	Caution: General warning symbol Caution: See documentation supplied with the product
	Date of manufacture

Intended use

The S⁴ Element Instruments are used for the implantation of implants of the S⁴ Element System.

Safe handling and preparation

Risk of injury caused by incorrect operation of the product!

- Adhere to instructions for use of all medical devices.

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, expertise or experience.
- Read, follow and keep safe the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product thoroughly, either by hand or mechanically, prior to initial sterilization.
- Store the new or unused product in a dry, clean and protected place.

- Prior to each use, visually inspect the product for: loose, bent, broken, cracked, worn, or fractured parts.
- Do not use the product if it is damaged or defective. Immediately set aside the product if it is damaged.
- Immediately replace damaged components with original spare parts.
- To avoid damage to the working end: Apply appropriate care when introducing the product through the working channel (e.g. trocar).

Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

- Bending or levering forces should not be made with the screwdrivers FW277R or FW276R.
- The threaded shaft of FW277R or FW276R should not be over tightened or tightened by mechanical means.

Risk of implant damages due to use of inappropriate screwdrivers!

- Only use monoaxial screwdriver FW276R with monoaxial screws.
- Only use polyaxial screwdriver FW277R with polyaxial screws.
- The screwdrivers FW277R or FW276R should not be used with torque wrenches or for interface with set screws.
- The tips of the screwdrivers FW277R and FW276R are fragile and should not be used as a lever.
- Should any part of an instrument become disassociated from the main instrument, it should be removed from the patient as it is not implantable.

Validated processing procedure

Note

Adhere to national law and regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Please note that successful processing of this medical product can be ensured only through a validated processing procedure.

Responsibility for the validation lies with the user/processor.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

To avoid unnecessary, excessive contamination of the entire instrument tray, care must be taken at the stage of application that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated surgical residues can make cleaning more difficult or ineffective, and can cause corrosion of stainless steel products. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehyde, alcohol) should be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions, in terms of visual reading and machine-readability, on stainless steel surfaces.

Residues containing chlorine or chlorides, e.g. surgical residues or residues from medicines, saline solutions and service water used for cleaning, disinfecting and sterilizing will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be adequately rinsed with fully desalinated water, followed by thorough drying.

Only tested and approved process chemicals (e.g. VAH/DGHM or FDA approval or CE mark) that are, according to the chemical manufacturers' recommendations, compatible with the product's materials may be used for processing the product. All application parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times must be observed under all circumstances. Failure to do so can have the following adverse effects:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing process chemicals (e.g. H₂O₂), as these can cause bleaching or layer loss.

S⁴ Element Instruments

- Do not use process chemicals that cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.
- Further detailed advice on hygienically safe and material-preserving/value-preserving reprocessing can be found at www.a-k-i.org
- Use suitable cleaning/disinfecting agents if the product is to be put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, thoroughly rinse the product with running water.

Preparations at the place of use

- Open jointed products.
- Remove sealing caps from Luer lock connectors.
- Open all valves/faucets.
- Rinse surfaces that cannot be visually inspected, e.g. on products with hidden gaps or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible residues as completely as possible, using a lint-free wet wipe.
- Have the product transferred, in dry condition and stored in a closed disposal container, for cleaning and disinfecting within 6 h.

Cleaning/Disinfecting

CAUTION

Damage to the product caused by inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning/disinfecting agents, according to the manufacturer's instructions,
 - that are approved for plastics,
 - do not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, prior to mechanical cleaning/disinfecting.
 - as an integrated mechanical supplement to mechanical cleaning/disinfecting.
 - as an additional process to remove residues left on products after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products in mechanical washers/disinfectors, provided the products can be securely fixed in the respective washing/disinfecting machine or on storage devices in such a way that they will be cleaned thoroughly.

Manual cleaning/disinfecting

- Clean products featuring movable links with the links open or being moved.
- Inspect visible surfaces for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfection (FW283R)

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	0.5	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make certain all accessible surfaces are moistened by the solution.
- Clean the product under running tap water with a suitable cleaning brush, until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden gaps, lumens or complex geometries, for at least 1 minute or until no more residues can be removed, using a suitable cleaning brush (TA-No. 011944/TE654202, GK469200).
- Mobilize non-rigid components such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning/disinfecting solution, using a disposable syringe (20 ml).

- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

S⁴ Element Instruments

Stage II

- Rinse (through) the entire product (all accessible surfaces) under running water.
- Allow residual water to drip off sufficiently.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain all accessible surfaces are moistened by the solution.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Make certain that all accessible surfaces are moistened.

Stage IV

- Rinse (through) the entire product (all accessible surfaces).
- Rinse lumens at least 5 times, using a disposable syringe (20 ml) and an appropriate rinsing adapter.
- Allow residual water to drip off sufficiently.

Stage V

- Dry the product with lint-free tissue or medical-quality compressed air.

Manual ultrasonic cleaning and immersion disinfection (FW276R, FW277R)

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	5	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Mount the jaws protection on the product.
- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened and acoustic shadows are avoided.
- Remove the jaws protection.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden gaps, lumens or complex geometries, for at least 1 minute or until no more residues can be removed, using a suitable cleaning brush (TA-No. 011944/TE654202, GK469200).
- Mobilize non-rigid components such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).

- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse (through) the entire product (all accessible surfaces) under running water.
- Allow residual water to drip off sufficiently.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain all accessible surfaces are moistened by the solution.
- Rinse lumens at least 5 times, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Make certain all accessible surfaces are moistened by the solution.

S⁴ Element Instruments

Stage IV

- Rinse (through) the entire product (all accessible surfaces) under running water.
- Rinse lumens at least 5 times, using a disposable syringe (20 ml) and an appropriate rinsing adapter.
- Allow residual water to drip off sufficiently.

Stage V

- Rinse (through) the entire product (all accessible surfaces) under running water.

Mechanical cleaning/disinfecting with manual pre-cleaning

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that A_0 is $>3\,000$ for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with brush (Countertorque Instrument FW283R)

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Rinse	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make certain all accessible surfaces are moistened by the solution.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden gaps, lumens or complex geometries, for at least 1 minute or until no more residues can be removed, using a suitable cleaning brush (TA-No. 011944/TE654202, GK469200).

- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse (through) the entire product (all accessible surfaces) under running water.

S⁴ Element Instruments

Manual pre-cleaning with ultrasound and brush (Screwdrivers FW276R/FW277R)

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Irrigation	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Clean the product in the ultrasonic cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened and acoustic shadows are avoided.
- Remove the jaws protection.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden gaps, lumens or complex geometries, for at least 1 minute or until no more residues can be removed, using a suitable cleaning brush (TA-No. 011944/TE654202, GK469200).
- Mobilize non-rigid components such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse (through) the entire product (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfector without ultrasound

- Place the product on a tray suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the special rinsing attachment of the injector carriage.
- Keep working ends open during cleaning.
- Position the product, with all links and joints open, on the tray.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER alkaline with tensides; application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfector program

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

S⁴ Element Instruments

Packaging

- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that any existing sharp edges are protected.
- Pack trays appropriately for the sterilization process (e.g. in Aesculap SterilContainer).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

- Check that the sterilizing agent will be in contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterilcontainer (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F– 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap Sterilcontainer cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Disposal

- Always adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Technical Service

Risk of injury and/or malfunction!

- **Do not modify the product.**

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.

3773 Corporate Parkway

Center Valley, PA 18034

USA

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1602

Fax: +49 7461 16-5621

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood, MO 63042

Aesculap Repair Hotline

Phone: +1 800 214-3392

Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

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Aesculap AG