



Aesculap Spine

GB	Instructions for use/Technical description
USA	S⁴ instruments
D	Gebrauchsanweisung/Technische Beschreibung S⁴-Instrumente
F	Mode d'emploi/Description technique Instruments S⁴
E	Instrucciones de manejo/Descripción técnica Instrumental S⁴
I	Istruzioni per l'uso/Descrizione tecnica Strumenti S⁴
P	Instruções de utilização/Descrição técnica Instrumentos S⁴
NL	Gebruiksaanwijzing/Technische beschrijving S⁴-instrumenten
CZ	Návod k použití/Technický popis S⁴-instrumenty


7

8

9

10

11

12

13


14

15

16

17

18

19

20

S⁴ instruments

Legend

- 1 Set screw starter FW177R
- 2 Counter FW178R (L handle)
- 3 Rod persuader FW208R
- 4 Reduction lever (consisting of threaded tube FW233R, screw driver FW234R and outer sleeve FW235R)
- 5 Counter FW236R
- 6 Wrench FW237R
- 7 Distractor (consisting of distractor FW238R and two distraction blades FW239R)
- 8 Rod insertion instruments FW240R
- 9 Spindle distractor FW241R (consisting of distractor and extension)
- 10 Rod length measuring instrument FW242R
- 11 Slotted hammer FW243R
- 12 Tissue protection sleeve FW244P
- 13 K-wire FW257S
- 14 K-wire aiming device FW258R
- 15 Screw driver monoaxial FW262R
- 16 Bone probe cannulated FW263R
- 17 Screw taps FW264R, FW265R, FW266R, FW267R, FW268R
- 18 Screwdriver polyaxial FW270R
- 19 Trocar FW271R
- 20 Tube sleeve FW272R

Symbols on product

	Sterilization with steam or dry heat
	Not for reuse in intended applications as defined by the manufacturer
	Attention, see instructions for use

Intended use

The S⁴ instruments are used for the implantation of the S⁴ implants with the following systems:

- S⁴ Spinal System
- S⁴ CS - Cannulated System
- S⁴ FRI

The S⁴ fracture reduction instruments are used for reduction and fracture stabilization with the S⁴ implants of the S⁴ Spinal System or the S⁴ CS - Cannulated System.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working tip: Exercise caution when passing products through the working channel (e.g. trocars).
- Read, follow and keep the surgical technics 026102/026402/026502/026702.

Safe operation

DANGER

Danger of death of the patient caused by pushing the guide wire into the aorta!

- When inserting the S⁴ implants or S⁴ instruments with the guide wire, always observe the markings on the guide wire.
- Hold the guide wire with a forceps.
- Remove the guide wire at the right moment.

WARNING

Trocar penetrating too deeply because of missing depth stop!

- Use trocar FW271R only in combination with aiming device FW258R.

S⁴ instruments

Disassembling


Fig. 1

Legend

A Button

- To dismount reduction lever **4**, see Fig. 1:
Firmly press button **A** and remove screwdriver FW234R.


Fig. 2

Legend

B Slide

- To dismount distractor **7** with the distraction blades, see Fig. 2:
Actuate slide **B** and remove distraction blade FW239R.


Fig. 3

- To dismount rod insertion instrument **8**, see Fig. 3:
unscrew the inner part.


Fig. 4

- Dismount spindle distractor **9**, see Fig. 4.


Fig. 5

- Dismount rod length measuring instrument **10**, see Fig. 5.

Assembling


Fig. 6

- Mount reduction lever **4**, see Fig. 6.


Fig. 7

- Mount distractor **7** with the distraction blades, see Fig. 7.


Fig. 8

- Mount rod insertion instrument **8**, see Fig. 8.

S⁴ instruments


Fig. 9

- Mount spindle distractor 9, see Fig. 9.


Fig. 10

- Mount rod length measuring instrument 10, see Fig. 10.

Validated processing procedure

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured through a validated processing procedure. The user/processor is responsible for the validation. Due to process tolerances, manufacturer's specifications can only serve as guide values for the processing procedures applied by the individual user.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures > 45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) should be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading of the laser engraving on stainless steel.

Residues containing chlorine or chlorides, as contained e.g. in surgical residues, tinctures, medicines, saline solutions, service water used for cleaning, and in cleaning/disinfecting agents, will cause corrosion damage (pitting, stress corrosion) and thereby destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals recommended by the chemicals' manufacturers as effective cleaning/disinfecting agents compatible with the materials used in the respective product. Strictly observe all application instructions, e.g. regarding temperatures, concentrations, exposure times, etc. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes or
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing process chemicals, which could cause bleaching/layer loss of the product.
- Do not use process chemicals that cause stress cracking of plastics such as PPSU, or attack and cause brittleness e.g. of silicone softeners.

- Clean the product immediately after use.
Encrusted debris on HF instruments is softened most effectively and conservatively by a 5 minutes immersion treatment in a 3 % H₂O₂ solution. Encrustations can be removed by hand, with a soft brush and/or in an ultrasound bath. After that treatment, proceed through the normal steps of processing.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications, Red Brochure – Proper maintenance of instruments.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open jointed products.
- Remove the sealing cap from the Luer lock connector.
- Remove visible residues as completely as possible, using a lint-free, single-use cleaning tissue.

Preparation prior to cleaning

- Have the product dry in a disposal container and ready for immediate cleaning and disinfecting within 30 min after use.

S⁴ instruments

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to manufacturer's instructions. The cleaning and disinfecting agent must
 - be approved for plastics and high-grade steel,
 - not attack softeners (e.g. silicone).
 - Observe specifications regarding concentration, temperature and exposure time.
 - Do not exceed the maximum allowable cleaning temperature of 93 °C.
-
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Manual cleaning/disinfecting

- Clean hinged or jointed products in open and closed positions.
- Check visible surface for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfecting and cleaning with brush

FW233R, FW234R, FW235R, FW238R, FW239R, FW241R, FW244P, FW258R, FW262R, FW263R, FW264R, FW265R, FW266R, FW267R, FW268R, FW270R, FW271R, FW272R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Final rinse	RT (cold)	0.5	-	FD-W	-
IV	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product under running tap water, if necessary with a suitable cleaning brush (TA no. 011944/TE654202, GK469200), until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives, which would damage the surfaces, for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage IV

- Dry the product completely with compressed air.

S⁴ instruments

Manual cleaning with immersion disinfecting and cleaning with brush and ultrasound

FW177R, FW208R, FW240R, FW242R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting ultrasound cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	1st Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
IV	2nd Intermediate rinse	RT	1	-	D-W	-
V	Intermediate drying	-	-	-	-	-
VI	Final rinse	RT (cold)	0.5	-	FD-W	-
VII	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush (TA no. 011944/TE654202, GK469200) until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.

- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use metal brushes or other abrasives for cleaning the product.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.

Stage V

- Dry the product with compressed air.

Stage VI

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage VII

- Dry the product completely with compressed air.

S⁴ instruments

Mechanical cleaning/disinfecting

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.
- Store any products with hinges or joints on the tray in such a way that the joints are kept open.

Mechanical alkaline cleaning and thermal disinfecting

FW234R, FW238R, FW239R, FW241R, FW271R

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	< 25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	- Concentrate, alkaline: pH = 10.9 < 5 % anionic tensides - 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD-W	- Concentrate, acid: pH = 2.6 Basis: Citric acid - 0.2 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD-W	-
V	Thermal disinfecting	94/201	10	FD-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

Mechanical cleaning/disinfecting with manual pre-cleaning

Manual pre-cleaning with brush

FW178R, FW208R, FW233R, FW235R, FW244P, FW258R, FW262R, FW263R, FW264R, FW265R, FW266R, FW267R, FW268R, FW270R, FW272R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush (TA no. 011944/TE654202, GK469200) until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives that would damage the surfaces for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

S⁴ instruments

Manual pre-cleaning with ultrasound

FW177R, FW240R, FW242R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting ultrasound cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Mount mouth protection on the product.
- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make certain that all accessible surfaces are moistened.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.
- Store products with hinges or joints on the tray in such a way that the joints are kept open.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	< 25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	- Concentrate, alkaline: pH = 10.9 < 5 % anionic tensides - 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD-W	- Concentrate, acid: pH = 2.6 Basis: Citric acid - 0.1 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD-W	-
V	Thermal disinfecting	94/201	10	FD-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

S⁴ instruments

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts (e.g. joints and links) with a maintenance oil suitable for the respective sterilization process (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Assemble the separable product, see Assembling.
- Mount reduction lever 4, see Fig. 6.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside the product if it is damaged.

Packaging

- Appropriately protect products with fine working tips.
- Secure product with retractor in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

Note

The product can be sterilized either in disassembled or in assembled condition.

Note

To prevent breakage due to stress crack corrosion, sterilize the instruments with the locking mechanism open, or locked on the first ratchet tooth.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilization acc. to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST 46-1993, validated acc. to EN ISO 17665 or EN 554/ISO 13683
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270–275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store processed products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

Technical Service

For service, maintenance and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of warranty rights and applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems Inc.
Attn: Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Designation	Art.-no.
Replacement nut (for threaded tube FW233R, distractor FW238R and spindle distractor FW241R)	FW233201

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG