

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-----|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| en  | <b>Instructions for use/Technical description</b><br>S <sup>4</sup> ® Long Tab Instruments                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| USA | Note for U.S. users<br>This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapImplantsystems.com">www.aesculapImplantsystems.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost. |
| de  | <b>Gebrauchsanweisung/Technische Beschreibung</b><br>S <sup>4</sup> ® Long Tab Instrumente                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| fr  | <b>Mode d'emploi/Description technique</b><br>Instruments S <sup>4</sup> ® Long Tab                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| es  | <b>Instrucciones de manejo/Descripción técnica</b><br>Instrumental S <sup>4</sup> ® Long Tab                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| it  | <b>Istruzioni per l'uso/Descrizione tecnica</b><br>Strumenti S <sup>4</sup> ® Long Tab                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| pt  | <b>Instruções de utilização/Descrição técnica</b><br>Instrumentos S <sup>4</sup> ® Long Tab                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| nl  | <b>Gebruiksaanwijzing/Technische beschrijving</b><br>S <sup>4</sup> ® Long Tab-instrumenten                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| sv  | <b>Bruksanvisning/Teknisk beskrivning</b><br>S <sup>4</sup> ® Long Tab-instrument                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ru  | <b>Инструкция по применению/Техническое описание</b><br>Инструменты S <sup>4</sup> ® Long Tab                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| cs  | <b>Návod k použití/Technický popis</b><br>Nástroje S <sup>4</sup> ® Long Tab                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| pl  | <b>Instrukcja użytkowania/Opis techniczny</b><br>Instrumenty S <sup>4</sup> ® Long Tab                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| sk  | <b>Návod na použitie/Technický opis</b><br>S <sup>4</sup> ® Long Tab nástroje                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| tr  | <b>Kullanım Kılavuzu/Teknik açıklama</b><br>S <sup>4</sup> ® Long Tab aletleri                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap® – a B. Braun brand

TA014587    2020-03    V6    Change No. 62188

# Aesculap®

## S<sup>4</sup>® Long Tab Instruments

### Legend

- 1 S<sup>4</sup>® Long Tab Counterholding instrument FW728R
- 2 S<sup>4</sup>® Long Tab Rod length measuring instrument FW729R
- 3 S<sup>4</sup>® Long Tab Screwdriver FW731R
- 4 S<sup>4</sup>® Long Tab Shaft for torque wrench FW732R
- 5 S<sup>4</sup>® Long Tab Insertion instrument for lock nut FW733R
- 6 S<sup>4</sup>® Long Tab FRI outer sleeve FW737R
- 7 S<sup>4</sup>® Long Tab Tab breaker FW738R
- 8 S<sup>4</sup>® Long Tab Tab protection ring FW739R
- A Pedicle screw in S<sup>4</sup>® Long Tab FRI outer sleeve FW737R
- a Movable measuring arms
- b Circumferential marking
- c Circumferential marking
- d Circumferential marking
- e Circumferential marking
- f Inspection window with marking
- g Scale
- h Body of the pedicle screw

### 1. About this document

*Note*  
General risk factors associated with surgical procedures are not described in these instructions for use.

#### 1.1 Scope

These instructions for use apply for the following products:

| Art. no. | Designation                                                 |
|----------|-------------------------------------------------------------|
| FW728R   | S <sup>4</sup> ® Long Tab Counterholding instrument         |
| FW729R   | 2 S <sup>4</sup> ® Long Tab Rod length measuring instrument |
| FW731R   | S <sup>4</sup> ® Long Tab Screwdriver                       |
| FW732R   | S <sup>4</sup> ® Long Tab Shaft for torque wrench           |
| FW733R   | S <sup>4</sup> ® Long Tab Insertion instrument for lock nut |
| FW737R   | S <sup>4</sup> ® Long Tab FRI outer sleeve                  |
| FW738R   | S <sup>4</sup> ® Long Tab Tab breaker                       |
| FW739R   | S <sup>4</sup> ® Long Tab Tab protection ring               |

*Note*  
The applicable CE mark for the product can be found on the label or packaging of the product.

- For article specific instructions for use and material compatibility and lifetime information, see B. Braun eIFU at eifu.bbraun.com

#### 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

**⚠ DANGER**  
Indicates a potential danger. If it is not prevented, death or severe injuries may result.

**⚠ WARNING**  
Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

**⚠ CAUTION**  
Indicates a possible threat of material damage. If not avoided, the product may be damaged.

### 2. Clinical use

#### 2.1 Areas of use and limitations of use

##### 2.1.1 Intended use

The S<sup>4</sup>® Long Tab instruments are used to complement the S<sup>4</sup>® Spinal System (instruments and implants) when implanting the S<sup>4</sup>® Long Tab screws.

##### 2.1.2 Indications

*Note*  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

##### 2.1.3 Contraindications

No known contraindications.

### 2.2 Safety information

#### 2.2.1 Clinical user

- General safety information**  
To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:
- Use the product only according to these instructions for use.
  - Follow the safety and maintenance instructions.
  - Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
  - Store any new or unused products in a dry, clean, and safe place.
  - Prior to use, check that the product is in good working order.
  - Keep the instructions for use accessible for the user.

*Note*  
The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

**Notes on surgical procedures**  
It is the user's responsibility to ensure that the surgical procedure is performed correctly. Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product. The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

#### 2.2.2 Product specific safety information

- Read, observe and store the surgical manual S<sup>4</sup>® Long Tab 082302 .
- Read, observe and store the operating instructions of the S<sup>4</sup>® Spinal System (instruments and implants), including:
  - S<sup>4</sup>® Spinal System – Lumbar/Deformity TA011187
  - S<sup>4</sup>® MIS Instruments TA014087
  - S<sup>4</sup>® Instruments TA012384
  - S<sup>4</sup>® MIS Percutaneous instruments TA013496

#### 2.2.3 Sterility

- The product is delivered in an unsterile condition.
- Clean the new product after removing its transport packaging and prior to its initial sterilization.

### 2.3 Application

**⚠ DANGER**  
Pushing the K wire forward into the aorta risks the life of the patient!  
► When inserting S<sup>4</sup>® implants or S<sup>4</sup>® instruments using the guidewire: Observe the markings on the guide wire.  
► Use X-ray control to check the correct positioning.  
► Use forceps to hold the guide wire.  
► Remove the guide wire in time.

**⚠ WARNING**  
Risk of injury and/or malfunction!  
► Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.  
► Always carry out a function test prior to each use of the product.

**⚠ WARNING**  
Danger of severe injuries from incorrectly positioning implants!  
The depth markings and scales on the instruments are for orientation only.  
► Observe the relevant reference points of the markings.  
► Check the correct position of the implants using x-ray.  
► Check the proper rod length using x-ray.

**⚠ WARNING**  
Risk of injury and/or malfunction!  
► When screwing the screws into the bone, apply only axial force. Avoid any transmission of lateral forces when turning the screw.  
► Do not use the insertion instrument for lock nut FW733R to loosen lock nuts.  
► To loosen lock nuts: Use the shaft of the torque wrench FW732R.  
► To tighten the lock nut: Always use the shaft of the torque wrench FW732R in conjunction with the intended torque wrench in combination with the corresponding counterholding instrument.

**⚠ CAUTION**  
The screw flanks may prematurely break off if the protection ring is removed!  
► Keep the tab protection ring FW739R mounted during all surgical steps until the screw flanks are broken off, see surgical manual S<sup>4</sup>® Long Tab.

- For easy handling, some S<sup>4</sup>® Long Tab instruments have markings or scales:
- S<sup>4</sup>® Long Tab Screwdriver FW731R: The circumferential marking **b** indicates the end of the protection ring if the screwdrivers, the protection ring FW739R and the S<sup>4</sup>® Long Tab screw are properly installed.
  - S<sup>4</sup>® Long Tab Shaft for torque wrench FW732R: The circumferential marking **c** indicates the end of the outer sleeve FW737R/counterholder FW728R if the lock nut has been tightened.
  - S<sup>4</sup>® Long Tab Insertion instrument for lock nut FW733R:
    - The circumferential marking **d** indicates the end of the outer sleeve FW737R/counterholder FW728R if the lock nut has been screwed in.
    - The circumferential marking **e** indicates the end of the protection ring FW739R, if the lock nut has been screwed in.
  - S<sup>4</sup>® Long Tab FRI outer sleeve FW737R:
    - The scale **g** on the outer sleeve refers to the body of the pedicle screw **h**, see Fig. A. When the edge of the screw is in line with the marking **f** in the inspection window, the screw is fully inserted in the outer sleeve.

### 3. Validated reprocessing procedure

#### 3.1 General safety information

*Note*  
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

*Note*  
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

*Note*  
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

*Note*  
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

*Note*  
If there is no final sterilization, then a virucidal disinfectant must be used.

*Note*  
For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com  
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

#### 3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

3.3 Reusable products

- No maximum number of applications and processing cycles has been specified for the product.
- The service life of the product is limited by damage, normal wear, the type and duration of usage, as well as the handling, storage and transportation of the product.
- Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product.

3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Preparing for cleaning

- Do not disassemble any anchoring screws/bolts that permanently connect the various components.
- Disassemble the product prior to cleaning, see Disassembly.

3.6 Disassembly

- Remove the movable measuring arm **a** from the rod length measuring instrument.

3.7 Cleaning/Disinfection

3.7.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning agents and disinfectants in accordance with the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable disinfection temperature of 95 °C.
- Use suitable cleaning/disinfectants for wet disposal. To prevent foam from forming and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly under running water

3.7.2 Validated cleaning and disinfection procedure

| Validated procedure                                                                                                                                 | Specific requirements                                                                                                                                                                                                                                                                                                                                     | Reference                                                                                                                                                                                          |
|-----------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection<br>■ FW728R to FW729R<br>■ FW731R to FW733R<br>■ FW737R to FW739R                                       | ■ Suitable cleaning brush<br>■ 20 ml disposable syringe<br>■ Clean products having movable hinges in the open position or while moving the joints.<br>■ Drying phase: Use a lint-free cloth or medical compressed air                                                                                                                                     | Chapter Manual cleaning/disinfection and subsection:<br>■ Chapter Manual cleaning with immersion disinfection                                                                                      |
| Alkaline machine cleaning and thermal disinfection<br>■ FW729R<br>■ FW732R<br>■ FW733R<br>■ FW738R                                                  | ■ Place product on a screen basket suitable for cleaning (make sure all areas will be reached by water jets).<br>■ Connect components with lumens and channels directly to the rinsing port of the injector carriage.<br>■ Open the links and joints before placing the product on the screen basket.                                                     | Chapter Machine cleaning/disinfection and subsections:<br>■ Chapter Mechanical alkaline cleaning and thermal disinfection                                                                          |
| Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection<br>■ FW728R<br>■ FW731R<br>■ FW737R<br>■ FW739R | ■ Suitable cleaning brush<br>■ 20 ml disposable syringe<br>■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).<br>■ Connect the lumens and channels of individual parts directly to the injector unit's special flushing connectors.<br>■ Open the links and joints before placing the product on the screen basket. | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection:<br>■ Chapter Manual pre-cleaning with a brush<br>■ Chapter Mechanical alkaline cleaning and thermal disinfecting |

3.8 Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfection process if necessary.

3.8.1 Manual cleaning with immersion disinfection

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|---------------|----------------------------------------------------------------|
| I     | Disinfecting cleaning | RT (cold)    | >15        | 2            | D-W           | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Intermediate rinse    | RT (cold)    | 1          | -            | D-W           | -                                                              |
| III   | Disinfection          | RT (cold)    | 5          | 2            | D-W           | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| IV    | Final rinse           | RT (cold)    | 1          | -            | FD-W          | -                                                              |
| V     | Drying                | RT           | -          | -            | -             | -                                                              |

D-W: Drinking water  
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
RT: Room temperature  
\*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.9 Machine cleaning/disinfection

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.9.1 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water quality | Chemical/Note                                                                                                |
|-------|----------------------|--------------|------------|---------------|--------------------------------------------------------------------------------------------------------------|
| I     | Pre-rinse            | <25/77       | 3          | D-W           | -                                                                                                            |
| II    | Cleaning             | 55/131       | 10         | FD-W          | ■ Concentrate, alkaline:<br>- pH = 13<br>- <5 % anionic surfactant<br>■ 0.5 % working solution<br>- pH = 11* |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W          | -                                                                                                            |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W          | -                                                                                                            |
| V     | Drying               | -            | -          | -             | According to the program for cleaning and disinfection device                                                |

D-W: Drinking water  
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.10 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.10.1 Manual pre-cleaning with a brush

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water quality | Chemical                                                       |
|-------|-----------------------|--------------|------------|--------------|---------------|----------------------------------------------------------------|
| I     | Disinfectant cleaning | RT (cold)    | >15        | 2            | D-W           | Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9* |
| II    | Rinsing               | RT (cold)    | 1          | -            | D-W           | -                                                              |

D-W: Drinking water  
RT: Room temperature  
\*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.10.2 Mechanical alkaline cleaning and thermal disinfecting  
Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                 | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                                                                                                                                               |
|-------|----------------------|--------------|------------|------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| I     | Prerinse             | <25/77       | 3          | D-W              | -                                                                                                                                                                                                                                      |
| II    | Cleaning             | 55/131       | 10         | FD-W             | <div><div>■</div>Concentrate, alkaline:<ul style="list-style-type: none"><li>- pH ~ 13</li><li>- &lt;5 % anionic surfactant</li></ul><div>■</div>working solution 0.5%<ul style="list-style-type: none"><li>- pH = 11*</li></ul></div> |
| III   | Intermediate rinse   | >10/50       | 1          | FD-W             | -                                                                                                                                                                                                                                      |
| IV    | Thermal disinfecting | 90/194       | 5          | FD-W             | -                                                                                                                                                                                                                                      |
| V     | Drying               | -            | -          | -                | In accordance with the program for the cleaning and disinfecting machine                                                                                                                                                               |

D-W: Drinking water  
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.11 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.11.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the surfaces for rough spots.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.11.2 Functional test

⚠ CAUTION

The product may become damaged (metal pitting/fretting corrosion) if insufficiently oiled!

- ▶ Oil any moving parts (e.g. joints, spool parts and threaded rods) prior to the functional test using oil suitable for the applied sterilization process (e.g. STERILIT® I oil spray JG600 for steam sterilization or STERILIT® I drip-feed lubricator JG598).
- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working property (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.12 Assembly

- ▶ Slide the movable measuring arm a onto the rod length measuring instrument.

3.13 Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.14 Steam sterilization

Note

FW729R: The product may only be sterilized in a disassembled state.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
  - Disassemble the product
  - Steam sterilization in a fractionated vacuum process
  - Steam sterilizer in accordance with DIN EN285 and validated in accordance with DIN EN ISO 17665
  - Sterilization in a fractionated vacuum process at 134°C, holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.15 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4. Technical service

⚠ CAUTION

Modifications to medical equipment may result in the voiding of any guarantee/warranty claims and any approvals.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses  
Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 7461 95-1601  
Fax: +49 7461 16-2887  
E-Mail: ats@aesculap.de  
Other service addresses can be obtained from the address indicated above.

5. Disposal

⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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