Aesculap Spine

GB
USA

Instructions for use

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www.aesculap.de

TA-Nr. 012865 05/10 Änd.-Nr. 41722

CE marking according to directive 93/42/EEC

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S⁴ Spinal System – augmentation screw

Note
Generally, the S⁴ Spinal System augmentation screw is covered by the instructions for use S⁴ Spinal System – Lumbar/Deformity TA011187! The present information on the augmentation screw is intended to complement the respective information contained in the instructions for use S⁴ Spinal System – Lumbar/Deformity.

Symbols on product

Symbol	Explanation
	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Adhere to instructions for use

Intended use

The S⁴ Spinal System implants are used for dorsal monosegmental and multi-segmental stabilization of the lumbar and thoracic spine.
The injection cannula is inserted in the S⁴ Spinal System augmentation screw to apply bone cement.

System components

- The S⁴ Spinal System augmentation screw comprises:
- Monoaxial/polyaxial screws (augmentation screw), supplied in sterile condition
 - Cement injection cannula (supplied in sterile condition), see TA013132

Note
There are special S⁴ instruments provided for the implantation of these system components and for the augmentation, distraction, compression and reduction of the lumbar and thoracic spine.

Materials

The materials used in the implants are listed on the packaging:

Materials
ISOTAN®: titanium forged alloy Ti6Al4V acc. to ISO 5832-3
ISOTAN®: pure titanium acc. to ISO 5832-2

The titanium implants are coated with a colored oxide. Slight color changes are possible, but do not affect the quality of the implant.
ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

- Surgically installed implants serve to support normal healing processes. They are not supposed to replace normal body structures or to support permanently loads occurring in cases where healing did not take place.
Use for:
- Fractures
 - Pseudarthrosis or delayed healing
 - Dislocation
 - Degenerative instability
 - Post-discectomy syndrome
 - Spondylolisthesis
 - Kyphosis
 - Degenerative scoliosis
 - Stenosis
 - Post-trauma instability
 - Tumors

Risk of screw breakage caused by inappropriate use of pedicle screws in cases of spondylolisthesis!
➤ Provide additional stabilization by intercorporeal fusion (ALIF, PLIF or TLIF).

Contraindications

- Do not use in the presence of:
- Fever
 - Acute or chronic vertebral infections in the spine, local or systemic
 - Severe osteoporosis or osteopenia:
 - The operating surgeon has to weigh possible risks affecting the success of treatment with the S⁴ Spinal System augmentation screw.
 - In cases of foreseeable reduced stability of the standard S⁴ Spinal System augmentation screw anchoring, additional stability can be achieved by fixation with bone cement.
 - Pregnancy
 - Infection
 - systemic
 - in the spine
 - local
 - Medical or surgical conditions that could preclude successful implantation
 - Mental illness
 - Dependency on pharmaceutical drugs, drug abuse, or alcoholism
 - Adiposity
 - Severely damaged bone structures that could prevent stable implantation of implant components
 - Neuromuscular disorders or illnesses
 - Bone tumors in the region of implant fixation
 - Wound healing disorders
 - Inadequate patient compliance
 - Foreign body sensitivity to the implant materials
 - Cases not listed under indications

Side effects and adverse interactions

- Implant failure due to overloading:
 - Bending
 - Loosening
 - Breakage
- Inadequate fixation
- Failed or delayed fusion
- Cement leakages (cement emerging from the vertebral body)
- Infection
- Fractured vertebral body or bodies
- Injuries to:
 - nerve roots
 - spinal cord
 - blood vessels
 - organs
- Changes of the normal spine lordosis
- Impairment of the gastrointestinal, urological and/or reproductive systems
- Pain or discomfort
- Bursitis
- Reduced bone density due to load avoidance
- Bone atrophy/fracture above or below the spine section treated with the implant
- Limited physical capacity
- Persistence of symptoms treated by the implantation
- Death

Safety notes

- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Only combine Aesculap modular implant components with each other.
- Do not use damaged or surgically excised components under any circumstances.
- Implants that have been used before must not be reused.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- Delayed healing can cause implant breakage due to metal fatigue.
- The physician in charge decides whether an implant component implanted should be removed,
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.
- Damage to the weight-bearing bone cement and/or bone structures can cause loosening of the components, fracture of the bone or implant and other grave complications.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
 - Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
 - Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
 - Do not use implant components that are past their expiry date or whose packaging is damaged.

Damage to implants caused by processing and resterilization!
➤ Do not process or resterilize the implants.

Application

- The operating surgeon draws up an operating plan that specifies and appropriately documents the following steps:
- Selection of the implant components and their dimensions
 - Positioning of the implant components in the bone
 - Location of intraoperative landmarks
- The following conditions must be fulfilled prior to application:
- All requisite implant components ready to hand.
 - Operating conditions highly aseptic.
 - Implantation instruments, including the special Aesculap implant system instruments, complete and in working condition.
 - Operating surgeon and operating team thoroughly conversant with the operating technique and the range of implants and instruments; information documents on these subjects complete and ready to hand.
 - The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific publications by medical authors.
 - The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the operating site.

- The intervention has been explained to the patient, whose consent concerning the following information has been documented:
- Delayed healing or incomplete fusion can cause the implant to fracture or loosen as a result of the extreme load to which it is subjected.
 - Smokers present an increased risk of fusion failure.
 - The lifespan of the implant depends on the patient's body weight.
 - Implant components must not be subjected to extreme load arising from intense exercise or strenuous physical work.
 - Revision surgery can become necessary in cases of implant loosening, fracture or loss of correction.
 - The patient must undergo regular medical follow-up examinations of the implant components.

- The implantation of the S⁴ Spinal System implant (augmentation screw) involves the following steps:
- Select the appropriate S⁴ Spinal System variant and anchoring implants according to indication, preoperative planning and intraoperative bone situation.
 - When choosing the S⁴ Spinal System implant, be careful to select the correct diameter, length and canal orientation.

Trauma to the spinal cord and nerve roots due to incorrect application!
➤ Align and screw in the instruments, the pedicle screw and the polyaxial screw only under X-ray control or using a navigation system.

- To avoid internal stress on, and weakening of, the implant: avoid scoring or scratching of the implant components.
- Only allow axial load on the screws during screw insertion. Avoid induction of lateral forces during screw insertion.

Risk of screw breakage caused by hard bone or contact with contralateral cortex!
➤ Prepare the vertebral body for the full length of the screw, using a thread cutter.

Damage to the hexagon socket in the screw head due to incorrect application of the Allen key!
➤ Make certain the Allen key is fully inserted in the hexagon socket of the screw head when turning the screw.

Note
Cementing is carried out only after all screws have been aligned and the connecting rods have been adapted, and before locking the S⁴ Spinal System augmentation screw!

Risk of breaking the injection cannula by inducing excessive shear forces!
➤ Do not apply excessive transversal load to the injection cannula.
➤ Do not bend the injection cannula.

Note

Once the S⁴ Spinal System augmentation screw has been cemented, it will be impossible to change the implant position, insert the reduction lever or apply corrections!

Note

With the S⁴ Spinal System augmentation screw finally cemented, explanting the screw can become extremely difficult, especially in cases of osteopenia.

- Prior to applying the bone cement, fully screw in the injection cannula at the thread in the bone screw head, without applying perceptible force.
- Check the bone cement for proper texture and viscosity before positioning the applicator on the injection cannula. Observe the manufacturer's advisories for the bone cement!
- Cement the S⁴ Spinal System augmentation screw only with cement that is suitable for screw augmentation.

Risk of screw breakage due to cement-free application!

- Insert the S⁴ Spinal System augmentation screw with cement.

Risk of injury from cement leakages!

- Always apply imaging techniques for control during the preparation and positioning of the implants.
- Make certain the slots of the S⁴ Spinal System augmentation screw are positioned fully inside the bone.
- Make sure that the pedicle cortex is not perforated.
- Do not penetrate the opposite cortex and ensure that the opposite cortex is undamaged.
- Do not insert or cement the S⁴ Spinal System augmentation screw in a fractured or partly resected vertebral body. Check by preoperative thin-layer CT.
- Screw the injection cannula into the screw head without applying any perceptible force. Make certain the axes of the pedicle screw and the injection cannula are aligned with each other.
- Apply cement only under high-quality radiographic control (max. volume of cement per polyaxial screw: 2 ml).

Risk of injury due to material incompatibility!

- Only use application systems that have been approved for use with bone cement.

Contamination of the locking thread/Insufficient fixation when removing the application system, due to incorrect viscosity of the cement!

- Make certain to be using cement of the correct viscosity.

- Do not allow the cement to cure completely before removing the injection cannula. Remove the cement application system incl. cannula as soon as the cement shows a doughy texture (allow for a certain amount of cement for reference purposes when preparing the cement).

Wound contamination due to insufficient curing of cement! Risk of caulking of the locking screw thread, resulting in system instability!

- Observe the correct curing period after removing the injection cannula.
- For monitoring the cement texture, prepare a reference volume of cement together with the cement to be applied.
- Make certain the system is not exposed to any load before the cement has reached its final strength.

Risk of connecting rod breakage!

- Do not bend back or bend excessively the connecting rods and cross connectors.

- If the screw head loosens during revision of the larger polyaxial screws (Ø 7.5 mm), remove the remaining screw component with the SW 3.5 screwdriver.

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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