- GB** **Instructions for use**
- USA** S⁴ Spinal System - dynamic connecting rod
- D** **Gebrauchsanweisung**
- S⁴ Spinal System - dynamischer Verbindungsstab
- F** **Mode d'emploi**
- S⁴ Spinal System - tige de liaison dynamique
- E** **Instrucciones de manejo**
- S⁴ Spinal System - barra de conexión dinámica
- I** **Istruzioni per l'uso**
- S⁴ Spinal System - barra di connessione dinamica
- P** **Instruções de utilização**
- S⁴ Spinal System - barra de união dinâmica
- NL** **Gebruiksaanwijzing**
- S⁴ Spinal System - dynamisch verbindingsstaafje
- S** **Bruksanvisning**
- S⁴ Spinal System - dynamisk förbindelsestav
- RUS** **Инструкция по применению**
- S⁴ Spinal System - динамический соединительный стержень
- CZ** **Návod k použití**
- S⁴ Spinal System - dynamická spojovací tyčka
- PL** **Instrukcja użytkowania**
- S⁴ Spinal System - dynamiczny pręt połączeniowy
- SK** **Návod na používanie**
- S⁴ Spinal System - dynamická spojovacia tyč
- TR** **Kullanım Kılavuzu**
- S⁴ Spinal System - dinamik bağlantı çubuğu

B | BRAUN
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 013452 11/11 V6

CE 0123 - DIR 93/42/EEC

Technical alterations reserved

Aesculap®

S⁴ Spinal System – dynamic connecting rod

Applicable to

- For item-specific instructions see also the Aesculap extranet at www.aesculap-extra.net

Intended use

The S⁴ Spinal System dynamic connecting rod is used for the dynamic monosegmental stabilization of the lumbar spine through the dorsal approach. It can be used to fuse one or two adjacent segments.

The S⁴ Spinal System dynamic connecting rod consists of a spring element, a short rod section on one side, and another rod section that varies in length on the other side.

S⁴ Spinal System dynamic connecting rods are supplied in a sterile condition and include:

- Straight dynamic connecting rods
- Pre-bent dynamic connecting rods

S⁴ Spinal System dynamic connecting rods are used with the following screws:

- S⁴ Spinal System cannulated polyaxial screws (supplied in an unsterile condition) or
- S⁴ Spinal System polyaxial augmentation screws (supplied in a sterile condition, observe TA012865)
- S⁴ Spinal System locking screws (supplied in an unsterile condition)

Note

Special instrumentation is required to implant these components.

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_T: Titanium forged alloy Ti6Al4V acc. to ISO5832-3
- ISOTAN®_P: Pure titanium acc. to ISO 5832-2

The titanium implants are coated with a colored oxide. Their color may change slightly, but this has no bearing on the quality of the implant.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

The S⁴ Spinal System dynamic connecting rod serves to support mobile segments.

Use the spring element of the S⁴ Spinal System dynamic connecting rod for dynamic stabilization in the event of:

- Signs of disk degeneration that have been confirmed via X-ray
- Central or lateral spinal canal stenosis in conjunction with a microsurgical decompression which has no effect on stability

The long rod section of the S⁴ Spinal System dynamic connecting rod for fusing one or two adjacent segments serves to support normal healing processes. It is not intended for use either as a replacement for natural body parts or to bear loads over the long term if healing does not occur.

Use the long rod section of the S⁴ Spinal System dynamic connecting rod for fusing in the event of:

- Fractures
- Pseudarthrosis or delayed healing
- Dislocation
- Degenerative instability
- Post-discectomy syndrome
- Spondylolisthesis
- Kyphosis
- Scoliosis
- Stenosis
- Post-trauma instability
- Tumors

WARNING

Risk of screw breakage caused by pedicle screws used for spondylolisthesis!
► Support the stabilization measure with interbody fusion (ALIF, PLIF or TLIF).

Contraindications

Do not use the S⁴ Spinal System dynamic connecting rod (spring element for dynamic stabilization and long rod section for fusing) in the event of:

- Fever
- Acute or chronic vertebral infections of a local or systemic nature
- Pregnancy
- Severe osteoporosis or osteopenia
- Medical or surgical conditions that could negatively affect the outcome of the implantation
- Mental illness
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Adiposity (BMI ≥ 30 kg/m²)
- Severely damaged bone structures that could prevent stable implantation of implant components
- Neuromuscular disorders or illnesses
- Bone tumors in the region of implant fixation
- Wound healing disorders
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials
- Treatment with immunosuppressives or long-term steroid treatment
- Severe muscular, nervous system or vascular system disease that impacts on the spine
- Any secondary conditions that would put the functionality or success of the implant at risk
- Vertebral body fractures
- Treatment of the thoracic or cervical spine
- Fusing of more than two segments
- Unilateral application of the implant
- Severe instabilities such as a laminectomy
- Facet joint arthrosis
- Cases not listed under indications

Furthermore, do not use the spring element for dynamic stabilization in the event of:

- Signs of an instability
- Support in the intervertebral space (e.g., bone chip, spacer or cages)
- An earlier failed attempt at fusion (pseudarthrosis), an earlier total facetectomy or trauma
- Degenerative scoliosis >10°

Side effects and interactions

The use or improper handling of the S⁴ Spinal System cannulated polyaxial screws or the S⁴ Spinal System polyaxial augmentation screws in conjunction with the S⁴ Spinal System dynamic connecting rod entails the following general risks:

- Implant failure due to overloading, wear or damage
 - Bending
 - Loosening
 - Breakage
- Inadequate fixation
- Infection
- Vertebral body fractures
- Tissue reaction to implant materials
- Hematomas and wound healing disorders
- Injuries to
 - nerve roots
 - spinal cord
 - blood vessels
 - organs
- Changes of the normal spine lordosis
- Impairment of the gastrointestinal, urological, and/or reproductive systems
- Pain or indisposition
- Bursitis
- Decreased bone density due to load avoidance
- Bone atrophy/fracture above or below the spine section provided for
- Limited physical ability
- Persistence of symptoms that were to be treated by the implantation
- Death

In addition, the following risks apply in the dynamically stabilized segment:

- Insufficient stabilization of the segment
- Advancing degeneration on, above or below the segment being treated

The following risk also applies in the fused segment:

- Failed or delayed fusion

Safety notes

CAUTION

Federal law restricts this device to purchase by, or on instruction by a physician!

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.

WARNING

Failure of the implant system due to incorrect screw selection!

- Only ever use cannulated polyaxial screws.
- Make sure they are correctly seated in the bone.
- Make sure the spring element is free of stress when seated.

WARNING

Never combine dynamic connecting rods and rods for fusion within the same segment!

- Use either only dynamic connecting rods or only rods for fusion within the same segment.

- The spring element must only be used for the dynamic stabilization of the treated segment.
- The disk and the facet joint must remain in the dynamically stabilized segment.
- The implant is extremely sensitive to damage; even small scratches can have a significant impact on its durability. It is essential to handle the implant with the utmost care.

WARNING

Risk of damage to the spring element during implantation!

- Prevent the spring element from coming into contact with other implants or instruments.
- Leave the insertion instrument on the spring element throughout the entire operation.

- Damage to the load-bearing structures of the segment being treated or the implant can lead to loosening or breakage of components, dislocation, migration, and other severe complications.
- Only the long rod section may be used to fuse an adjacent segment.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

- Delayed healing in the fused segment can cause implant breakage due to metal fatigue.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.
- The attending physician shall make any decision with regard to the removal of implant components that have been used,

Sterility

- S⁴ Spinal System dynamic connecting rods are supplied in a sterile condition.
- S⁴ Spinal System cannulated polyaxial screws are supplied in an unsterile condition.
- S⁴ Spinal System locking screws are supplied in an unsterile condition.

Note
Different reprocessing instructions apply to the S⁴ Spinal System dynamic connecting rods (sterile), the S⁴ Spinal System cannulated polyaxial screws (unsterile), and the S⁴ Spinal System locking screws (unsterile).

- The following applies to S⁴ Spinal System dynamic connecting rods:
- The implant components come individually packed in protective packaging that is labeled according to its contents.
 - The implant components are gamma-sterilized.
 - Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
 - Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
 - Do not use implant components that are past their expiration date or whose packaging is damaged

WARNING

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants

- The following applies to S⁴ Spinal System cannulated polyaxial screws and S⁴ Spinal System locking screws:
- The implant components are supplied in an unsterile condition.
 - The implant components are packaged individually.
 - Store the implant components in their original packaging and only remove them from their original and protective packaging immediately prior to application.
 - Use the implant system storage devices for processing, sterilization and sterile setup.
 - Use a suitable tray for cleaning/disinfection.
 - Use the system storage device only for sterilization and sterile provision.
 - Ensure that the implant components in their implant system storage devices do not come into contact with each other or with instruments.
 - Ensure that the implant components are not damaged in any way.

Prior to initial sterilization and subsequent resterilization, the implant components must be cleaned using the following validated reprocessing procedure:

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.aesculap-extra.net
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Validated reprocessing procedure

Note
The reprocessing procedure that follows applies only to S⁴ Spinal System cannulated polyaxial screws and S⁴ Spinal System locking screws.

CAUTION

Risk of damage to implants during reprocessing!

- Do not reprocess S⁴ Spinal System dynamic connecting rods.

Cleaning/Disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Risk of damage to the product due to inappropriate cleaning/disinfecting agents!

- Following the manufacturer's instructions, use cleaning and disinfecting agents which
 - are approved for the material in question (e.g., aluminum, plastics, high-grade steel),
 - do not attack softeners (e.g., in silicone).
- Observe specifications regarding concentration, temperature, and exposure time.

- Carry out ultrasound cleaning:
 - As an effective mechanical supplement to manual cleaning/disinfection

Validated cleaning and disinfecting procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with ultrasound and immersion disinfection	■ 20 ml disposable syringe ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air.	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with ultrasound and immersion disinfection

Manual cleaning/disinfection

- Prior to manual disinfection, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, visually check visible surfaces for residues.
- Repeat the cleaning/disinfection process if necessary.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Ultrasound cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	–	FD–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended:BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfecting procedure.

Phase I

- Clean the product in an ultrasound cleaning bath (frequency of 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Do not use brushes or other abrasives that would damage the product surfaces, as this would cause corrosion.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 minute.
- Mobilize non-rigid components, such as set screws and hinges, during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a suitable disposable syringe.
- To prevent internal stresses and weakening of the implants, avoid scoring or scratching any components.

Phase II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws and hinges, during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the instrument in the disinfectant solution.
- Mobilize non-rigid components, such as set screws and hinges, during disinfection.
- Rinse lumens at least five times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws and hinges, during final rinsing.
- Rinse lumens at least five times using an appropriate disposable syringe.
- Drain any remaining water fully.

Phase V

- In the drying phase, dry the product via appropriate means (e.g., wipes, compressed air), see Validated cleaning and disinfecting procedure.

For implant components that are to be resterilized:

WARNING

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- Handle the implants with new gloves only.
- Keep the implant system storage devices covered or closed.
- Process implant system storage devices separately from instrument trays.
- Clean implants must not be processed together with contaminated implants.
- Process the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.
- Manually clean and disinfect the implant components.
- Do not reuse surgically contaminated implants!

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- Do not reprocess implants that have been directly or indirectly contaminated with blood.

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- Immediately sort out damaged or inoperative products.

Packaging

- Place the product in its holder or on a suitable tray.
- Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Make sure that the packaging will prevent a recontamination of the product whilst in storage.

Sterilization

- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

WARNING

Risk of injury caused by incorrect operation of the product!

- ▶ Attend appropriate product training before using the product.
- ▶ For information about product training, please contact your national B. Braun/Aesculap agency.

- The operating surgeon shall devise an operation plan that specifies and accurately documents the following:
- Selection of the implant components and their dimensions
 - Positioning of the implant components in the bone
 - Location of intraoperative landmarks
- The following conditions must be fulfilled prior to application:
- All requisite implant components are ready to hand.
 - Operating conditions are highly aseptic.
 - All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
 - The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
 - The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
 - Information has been obtained from the manufacturer in the event of an ambiguous preoperative situation and if implants are present in the area to be treated

- The intervention has been explained to the patient, whose consent concerning the following information has been documented:
- In the case of delayed or incomplete fusion in the fused segment, the implants can break and loosen due to high loads
 - The life-span of the implant depends on the patient's body weight.
 - The implant components must not be overloaded by extreme strains, hard physical labor or sports.
 - Corrective surgery may be necessitated by implant loosening, fracture or loss of correction.
 - Smokers present an increased risk of bone fusion failure.
 - The patient must undergo medical check-ups of the implant components at regular intervals.
- The following steps need to be carried out in order to implant the S⁴ Spinal System implant:
- ▶ Select the appropriate S⁴ Spinal System variant and anchoring implants according to indication, preoperative planning, and intraoperative bone situation.

WARNING

Pedicle screw clearance too small for spring element!

- ▶ Estimate the pedicle screw clearance prior to surgery.
- ▶ Use K-wire to determine the insertion point of the pedicle screws.
- ▶ Use a test implant to establish whether the pedicle screw clearance is sufficient as soon as the K-wires have been put in place.

- ▶ Only ever use S⁴ Spinal System cannulated polyaxial screws or S⁴ Spinal System polyaxial augmentation screws (TA012865).
- ▶ When selecting S⁴ Spinal System screw types, make certain to select the correct diameter, length, and canal orientation.

WARNING

Incorrect application can result in trauma to the spinal cord and nerve roots!

- ▶ Align and screw in the instruments, the pedicle screw and the polyaxial screw only under X-ray control or using a navigation system.

- ▶ Use a test implant to establish and clear the space required for the spring element.

WARNING

Risk of the spring element breaking due to contact with osseous structures!

- ▶ Prepare a space for the spring element.
- ▶ Use a test implant.

- ▶ Avoid any height difference between the cranial and caudal pedicle screws of the spring element.
- ▶ To avoid internal stress on, and weakening of, the implant: avoid scoring or scratching of the implant components.
- ▶ Do not bend metal implants, except for the long rod sections of the S⁴ Spinal System dynamic connector rod.

WARNING

Risk of damage to the spring element when bending the long rod sections!

- ▶ Only ever bend the long rod sections.
- ▶ Ensure that adequate clearance (min. 10 mm) is maintained between the bending instrument and the spring element during the bending procedure.
- ▶ Prevent the bending instrument from coming into contact with the spring element.

- ▶ Do not bend the long rod sections back or bend them excessively.
- ▶ Only use bending instruments from the S⁴ Spinal System instrument set for bending the long rod sections.

WARNING

Risk of damage to the spring element when cutting the long rod sections!

- ▶ Only ever cut the long rod sections.
- ▶ Ensure that adequate clearance (min. 10 mm) is maintained between the rod cutter and the spring element during the cutting procedure.
- ▶ Prevent the rod cutter from coming into contact with the spring element.

- ▶ Only use rod cutters from the S⁴ Spinal System instrument set for cutting the long rod sections.
- ▶ Always use the insertion forceps to put the spring element in position. Avoid making contact with other implant components.
- ▶ When putting the spring element in position, make sure that the three loops on the spring element are facing the dorsal direction and the two loops are facing the ventral direction.
- ▶ To position the locking screw, always use the locking screw insertion instrument intended for this purpose.
- ▶ Do not carry out any derotation maneuvers.
- ▶ For tightening and loosening the locking screw, always use the screwdriver and countering instrument intended for this purpose. If there is not much space between the screw and the spring element, use the sleeve for dynamic rods.
- ▶ Always tighten the locking screws with the torque wrench intended for this purpose.

WARNING

Risk of injury due to improper mounting of the locking screw!

- ▶ Position the locking screw correctly.
- ▶ Make certain that the rods are correctly positioned on the floor of the groove.
- ▶ Fully tighten the locking screw with the torque wrench.
- ▶ Break off the flanks only once the locking screw has been fully tightened.

WARNING

Correction loss due to insufficient fixation of the implant head!

- ▶ Never loosen the connection of the implant head after it has been tightened.
- ▶ Tighten the locking screw only after every necessary corrective measure has been performed.

WARNING

Risk of damage to the implant due to overtightening of the locking screw!

- ▶ Tighten the locking screw only with the exact torque specified for this purpose.

WARNING

Expansion of the implant head caused by tightening of the locking screw without application of the countering instrument!

- ▶ Always use the countering instrument when tightening the locking screw.
- ▶ If there is not much space between the screw and the spring element, use the sleeve for dynamic rods.

WARNING

Bending or twisting of the spring element while tightening the locking screws!

- ▶ Use the countering instrument with care
- ▶ Make certain that the heads of the pedicle screws are in parallel.
- ▶ Clamp the spring element during tightening with the insertion forceps.

- ▶ Check that the spring element is free of stress when seated before and after tightening the locking screws.

WARNING

Damage to the implant caused by applying the reduction instruments too high at the implant!

- ▶ Always apply reduction instruments (e.g. distraction and compression forceps) below the rod at the implant.

WARNING

Bracing/twisting of the spring element via repositioning instruments!

- ▶ Do not perform any manipulations (using distraction and compression forceps, for example) in the dynamically stabilized segment.

WARNING

Inadequate fixation due to incorrect positioning of insert!

- ▶ When placing the rod, make certain that the insert is correctly positioned.
- ▶ After manipulating the rod, make certain that the insert is correctly positioned.
- ▶ After final locking, check that the insert is correctly positioned.

WARNING

Inadequate fixation due to incorrect positioning of the rod in the screw head!

- ▶ Adjust the curvature of the rod to the anatomical situation.
- ▶ Do not bend the rod in the region of the screw head.
- ▶ Make certain and check that the screw head is at right angles to the rod.
- ▶ Make certain and check that the rod is positioned securely in the groove base.

WARNING

Risk of damage to the implant caused by spondylolisthesis repositioning via the locking screw!

- ▶ Always use the rod persuader for spondylolisthesis reposition.
- ▶ Spondylolisthesis repositioning is only permitted in the fused level.

- ▶ Always use the special flank breaking forceps for breaking off the flanks.
- ▶ Once the locking screws have been locked, remove the insertion forceps.
- ▶ If the screw head comes loose during revision of the larger polyaxial screws (Ø 7.5 mm or larger), remove the remaining screw component with the SW 3.5 screwdriver.
- ▶ For tightening the locking screws of the cross connector, always use the torque wrench and countering instrument intended for this purpose.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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