Aesculap Spine

GB

Instructions for use

USA

S⁴ Spinal System - injection cannula

D

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S⁴ Spinal System - Injektionskanüle

F

Mode d'emploi

S⁴ Spinal System - Canule d'injection

E

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TA-Nr. 013132 07/09 V3

0123

CE marking according to directive 93/42/EEC

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Symbols on product and packages

Symbol	Explanation
	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution: See documentation supplied with the product

Intended use

The S⁴ Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar and thoracic spine.

The S⁴ Spinal System - injection cannula is used for applying bone cement through the S⁴ Spinal System - augmentation screw.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Note

Generally, the instructions for use of the S⁴ Spinal System - Lumbar/Deformity TA011187 and the S⁴ Spinal System augmentation screw TA012865 apply to the S⁴ Spinal System - augmentation screw!

The present information on the S⁴ Spinal System - injection cannula complements the respective descriptions in the instructions for use of the S⁴ Spinal System - Lumbar/Deformity and the augmentation screw.

Note

The S⁴ Spinal System - injection cannula is intended only for single use in connection with the S⁴ Spinal System - augmentation screw.

- The product is gamma-sterilized and ships in sterile packaging.
- The product must not be reused.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
 - Read, follow and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Do not use products from open or damaged sterile packaging.
 - Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Do not re-sterilize the product.
 - Do not use the product after its use-by date.

Safe operation

- Risk of injury from cement leakages!**
- Insert the injection cannula into the screw head by turning it gently, without applying palpable force.
 - Make certain the injection cannula is in axial alignment with the pedicle screw.
 - Cement application only under high-quality radiographic control (max. volume of cement per pedicle screw: 2 ml).

- Risk of injury and/or malfunction!**
- Always carry out a function check before using the product.

- Risk of breaking the injection cannula by inducing excessive shear forces!**
- Do not subject the injection cannula to transverse loads; do not bend the injection cannula.

- Contamination of the locking thread/ball head of the S⁴ Spinal System - augmentation screw!**
- Inadequate fixation during removal of the application system due to incorrect viscosity of the cement!**
- Be sure to use cement of the correct viscosity.
 - Make certain the locking thread/ball head are not contaminated with cement.

Single-use products

- Patients at risk from damaged, broken or contaminated products!**
- These products are intended for single use only.**
- Do not reprocess these products.

Storage

- Store single-use products in sterile packaging in dust-protected conditions in a dry, dark and temperature-controlled room.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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