AESCULAP®

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Legend

- 1 Outer lid
- 2 Outer lid lock
- 3 Inner lid
- 4 Inner lid lock
- 5 Bottom
- 6 Indicator seal retainer
- 7 Indicator seal
- 8 Handle
- 9 Slot (for the indicator seal on the inner lid lock)
- 10 Universal filter retainer
- 11 Cap
- 12 Push buttons
- 13 Lid latch (for optional outer lid of BASIS container)
- 14 Plastic cover
- 15 Perforation field cover
- 16 Germ retention system
- 17 Adapter frame
- 18 Finger tray
- 19 Handling pin
- 20 Ribbed cover grid
- 21 Mounting handle
- 22 Plastic cover lock

Symbols on the product

	Attention, see instructions for use
---	--

Contents

1.	Safe handling	8
2.	Product description.....	9
2.1	Intended use	9
2.2	Mode of operation	9
3.	Preparation and setup	9
3.1	First use.....	9
4.	Working with the sterile container.....	9
4.1	System set-up.....	9
4.2	Function check	10
4.3	Safe operation.....	10
5.	Maintenance	11
6.	Troubleshooting list	12
7.	Technical Service	13
8.	Accessories/Spare parts	13
9.	Technical specifications.....	13
10.	Extracts from relevant standards	13
10.1	Standards cited.....	13
10.2	Storage periods for sterile materials: extract from DIN 58953-9, sterile materials supply	13

1. Safe handling

- Prior to use, check for proper condition and functioning of the product.
- In order to avoid damaging the device through improper assembly/operation and mitigating the manufacturer's liability:
 - Use the product only according to the instructions.
 - Follow the safety instructions and maintenance advisories.
 - Only combine Aesculap products with each other.
 - Never use damaged or faulty sterile containers. Replace any damaged components immediately with original spare parts.
 - If sterile container components are repaired in a way that could affect the container's germproof qualities: Inspect the sterile container thoroughly before use.
- Make sure that the product and accessories are operated and used only by qualified or trained and experienced personnel.
- Keep the instructions for use accessible to all staff.
- Adhere to applicable norms.
- Follow the general guidelines and aseptic principles when handling contaminated items that have undergone or are to undergo sterilization.

2. Product description

2.1 Intended use

The Aesculap sterile container system is intended for use as a sterile packaging for instruments and textiles due for steam sterilization according to a validated steam sterilization process (e.g. in a sterilizer complying with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated according to EN 554/ISO 13683).

After being sterilized, the sterile materials are stored in the sterile container until they are used.

Note

Please contact your Aesculap representative if your Aesculap sterile containers are to be used in any other steam sterilization process.

2.2 Mode of operation

The Aesculap sterile container system complies with the stipulations of DIN EN 868-1, DIN EN 868-8, DIN 58953-9 and ISO 11607.

Sterile containers with perforated lids and closed bottoms are suitable for:

- fractional flow processes
- fractional vacuum processes
- steam injection processes

Sterile containers with perforated lids and perforated bottoms are suitable for:

- gravitational procedures
- pre-vacuum procedures

3. Preparation and setup

3.1 First use

- Prior to the first sterilization, clean the new sterile container either manually or mechanically.
- After cleaning, use a suitable filter, see Changing filters.

Sterile container system PrimeLine:

The germ retention system 16 is integrated with the system.

4. Working with the sterile container

4.1 System set-up

Removing the outer lid/inner lid/plastic cover

If an outer lid is used, it can be removed to clean the sterile container and, if it is soiled, it can be separated from the inner lid.

VARIO container (by default with outer and inner lid):

- Push down the outer lid lock 2 and remove the outer lid 1 and the inner lid 3.

BASIS container (with retrofitted outer lid):

- Dismount the combined outer lid 1 and inner lid 3 from the bottom 5.
- Loosen the lid latch 13 and remove the outer lid 1.

Sterile container system PrimeLine:

- Remove the plastic cover 14 from the bottom 5.

Disinfecting and cleaning sterile container

Damage to the plastic cover caused by Acetone!

- Clean the plastic cover and filter parts with pH-neutral detergents only.

Note

Literature on disinfecting and cleaning:

"Test series and statements; Working group instrument preparation; C 759", available through Aesculap.

- For disinfecting and cleaning the sterile containers check with the manufacturers of the washers/disinfectors, water preparation, and chemical cleaning solutions.
- Use AESCULAP Eloxal Cleaner JG 601 to remove adhesive residues, inscriptions and similar serious contaminations on Eloxal surfaces.

Manual cleaning and disinfecting

- Wipe surfaces with a clean, lint-free cloth and a pH-neutral cleaning and disinfecting solution.
- Remove encrusted materials with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Carry out the final rinse in distilled or demineralized water.

Mechanical cleaning and disinfecting

Note

Repeated mechanical cleaning may lead to a discoloring (bleaching, staining) of the aluminum lids. This has no effect whatsoever on their performance or efficiency.

For colored aluminum lids, apply, if possible, a process involving a neutral detergent and fully demineralized water.

- Use suitable storage racks for mechanical cleaning, according to the manufacturers of the respective devices.
- Carry out the final rinse in distilled or demineralized water.

Changing filters

Change filters at the following intervals, depending upon filter type:

- Disposable filter: prior to each sterilization cycle
- Textile filter: every 50 sterilization cycles
- Reusable filter: every 1000 sterilization cycles
- PrimeLine germ retention system: every 5000 sterilization cycles

VARIO container and BASIS container:

- Press simultaneously on both push buttons **12** on the universal filter retainer **10**.
- Remove the universal filter retainer **10**.
- Insert a new filter and remount the universal filter retainer **10**.
- Push down the cap **11** on the universal filter retainer **10** until you hear it click into place.

Sterile container system PrimeLine:

- Turn the ribbed cover grid **20** anticlockwise until it is unlocked.
- Remove the ribbed cover grid **20** of the germ retention system **16**.
- Turn the germ retention system **16** anticlockwise at the finger tray **18** or the mounting handle **21** until it disengages from the adapter frame **17**.
- Lift the germ retention system **16**, gripping it at the handling pin **19**, and remove it.
- For installing the germ retention system **16**, retrace the above steps.

4.2 Function check

- Visually inspect all components of the sterile container for damage and correct functioning.
- Use sterile containers only if they are in mint condition. Replace any damaged components immediately with original spare parts or have the parts repaired.

4.3 Safe operation

Loading sterile container

Instrument system

According to DIN EN 868-8 and DIN 58953-9, the following maximum load of the container is recommended:

- Standard container: 10 kg
- Half container: 5 kg
- 3/4-container: 7 kg

- Store instruments in wire baskets in suitable supports.
In doing this, lay hollow components, bowls, plates, etc. face down and on a slant.
- For double-packaging of sterile materials in accordance with DIN 58953-9 wrap the loaded wire basket separately.

VARIO container and BASIS container:

- Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed (VARIO and BASIS container).
- Use the inner lid lock **4** to closure the inner lid **3** to the container bottom **5**.
The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that the ribbed cover grid **20** in the plastic cover **14** remains unobstructed.
- Use the plastic cover lock **22** to lock the plastic cover **14** on the bottom **5**.
The plastic cover lock **22** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Textiles

- Pack the textiles together in such a way that they are upright in the sterile container.
- Make certain that when the sterile container is fully loaded, it is still possible to insert an outstretched hand between the individual items without difficulty.
- For double-packaging of sterile materials in accordance with DIN 58953-9, wrap the textile items for sterilization separately.
- Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed (with VARIO and BASIS containers).
- Use the inner lid lock **4** to closure the inner lid **3** to the container bottom **5**.
The inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that the ribbed cover grid **20** in the plastic cover **14** remains unobstructed.
- Use the plastic cover lock **22** to lock the plastic cover **14** on the bottom **5**.
The plastic cover lock **22** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Inserting the indicator seal

- Having loaded the sterile container, make notes of the following data on the indicator seal **7**: Date of sterilization, sterilization number, expiry date, name and signature.
- Slide the indicator seal **7** from the outer side into the indicator seal retainer **6** in such a way that the indicator segment engages in the slot on the lid lock and seals the closure.

Loading the sterilizer

- Prepare the sterile container and the sterilizer for loading of the latter in the following way:
 - Insofar as possible, load the sterilizer with a mixture of textiles and instrument containers.
 - Do not use outer packaging for the sterile containers.
 - Always place heavy sterile containers at the bottom of the sterilizer.
 - Do not exceed the maximum height of 60 cm for stacks of containers.
 - For all container versions, never allow the airtight closure (e.g. by packaging foil) of the media passages in the bottom and inner lid.
 - Sterilize VARIO containers with their outer lids attached.
- Always carry the sterile containers by their handles.
- Transport stacks of sterile containers in such a way that the stacks do not topple over.
- Follow the sterilizer manufacturer's instructions.
- Sterilize with steam, taking note of the following:

The sterilization has to be done according to a validated steam sterilization procedure (e.g. in a sterilizer according to EN285/ANSI/AAMI/ISO11134-1993, ANSI/AAMI ST46-1993 and validated according to EN 554/ISO 13693). In case of application of the fractionated vacuum procedure the sterilization has to be done using the 134 °C/2 bar program with a minimal holding time of 5 minutes.

Unloading the sterilizer

Risk of burns due to a hot sterile container after sterilization!

- Always wear gloves when handling a container.

- With heavy sterilizer loads, allow the sterile containers to cool to room temperature by leaving the sterilizer door slightly ajar.
- Always carry the sterile containers by their handles.
- Transport stacks of sterile containers in such a way that the stacks do not topple over.
- Store sterile containers during cooling down in air-conditioned rooms in accordance with DIN 58953-9 standards.

To prevent condensate from forming: Never allow sterile containers to cool down on a cold floor or in a drafty place.

Checking and providing sterile materials

The contents of a sterile container are only deemed sterile if the container has been properly sterilized.

If this is not the case, the sterile materials must be reprocessed again.

Risk of contamination from improperly sterilized materials!

- Before providing sterile materials, check to ensure that the sterilization was successful.

- Before opening the sterile container, make certain that its indicator seal is not damaged.
- Make certain that the color of the indicator point has changed.
- Make certain that all container components, particularly the lid locks, are intact.
- Make certain that there is no condensate in the container.

If any condensate is found, eliminate the cause, see Troubleshooting list.

Storing sterile containers

- Set the sterile containers in place one at a time.
 - or –
 - Stack the sterile containers (maximum stack height 60 cm).
- Store sterile containers in a dry, clean and protected place.
- Comply with DIN 58953-9 standards for storage periods and conditions, see Storage periods for sterile materials: extract from DIN 58953-9, sterile materials supply.

Recommended storage periods according to DIN 58953-9 guidelines: 6 months.

5. Maintenance

The seals in the container lids have been tested, in conformance with DIN EN 868-8, Appendix G for a minimum service life of 5000 cycles.

- In case of visible damage, change the seal immediately.
- Have the lid sent to Aesculap, see Technical Service.
 - or –
 - Order seal and sticker, see Accessories/Spare parts and replace seal.

6. Troubleshooting list

Error	Cause	Remedy
Too much condensate inside the sterile container	Temperature of sterile materials too low prior to sterilization	Allow sterile materials to come to room temperature (approx. 20 °C) before sterilization
	Textiles too damp	Sterilize dry textiles only
	Sterile container too heavy	Standard container with instruments: max. load 10 kg
		Half container with instruments: max. load 5 kg
		Standard container with textiles: max. load 8 kg
		3/4-container with instruments: max. load 7 kg
	Sterile materials improperly positioned	Lay hollow components, bowls, plates etc. face down and on a slant
		Arrange textiles in loose vertical piles, do not press them together
	No inner wrapping used	Use suitable inner wrapping
	Sterile container incorrectly positioned in sterilizer	Always place heavy sterile containers at the bottom
	Sterile containers supplied for use immediately after sterilization	Allow sterile containers to cool down to room temperature prior to processing
	Sterile containers improperly positioned during cooling phase	Do not store sterile containers on a floor or in a drafty place Store sterile containers in air-conditioned areas according to DIN 58953-9 standards
	Sterilizer properties do not meet DIN EN 285 standards	Have sterilizer serviced regularly Check drying vacuum Check drying time
		Check steam quality and upgrade if necessary
Condensate in groove of (inner) lid	Empty-cycle and vacuum test not run daily before sterilization begins	Run empty-cycle and vacuum test daily before beginning sterilization
	Unsuitable sterilizer cycle selected	Select a program according to the load
	Sterilizer door left open too long, sterilizer cooled down	Load and unload sterilizer quickly
Indicator seal color unchanged	Sterilization improperly performed, sterilizer faulty	Have sterilizer repaired by manufacturer
BASIS/VARIO sterile container or PrimeLine sterile container system deformed	Perforations covered during sterilization	Never cover perforations
Lid lock faulty	Outer lid used to carry sterile container	Always carry the sterile containers by their handles
Textile filter has brown discoloration	Textile filter no longer useable or its shelf life has expired	Check steam quality and upgrade if necessary Change filters

7. Technical Service

For service, maintenance or repairs, contact your national Aesculap/ B. Braun representative.

All guarantees and warranty rights as well as licenses will be void if modifications on the device are performed.

Service addresses

AESCU LAP Technischer Service

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78532 Tuttlingen/Germany

Phone: +49 (7461) 95 27 00

Fax: +49 (7461) 16 28 87

E-mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

After repairs have been performed

- If sterile container components are repaired in a way that could affect the container's germproof qualities: Inspect the sterile container thoroughly before use.

8. Accessories/Spare parts

Accessories and consumables can be found in brochure no. C 404 01.

9. Technical specifications

The various models and sizes of the sterile containers can be found in brochure no. C 404 01.

10. Extracts from relevant standards

10.1 Standards cited

The following standards are cited in connection with the sterile containers:

- DIN 58953-9
- DIN 58946-1, DIN EN 285, DIN EN 866-3
- EN 868-1, EN 868-8
- EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, with applicable EN 554/ISO 13683
- ISO 11607

10.2 Storage periods for sterile materials: extract from DIN 58953-9, sterile materials supply

Loss of sterility is not so much connected to the storage periods as to outside influences and the effects of storage, transport and handling. Therefore, blanket statements cannot be made regarding appropriate storage periods. The appropriate length of storage periods is determined by the hygiene control committee. Responsibility for storage conditions and duration of storage periods lies with the hospital operator or its medical director.

Storage periods should be defined differently (where appropriate) for each department, and in writing. In doing this, the type of packaging and the storage conditions should be taken into account (see Table).

Recommended storage periods for sterile materials in sterile containers

Packaging for sterile materials	Type of packaging	Storage periods
Sterile container according to DIN EN 868-1 or DIN EN 868-8 standards	Single-packaging for sterile materials Double-packaging for sterile materials	6 months
NOTE 1: The information contained in this table are recommendations which, if followed, will decrease the risk of contamination when sterile containers are transported and opened. Storage for longer periods can lead to contamination of outer surfaces, depending upon storage conditions. Although such longer periods do not necessarily lead to recontamination of the contents of the package during storage, they do increase the risk of contamination when sterile containers are transported or opened. The storage recommendations contained in this table are made on the assumption that the materials will be stored in a dry area that is protected against dust. Longer storage periods are feasible, depending upon conditions at the storage site.		
NOTE 2: The use of an inner wrapping for double-packaging sterile materials improves aseptic presentation and is therefore recommended.		

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AESCULAP AG & CO. KG