GB	Instructions for use TARGON® femur and tibia nail system
D	Gebruuchsanweisung TARGON®-Femur- und -Tibia-Nagelsystem
F	Mode d'emploi Système de clous TARGON® pour fémur et tibia
E	Instrucciones de manejo Sistema de clavos para fémur y tibia TARGON®
I	Istruzioni per l'uso Sistema di chiodi per femore e tibia TARGON®
P	Instruções de utilização Sistema de pregos femorais e tibiais TARGON®
NL	Gebruiksaanwijzing TARGON®-femur- en -tibia-nagelsysteem
RUS	Инструкция по применению Система гвоздей TARGON® для бедренной и большеберцовой кости
CZ	Návod k použití TARGON® - systém femorálních a tibialních hřebů
PL	Instrukcja użytkowania System gwoździ udowych i piszczelowych TARGON®
SK	Návod na používanie Femorálny a tibialný klincový systém TARGON®
CN	使用说明 TARGON® 带锁髓内钉系统

B | BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 010451 01/11 V4 Änd.-Nr. 32365/32607

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Simbolo CE, em conformidade com a Directiva 93/42/CEE
CE-certificering conform richtlijn 93/42/EEG
Маркировка CE согласно директиве 93/42/ЕЭС
CE-označení podle směrnice 93/42/EHS
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
Označenie CE podľa smernice 93/42/EHS
根据 93/42/EWG 指令的 CE 标志

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Technische wijzigingen voorbehouden
Мы оставляем за собой право вносить технические изменения
Technické zmeny vyhradeny
Zmiany techniczne zastrzeżone
Technické zmeny vyhradené
保留技术变动的权力

Intended use

The TARGON® femoral and tibial nail system is used for splinting of fractures of long, tubular bones (femur and tibia).

System configurations

- TARGON® femoral universal nail various diameters and lengths
- TARGON® tibial universal nail various diameters and lengths
- TARGON® tibial solid titanium nail various diameters and lengths
- TARGON® femoral solid titanium nail various diameters and lengths
- For further information regarding diameters and length, see operating manual O 113 and brochure O 111.
- TARGON® retrograde femoral nail various diameters and lengths
- For further information regarding diameters and lengths, see brochure O 190.

Associated implants:

- Fixation screws in various diameters and lengths
- Threaded sleeves for fixation screws
- Fixation screws for the end of the nail on the connecting side

Materials

The materials used in the implant are listed on the package.

- ISOTAN®_F titanium forged alloy Ti6Al4V according to ISO 5832-3
- Implant steel acc. to ISO 5832-1, cold-worked for high strength, with chemically passive surface

The titanium implants are coated with a colored oxide. Their color may change slightly, but this has no bearing on the quality of the implant.

ISOTAN®_F is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use the TARGON® femoral universal nail, the TARGON® tibial universal nail, the TARGON® tibial solid titanium nail and the TARGON® femoral solid titanium nail in cases of:

- Open or closed fractures of the femur shaft or tibia shaft
- Pseudoarthrosis and transposition osteotomies in the shaft region
- Reconstruction following tumor resection in the shaft region
- Osteotomies in the shaft region to lengthen or shorten the limb
- Pseudoarthrosis
- Knee arthrodesis (TARGON® F special length)

Use the TARGON® retrograde femoral nail for:

- Femur shaft fractures
- Extraarticular or intraarticular fractures of the distal femur (these fractures also with ipsilateral femur shaft fractures)
- Patients with ipsilateral knee prosthesis, hip prosthesis or other proximal implants, e.g. a dynamic hip screw
- Shaft fractures for which additional proximal implants are necessary, e.g. ipsilateral femur fractures
- Shaft fracture and additional nailing of the ipsilateral tibia (identical, extended approach incision)
- Pathological fractures
- Pseudoarthroses
- Ankle arthrodesis

The operating surgeon is responsible for establishing these indications, as well as other indications not listed above, taking into account the individual clinical, biological and biomechanical situation.

Contraindications

Generally, do not use in cases of:

- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Do not use the TARGON® femoral universal nail, the TARGON® tibial universal nail, the TARGON® tibial solid titanium nail and the TARGON® femoral solid titanium nail in cases of:

- Open or closed fractures in joint regions
- Trochanteric fractures of the femur

Do not use the TARGON® retrograde femoral nail in cases of:

- Septic knee arthritis

Side-effects and adverse interactions

- Position change, loosening, bending or fracture of the implant
- Position change and loosening of fragments
- Primary and secondary infections
- Venous thrombosis, lung embolism, and cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Delayed healing or nonunion
- Pseudoarthrosis
- Compartment syndrome
- Pain in the region of the nail insertion point and in the region of the locking components

Safety notes

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- The operating surgeon must have a thorough command of all hands-on and conceptual aspects of the established operating techniques.
- Read and keep the instructions for use; follow the instructions in the OP Manual for the system variants.
- General risks entailed by surgical procedures are not described in these instructions for use.
- Follow the instructions for use of the respective Aesculap implant components.
- It is the operating surgeon's responsibility to ensure that the correct combination of implant components is selected and implanted properly.
- TARGON® implants may only be implanted using the Aesculap implantation instruments specially designed for this purpose.

- Use implant components only in the nail holes intended for this purpose.
- The operating surgeon must be thoroughly familiar with the bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- Under no circumstances may implant components from different suppliers be combined.
- Implants that have been used once for their designated purpose may not be reused.
- Never use damaged implant components or implant components that have been removed surgically.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- To allow the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals through appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- Patients must be provided with information regarding dynamic load limitations, specific behavior guidelines and the risk (e.g., implant failure) of not adhering to these parameters. The implants act as a splint for the bone until the healing process is completed. The implant cannot replace the function, nor bear the load, of a healthy bone. The patient must be very careful with the injured limb.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudoarthrosis develops, or if the implant is overloaded for an excessively long period. One way to achieve this is through mobilization.
- Mobilization entails the risk that the nail in the knee, hip or ankle joint will migrate. If mobilization is planned, make certain to leave sufficient space in front of the end of the nail.
- Always apply TARGON® implants in such a way that only mild forces have to be transferred through the implant and that any such forces will be born by the bone at an early stage.
- If fixation screws are positioned in the same plane as, or in the region of, the fracture line, increased lever action will lead to a considerably higher load being placed on the nail. This can lead to implant failure. The postoperative load on the nail must be kept appropriately low: the implant should be fully loaded only after a callus has developed.
- After fixing the nail, remove any transmedullary bone or fragment correction screws; or, within the patient's allowable limits, make certain that the nail can tolerate a high dynamic load engendered by unpropitious lever action.
- Never use damaged implant components or implant components that have been removed surgically. A damaged or scratched implant surface, particularly in the region of the nail drill hole, can lead to the development of cracks and early implant failure.
- Do not apply excessive force when implanting and explanting TARGON® locking nails and the locking screws and implant components and when working with the instruments. Should a problem arise, check the position and orientation of the implants, instruments and fragments and analyze the origin of the problem. If necessary, repeat one or more steps of the procedure and inspect the instruments (e.g., to see if the drill helixes are clogged).
- Only use cutting instruments, e.g. drill bits, guide pins, etc., that are sharp.
- Increase the load on the fracture and implant in accordance with the progress of healing.
- Delayed healing can cause implant breakage due to metal fatigue.
- At the time of follow-up exams, the scope of postoperative load and aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. The dimensions of the implant should also be taken into account.

Recommended times of follow-up exams:

- prior to the patient's discharge from hospital
- 10-12 weeks post operationem
- 6 months post operationem
- 12 months post operationem

For implanting the nail via a guide pin:

- Use a fitting guide pin
- Avoid jamming the implant or guide pin
- Monitor the guide pin end by X-ray control
- Remove the guide pin prior to locking

Aesculap and the Aesculap Academy regularly hold general and product-specific training events, offering users helpful application advice, especially for complicated or critical situations.

Sterility

- The implants are delivered unsterile.
- The implant components are delivered each in an individual packaging.

For implant components in original packaging:

- Store the implant components in their original packaging and only remove them from their original and protective packaging immediately prior to application.
- Use the implant system storage devices for processing, sterilization and sterile setup.
- Make certain that the implant components in their implant system storage devices will not get into contact with each other or with instruments.
- Make certain that the implant components will not be damaged under any circumstances.

For implant components that are to be resterilized:

Intraoperative contamination with blood, secretions and other fluids may render the affected component unsuitable for resterilization!

- Handle the implants with new gloves only.
- Keep the implant system storage devices covered or closed.
- Put away implant system storage devices separately from instrument trays.
- Process the implant components individually and separately if no implant system storage devices are available. When doing this, make certain that the implant components will not be damaged.
- Carry out manual and, if necessary, additional ultrasound pre-cleaning of soiled implant components.
- Mechanically clean and disinfect the implant components.
- Carry out the final rinse with distilled, demineralized or fully desalinated water.
- Always observe the applicable hospital guidelines concerning the supply of sterile materials.

Risk of contamination of the patient!

- Do not resterilize or reuse implants contaminated e.g. with blood or secretions, under any circumstances.

TARGON® implants to which above-mentioned restrictions do not apply can be resterilized after proper processing.

Sterilization method and parameters

Note
Observe all relevant national regulations and standards with regard to processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note
Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

- Sterilize with steam, always observing the following rule:
Sterilization must be carried out through a validated steam sterilization process (e.g. in a sterilizer compliant with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, validated according to EN 554/ISO 13683).
For the fractionated vacuum process: Sterilization has to be carried out through the 134 °C/2 bar program with a minimum holding time of 5 minutes.

Validated processing procedure

Cleaning/Disinfecting

Manual cleaning with immersion disinfecting and cleaning with brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Final rinse	RT (cold)	0.5	–	FD–W	–
IV	Drying	RT	–	–	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized)
RT: Room Temperature

- Stage I**
- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
 - Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surface.
 - Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed.
 - After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
 - To avoid the risk of corrosion, do not use a metal brush or other abrasives, which would damage the surfaces, for cleaning the product.

- Stage II**
- Rinse the product completely (all accessible surfaces) under running water.

- Stage III**
- Carry out a full rinse of the product (all accessible surfaces) under running water.
 - Allow water to drip off for a sufficient length of time.

- Stage IV**
- Dry the product completely with compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

Manual pre-cleaning with brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room Temperature

- Stage I**
- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
 - Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
 - Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed.
 - After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
 - To avoid the risk of corrosion, do not use a metal brush or other abrasives that would damage the surfaces for cleaning the product.

- Stage II**
- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	– Concentrate, alkaline: pH = 10.9 <5 % anionic tensides – 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD–W	– Concentrate, acid: pH = 2.6 Basis: Citric acid – 1 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD–W	–
V	Thermal disinfecting	94/201	10	FD–W	–
VI	Drying	90/194	40	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized)

Application

Detailed operating instructions for the TARGON® nail systems are available at any time from the relevant Aesculap agency.

The operating surgeon draws up an operating plan that specifies and appropriately documents the following steps:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the Aesculap implant system special instruments, must be complete and operational
- The operating surgeon and operating room team are thoroughly conversant with the operating technique, as well as the range of implants and implant instruments to be applied; complete information materials on these subjects are readily available at the workplace
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors
- The manufacturer has been consulted if the preoperative situation was ambiguous or if implants were found in the operating site.

The intervention has been explained to the patient, whose consent concerning the following information has been documented:

- Under certain circumstances it can happen that after a fracture is treated with an implant, limb anatomy will not be fully restored.
- Following the fracture treatment, the functionality of adjacent joints may be restricted.
- Pains may occur following the fracture treatment.
- Implants must not be subjected to extreme load arising from sports activities or strenuous physical work. The implant may loosen or suffer material fracture if subjected to excessive load.
- Overloading the implant very briefly as the consequence of an event such as a fall can cause the implant to fracture.
- The patient must be informed of the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The risks of transgressing these rules must be explained to the patient.
- Implant loosening, implant breakage, pseudarthrosis or other complications can necessitate revision surgery.
- The patient must undergo regular medical follow-up examinations.
- The attending physician decides when to remove, in part or completely, the osteosynthesis implants through surgical intervention.

WARNING

Complications due to improper application and/or incorrect fracture assessment/treatment!

- Only apply nails of appropriate length and diameter.
- Make certain that implants are aligned correctly.

To check that the targeting device with the targeting attachment is correctly mounted on the retrograde TARGON® femoral nail:

- Pass the centering pin through one of the proximal drill holes in the targeting attachment, and through the nail.

To check that the targeting device without targeting attachment is correctly mounted on the TARGON® universal and titanium locking nails:

- Use the drill sleeve and the core-hole drill.

WARNING

Risk of shattered bones and jammed or bent nails due to improper use!

- Do not use excessive force, under any circumstances, when handling/applying the instruments and implants.
- Only apply mild hammer taps to TARGON® F/T and, if necessary, RF nails to insert them into the medullary cavity.
- Check the nail insertion site, and correct if necessary.
- Should a nail be difficult to insert, drill open the medullary cavity to a smaller (for solid nails) or larger (for drilled nailing) extent.

Depending on the indication, the surgeon decides:

- if implant locking is necessary
- if the implant is to be locked, whether it will be locked statically or dynamically.

CAUTION

Drill targets missed, drill bit breakage and guide pin breakage as well as early implant failure due to crack formation caused by surface damage in the drill hole regions!

- Use a drill bit of the correct diameter.
- Level the drill insertion site with a flat counterbore or use a center punch.
- Only use sharp drill bits and guide pins and drill carefully and without applying excessive force. Avoid any deviations from the drilling direction.
- Do not exert any force on the targeting device during the drilling process. Do not bend or twist the targeting device.

Note
For locking an implant without a targeting attachment, Aesculap recommends free-hand locking supported by a C-bow and laser pilot beam, or with a radiolucent angle transmission!

Implant removal:

- When removing the TARGON® nails with the extractor instrument and the slotted hammer, only apply a quick succession of gentle taps, which is more effective than more powerful taps.
- Make certain that all locking screws have been explanted before removing the TARGON® nails.

Note

During implant removal, complications can arise due to jammed implants, grown-in bone tissue, etc. Under certain circumstances, the implants and/or the instruments could get damaged. For such cases we recommend using a set of special instruments, which can be ordered from Aesculap. This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Follow the instructions for use for such cases!

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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