- GB** Instructions for use/Technical description
Targon® FN
- D** Gebrauchsanweisung/Technische Beschreibung
Targon® FN
- F** Mode d'emploi/Description technique
Targon® FN
- E** Instrucciones de manejo/Descripción técnica
Targon® FN
- I** Istruzioni per l'uso/Descrizione tecnica
Targon® FN
- P** Instruções de utilização/Descrição técnica
Targon® FN
- NL** Gebruiksaanwijzing/Technische beschrijving
Targon® FN
- S** Bruksanvisning/Teknisk beskrivning
Targon® FN
- RUS** Инструкция по применению/Техническое описание
Targon® FN
- CZ** Návod k použití/Technický popis
Targon® FN
- PL** Instrukcja użytkowania/Opis techniczny
Targon® FN
- SK** Návod na používanie/Technický opis
Targon® FN
- TR** Kullanım Kılavuzu/Teknik açıklama
Targon® FN

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 012010 01/13 V6 Änd.-Nr. 45787

Technical alterations reserved

Aesculap® Targon®
Targon® FN

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bb.raun.com

Intended use

The Targon® FN system is used for splinting, stabilization and fixation of the proximal femur.

System configurations:

- Targon® FN femoral plate
- Targon® FN Telescrew support screw system in various lengths

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®: Titanium forged alloy Ti6Al4V acc. to ISO5832-3

The titanium implants are anodized with a colored oxide layer. Slight changes in coloration may occur, but do not affect the implant quality.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for:

- Medial femoral neck fractures

The operating surgeon is responsible for establishing the above indications and indications beyond those, taking into account the specific clinical, biological, and biomechanical situation.

Contraindications

Do not use in the presence of:

- Trochanteric femur fractures
- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side effects and interactions

- Changes in position, loosening, wear and tear or fracture of implant components
- Position change and loosening of fragments
- Delayed healing, or non-union and development of pseudarthrosis
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons, or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Pain in the implant entry region and in the region of locking components

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present instructions for use.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Targon® implants may only be implanted using the Aesculap implantation instruments specially designed for this purpose.
- Aesculap is not responsible for any complications arising from wrong indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The instructions for use for individual Aesculap implant components must be followed.
- Follow the operating manual.
- Use implant components only in the appropriate plate drill holes.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically removed components.
- Implants that have been used before must not be reused.
- Mobilization, dynamic locking or use of compression involve the risk of migration of the nail and/or the locking component into the joint region or out of the bone. Therefore the migration paths of the implants and those of the fragments must be considered in the implantation and taken into account when selecting the implant length.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete callus formation.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudarthrosis develops, or if the implant is overloaded for an excessively long period. One way to achieve this is through mobilization.
- Always insert the Targon® implants in such a way that forces to be transmitted are slight and can be transmitted by the bone at an early stage.
- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, instruments, and fragments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g. for blocked drill helix).

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant bearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.
- Increase the load on the fracture and implant in accordance with the progress of healing.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- At the time of follow-up exams, the scope of postoperative load and aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. Also, the dimensions of the implant should be taken into account.
- Recommended follow-up intervals
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.
- Use the implant system storage devices for processing, sterilization and sterile setup.
- Ensure that the implant components in their implant system storage devices do not come into contact with each other or with instruments.
- Ensure that the implant components are not damaged in any way.

Prior to resterilization, the implant components must be cleaned using the following validated reprocessing procedure:

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bb.raun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Validated reprocessing procedure

Mechanical alkaline cleaning and thermal disinfection

- Place the implants on a tray that is suitable for cleaning (avoiding rinsing blind spots).

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - < 5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

For implant components that are to be resterilized:

WARNING

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- Handle the implants with new gloves only.
- Keep the implant system storage devices covered or closed.
- Process implant system storage devices separately from instrument trays.
- Clean implants must not be processed together with contaminated implants.
- Process the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.
- Mechanically clean and disinfect the implant components.
- Do not reuse surgically contaminated implants!

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- Do not reprocess implants that have been directly or indirectly contaminated with blood.

Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Immediately set aside damaged or inoperative products.

Packaging

- ▶ Place the product in its holder or on a suitable tray.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Make sure that the packaging will prevent a recontamination of the product during storage.

Sterilization

- ▶ Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly familiar with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- Under certain circumstances it can happen that, when a fracture is treated with an implant, humerus anatomy will not be fully restored.
 - Following treatment of a fracture, shoulder and elbow joint function may be impaired.
 - Pains may occur following the fracture treatment.
 - The implants must not be exposed to overload due to excessive strains, hard physical work or sports activities. Overstraining entails the risk of implant loosening or breakage.
 - The patient must be informed about the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The risks of transgressing these rules must be explained to the patient.
 - Revision surgery may become necessary in cases of implant loosening or pseudarthrosis.
 - The patient must undergo regular medical follow-up examinations.
 - The attending physician decides when to remove, in part or completely, the osteosynthesis implants through another surgical procedure.
- ▶ Check for correct mounting of the plate on the aiming device: Insert the drill through the drill sleeve and the aiming hole and guide it through the plate hole.

WARNING

Complications due to incorrect use and/or incorrect assessment/treatment of the fracture!

- ▶ Use support screws of the appropriate length.
- ▶ Make certain that implants are aligned correctly.
- ▶ Handle the implants with appropriate care.

Explantation

Note
Implant removal can give rise to complications due to implants that are locked in place, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments. For such cases a set of special instruments is recommended, which can be ordered from Aesculap. This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such cases!

Further information on Aesculap implant systems is available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

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