Targon®

GB

Instructions for use

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Intended use

Targon® locking nails are used for intramedullary splinting, stabilization and fixation of fractures of long tubular bones.

Systems:

- Targon® PH/H
 - Proximal humerus nail
 - Humerus shaft nail
- Targon® PF, Targon® PFT
 - Proximal femur nail
- Targon® TX
 - Tibia nail

Each system consists of locking nails, locking components (e.g. screws, sleeves, pins, etc.) and associated special implantation instruments.

Materials

The following materials are used for Targon® implants.

- ISOTAN®F titanium forged alloy Ti6Al4V acc. to ISO 5832-3

For the Targon® PH system:

- Nails with additional screw brake made of PEEK-OPTIMA®
- Fixation washer made of ISOTAN®P pure titanium acc. to ISO 5832-2

The materials used in the respective implants are listed on the individual packages.

PEEK-OPTIMA® is a registered trademark of Invibio Ltd, Lancashire FY5 4QD, UK.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

The titanium implants are coated with a colored oxide. Slight color changes can occur, but will not affect the quality of the implant.

Indications

Use Targon® PH/H in cases of:

- Multi-fragment fractures of the humerus head, up to 4-fragment fracture (TARGON® PH)
- Combined humerus head and shaft fractures (TARGON® PH long)
- Stable or unstable humerus shaft fractures (TARGON® H)
- Pathological fractures
- Pseudarthroses

Use Targon® PF/Targon® PFT in cases of:

- Per-, inter- and subtrochanteric femur fractures
- Above-mentioned fractures combined with ipsilateral femoral shaft fractures (Targon® PF long)
- Lateral femoral neck fractures combined with ipsilateral femoral shaft fractures (Targon® PF long)
- Peritrochanteric fractures (lateral femoral neck fractures extending into the trochanter zones)
- Pathological fractures
- Pseudarthroses

Use Targon® TX in cases of:

- Open and closed tibial shaft fractures
- Repositioning osteotomies in the shaft region
- Reconstruction following tumor resection in the shaft region
- Osteotomies in the shaft region to lengthen or shorten the limb
- Pathological fractures
- Pseudarthroses
- Metaphyseal fractures of the tibia without articular involvement

The operating surgeon is responsible for establishing the above indications and indications beyond those, taking into account the specific clinical, biological and biomechanical situation.

Contraindications

Do not use Targon® PH/H in cases for:

- Fractures up to 7 cm proximal of the olecranon
- Crushed head cap of the humerus head (C3 fractures)
- Isolated avulsion fracture of a tuberosity

Do not use Targon® PF/Targon® PFT in cases of:

- Medial intermediary femoral neck fractures

Do not use Targon® TX as the sole implant system in cases of :

- Fractures with articular involvement

Generally, do not use in cases of:

- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side-effects and adverse interactions

- Dislocation, loosening, wear, bending and fracture of implant components
- Position change and loosening of fragments
- Delayed healing, or non-union and development of pseudarthrosis
- Primary and secondary infections
- Venous thrombosis, lung embolism, and cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Compartment syndrome
- Pain in the region of the nail insertion point and in the region of locking components

Safety notes

- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- General risk factors associated with surgical procedures are not described in the present documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Read and keep the instructions for use; observe the OP Manual for the system variants.
- Use implant components only in the nail holes intended for this purpose.
- The operating surgeon must be thoroughly familiar with the bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Prior to the operation, the surgeon must make themselves familiar with the Targon® implants and the operating technique to be used.

- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- The Targon® implants may only be implanted with the Aesculap implantation instruments intended for this purpose.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- Under no circumstances may implant components from different suppliers be combined.
- Do not use damaged or surgically excised components under any circumstances.
- Implants that have been used before must not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete callus formation.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudarthrosis develops, or if the implant is overloaded for an excessively long period. This can be achieved e.g. through dynamizing.
- Once the nail is locked, remove transmedullary bone- and fragment correction screws, or take into account, when determining the allowable patient load, that adverse lever effects can expose the nail to increased load.
- Always insert the TARGON® implants in such a way that forces to be transferred are slight and can be transmitted by the bone at an early stage.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- Never use damaged implant components or implant components that have been removed surgically.
- Always use cutting instruments that are sharp, e.g. drills, guide pins, etc.
- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, instruments and fragments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g. for blocked drill helix).
- Dynamizing, dynamic locking or use of compression involve the risk of migration of the nail and/or the locking component into the joint region or out of the bone. Therefore the migration paths of the implants and those of the fragments must be considered in the implantation and taken into account when selecting the implant length.
- Increase the load on the fracture and implant in accordance with the progress of healing.
- Aesculap and the Aesculap Academy offer regular general and product-specific training courses, where users receive helpful advice especially for problem situations. Do not hesitate to use such events.
- At the time of follow-up exams, the scope of postoperative load and aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. Also, the dimensions of the implant should be taken into account.
- Recommended follow-up intervals:
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem
- Implanting the nail via a guide pin:
 - Use an appropriate guide pin
 - Avoid jamming the implant or guide pin
 - Monitor the guide pin end by X-ray control
 - Remove the guide pin before locking the implant

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.

Risk of contamination of the patient!
➤ **Do not resterilize or reuse implants contaminated e.g. with blood or secretions, under any circumstances.**

Targon® implants that are not subject to above restrictions can be resterilized after proper processing.

Prior to the initial sterilization and subsequent resterilization, the implant components must be cleaned through the following validated processing procedure:

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, the manufacturer's specifications can only serve as an approximate guide for assessing the processing procedures applied by the individual operator/processors.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Validated processing procedure

Mechanical alkaline cleaning and thermal disinfecting

- Place the implants on a tray that is suitable for cleaning (avoiding rinsing blind spots).

Machine type: Single-chamber washer/disinfector without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	BBRAUN HELIMATIC CLEANER alkaline with tensides, application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to disinfecter program

D–W: Drinking water
FD–W: Fully desalinated water (demineralized)

WARNING

Intraoperative contamination with blood, secretions and other fluids may render the affected component unsuitable for resterilization!

- Handle the implants with new gloves only.
- Keep the implant system storage devices covered or closed.
- Put away implant system storage devices separately from instrument trays.
- Clean implants must not be processed together with contaminated implants.
- Process the implant components individually and separately if no implant system storage devices are available. When doing this, make certain that the implant components will not be damaged.
- Mechanically clean and disinfect the implant components.
- Do not reuse intraoperatively contaminated implants!

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- Do not reprocess implants that were directly or indirectly contaminated with blood.

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational and undamaged.
- Immediately sort out damaged or inoperative products.

Packaging

- Sort the product into its appropriate storage device or put it on a suitable tray.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical intervention has been explained to the patient and their agreement concerning the following information has been recorded:

- Treatment of a fracture with implants may not always be successful in completely restoring the original bone anatomy.
- Following the fracture treatment, the functionality of adjacent joints may be restricted.
- Pains may occur following the fracture treatment.
- The implants must not be subjected to overload through extreme strains, hard physical work or sports activities. The implant may loosen or suffer material fracture due to excessive strain e.g. suffered in a fall.
- The patient must be informed of the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The implants will not sustain the load and functionality of a healthy bone. The risks of transgressing these rules must be explained to the patient.
- Corrective surgery may become necessary in cases of implant loosening, fragment migration, pseudarthroses or other complications.
- The patient must undergo regular medical follow-up examinations.
- The attending physician decides when to remove, in part or completely, the osteosynthesis implants through surgical intervention.

WARNING

Complications due to incorrect use and/or incorrect assessment/treatment of the fracture!

- Make certain that nails of a suitable length and diameter are used.
- Make certain that implants are aligned correctly.

- Check for correct mounting of the nail on targeting device: Insert the drill through the drill sleeve and the targeting hole and guide it through the nail hole.

WARNING

Risk of nail seizure, nail bending or bone rupture!

- Do not apply excessive force on instruments or implants.
- Apply only mild hammer taps to insert Targon® TX and, if necessary, PF/PFT nails into the medullary cavity.
- To avoid bone ruptures, do not impact Targon® PH/H nails under any circumstances.
- Check and, if necessary, correct the nail insertion point.
- If a nail proves difficult to insert, slightly drill open the medullary cavity.

Depending on the indication, the surgeon decides:

- whether implant locking is necessary
- if the implant is to be locked, whether it will be locked statically or dynamically.

Note

For locking an implant without a targeting attachment, Aesculap recommends free-hand locking supported by a C-bow and laser pilot beam, or with a radiolucent angle transmission!

WARNING

Avoid drill mistargeting, drill breakage and guide pin breakage as well as premature implant failure due to crack formation as a result of damaged surfaces in the drill hole regions!

- Use a drill bit of the correct diameter.
- Level the drill insertion site with a flat counterbore or use a center punch.
- Only use sharp drill bits and guide pins and drill carefully and without applying excessive force. Avoid direction deviations.
- During center punching, flat countersinking and drilling, do not apply force to the targeting device; do not bend or twist the targeting device.

Implant removal

- Prior to removing the nail, remove all screws.
- Remove the nail with quick, light taps, and with the aid of the impactor, extraction pin and slotted hammer.

Note

Implant removal can give rise to complications due to implants that are stuck, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments. For such cases a set of special instruments is recommended, which can be ordered from Aesculap. This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such case!

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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