Aesculap®

TSPACE® PEEK/TSPACE® XP/TSPACE® 3D

Aesculap Spine

en Instructions for use/Technical description
Instruments

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapImplantsystems.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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



Instruments

Legend

- A Insertion instrument SN305R (assembled)
- B Insertion instrument SN305R (disassembled):
 - 1 Tube
 - 2 Switch
 - 3 Rotary nut
 - 4 Irrigation connector
 - 5 Pin
 - 6 Rod
 - 7 Dihedron
- C Impact cylinder SN320R
- D Trial implant
- E Filling holder SN304R

Contents

1.	About this document.....	2
1.1	Scope	3
1.2	Safety messages.....	4
2.	Clinical use.....	4
2.1	Areas of use and limitations of use	4
2.1.1	Intended use.....	4
2.1.2	Indications	4
2.1.3	Contraindications	4
2.2	Safety information.....	4
2.2.1	Clinical user	4
2.2.2	Product.....	4
2.2.3	Sterility.....	4
2.3	Application.....	5
2.3.1	Preparing the implant bed	5
2.3.2	Determining the implant height with the trial implant.....	5
2.3.3	Inserting the implant in the implant bed ..	7
3.	Validated reprocessing procedure	8
3.1	General safety information.....	8
3.2	General information	9
3.3	Reusable products	9
3.4	Preparations at the place of use.....	9
3.5	Preparing for cleaning.....	9
3.6	Disassembly	9
3.7	Cleaning/Disinfection	10

3.7.1	Product-specific safety information on the reprocessing method	10
3.7.2	Validated cleaning and disinfection procedure	10
3.8	Manual cleaning/disinfection	11
3.8.1	Manual cleaning with immersion disinfection.....	12
3.8.2	Manual cleaning with ultrasound and immersion disinfection	13
3.9	Mechanical cleaning/disinfection.....	14
3.9.1	Mechanical alkaline cleaning and thermal disinfection.....	14
3.10	Mechanical cleaning/disinfection with manual pre-cleaning	15
3.10.1	Manual pre-cleaning with a brush.....	15
3.10.2	Manual pre-cleaning with ultrasound and brush.....	16
3.10.3	Mechanical alkaline cleaning and thermal disinfecting.....	17
3.11	Inspection	17
3.11.1	Visual inspection	17
3.11.2	Functional test.....	18
3.12	Assembly.....	18
3.13	Packaging	18
3.14	Steam sterilization.....	18
3.15	Storage	18
4.	Technical service	18
5.	Disposal.....	19
6.	Symbols on product and packaging	19

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply for the following products:

Art. no.	Designation
SN304R	TSPACE PEEK/TSPACE ^{XP} filling holder
SN305R	TSPACE PEEK/TSPACE ^{XP} insertion instrument
SN320R	TSPACE PEEK/TSPACE ^{XP} impact cylinder
SN704R	TSPACE 3D filling holder

TSPACE PEEK/TSPACE^{XP}/TSPACE 3D trial implants 26 mm

Art. no.	Designation
SN322R	Trial implant 5° 26 x 11.5 x 7 mm
SN323R	Trial implant 5° 26 x 11.5 x 8 mm
SN324R	Trial implant 5° 26 x 11.5 x 9 mm
SN325R	Trial implant 5° 26 x 11.5 x 10 mm
SN326R	Trial implant 5° 26 x 11.5 x 11 mm
SN327R	Trial implant 5° 26 x 11.5 x 12 mm
SN328R	Trial implant 5° 26 x 11.5 x 13 mm
SN329R	Trial implant 5° 26 x 11.5 x 14 mm
SN330R	Trial implant 5° 26 x 11.5 x 15 mm
SN332R	Trial implant 5° 26 x 11.5 x 17 mm

TSPACE PEEK/TSPACE^{XP}/TSPACE 3D trial implants 30 mm

Art. no.	Designation
SN352R	Trial implant 5° 30 x 11.5 x 7 mm
SN353R	Trial implant 5° 30 x 11.5 x 8 mm
SN354R	Trial implant 5° 30 x 11.5 x 9 mm
SN355R	Trial implant 5° 30 x 11.5 x 10 mm
SN356R	Trial implant 5° 30 x 11.5 x 11 mm
SN357R	Trial implant 5° 30 x 11.5 x 12 mm

Art. no.	Designation
SN358R	Trial implant 5° 30 x 11.5 x 13 mm
SN359R	Trial implant 5° 30 x 11.5 x 14 mm
SN360R	Trial implant 5° 30 x 11.5 x 15 mm
SN362R	Trial implant 5° 30 x 11.5 x 17 mm

TSPACE PEEK/TSPACE^{XP}/TSPACE 3D trial implants 34 mm

Art. no.	Designation
SN382R	Trial implant 5° 34 x 11.5 x 7 mm
SN383R	Trial implant 5° 34 x 11.5 x 8 mm
SN384R	Trial implant 5° 34 x 11.5 x 9 mm
SN385R	Trial implant 5° 34 x 11.5 x 10 mm
SN386R	Trial implant 5° 34 x 11.5 x 11 mm
SN387R	Trial implant 5° 34 x 11.5 x 12 mm
SN388R	Trial implant 5° 34 x 11.5 x 13 mm
SN389R	Trial implant 5° 34 x 11.5 x 14 mm
SN390R	Trial implant 5° 34 x 11.5 x 15 mm
SN392R	Trial implant 5° 34 x 11.5 x 17 mm

TSPACE 3D trial implants 38 mm

Art. no.	Designation
SN783R	Trial implant 5° 8 x 11.5 x 38 mm
SN784R	Trial implant 5° 9 x 11.5 x 38 mm
SN785R	Trial implant 5° 10 x 11.5 x 38 mm
SN786R	Trial implant 5° 11 x 11.5 x 38 mm
SN787R	Trial implant 5° 12 x 11.5 x 38 mm
SN788R	Trial implant 5° 13 x 11.5 x 38 mm
SN790R	Trial implant 5° 15 x 11.5 x 38 mm

Note

The applicable CE mark for the product can be seen on the label or packaging of the product.

- For article specific instructions for use and material compatibility and lifetime information, see B. Braun eFU at eifu.bbraun.com

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The TSPACE PEEK/TSPACE^{XP}/TSPACE 3D instruments are used to insert TSPACE PEEK (third generation), TSPACE^{XP} and TSPACE 3D trial implants and implants into the thoracolumbar spine.

The SN305R insertion instrument is used for implanting implants and for determining the implant size by inserting trial implants into the intervertebral space.

The impact cylinder SN320R is coupled to the insertion instrument. The impact cylinder is used to expel the trial implants from the intervertebral space or to revise the implants.

The filling holders SN304R (for TSPACE PEEK/TSPACE^{XP}) and SN704R (for TSPACE 3D) are used to hold the implants on the surgical table while filling the inner window with bone material.

2.1.2 Indications

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

For indications, see Intended use.

2.1.3 Contraindications

No known contraindications.

2.2 Safety information

2.2.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- Use the product only according to these instructions for use.
- Follow the safety and maintenance instructions.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to use, check that the product is in good working order.
- Keep the instructions for use accessible for the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.2.2 Product

Product-specific safety information

- Only use the product with visual control.

2.2.3 Sterility

The product is delivered in an unsterile condition.

- Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.3 Application

⚠ WARNING

Risk of injury and/or malfunction!

- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Always carry out a function test prior to each use of the product.

2.3.1 Preparing the implant bed

⚠ WARNING

Risk of injury to nerves, blood vessels, and tissue caused by sharp instruments!

- To avoid injuries, proceed with extreme care during the operation.
- Resect the facet joint.
- Open the intervertebral disk space.
- Restore the height of the intervertebral disk using distractors.
- Clear the intervertebral disk space using rongeurs, sharp spoons and curettes.
- If necessary, apply bone material.

2.3.2 Determining the implant height with the trial implant

Pick up the trial implant with the insertion instrument

- Turn the rotary nut **3** on the insertion instrument **A** counterclockwise ("loosen" direction) to push the rod **6** forward, see Fig. 1.


Fig. 1 Close-up view of the insertion instrument, switch open

- Couple the trial implant **D** to the dihedron **7** on the insertion instrument, **A** see Fig. 2. When doing so, pay attention to the orientation of the trial implant (orientation according to the marking on the insertion instrument).


Fig. 2 Coupling the trial implant to the insertion instrument

Note

The switch **2** can only be turned when the rod is in the front position and the pin **5** is fully visible in the switch.

- Turn the switch **2** to the right (direction ) , see Fig. 3.


Fig. 3 Close-up view of the insertion instrument, switch closed

- Hold the trial implant **D** at 0° position on the insertion instrument **A**, see Fig. 8.
- Turn the rotary nut **3** on the insertion instrument **A** clockwise ("tighten" direction) until the trial implant is clamped, see Fig. 4.


Fig. 4 Close-up view of the insertion instrument, trial implant/implant clamped

- Check manually that the trial implant is securely in place.

Inserting the trial implant

- Insert the trial implant **D** with the insertion instrument **1** straight (0° position) into the prepared implant bed. To do so, carefully hit with a hammer on the impact surface of the rotary nut **3**.
- To turn the trial **D** implant:
 - Gently turn the rotary nut **3** counterclockwise ("loosen" direction). When doing so, make certain that the pin **5** remains outside of the switch **2**.
 - Impact the trial implant **D** further with the hammer until the required end position is reached.
- Check the position of the trial implant **D** using intraoperative X-ray control.

Removing the trial implant

- Mounting the impact cylinder **C** onto the insertion instrument **A**, see Fig. 5.


Fig. 5 Mounting the impact cylinder onto the insertion instrument

- Carefully expel the trial implant **D** using the impact cylinder **C** out of the intervertebral space, see Fig. 6.


Fig. 6 Expelling the trial implant/implant with the impact cylinder

- Uncoupling the impact cylinder **C** from the insertion instrument **A**.

Uncoupling the trial implant from the insertion instrument

- Turn the rotary nut **3** counterclockwise ("loosen" direction), until the pin **5** is fully in the switch **2**, see Fig. 3. When doing so, ensure that the pin **5** does not jam the switch **2**.
- Turn the switch **2** to the left (direction ) , see Fig. 1.
- Remove the trial implant.
- If necessary, repeat the process with the next biggest trial implant, until the correct implant size is determined.

2.3.3 Inserting the implant in the implant bed

Picking up the implant with the insertion instrument

- ▶ Select an implant according to the size of the trial implant.
- ▶ Remove the implant with protective packaging from the sterile packaging.
- ▶ Mount the implant **a** in the protective packaging onto the dihedron **7** of the insertion instrument **A**, see Fig. 7. When doing so, pay attention to the orientation of the implant (orientation according to the marking on the insertion instrument).


Fig. 7 Coupling the implant to the insertion instrument

Legend

a Implant

- ▶ Turn the switch **2** to the right (direction ) , see Fig. 3.
- ▶ Hold the implant at 0° position, see Fig. 8.


Fig. 8 Implant at 0° position

- ▶ Turn the rotary nut **3** on the insertion instrument **A** clockwise (:"tighten" direction) until the trial implant is clamped, see Fig. 4.
- ▶ Remove the protective packaging from the implant.
- ▶ Check manually that the implant is securely in place.

Filling the implant with bone material (optional)

- ▶ Place implant **a** in filling holder **E**, see Fig. 9. When doing so, observe the implant length (marking on the filling holder: L26/L30/L34).


Fig. 9 Implant in filling holder

- ▶ Compress the bone grafts or bone replacement material in the inner window of the implant **a** using the pestle provided in the set. When doing so, ensure that the pestle does not damage the surface of the implant.

Inserting the implant

- ▶ Insert the implant **a** with the insertion instrument **A** straight (0° position) into the prepared implant bed. To do so, carefully hit with a hammer on the impact surface of the rotary nut **3**.


Fig. 10 Inserting the implant into the implant bed

- ▶ To turn the implant **a**:
 - Gently turn the rotary nut **3** counterclockwise ("loosen" direction). When doing so, make certain that the pin **5** remains outside of the switch **2**.
 - Impact the implant **a** further with the hammer until the required end position is reached.
- ▶ Check the position of the implant **a** using intraoperative X-ray control.
- ▶ Once the end position (see Fig. 11) is reached: turn the rotary nut **3** counterclockwise ("loosen" direction), until the pin **5** is fully in the switch **2**, see Fig. 3. When doing so, ensure that the pin **5** does not jam the switch **2**.


Fig. 11 Implant in end position

- ▶ Turn the switch **2** to the left (direction , see Fig. 1.
- ▶ Carefully remove the insertion instrument **A** from the implant **a**. Avoid any tilting or canting of the instrument when doing so.
- ▶ Place bone or bone replacement material around the implant **a**.

Removing the implant intraoperatively

If, in the opinion of the surgeon, it becomes necessary to remove the implant intraoperatively, carry out the following steps:

- ▶ Check the position for coupling the insertion instrument **A**, see Fig. 1.
- ▶ Couple insertion instrument **A** to implant **a** in situ, see Fig. 7.
- ▶ Turn the switch **2** to the right (direction , see Fig. 3.
- ▶ Mounting the impact cylinder **B** onto the insertion instrument **A**, see Fig. 5.
- ▶ Carefully expel the implant **a** using the impact cylinder **C** out of the intervertebral space, see Fig. 6.

3. Validated reprocessing procedure

3.1 General safety information

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see B. Braun eIFU at eifu.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures $>45^{\circ}\text{C}$ nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of $\text{pH} > 8$ to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to "AKI-Brochures", "Red brochure".

3.3 Reusable products

- There is no set maximum number of uses and preparation cycles for the product.
- The life of the product is limited by damage, normal wear and tear, type and duration of use, as well as handling, storage and transport of the product.
- A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional.

3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Preparing for cleaning

- Disassemble the product prior to cleaning, see Disassembly.
- For mechanical cleaning: mount the irrigation connector 4 to the rear end of the rotary nut 3.

3.6 Disassembly

SN305R insertion instrument

- Turn the rotary nut 3 counterclockwise ("loosen" direction).
- Turn the switch 2 to the left (direction , see Fig. 1).
- Turn the rotary nut 3 clockwise ("tighten" direction), until the rod 6 is out of the thread.
- Remove the rod 6 from the tube 1, see Fig. 12.


Fig. 12 Insertion instrument disassembled

3.7 Cleaning/Disinfection

3.7.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents
 - that are approved for plastic material and high-grade steel.
 - that do not attack softeners (e.g., in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable disinfection temperature of 95 °C.
- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water

- ▶ SN304R filling holder: manually pre-clean the product (with a cleansing brush) if there are residues of bone, tissue or ancillary materials (such as plaster or bone cement).

3.7.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ SN304R ■ Trial implants	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Drying phase: use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ SN305R ■ SN320R	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ SN304R	■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfection

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ SN320R ■ Trial implants	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ SN305R	■ Suitable cleaning brush ■ 20 ml disposable syringe ■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

3.8 Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfection process if necessary.

3.8.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.8.2 Manual cleaning with ultrasound and immersion disinfection

Phase	Step	D [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	5	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.9 Mechanical cleaning/disinfection

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.9.1 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfecting machine without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant 0.5 % working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.10 Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.10.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.10.2 Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
 RT: Room temperature
 *Recommended: BBraun Stabimed fresh

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

3.10.3 Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 <5 % anionic surfactant 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.
- ▶ Check the product for burrs that could damage tissue or surgical gloves.
- ▶ Check the product for loose or missing parts.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.11 Inspection

- ▶ Allow the product to cool down to room temperature.
- ▶ Dry the product if it is wet or damp.

3.11.1 Visual inspection

- ▶ Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- ▶ If the product is dirty: repeat the cleaning and disinfection process.
- ▶ Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- ▶ Check the product for missing or faded labels.
- ▶ Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- ▶ Check the surfaces for rough spots.

3.11.2 Functional test

CAUTION

Damage (metal cold welding /friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).
- ▶ Assemble disassembled products, see Assembly.
- ▶ Check that the product functions correctly.
- ▶ Check that all moving parts are working property (e.g. hinges, locks/latches, sliding parts etc.).
- ▶ Check for compatibility with associated products.
- ▶ Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.12 Assembly

SN305R insertion instrument

- ▶ Turn the switch **2** to the left (direction .
- ▶ Insert the rod **6** into the tube **1**, see Fig. 12.
- ▶ Turn the rotary nut **3** counterclockwise ("loosen" direction) to push the rod **6** forward, see Fig. 1.

3.13 Packaging

- ▶ Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
- ▶ Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.14 Steam sterilization

Note

SN305R insertion instrument may only be sterilized in disassembled state.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the product
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 C/holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.15 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4. Technical service

CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ **Do not modify the product.**
- ▶ **For service and repairs, please contact your national B. Braun/Aesculap agency.**

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

5. Disposal

WARNING

Risk of infection due to contaminated products!

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany