Aesculap Sterile Technology

GB

Instructions for use

Important advice on the sterile container system

USA

Instructions for use

Note for U.S. users

U.S. cleared device instructions for use can be found at

www.aesculapusa.com "Products & Services".

If you do not have internet access, please call your local Aesculap representative or customer service at 1-800-258-1946.

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www.aesculap.de

TA-Nr. 013139 11/10 Änd.-Nr. 42641

CE marking according to directive 93/42/EEC

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Sterrad® 100 S, Sterrad® 200, Sterrad® NX and Sterrad® 100 NX are registered trademarks of ASP.

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Sterile container system

Note

These brief instructions for use apply to standard mini, dental and optics containers.

Note

Observance of instructions for use C48381 (Revision 07/2009) of the standard container and C29581 (Revision 11/2010) of the mini, dental and optics container is absolutely imperative to ensure safe use of the respective containers. If the most recent edition of the instructions for use is not available to you, please contact Aesculap.

Intended use

The Aesculap sterile container system is a reusable sterile barrier system that preserves the sterility of medical products until they are used or reach their use-by date.

Medical products can be sterilized, stored and transported in the Aesculap sterile container system. It can also be used for transporting medical products back to sterile reprocessing after use.

The Aesculap sterile container system is suitable for steam sterilization. The optics container is also suitable for sterilization in ethylene oxide.

Specially marked containers (Sterilcontainer S) can also be sterilized in Sterrad® 100 S, Sterrad® 200, Sterrad® NX, Sterrad® 100 NX, V-Pro® 1 and V-Pro® 1 Plus.

WARNING

Only combine Aesculap sterile container products with each other!

Combining Aesculap sterile container products with components from other manufacturers will lead to the risk of container leaks.

The functionality of Aesculap sterile containers as germ barriers during transport and storage of sterile materials cannot be guaranteed if Aesculap sterile container products are combined with components from other manufacturers.

Function check

- Visually inspect all components for any damage and correct functioning (seals intact, plastic parts free from cracks).

Cleaning and disinfection

WARNING

Risk of corrosion and tension cracks due to inappropriate cleaning and disinfection!

- When cleaning and disinfecting, always observe the specifications given by the manufacturers of the cleaning agents and disinfectants.
- Exclusively use pH-neutral or mildly alkaline cleaning agents.

WARNING

Additional instructions for PrimeLine:

- Do not use solvents such as acetone or benzene.
- Exclusively use fully demineralized water for the final rinse.
- Make certain that no residues from the cleaning process remain on the lid.

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V-Pro® 1 and V-Pro® 1 Plus are registered trademarks of Steris.