- GB

Instructions for use/Technical description
- USA

Reusable electrode handles (monopolar)
- D

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Technical alterations reserved

Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Intended use

The reusable electrode handles (monopolar), which are fitted with a fixed cable, are used in open surgical procedures. The reusable electrode handles (monopolar) are used to conduct the required HF current from the HF device to the operating site, to hold the required working electrode and, if applicable, to activate the cutting or coagulating current from the HF device (handles with activation keys).

Product variants

The reusable electrode handles (monopolar) are available with various handle shapes and connectors and in different sizes (not all combinations available).

Handle shapes

Fig. 1 Narrow electrode handle without activation keys

Fig. 2 Broad electrode handle with activation keys

Fig. 3 Narrow electrode handle with activation keys

Pull out plug

- Check the plug for compatibility with the socket of the HF device.

Designation	Illustration
3-pin plug (e.g. for Aesculap GN300 and GN640 or Valleylab devices)	
Koax plug (e.g. for Aesculap GN300 and GN640)	
Metal plug (e.g. for Aesculap GK160, GK170 or GK450)	
5 mm plug ERBE (e.g. for ERBE-ICC units)	
8 mm plug (e.g. Valleylab units)	

Electrode holders


Fig. 4 Hexagon anti-twist guard

The electrode holders have a hexagon anti-twist guard. The electrode handles can hold electrodes of the following diameters (depending on variant):

- 1.6 mm
- 1.7 mm
- 2.4 mm
- 4.0 mm

Note

Further information about reusable electrode handles (monopolar) can be found in Aesculap brochure C-304-81.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).
- Disconnect the plug connection only by pulling at the plug.

WARNING

Risk of injury from ignition or explosion of flammable gases!
Sparks may occur when using the HF device as directed.

- Observe the safety guidelines in the instructions for use of the HF device.

WARNING

Thermal injuries to patients/users due to insufficient insulation of leads in active accessories!

- Adjust the HF device to an appropriate setting to ensure that the peak output voltage does match or not exceed the accessory voltage rating specified for the product.

The manufacturer has tested the product and verified that its insulation can withstand 100 reprocessing cycles. In clinical practice, the service life will depend on the individual intraoperative usage and the hospital's specific reprocessing conditions.

- Prior to each use, inspect products for: damage or surface changes to the insulation.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Adjust the HF power output to the intended surgical intervention. Take into account clinical experience or reference values.
- Select the lowest possible HF power output.
- Keep the product's contact surfaces clean during surgery. Remove encrusted tissue residues or body fluids with a moistened swab.

The product is furnished with the following connection on the plug side: see Product variants.

The accessory voltage rating of the product is 6 000 Vp.

The accessory voltage rating must exceed or match the peak output voltage with which the product is operated in combination with a suitable HF device at an appropriate operating mode/setting (see IEC/DIN EN 60601-2-2).

To avoid HF burns:

- Guide the handle cable to the operating site in such a way that it is not in contact with the patient or with other leads.
- Always keep the working end of the product in the user's field of vision whenever the HF power is activated.
- Prior to activating the HF device, check that the working end of the product is not touching any electrically conductive accessories.
- Never place the product on or next to the patient.
- When using accessories for endoscopy or laparoscopy, deactivate the automatic switch-on mode of the HF device.
- Follow the instructions for use of the HF device.

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

Connecting the working electrode

- Insert the working electrode in the handle so that the hexagon anti-twist guard is effective.

Fig. 5 Inserting the working electrode in the handle

Connecting the handle at the HF device

- ▶ Insert the plug in the appropriate monopolar output socket of the HF device.


Fig. 6 E.g. inserting the coax plug in the monopolar output socket (GN640)

Activate HF

When using a handle with activation keys, the HF output can be activated by pressing the keys.

- ▶ Press the yellow key: activation of the "Cutting" operating mode.
- ▶ Press the blue key: activation of the "Coagulating" operating mode.

Note

For correct functioning of the HF device, configure the device for handle activation according to the instructions for use of the HF device!

When using a handle without keys, the HF output of the HF device is controlled by a foot switch connected to the HF device.

Validated reprocessing procedure

Note

Note the validated processing procedure brochure (AVA-V6), available separately, see Validated cleaning and disinfection procedure.

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- ▶ Remove the working electrodes from the handle.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- ▶ Disassemble the product prior to cleaning.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for use, for example, on aluminum, plastic materials, and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- ▶ Do not exceed the maximum permitted cleaning temperature of 93 °C.

Immersion treatment in a 3 % H₂O₂ solution for approx. 5 minutes is a particularly effective and gentle method to dissolve encrustations from HF instruments. The debris can be removed manually with a soft brush. This is followed by the conventional reprocessing steps.

- ▶ Do not clean the product in an ultrasound bath.
- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Mechanical alkaline cleaning and thermal disinfection	■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).	Chapter Mechanical cleaning/disinfecting and sub-chapter: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	■ Concentrate, alkaline: – pH = 13 – <5 % anionic surfactant ■ 0.5 % working solution – pH = 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disin- fecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Place the product in its holder or on a suitable tray.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the instrument
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
 - Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.
- To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

Risk of injury and/or malfunction!

► Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Note
Information about accessories and spare/replacement parts can be found in Aesculap brochure C-304-81.

Technical Data

Classification acc. to Directive 93/42/EEC

Designation	Class
Reusable electrode handles (monopolar) without keys	IIa
Reusable electrode handles (monopolar) with keys	IIb

Conforming to standard

IEC/DIN EN 60601-2-2

Ambient conditions

	Storage and transport conditions
Relative humidity	0 % to 75 %, no condensation
Ambient temperature	-20 °C to +50 °C
Atmospheric pressure	500 hPa to 1 600 hPa
Direct sunlight	Protect against direct sunlight

Disposal

Note
The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Always adhere to national regulations when disposing of or recycling the product or its components!

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass contains instructions on disassembly of the product, as well as information for proper disposal of environmentally hazardous components)

Any product carrying this symbol must be disposed of separately through dedicated electrical and electronic devices recycling. Within the EU, such disposal is taken care of by the manufacturer free of charge.

Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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