





Aesculap® Acculan 4

- (GB)** Instructions for use/Technical description
- (USA)** Attachments
- (D)** Gebrauchsanweisung/Technische Beschreibung
Aufsätze
- (F)** Mode d'emploi/Description technique
Embouts
- (E)** Instrucciones de manejo/Descripción técnica
Cabezales
- (I)** Istruzioni per l'uso/Descrizione tecnica
Terminali
- (P)** Instruções de utilização/Descrição técnica
Cabeçotes
- (NL)** Gebruiksaanwijzing/Technische beschrijving
Opzetstukken

- (S)** Bruksanvisning/Teknisk beskrivning
Tillsatser
- (RUS)** Инструкция по применению/Техническое описание
Насадки
- (CZ)** Návod k použití/Technický popis
Násadce
- (PL)** Instrukcja użytkowania/Opis techniczny
Nasadki
- (SK)** Návod na použitie/Technický opis
Nadstavec
- (TR)** Kullanım Kılavuzu/Teknik açıklama
Başlıklar
- (KR)** 사용 설명서 / 기술 설명
부착 장치

Aesculap® Acculan 4

Attachments

Legend

- 1 Drill
- 2 Knob (for speed control)
- 3 Safety catch
- 4 Rotating sleeve
- 5 Arrow
- 6 Reamer attachment with quick-action chuck adapter
- 7 Unlocking sleeve
- 8 Tool coupling
- 9 Drill attachment with quick-action chuck adapter
- 10 Drill attachment keyless Jacobs chuck
- 11 Locking sleeve
- 12 Attachment with Jacobs chuck
- 13 Tightening key
- 14 Zimmer chuck with Hudson/Zimmer coupling
- 15 Unlocking sleeve
- 16 Adapter AO large to DHS/DCS step drill
- 17 Sagittal saw attachment
- 18 Tool attachment with L sagittal tool coupling
- 19 Knob for tool release
- 20 Drill attachment for Synthes radiolucent angle transmission
- 21 Kirschner wire attachment
- 22 Clamping lever
- 23 Adjustment sleeve
- 24 Kirschner wire chuck
- 25 L sagittal saw blade
- 26 Oil spray adapter

The diagrams are only schematic.

Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
	Machine-readable, two-dimensional code The code contains a unique serial number which can be used for electronic tracking of the individual instrument. The serial number is based on the global standard sGTIN (GS1).
	Date of manufacture
	Manufacturer's batch designation
	Manufacturer's serial number
	Manufacturer's article number
	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage

Contents

1.	Applicable to.....	3
2.	General information.....	3
2.1	Intended use.....	3
2.2	Main functions and design characteristics	4
2.3	Indications	4
2.4	Absolute contraindications	4
2.5	Relative contraindications.....	4
3.	Safe handling	4
4.	Product description	5
4.1	Scope of supply.....	5
4.2	Components required for operation	5
4.3	Operating principle.....	5
5.	Preparation	5
6.	Working with the device	5
6.1	System set-up.....	6
6.1.1	Connecting the accessories	6
6.2	Protection against inadvertent activation	6
6.3	Coupling and uncoupling of attachment to drill and reamer..	6
6.4	Attaching and detaching the tool from the attachment	6
6.4.1	Attachments with quick-action chuck GB623R/GB657R/ GB663R/GB664R/GB665R/GB668R/GB669R/GB670R	6
6.4.2	Drill attachment keyless Jacobs chuck GB620R	7
6.4.3	Attachments with three-jaw chuck GB621R/GB639R/GB667R	7
6.4.4	Zimmer chuck with Hudson/Zimmer shaft GB176R	7
6.4.5	Harris chuck with AO large shaft GB184R and DIN chuck with AO large shaft XF457R.....	7
6.4.6	Adapter AO large to DHS/DCS step drill GB628R	7
6.4.7	Sagittal saw attachment GB660R.....	7
6.4.8	Drill attachment for Synthes radiolucent angle transmission GB645R.....	8
6.4.9	Insert and clamp Kirschner wire attachment GB641R	8
6.5	Function checks	8
6.6	Safe operation	9
7.	Validated reprocessing procedure.....	10
7.1	General safety notes	10
7.2	General information.....	10
7.3	Preparations at the place of use.....	10
7.4	Preparation before cleaning	11
7.5	Cleaning/disinfection.....	11
7.5.1	Product-specific safety instructions for the reprocessing procedure	11
7.6	Manual cleaning with wipe disinfection	12
7.7	Automatic cleaning/disinfection with manual pre-cleaning ..	13
7.7.1	Manual pre-cleaning with a brush	13
7.7.2	Mechanical alkaline cleaning and thermal disinfection	14
7.8	Inspection, maintenance and checks	14

7.9	Packaging	14
7.10	Steam sterilization.....	14
7.11	Sterilization for the US market.....	15
7.12	Storage	15
8.	Maintenance	15
9.	Troubleshooting list.....	16
10.	Technical Service	18
11.	Accessories/Spare parts.....	18
12.	Technical data	18
12.1	Classification acc. to Directive 93/42/EEC.....	18
12.2	Performance data, information about standards	19
12.3	Operating mode	21
12.4	Environmental conditions	21
13.	Disposal	21
14.	Distributor in the US/Contact in Canada for product information and complaints	21

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

2. General information

2.1 Intended use

Task/Function

The drill and reamer GA330, combined with the relevant attachment and tool is used to process hard tissue, cartilage and the like, as well as bone replacement material, to screw bone pins in and out and to set Kirschner wires.

Application Environment

The product is used in operating rooms in sterile zones outside of the explosion risk zone (such as areas with pure oxygen or anesthesia gases).

Aesculap® Acculan 4

Attachments

2.2 Main functions and design characteristics

Speed	GB667R GB668R GB669R GB670R	0 min ⁻¹ to max. 250 min ⁻¹
	GB657R	0 min ⁻¹ to max. 370 min ⁻¹
	GB645R	0 min ⁻¹ to max. 900 min ⁻¹
	GB620R GB621R GB663R GB664R GB665R	0 min ⁻¹ to max. 1 000 min ⁻¹
	GB623R GB639R GB641R	0 min ⁻¹ to max. 1 250 min ⁻¹
	GB660R	0 min ⁻¹ to max. 17 000 min ⁻¹
Rotational direction	Right and left rotation, oscillation	
Operating mode	Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1) Right/left rotation: ■ 60 second application, 60 second pause ■ 20 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C Milling (clockwise/counterclockwise): ■ 30 second application, 30 second pause ■ 8 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C Drilling (oscillation): ■ 60 second application, 60 second pause ■ 4 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C Saw mode with GB660R: ■ 30 second application, 60 second pause ■ 3 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C	

2.3 Indications

The type and area of application depend on the tool selected.

2.4 Absolute contraindications

The product is not licensed for use on the central nervous system or central circulatory system.

2.5 Relative contraindications

The safe and effective use of the product greatly depends on influences which can only be controlled by the user. Therefore the specifications provided represent framework conditions only.

Clinically successful use of the product is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

3. Safe handling

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

WARNING

Risk of injury and material damage due to inappropriate use of the product!

► Use the product only in accordance with the intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

► Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
 - Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
 - Prior to use, check that the product is in good working order.
 - To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.

- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.
- Ensure that the electrical installation of the room is consistent with the requirements of IEC/DIN EN.
- Do not operate the product in explosion-hazard areas.
- Sterilize product before use.
- When handling, observe instructions for use that are relevant to Aesculap holder systems TA009721, see Aesculap Extranet at <https://extranet.bbraun.com>

4. Product description

4.1 Scope of supply

Art. no.	Designation
GB176R	One of the following attachments/adapters:
GB184R	Zimmer chuck with Hudson/Zimmer shaft
GB620R	Harris chuck with AO large shaft
GB621R	Drill attachment keyless Jacobs chuck
GB623R	Drill attachment large Jacobs chuck
GB628R	Drill attachment AO small
GB639R	Adapter AO large to DHS/DCS step drill
GB641R	Drill attachment small Jacobs chuck
GB645R	Kirschner wire attachment
GB657R	Drill attachment for Synthes radiolucent angle transmission
GB660R	Medullary reamer attachment AO large
GB663R	Sagittal saw attachment
GB664R	Drill attachment Trinkle
GB665R	Drill attachment Aesculap hexagon
GB667R	Drill attachment Hudson/Zimmer
GB668R	Reamer attachment large Jacobs chuck
GB669R	Reamer attachment AO large
GB670R	Reamer attachment Hudson/Zimmer
XF457R	Reamer attachment Harris
GA031R	DIN chuck with AO large shaft
GA062R	Chuck key for large Jacobs chuck (GB620R, GB621R and GB667R)
TA014541	Key for small Jacobs chuck (GB639)
	Instructions for use for Acculan 4 attachments (flyer)

4.2 Components required for operation

- Drill and reamer GA330 (ready)

4.3 Operating principle

The attachment can be coupled to the drill and reamer in three different positions, each offset by 120°.

The attachment locks automatically when attached to the drill and reamer. The attachment can be removed by operating a rotating sleeve on the drill and reamer.

The attachments have various integrated couplings on the working end, to enable attachment of various tools, Kirschner wires or adapters.

5. Preparation

Aesculap assumes no liability if the following rules are not followed:

- Do not use products from open or damaged sterile packaging.
- Prior to use, inspect the product and its accessories for any visible damage.
- Use the products and their accessories only if they are in perfect technical condition.

6. Working with the device

WARNING

Risk of infection and contamination!

Product is delivered unsterilized!

- Sterilize the product before use pursuant to the operating instructions.

WARNING

Risk of injury and material damage due to inappropriate use of tools!

- Always follow the safety advice and information given in the instructions for use.
- When coupling/uncoupling, handle tools with cutting edges with care.

WARNING

Damage to the product if dropped!

- Use the products only if they are in perfect technical condition, see Function check

WARNING

Risk of burns to skin and tissue through blunt tools or if product has not been maintained properly!

- Use tools only if they are in perfect condition.
- Replace blunt tools.
- Maintain the product properly, see maintenance guide.

Aesculap® Acculan 4

Attachments

6.1 System set-up

6.1.1 Connecting the accessories

DANGER

Risk of injury due to unapproved configuration using additional components!

- Ensure that the classification of all components used is consistent with the classification of the product (such as type BF or type CF).

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

All configurations must fulfill the fundamental standard IEC/DIN EN 60601-1. The person connecting the devices with each other is responsible for the configuration and must ensure that the fundamental standard IEC/DIN EN 60601-1 or relevant national standards.

- Follow the instructions for use of individual accessories.
- Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

6.2 Protection against inadvertent activation

The rotational speed control button can be locked to prevent inadvertent activation of the drill and reamer during a tool/attachment change.

Locking the knob for speed control 2:

- Twist the safety catch 3 to position OFF.
The speed control knob 2 is blocked and the product 1 cannot be operated.

Unlocking the knob for speed control 2:

- Twist the safety catch 3 to position ON.
The speed control knob 2 is unlocked and the product 1 can be operated.

Note

For further information on the drill and reamer, see TA014538 or TA014539 (flyer).

6.3 Coupling and uncoupling of attachment to drill and reamer

WARNING

Risk of injury when attaching/removing attachments/tools in the ON position through inadvertent activation of the product!

- Only attach/remove attachments/tools in the OFF position.

- Secure drill and reamer 1 against inadvertent activation with the safety catch 3, see Protection against inadvertent activation.

Coupling

- Push the attachment onto drill and reamer 1 until it engages with a click.
- Pull at attachment to check that it is securely coupled.

Uncoupling

- Turn the sleeve 4 in the direction of the arrow 5 and simultaneously remove the attachment from the drill and reamer 1.

6.4 Attaching and detaching the tool from the attachment

WARNING

Risk of injury when attaching/removing attachments/tools in the ON position through inadvertent activation of the product!

- Only attach/remove attachments/tools in the OFF position.

- Ensure that the tool connector and attachment type are compatible.

6.4.1 Attachments with quick-action chuck GB623R/GB657R/GB663R/GB664R/GB665R/GB668R/GB669R/GB670R

Coupling

- Pull back on the unlocking sleeve 7.
- Insert the tool shaft in the correct position up to the stop in the tool holder 8 of the attachment 6/9.
- Release the unlocking sleeve 7.
Tool is attached.
- Pull at tool to check that it is securely coupled.

Uncoupling

- Pull back the unlocking sleeve 7.
- Remove tool.

6.4.2 Drill attachment keyless Jacobs chuck GB620R

Note

The drill attachment keyless Jacobs chuck GB620R may open during left-hand rotation of tools (e.g. when pulling out tools).

To operate tools in left-hand rotation, use the Jacobs chuck with key GB621R, GB639R or GB667R.

Note

If necessary, the key GA031R can be used for coupling and decoupling.

Coupling

- ▶ Rotate the locking sleeve **11** clockwise and open Jacobs chuck.
- ▶ Insert the tool shaft in the correct position up to the stop in the tool holder **8** of the attachment **10**.
- ▶ Rotate locking sleeve **11** counterclockwise and close Jacobs chuck.
- ▶ Pull at tool to check that it is securely coupled.

Uncoupling

- ▶ Rotate the locking sleeve **11** clockwise and open Jacobs chuck.
- ▶ Remove tool.
- ▶ Rotate the locking sleeve **11** counterclockwise and close Jacobs chuck.

6.4.3 Attachments with three-jaw chuck GB621R/GB639R/GB667R

Coupling

- ▶ Open the Jacobs chuck with the key **13**.
- ▶ Insert the tool shaft in the correct position up to the stop in the tool holder **8** of the attachment **12**.
- ▶ Close the Jacobs chuck with the key **13** and pull tight.
- ▶ Pull at tool to check that it is securely coupled.

Uncoupling

- ▶ Open the Jacobs chuck with the key **13**.
- ▶ Remove tool.

6.4.4 Zimmer chuck with Hudson/Zimmer shaft GB176R

Coupling

- ▶ Push unlocking sleeve forward **15** and hold.
- ▶ Insert the tool shaft in the correct position up to the stop in the tool holder **8** of the attachment **14**.
- ▶ Release the unlocking sleeve **15**.
- ▶ Pull at tool to check that it is securely coupled.

Uncoupling

- ▶ Push unlocking sleeve forward **15** and hold.
- ▶ Remove tool.

6.4.5 Harris chuck with AO large shaft GB184R and DIN chuck with AO large shaft XF457R

Coupling

- ▶ Push unlocking sleeve back **15** and hold.
- ▶ Insert the tool shaft in the correct position up to the stop in the tool holder **8** of the attachment **14**.
- ▶ Release the unlocking sleeve **15**.
- ▶ Pull at tool to check that it is securely coupled.

Uncoupling

- ▶ Push unlocking sleeve back **15** and hold.
- ▶ Remove tool.

6.4.6 Adapter AO large to DHS/DCS step drill GB628R

Coupling

- ▶ Push unlocking sleeve back **15** and hold.
- ▶ Insert the tool shaft in the correct position up to the stop in the tool holder **8** of the attachment **16**.
- ▶ Release the unlocking sleeve **15**.
- ▶ Pull at tool to check that it is securely coupled.

Uncoupling

- ▶ Push unlocking sleeve back **15** and hold.
- ▶ Remove tool.

6.4.7 Sagittal saw attachment GB660R

Coupling

- ▶ Insert saw blade **25** into the slot of the tool holder **18**, see Fig. A. Ensure that the side stops of the saw blade are aligned with the tool holder.
The tool engages.
- ▶ Pull on the attachment **22** to ensure a secure fit.

Uncoupling

- ▶ Press knob for tool release **19** and pull saw blade **25** out of the tool holder **18**.

Aesculap® Acculan 4

Attachments

6.4.8 Drill attachment for Synthes radiolucent angle transmission GB645R

The drill attachment allows connecting a radiolucent angle transmission for applying drill holes in the bone under radiographic control.

The drill attachment is only suitable for use of the radiolucent angle transmission 511.300 from the Synthes company.

- Follow the instructions for the radiolucent angle transmission 511.300 in the operating instructions of the Synthes company.

Coupling the angle transmission

- Push the radiolucent angle transmission up to the stop on holder **8** of the drill attachment **20**.
- If necessary, slightly move the radiolucent unit.

Decouple angled handpiece

- Forcefully pull out the radiolucent angle transmission from the drill attachment **20**.

Couple/decouple the tool

- Follow the instructions for the radiolucent angle transmission 511.300 in the operating instructions of the Synthes company.

6.4.9 Insert and clamp Kirschner wire attachment GB641R

Risk of injury when using long guide wires!

- Use Kirschner wire protection sleeve when using long guide wires.

Note

The special Kirschner wire chuck is recommended for the placement of guide wires. This quick-action chuck adapter allows quick and easy tensioning of drilling wires.

For further information see TA014535 or TA014539 (flyer).

The following diameters can be set on the Kirschner wire attachment:

- 0.6 mm to 1.8 mm
- 1.8 mm to 3.0 mm
- 3.0 mm to 4.0 mm

Insert Kirschner wire

- Ensure that the clamping lever **22** is in the starting position (deactivated status).
- Set the adjusting sleeve **23** of the Kirschner wire attachment to the desired diameter range:
 - Push back the adjusting sleeve **23** and turn until the desired diameter range has been reached.
 - Release adjustment sleeve **23**, making certain that adjustment sleeve **23** clicks into position.
- Insert the Kirschner wire into the Kirschner wire chuck **24** until the desired exposed length is achieved.

The Kirschner wire is held in the intended position in the attachment by an automatic clamping mechanism.

Clamping the Kirschner wire

- Pull on the clamping lever **22** and hold in the desired position. The further the clamping lever is pulled back, the higher the tension of the Kirschner wire.

Note

The Kirschner wire will only remain tensioned with the clamping lever pulled back. If the clamping lever is released, it will move back into the starting position and the Kirschner wire is freely mobile.

6.5 Function checks

The function test must be performed before each use and after each intra-operative change of attachment and tool.

- Check that the attachment is coupled properly. To test this, pull on the attachment.
- Check that the tool is coupled properly. To do this, pull on the tool.
- Ensure that the tool blades are not mechanically damaged.
- Release drill and reamer for use (ON position).
- Briefly run the drill and reamer at maximum speed in clockwise and counterclockwise mode.
- Make certain that the rotational direction is correct in each case.
- Check the product for damage, abnormal operation noises, heavy vibration and excessive heat.
- Do not use the product if it is damaged or defective.
- Set aside the product if it is damaged.

6.6 Safe operation

WARNING

Coagulation of patient tissue or risk of burns for patients and user through hot product!

- ▶ Do not use drill attachments to operate reamer tools such as the acetabulum reamer and the medullary reamer.
- ▶ Operate any reamer tool only with a reamer attachment.
- ▶ Cool the tool during operation.
- ▶ Keep product/tool out of the reach of patients.
- ▶ Let the product/tool cool down.
- ▶ Use a cloth to protect against burns when changing the tool.

WARNING

Risk of infection or injury through aerosol formation or particles ejected from the tool!

- ▶ Use suitable protection (such as waterproof protective clothing, face mask, safety glasses, suction).

WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!

- ▶ Apply the product only under visual control.

WARNING

Risk of injury and damage to the tool/system!

The rotating tool may get caught in drapes (such as textiles).

- ▶ Do not let the tool come into contact with drapes (such as textiles) during operation.

Operating the product

Note

For further information on the drill and reamer, see TA014538 or TA014539 (flyer).

- ▶ Start the drill and reamer at moderate speed.
- ▶ Apply moderate pressure to prevent slipping.
- ▶ Do not bend the tool, otherwise there is a danger of breakage.
- ▶ When screwing bone screws and bone pins in and out: ensure correct torque setting.

Use the following attachments for drilling:

- GB620R
- GB621R
- GB623R
- GB639R
- GB663R
- GB664R
- GB665R

Use the following attachments for milling (medullary reaming):

- GB667R
- GB668R
- GB669R
- GB670R

Use the following attachment for medullary reaming:

- GB657R

Use the following attachment for the placement of Kirschner wires:

- GB641R

Use only the following attachment for using the radiolucent angle transmission 511.300 from the Synthes company:

- GB645R

Use the following attachment for drilling:

- GB660R

Aesculap® Acculan 4

Attachments

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this. The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

*For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.*

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

On stainless steel, residues containing chlorine or chloride (such as surgical residues, drugs, saline solutions in water for cleaning, disinfection and sterilization) may lead to corrosion (pitting corrosion, tensile corrosion) and thus to the destruction of the product. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Visual material changes (such as fading or color changes in titanium or aluminum). For aluminum, the application/process solution only needs to be pH >8 to cause visible surface changes.
- Material damage (such as corrosion, cracks, breaks, premature aging or swelling).
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

7.3 Preparations at the place of use

- ▶ Remove all attached components from the product (tool and accessories).
- ▶ Remove the attachment from drill and reamer.
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.4 Preparation before cleaning

- ▶ Before the first automatic cleaning/disinfection: attach ECCOS fixations in a suitable sieve basket.
- ▶ Insert the products in the correct position into the ECCOS fixations, see Fig. B.

Kirschner wire attachment GB641R

- ▶ Turn adjustment sleeve **20** to setting for largest Kirschner wire diameter.

7.5 Cleaning/disinfection

7.5.1 Product-specific safety instructions for the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for plastic material and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C during chemical cleaning and/or disinfection.
- ▶ Do not exceed the maximum temperature of 96 °C during thermal disinfection with FD water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

Note

The indicated drying time is a guide time only. It must be checked taking into account the specific conditions (e.g. load) and if applicable adjusted.

Aesculap® Acculan 4

Attachments

7.6 Manual cleaning with wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	≥2	-	D-W	Until visually clean
II	Cleaning with enzyme solution	RT (cold)	≥2	0.8	D-W	pH-neutral*
III	Intermediate rinse	RT	≥5	-	D-W	-
IV	Drying	RT	-	-	-	-
V	Wipe disinfection	-	>1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol
VI	Final rinse	RT (cold)	0.5	-	FD-W	-
VII	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

* Validated with enzyme cleaner "Cidezyme Johnson & Johnson"

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

Phase I

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Clean the product under running water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- ▶ Brush cannulation with cleaning brush TA011944 and difficult to access surfaces with a suitable cleaning brush made of plastic for at least 1 min.

Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available on Aesculap Extranet at <https://extranet.bbraun.com>).

Phase II

- ▶ Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- ▶ Spray products with a pH neutral enzyme solution, let soak in for at least 2 minutes and then wipe off.

Phase III

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Rinse product under running tap water for at least 5 minutes.
- ▶ Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- ▶ Contamination should be removed with a lint-free cloth or soft brush moistened with enzyme cleaner.
- ▶ Flexible components (such as sleeves) and cannulations should be rinsed for 20 seconds with the water pistol (cold water, at least 2.5 bar).
- ▶ After manual cleaning, check visible surfaces and areas of flexible components for residues.
- ▶ If necessary, repeat the cleaning process (phase 1 to 3).

Phase IV

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

Phase V

- ▶ Wipe all surfaces of the product with a single-use disinfecting wipe.

Phase VI

- ▶ Rinse disinfected surfaces after the prescribed reaction time for at least 1 minute under running demineralized water.
- ▶ Drain any remaining water fully.

Stage VII

- ▶ Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

7.7 Automatic cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must fundamentally have a tested efficacy (such as FDA approval or CE label pursuant to DIN EN ISO 15883).

Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

7.7.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry/Note
I	Rinsing	RT (cold)	-	-	D-W	Until visually clean
II	Brushes	RT (cold)	-	-	D-W	Until visually clean

D-W: Drinking water

RT: Room temperature

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.
- ▶ For Kirschner wire attachment GB641R: Set the adjustment sleeve **20** to the greatest Kirschner wire diameter.

Phase I

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Thoroughly clean the product under running water.

Phase II

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Brush cannulation with cleaning brush TA011944 and difficult to access surfaces with a suitable cleaning brush made of plastic for at least 1 min.
- ▶ After manual preparation, check visible surfaces for residue and repeat the pre-cleaning process as needed.

Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available on Aesculap Extranet at <https://extranet.bbraun.com>).

Aesculap® Acculan 4

Attachments

7.7.2 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemistry/Note
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH ~ 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	mind. 10 min at max. 120 °C

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ For Kirschner wire attachment GB641R: Set the adjustment sleeve 20 to the greatest Kirschner wire diameter.
- ▶ Insert the product in its proper position in the ECCOS fixation.
- ▶ Then connect the interior rinse device to the ECCOS fixation and connect to the rinse connector of the cleaning/disinfection device/rinse cart.
- ▶ After automatic cleaning/disinfection, check visible surfaces for residues and repeat the cleaning/disinfection process as needed.

7.8 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Spray product after every cleaning and disinfection with Acculan 4 oil spray adapter 23 GB600480 (black) for approx. 2 sec. with Aesculap STERILIT Power Systems oil spray GB600, see Fig. C.

Note

Aesculap recommends occasional spraying of moving parts (such as knobs, couplings, cover caps) with the Aesculap STERILIT-Power-Systems oil spray.

- ▶ Check the product for any damage, abnormal running noise, overheating or excessive vibration.
- ▶ Set aside the product if it is damaged.

7.9 Packaging

- ▶ Always follow the instructions for use of the respective packaging and storage devices (e.g. instructions for use TA009721 for Aesculap-ECCOS holder system).
- ▶ Insert the products in the correct position into the ECCOS fixations, see Fig. B.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

7.10 Steam sterilization

Note

Remove all attached components from the product (tools, accessories) before sterilization.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Use a validated sterilization method:
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer DIN EN 285 and validated pursuant to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min

When sterilizing multiple products in one steam sterilizer:

- ▶ Ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

7.11 Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

7.12 Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

8. Maintenance

To ensure reliable operation, the product must be maintained at least once a year.

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

Aesculap® Acculan 4

Attachments

9. Troubleshooting list

► Have defective products repaired by Aesculap Technical Service, see Technical Service.

Issue	Cause	Detection	Remedy
Product becomes too hot	Excessive use	Attachment very hot	Follow operating instructions (operating mode).
	Defective attachment drive/ball bearings	Attachment very hot	Follow operating instructions (preparation, care). Preventative: lubricate attachment before every sterilization. Have the manufacturer repair the attachment.
	Fall damage, product defective	Attachment very hot	Have the attachment repaired by the manufacturer
	Blunt tool	Heating of tool and product	Change tool.
Insufficient power	Product is operating in left rotation mode	Serrated tool run in counterclockwise mode	Use serrated tools only in clockwise mode
	Attachment defective	Attachment very hot	Follow operating instructions (preparation, care). Preventative: lubricate attachment before every sterilization. Follow operating instructions (operating mode). Have the attachment repaired by the manufacturer
	Blunt tool	Tool blades worn out	Change tool.
Loud running noise	Defective attachment drive/ball bearings	Loud, perceptible noise during operation	Follow operating instructions (preparation, care). Preventative: lubricate attachment before every sterilization. Have the manufacturer repair the attachment.
Attachment does not attach or detach.	Attachment not compatible	Attachment doesn't engage	Use the compatible attachment for drill.
	Connection on attachment deformed/defective	Attachment is difficult to couple or cannot be coupled	Have the manufacturer repair the attachment.
	Coupling on drill deformed/defective	Attachment is difficult to couple or cannot be coupled	Have the manufacturer repair the drill.
	Sleeve on drill hard to move	Attachment doesn't engage	Turn and hold the drill sleeve, then attach the attachment. Turn the sleeve back if necessary. Preventative: lubricate the sleeve on the drill before each sterilization.

Issue	Cause	Detection	Remedy
Tool does not attach or detach.	Incompatible tool	Tool does not click into place	Use the tool compatible with the attachment.
	Tool connector deformed/defective	Tool is difficult to couple or cannot be coupled	Use a new tool.
	Coupling on attachment deformed/defective	Tool is difficult to couple or cannot be coupled	Have the manufacturer repair the attachment.
	Tool coupling dirty	Tool is difficult to couple or cannot be coupled	Clean tool or use a new tool. Clean the attachment.
For GB641R: Kirschner wire does not insert into the Kirschner wire chuck	Kirschner wire chuck not set correctly	Diameter of the Kirschner wire is not consistent with the scale setting	Set the adjustment sleeve to the right Kirschner wire diameter.
For GB641R: Kirschner wire not turning	Kirschner wire chuck not set correctly	Diameter of the Kirschner wire is not consistent with the scale setting	Set the adjustment sleeve to the right Kirschner wire diameter.
	Open the clamping lever	Open the clamping lever	Close the clamping lever.
	Clamping lever not pulled back with sufficient force	Kirschner wire not turning	Pull back the clamping lever forcefully and hold it.
Tool not moving	Attachment not fully attached to the drill	Attachment can be pulled out of the drill	Correctly attach the attachment and perform a function test.
	Tool not fully attached to the attachment	Tool can be pulled out of the tool coupling	Correctly attach the tool and perform a function test.
	Attachment defective	Drill turning, but attachment is not turning.	Have the manufacturer repair the attachment.
	Drill defective	Drill not turning	Have the manufacturer repair the drill.
	The safety catch of the drill is in the OFF position	Safety catch is in OFF position	Turn the safety catch into the ON position.
Oil spray adapter not attachable	Oil spray adapter not compatible	Oil spray adapter not attachable	Use the compatible oil spray adapter for GA344.

Aesculap® Acculan 4

Attachments

10. Technical Service

WARNING

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- ▶ **Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.**
- ▶ **Do not modify the product.**

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95 -1601
Fax: +49 (7461) 14 -939
E-Mail: ats@aesculap.de
Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042 USA
Aesculap Repair Hotline
Phone: +1 (800) 214 -3392
Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

11. Accessories/Spare parts

Art. no.	Designation
GB243R	ECCOS set for 2 large machines
GB600	STERILIT Power Systems oil spray
GB600840	Oil spray adapter
TA011944	Cleaning brush
GB496R	ECCOS fixation for 3 attachments
GB497R	ECCOS fixation for one attachment
GA031R	Key for large Jacobs chuck
GA062R	Key for small Jacobs chuck
TA014541	Instructions for use for Acculan 4 attachments (flyer)

12. Technical data

12.1 Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GB176R	Zimmer chuck with Hudson/Zimmer shaft	IIa
GB184R	Harris chuck with AO large shaft	
GB620R	Drill attachment keyless Jacobs chuck	
GB621R	Drill attachment large Jacobs chuck	
GB623R	Drill attachment AO small	
GB628R	Adapter AO large to DHS/DCS step drill	
GB639R	Drill attachment small Jacobs chuck	
GB641R	Kirschner wire attachment	
GB645R	Drill attachment for Synthes radiolucent angle transmission	
GB657R	Medullary reamer attachment AO large	
GB660R	Sagittal saw attachment	
GB663R	Drill attachment Trinkle	
GB664R	Drill attachment Aesculap hexagon	
GB665R	Drill attachment Hudson/Zimmer	
GB667R	Reamer attachment large Jacobs chuck	
GB668R	Reamer attachment AO large	
GB669R	Reamer attachment Hudson/Zimmer	
GB670R	Reamer attachment Harris	
XF457R	DIN chuck with AO large shaft	

12.2 Performance data, information about standards

Rotational direction	Right and left rotation, oscillation
Conforming to standard	IEC/DIN EN 60601-1

The product was subject to 500 preparation cycles as a test by the manufacturer and has passed.

The performance data for the individual attachments can be found in the table below.

Drill attachments

Attachment	Connector	Dimensions Length x Ø [mm] ± 5 %	Weight [g] ± 10 %	Max. speed [min ⁻¹]	Max. torque [Nm]	Cannulation [mm]
GB620R	Jacobs chuck keyless Ø 0.3 mm to 7.4 mm	125 x 35	397	1 000	5	4
GB621R	Jacobs chuck with key Ø 0.5 mm to 7.4 mm	110 x 35	300	1 000	5	4
GB623R	AO small	90 x 35	210	1 250	4	4
GB639R	Jacobs chuck with key Ø 0.5 mm to 4.0 mm	105 x 35	243	1 250	4	4
GB645R	Synthes radiolucent angle transmission	110 x 35	225	900	5	4
GB663R	Trinkle	93 x 35	217	1 000	5	4
GB664R	Aesculap hex	86 x 35	225	1 000	5	4
GB665R	Hudson/Zimmer	100 x 35	230	1 000	5	4

Milling attachments

Attachment	Connector	Dimensions Length x Ø [mm] ± 5 %	Weight [g] ± 10 %	Max. speed [min ⁻¹]	Max. torque [Nm]	Cannulation [mm]
GB657R	AO large	108 x 35	300	370	12	4
GB667R	Jacobs chuck with key Ø 0.5 mm to 7.4 mm	133 x 35	280	250	19	4
GB668R	AO large	108 x 35	300	250	19	4
GB669R	Hudson/Zimmer	118 x 35	300	250	19	4
GB670R	Harris	106 x 35	290	250	19	4

Aesculap® Acculan 4

Attachments

Saw Attachments

Attachment	Connector	Dimensions L x W x H [mm] ± 5 %	Weight [g] ± 10 %	Max. oscillation frequency [min ⁻¹]	Cannulation [mm]
GB660R	Rapid Action	89 x 36 x 44	270	17 000	–

Kirschner wire attachment

Attachment	Connector	Dimensions L x W x H [mm] ± 5 %	Weight [g] ± 10 %	Max. speed [min ⁻¹]	Cannulation [mm]
GB641R	Kirschner wire Ø 0.6 mm to 4 mm	93 x 36 x 137	285	1 250	4

Adapter

Attachment	Connector	Dimensions Length x Ø [mm] ± 5 %	Weight [g] ± 10 %	Max. speed [min ⁻¹]	Max. torque [Nm]	Cannulation [mm]
GB184R	Harris	88 x 22	85	With GB668R: 250 With GB657R: 370	With GB668R: 19 With GB657R: 12	–
GB176R	Zimmer	50 x 22	55	With GB669R: 250 With GB665R: 1 000	With GB669R: 19 With GB665R: 5	3
GB628R	DHS/DCS	61 x 22	70	With GB668R: 250 With GB657R: 370	With GB668R: 19 With GB657R: 12	4
XF457R	DIN	75 x 25	100	With GB668R: 250 With GB657R: 370	With GB668R: 19 With GB657R: 12	–

12.3 Operating mode

Operating mode	Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1) Right/left rotation: ■ 60 second application, 60 second pause ■ 20 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C Milling (clockwise/counterclockwise): ■ 30 second application, 30 second pause ■ 8 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C Drilling (oscillation): ■ 60 second application, 60 second pause ■ 4 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C Saw mode with GB660R: ■ 30 second application, 60 second pause ■ 3 repetitions ■ 30 min cooling time ■ Max. Temperature 48 °C
----------------	--

12.4 Environmental conditions

	Operation	Storage and transport
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

14. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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