
Aesculap Spine

-
- | | |
|------------|----------------------------------|
| GB | Instructions for use |
| USA | activ C insertion instrument |
| D | Gebrauchsanweisung |
| | activ C-Einsetzinstrument |
| F | Mode d'emploi |
| | Instrument d'insertion activ C |
| E | Instrucciones de manejo |
| | Instrumento de inserción activ C |
| I | Istruzioni per l'uso |
| | Strumento inseritore activ C |
| P | Instruções de utilização |
| | Instrumento introdutor activ C |
| NL | Gebruiksaanwijzing |
| | activ C-inzetinstrument |
| CZ | Návod k použití |
| | activ C-nasazovací instrument |

activ C insertion instrument

Legend

- 1 Insertion instrument
- 2 Spacer
- 3 Key (for insertion instrument)
- 4 Repositioning instrument
- 5 Irrigation connector

Symbols on product

	Attention, see instructions for use
---	-------------------------------------

Intended use

The insertion instrument is used for picking up the activ C intervertebral disk prosthesis from its sterile packaging, inserting it in the prepared intervertebral space, and carrying out position corrections of the prosthesis.

Available sizes

Art. no.	Designation
FW857R	Insertion instrument with irrigation connector H6/H7
FW863R	Spacer for implant height 5 mm
FW864R	Spacer for implant height 6 mm
FW865R	Spacer for implant height 7 mm
FW866R	Insertion instrument with irrigation connector H5
FW867R	Repositioning instrument
FW945R	Key for insertion instrument

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working tip: Introduce the product carefully through the working channel (e.g. trocar).

Safe operation

Note

Follow OP-Manual no. 031302 and instructions for use "activ C intervertebral disk prosthesis" TA011995.

Assembling the instrument

Note

Mind the implant height when selecting the insertion instrument!

- Remove irrigation connector 5 from the handle end of insertion instrument 1.
- Verify that clamping sleeve 7 on insertion instrument 1 has been loosened.
- Turn clamping sleeve 7 counterclockwise until it is not in connection anymore with connecting piece 6 of insertion instrument 1.
- Connect spacer 2 for the intraoperatively selected implant height with insertion instrument 1.
- Insert spacer 2 in the hole of the connecting piece, observing the "CRANIAL" and "CAUDAL" markings on spacer 2 and insertion instrument 1. At the same time keep button 9, "push to clean", at the handle of insertion instrument 1 pressed down.
- Insert spacer 2 down to the stop in insertion instrument 1.

- Release "push to clean" button 9. By axially pulling at spacer 2, check that the spacer is properly locked in insertion instrument 1.
- Turn adjusting wheel 8 counterclockwise to set the depth stop to minimum insertion depth.

Coupling the implant

WARNING

Damage to the implant due to improper handling!

- Observe the instructions for use of activ C intervertebral disk prosthesis (TA011995).
- Mount the implant on the insertion instrument immediately after removing it from the implant packaging.
- Handle with caution the friction surfaces of the prosthesis plates and of the polyethylene inlay.

- Open the implant packaging.

Note

Do not remove the implant from the packaging!

- Mount insertion instrument 1 on the disk prosthesis with the latter still in its packaging. Observe the "CRANIAL" and "CAUDAL" markings. The superior prosthesis plate (the one with the spikes) must abut the part of insertion instrument 1 marked "CRANIAL", while the inferior prosthesis plate (with the fin) has to abut the part of insertion instrument 1 marked "CAUDAL".
- Turn clamping sleeve 7 clockwise to fixate the intervertebral disk prosthesis on insertion instrument 1.
- Remove disk prosthesis 1 from its packaging, holding it with the insertion instrument.


Fig. 1

activ C insertion instrument

Inserting the implant

WARNING

Damage to vertebral body end plates when the intervertebral disk prosthesis is tapped in!

- **Apply proper care when tapping in the intervertebral disk prosthesis.**
- Carefully, under X-ray control, insert the disk prosthesis into the intervertebral disk space. The superior prosthesis plate (plate with the spikes) must be oriented towards cranial, the inferior prosthesis plate (the plate with the fin) towards caudal.
- If necessary impel the implant with mild hammer taps on the proximal end of insertion instrument 1.

WARNING

Compression of the spinal canal and of other posterior elements caused by the intervertebral disk prosthesis inserted too deeply!

- **Prior to inserting the intervertebral disk prosthesis, adjust the depth guide to minimum insertion depth.**
- **Introduce the disk prosthesis into the intervertebral space under X-ray control.**
- If necessary, correct the depth stop position with adjusting wheel 8 until the intended position is reached.
- Should repositioning become necessary, screw repositioning instrument 4 onto the proximal end of insertion instrument 1 and carefully, with a slotted hammer, move the implant to its intended position.

Decoupling the implant

Note

To avoid damage to the implant, do not couple again an insertion instrument that already has been decoupled once.

- Loosen clamping sleeve 7 on insertion instrument 1 by turning it counterclockwise, and remove insertion instrument 1.
- If clamping sleeve 7 is difficult to loosen, insert key 3 in one of the lateral drill holes of clamping sleeve 7 and loosen clamping sleeve 7 in this way.
- On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the intervertebral disk implant.

Disassembling

- Uncouple spacer 2 from insertion instrument 1 by pressing "push to clean" button 9 at the handle of insertion instrument 1 and simultaneously pulling it out, axially, towards the working end.
- If applicable, unscrew repositioning instrument 4 from insertion instrument 1 by turning it counterclockwise.
- Prior to processing make sure that clamping sleeve 7 sits loosely on the instrument shaft.

Validated processing procedure

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, manufacturer's specifications can only serve as guide values for the processing procedures applied by the individual user.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of "stainless" steels. Consequently the time interval between application and processing should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures > 45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) should be used.

Excessive doses of neutralizers or basic detergents can cause chemical deterioration and/or fading of the laser engraving on stainless steel.

Residues containing chlorine or chlorides, as contained e.g. in surgical residues, tinctures, medicines, saline solutions, service water used for cleaning, and in cleaning/disinfecting agents, cause corrosion damage (pitting, stress corrosion) and thereby destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals recommended by the chemicals' manufacturers as effective cleaning/disinfecting agents compatible with the materials used in the respective product. Strictly observe all application instructions, e.g. regarding temperatures, concentrations, exposure times, etc. Failure to do so can result in the following unwelcome effects:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes or
- material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use process chemicals that cause stress cracking of plastics such as PPSU, or attack e.g. silicone softeners and cause brittleness.

activ C insertion instrument

- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications, Red Brochure – Proper maintenance of instruments.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove visible residues as completely as possible, using a lint-free, single-use cleaning tissue.

Preparations prior to cleaning

- Have the product dry in a disposal container and ready for immediate cleaning and disinfecting within 30 min after use.

Note

Clamping sleeve 7 must sit loosely on the instrument shaft during cleaning.

Note

To ensure proper cleaning of spacers 2, the PEEK forks must sit loosely (so that they can be turned in both directions) on the thread.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use a cleaning and disinfecting agent according to manufacturer's instructions that
 - must be approved for the product materials (e.g. aluminum, plastics, high-grade steel)
 - must not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 60 °C.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Manual cleaning/disinfecting

FW857R, FW863R–FW867R, FW945R

- Clean hinged or jointed products in open and closed positions.
- Check visible surfaces for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfecting and cleaning with brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Final rinse	RT (cold)	0.5	-	FD-W	-
IV	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives, which would damage the surfaces, for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage IV

- Dry the product completely with compressed air.

activ C insertion instrument

Mechanical cleaning/disinfecting with manual pre-cleaning

FW857R, FW863R–FW867R, FW945R

Note

For mechanical cleaning: Install the irrigation connector on the handle end of the insertion instrument.

Manual pre-cleaning with brush

FW857R/FW866R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	–	D–W	–

D–W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives that would damage the surfaces for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.
- Store any products with hinges or joints on the tray in such a way that the joints are kept open.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	- concentrate, alkaline: pH = 10.9 <5 % anionic tensides - 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD-W	- concentrate, alkaline: pH = 2.6 Basis: Citric acid - 0.2 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD-W	-
V	Thermal disinfecting	94/201	10	FD-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

activ C insertion instrument

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts (e.g. joints and links) with a maintenance oil suitable for the respective sterilization process (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Reassemble separable products.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside the product if it is damaged.

Packaging

- Appropriately protect products with fine working tips.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that any cutting edges are protected.
- Package baskets appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

Note

The product may only be sterilized in disassembled condition.

Note

The irrigation connector must be removed from the handle end of the insertion instrument prior to sterilization.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilization acc. to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST 46-1993, validated acc. to EN ISO 17665 or EN 554/ISO 13683
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270— 275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store processed products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

Technical Service

For service, maintenance or repairs, contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG