

Aesculap Spine

-
- | | |
|------------|--|
| GB | Instructions for use/Technical description |
| USA | activ C handle, reamer guide, chisel guide, trial implants |
| D | Gebrauchsanweisung/Technische Beschreibung |
| | activ C Handgriff, Fräseführung, Meißelführung, Probeimplantate |
| F | Mode d'emploi/Description technique |
| | Poignée activ C, guide-fraise, guidage de burin, implants d'essai |
| E | Instrucciones de manejo/Descripción técnica |
| | Mango activ C, guía de fresado, guía de cincel, implantes de prueba |
| I | Istruzioni per l'uso/Descrizione tecnica |
| | Impugnatura, guida per la fresa, guida per lo scalpello, impianti di prova activ C |
| P | Instruções de utilização/Descrição técnica |
| | Cabo activ C, guia de fresa, guia de cincel, implantes de ensaio |
| NL | Gebruiksaanwijzing/Technische beschrijving |
| | activ C-handgreep, freesgeleider, beitelgeleider, proefimplantaten |
| CZ | Návod k použití/Technický popis |
| | Rukojeť, vedení frézy, vedení dláta, zkušební implantáty activ C |

FW647R

FW645R

activ C handle, reamer guide, chisel guide, trial implants

Legend

- 1 Trial implant
- 2 Handle
- 3 Reamer guide
- 4 Mounting block (for reamer guide)
- 5 Reamer
- 6 Reamer handpiece
- 7 HiLAN pneumatic motor
- 8 active C chisel insert
- 9 Chisel holder
- 10 Chisel guide

Symbols on product and packages

Symbol	Explanation
	Adhere to instructions for use
	Not for reuse in intended applications as defined by the manufacturer
	Position "closed"
	Position "open"
	Ready

Intended use

The trial implants are used in conjunction with the handpiece to verify the size, height and shape of the activC intervertebral disk prosthesis to be implanted, and in conjunction with the reamer guide or the chisel guide to prepare the fin bed for the inferior implant plate of the activC intervertebral disk prosthesis.

Available sizes

Art. no.	Designation
FW789R ¹	Trial implant flat size XS, H5
FW790R ¹	Trial implant flat size XS, H6
FW791R ¹	Trial implant flat size S, H5
FW792R ¹	Trial implant flat size S, H6
FW793R ¹	Trial implant flat size M, H5
FW794R ¹	Trial implant flat size M, H6
FW790R ¹	Trial implant flat size L, H5
FW790R ¹	Trial implant flat size L, H6
FW797R ¹	Trial implant flat size XL, H5
FW798R ¹	Trial implant flat size XL, H6
FW799R ¹	Trial implant flat size XXL, H5
FW800R ¹	Trial implant flat size XXL, H6
FW8701 ¹	Handle
FW871R ²	Reamer guide (with mounting block)
GB771R ²	Hi-Line XS handpiece angled XL-I
GE700SU ⁴	Hi-Line XS activ C drill XL-I D1, 5
FW784R ³	activ C chisel guide
FW786R ³	active C chisel holder
FW787R ⁴	active C chisel

¹ Trays FW645R and FW647R

² only tray FW647R

³ only tray FW645R

⁴ single-use products

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

WARNING

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information regarding such training, please contact your national B. Braun/Aesculap agency.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only for its intended purpose, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Exercise caution when passing products through the working channel (e.g. trocar).

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

Note

Follow OP-Manual no. O31302 and instructions for use "activ C intervertebral disk prosthesis", see TA011995.

Mounting the trial implant of the selected size on the handle


Fig. 1 Handle for trial implants in open (left) and closed (right) position

- Turn locking pipes 11 of handle 2 to position "open"

activ C handle, reamer guide, chisel guide, trial implants

- Mount the trial implant with the selected size, height and shape onto the handle **2**. Note the "CRANIAL" and "CAUDAL" markings on the handle and on the trial implant.
- Turn locking pipes **11** of handle **2** to position "closed" .
- Check the trial implant for firm seating.

WARNING

Compression of the spinal canal and other posterior elements caused by the trial implant inserted too deeply!

- Before inserting the trial implant, set the depth stop to the minimum insertion depth.
- Insert the trial implant into the intervertebral space under X-ray control.

- Use adjusting wheel **12** to set the depth stop of the trial implant to the minimum insertion depth.
- Insert the trial implant in the intervertebral space under X-ray control. Ensure correct cranial and caudal alignment of the trial implant.
- Align the markings of the trial implant center to the midline marking.
- If necessary, correct the position of the depth stop with adjusting wheel **12** until the intended position is reached. Ensure that depth stop abuts to the anterior edge of the vertebral body.

WARNING

Increased risk of migration if the selected implant size is too small!

- Select the trial implant or implant size, respectively, to achieve maximum coverage of the vertebral body endplates.
- On the X-ray, in AP and lateral view, check to ensure that the size, position and alignment of the trial implant are correct.

WARNING

Irritation of the facet joints and strain on the spinal cord and nerve roots due to excessive or insufficient distraction!

- Select the correct trial implant/ implant height.
- Check for correct height of the trial implant in a lateral X-ray.

- On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the trial implant.
- With the trial implant properly positioned, turn locking pipes **11** to position "open"  and pull off handle **2** carefully in an axial direction from the trial implant.
The trial implant remains in the intervertebral space.

Note

The fin bed can either be prepared by reaming or chiseling. This can be carried out either with the reamer guide with handpiece or the chisel guide with chisel holder, depending on which set is used.

Reaming the fin bed/ Mounting the reamer guide

Note

The handle and the reamer guide are fitted with identical locking mechanisms.

- Turn locking pipes **11** of reamer guide **3** to position "open" .
- Mount the reamer guide **3** onto the trial implant in the intervertebral space. Note the "CRANIAL" and "CAUDAL" markings on the reamer guide and on the trial implant.
- Turn locking pipes **11** of reamer guide **3** to position "closed" .
- Reamer guide **3** is fixed to the trial implant.
- Check again, by X-ray imaging, the position of the trial implant.

Reaming the fin bed/Mounting the reamer handpiece

WARNING

Damage to nerves, spinal cord and other tissue because of fin bed reamed too deeply due to using too long a reamer handpiece or reamer.

- Only use reamer handpiece Hi-Line XS XL-1 (GB771R).
- Only use activC reamer (GE700SU).

- Prior to mounting, verify that the safety stop of the pneumatic motor is in position "0".
- Assemble the compressed-air tube, the HiLAN pneumatic motor and reamer handpiece 6 as described and illustrated in the respective instructions for use supplied with the products.


Fig. 2 Assembling

- Keep pressed locking clip 13 of mounting block 4 and push the clip over reamer handpiece 6, see Fig. 3.

Fig. 3 Mounting block with pressed locking clip

- Slide on mounting block 4 as far as it will go and release locking clip 13, see Fig. 4.

Fig. 4 Mounting block with locking clip engaged

- Pull at the mounting block 4 to check whether it is engaged properly. The arrow marking on the locking clip 13 must point to the "locked"  pictogram.
- Mount reamer 5 as described in the instructions for use of reamer handpiece 6.

WARNING

Risk of injury to blood vessels, nerves, spinal cord or other tissue if wrong reamer is used!

- Only use activ C reamer GE700SU.
- Check for technical faults prior to use.
- Do not use bent reamers.
- Do not reuse reamers after first use.

activ C handle, reamer guide, chisel guide, trial implants

DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

Note

The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk. The fin bed must only be prepared for the inferior implant plate.

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- Ensure that the position of the trial implant does not change during reaming.

WARNING

Risk of injury from sharp reamer!

- To avoid injuries to critical blood vessels, nerves and other tissue, as well as injuries to the operating surgeon, always work with appropriate care when a reamer is clamped in the machine.
- Do not activate the motor before the handpiece with the coupled reamer is properly positioned in the guide hole and the guiding slot of the reamer guide.
- **A:** Position reamer handpiece **6** with the coupled reamer through the frontal guide hole of reamer guide **3**, see Fig. 5.
- **B:** Position the laterally protruding pins of mounting block **4** through the rear guiding slots of reamer guide **3**, see Fig. 5.


Fig. 5 Sliding the reamer handpiece into the reamer guide

- With the hand control inactivated, push the safety stop to . The system is ready for operation.

Note

The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk.

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- Ensure that the position of the trial implant does not change during reaming.

WARNING

Risk of injury to nerves and spinal cord caused by too deep insertion of the reamer!

- Only use activC instruments.
- Use the instruments according to their intended use.
- For reaming the holes, apply the reamer only in combination with a correctly selected and positioned trial implant and with the reamer guide coupled.
- Ensure that the depth stop of the trial implant abuts to the anterior edge of the vertebral body.
- Always maintain lateral X-ray control while reaming the holes.

- To ream the first hole under lateral X-ray control: With the motor running, push reamer handpiece **6** at an even speed and with slight forward pressure into the first guiding slot down to the stop.
- Pull back reamer handpiece **6** and, after moving the pins of mounting block **4** to the second guiding slot, repeat the reaming procedure under lateral X-ray control.

The reaming depth will be limited by the stop at reamer guide **3**.

- With the hand control inactivated, shift the safety stop on the pneumatic motor to position "0".
- Cautiously pull out reamer handpiece **6** first from the rear and then from the frontal reamer guide **3**.
- Remove the reamer guide **3** and trial implant from the intervertebral disk space. Clean the fin bed with nerve hooks from the discectomy set.

Chiseling the fin bed/Mounting the reamer guide

Note

The handle and the chisel guide are fitted with identical locking mechanisms.

- Turn the locking pipes **11** of the chisel guide **10** to the "open" position .
- Mount the chisel guide **10** onto the trial implant in the intervertebral space. Note the "CRANIAL" and "CAUDAL" markings on the handle and on the trial implant.
- Turn locking pipes **11** of chisel guide **10** to position "closed" .
- Chisel guide **10** is fixed to the trial implant.
- Check again, by X-ray imaging, the position of the trial implant.

WARNING

Damage to nerves, spinal cord and other tissue because of fin bed reamed too deeply due to using too long a chisel.

- Only use chisel holder FW786R.
- Only use chisel FW787R.

- Mount chisel onto the chisel holder as shown.
- To do this, attach the chisel on the side of the chisel holder **9** and push back until it engages.
- Pull to check if the chisel is engaged properly.

activ C handle, reamer guide, chisel guide, trial implants

Fig. 6 Chisel and chisel holder

WARNING

Risk of injury to blood vessels, nerves, spinal cord or other tissue if wrong chisel is used!

- Only use chisel FW787R.
- Do not use bent chisels.
- Do not reuse reamers after first use.

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- Make certain that the position of the trial implant does not change during reaming.

WARNING

Risk of injury caused by sharp chisel!

- To avoid injuries to critical blood vessels, nerves and other tissue, as well as injuries to the operating surgeon, always work with appropriate care when a chisel is clamped in the chisel holder.
- A: Carefully push the chisel holder 9 with the chisel through the rear chisel guide.
- B: Carefully push the chisel holder 9 with the chisel through the front chisel guide.

Note

The chiseling of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk. The fin bed must only be prepared for the inferior implant plate.


Fig. 7 Chisel holder and chisel guide

Note

The chiseling of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk.

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

- Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
- Make certain that the position of the trial implant does not change during chiseling.

WARNING

Risk of injury to nerves and spinal cord caused by too deep insertion of the reamer!

- Only use activ C instruments.
- Use the instruments according to their intended use.
- Only use the fin in combination with a correctly selected and positioned trial implant and with the chisel guide coupled.
- Take care that the depth stop of the trial implant abuts the anterior edge of the vertebral body.
- Always maintain lateral X-ray control while chiseling the fin.

- Carefully mount the chisel guide with chisel on the vertebral body.
- With slight pressure and at an even speed, chisel under lateral X-ray control. The chiseling depth is limited by the stop on the chisel guide 10.

- Carefully pull the chisel backwards out of the chisel guide.
- Remove the chisel guide 10 and the trial implant from the vertebral disk space.
- Clean the fin bed with nerve hooks from the discectomy set.

Disassembly

Uncoupling the reamer guide/handle from the trial implant

- Turn locking pipes 11 of reamer guide 3 or handle 2, respectively, to position "open" .
- Pull off the trial implant in axial direction.
- Decouple reamer 5 from reamer handpiece 6 as described in instructions for use TA011267.

Detaching the mounting block from the reamer handpiece

- Hold locking clip 13 of mounting block 4 pressed down (arrow marking points to position "open" ) and pull off the mounting block from reamer handpiece 6.

Uncoupling the chisel guide/handle from the trial implant

- Turn locking pipes 11  of chisel guide 10 or handle 2, respectively, to position "open".
- Pull off the trial implant in axial direction.
- Raise the chisel by the tip on the instrument and push out forwards.

activ C handle, reamer guide, chisel guide, trial implants

Validated reprocessing procedure

Note

National laws, national and international standards and directives, and product-specific hygiene regulations for processing must be observed.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful processing of this medical device can only be guaranteed following prior validation of the processing method. The operator/sterile processing technician is responsible for this.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Single-use products

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

➤ **Do not reprocess the product.**

General information

To prevent increased contamination of loaded instrument trays during use, please ensure that contaminated instruments are collected separately and not returned to the instrument tray.

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion of stainless steel. You should therefore leave no more than 6 hours between use and cleaning, not use pre-clean fixing temperatures >45°C, and ensure that non-fixing disinfectants (aldehyde-free/alcohol-free) are used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chloride - e.g. in surgical residues, drugs, saline solutions and water for cleaning, disinfection and sterilization - may cause corrosion damage to stainless steel (pinholing, stress corrosion), thus rendering the products unusable. These must be removed by rinsing thoroughly with demineralised water and then drying.

Only use process chemicals which have been tested and approved (e.g. VAH/DGHM or FDA-certified or CE marking) and recommended by the chemical manufacturer as being compatible with the materials. All the chemical manufacturer's application specifications regarding temperature, concentration and contact time should be strictly observed. Failure to do so can result in the following problems:

- Visible changes to materials e.g. fading or changes in the color of titanium or aluminum. As regards aluminum, visible changes to the surface may already occur at a pH level of >8 in the application/working solution.
- Material damage such as corrosion, cracks, fractures, premature deterioration or swelling.
- Do not use process chemicals which cause stress cracks or brittleness in plastics.
- Clean the product immediately after use.

Please see www.a-k-i.org for more detailed information on hygienically safe reprocessing methods which serve to protect the materials as well as retain their value.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Rinse non-visible surfaces such as those in instruments with concealed crevices, lumens or complex geometries, preferably with distilled water, using a disposable syringe for instance.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Disassemble the product prior to cleaning, see Disassembly.

Cleaning/disinfection

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for use e.g. on aluminum, plastic materials and stainless steel,
 - do not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 60°C.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfection.

activ C handle, reamer guide, chisel guide, trial implants

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.

Note

Clean reamer guide **3**, chisel guide **10** and handle **2** with locking pipes **11** in closed position () and do not move them during cleaning.

- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Where necessary, repeat the cleaning process.

Manual cleaning with immersion disinfection

FW789R–FW800R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	>15	2	D–W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfection	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
IV	Final rinse	RT (cold)	1	–	DM–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution for at least 15 minutes. Ensure that all accessible surfaces are moistened.
- Clean the product under running tap water with a suitable cleaning brush where necessary for as long as it takes to remove all discernable residues.
- Use a suitable cleaning brush (TA011944) to clean all surfaces which are not visible, e.g. in instruments with concealed crevices, lumens or complex geometries, for at least 1 min or until no further residues can be removed.
- Mobilize non-rigid components, such as set screws, joints, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Stage III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during disinfection.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Ensure that all accessible surfaces are moistened.

Stage IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse the lumens at least five times using a disposable syringe (20 ml) and a suitable irrigation adapter.
- Drain any remaining water fully.

Stage V

- Dry the product with a lint-free cloth or medical compressed air.

activ C handle, reamer guide, chisel guide, trial implants

Manual cleaning with ultrasound and immersion disinfection

FW870/FW871R/FW784R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	BBraun Stabimed; phenol-free, aldehyde-free and QAV-free
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	BBraun Stabimed; phenol-free, aldehyde-free and QAV-free
IV	Final rinse	RT (cold)	1	-	DM-W	-
V	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

Stage I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all discernable residues have been removed.
- Use a suitable cleaning brush (TA011944) to clean all surfaces which are not visible, e.g. in instruments with concealed crevices, lumens or complex geometries, for at least 1 min or until no further residues can be removed.
- Mobilize non-rigid components, such as set screws, joints, etc. during cleaning.

- After cleaning, use a 20-ml single-use syringe to rinse thoroughly, for at least 5 times, these difficult to access parts of the product.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Drain any remaining water fully.

Stage III

- Fully immerse the product in the disinfectant solution.
- Rinse lumens at least 5 times at the beginning of the exposure time, using a disposable syringe (20 ml) and an appropriate rinsing adapter. Ensure that all accessible surfaces are moistened.

Stage IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.

- Rinse the lumens at least five times using a disposable syringe (20 ml) and a suitable irrigation adapter.
- Drain any remaining water fully.

Stage V

- Dry the product thoroughly with a lint-free cloth or medical compressed air.

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note

For thermal disinfection, always use fully desalinated water (demineralized, low microbial contamination: drinking water quality at least) and ensure that Ao is >3 000 for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

activ C handle, reamer guide, chisel guide, trial implants

Manual pre-cleaning with a brush

FW789R–FW800R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D–W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Rinsing	RT (cold)	1	–	D–W	–

D–W: Drinking water

RT: Room temperature

Stage I

- Fully immerse the product in the cleaning/ disinfecting solution for at least 15 minutes. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all discernable residues have been removed.
- Use a suitable cleaning brush (TA011944) to clean all surfaces which are not visible, e.g. in instruments with concealed crevices, lumens or complex geometries, for at least 1 min or until no further residues can be removed.
- Mobilize non-rigid components, such as set screws, joints, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual Pre-Cleaning with ultrasound and brush

FW870/FW871R/FW784R

Note

Clean reamer guide **3**, chisel guide **10** and handle **2** with locking pipes **11** in closed position () and do not move them during cleaning.

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

Stage I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all discernable residues have been removed.
- Use a suitable cleaning brush (TA011944) to clean all surfaces which are not visible, e.g. in instruments with concealed crevices, lumens or complex geometries, for at least 1 min or until no further residues can be removed.
- Mobilize non-rigid components, such as set screws, joints, etc. during cleaning.
- After cleaning, use a 20-ml single-use syringe to rinse thoroughly, for at least 5 times, these difficult to access parts of the product.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.

activ C handle, reamer guide, chisel guide, trial implants

Mechanical alkaline cleaning and thermal disinfection

FW789R–FW800R, FW870, FW871R, FW784R

Machine type: single-chamber cleaning/disinfection device without ultrasound

- Place the instrument in a tray suitable for washing (to ensure that all parts are cleaned).
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Keep working ends open for cleaning.
- Place instruments in the tray with their hinges open.

Note

Clean reamer guide **3**, chisel guide **10** and handle **2** with locking pipes **11** in closed position () and do not move them during cleaning.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	BBRAUN HELIMATIC CLEANER alkaline with tensides; application solution 0.5%
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfection	90/194	5	FD–W	–
V	Drying	–	–	–	According to disinfectant program

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts (e.g. joints and links) with a maintenance oil suitable for the respective sterilization process (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble separable products.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Ensure that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization

Note

The product may only be sterilized in disassembled condition.

Note

To avoid breakage due to stress crack corrosion, sterilize the instruments with the lock fully open or locked no further than on the first ratchet tooth.

- Ensure that the sterilant can reach all external and internal surfaces (e.g. by opening valves and stopcocks).
- Validated sterilization process
 - Disassemble the product.
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer specified in DIN EN 285 and validated in accordance with DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer: ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

activ C handle, reamer guide, chisel guide, trial implants

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F– 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Technical Service

WARNING

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of warranty rights and applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems, LLC.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems, LLC.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG