

- (GB)** Instructions for use/Technical description
- (USA)** activC handle, reamer guide, chisel guide, trial implants
- (D)** Gebrauchsanweisung/Technische Beschreibung
- activC Handgriff, Fräserführung, Meißelführung, Probeimplantate
- (F)** Mode d'emploi/Description technique
- Poignée, guidage de fraise, guidage de burin, implants d'essai activC
- (E)** Instrucciones de manejo/Descripción técnica
- Mango, guía de fresado, guía de cinzel, implantes de prueba activC
- (I)** Istruzioni per l'uso/Descrizione tecnica
- Impugnatura, guida per la fresa, guida per lo scalpello, impianti di prova activC
- (P)** Instruções de utilização/Descrição técnica
- Cabo activC, guia de fresa, guia de cinzel, implantes de ensaio
- (NL)** Gebruiksaanwijzing/Technische beschrijving
- activC handgrip, freesgeleider, beitjelgeleider, proefimplantaten
- (S)** Bruksanvisning/Teknisk beskrivning
- activC handtag, frässtyrning, mejselstyrning, provimplantat
- (RUS)** Инструкция по применению/Техническое описание
- Ручка activC, направляющая фрезы, направляющая резца, тестовые имплантаты
- (CZ)** Návod k použití/Technický popis
- Rukojeť, vedení frézy, vedení dláta, zkušební implantáty activC
- (PL)** Instrukcja użytkowania/Opis techniczny
- Rękojeść, prowadnica freza, prowadnica dłuta, implanty próbne activC
- (SK)** Technical alterations reserved
- Rukoväť, frézové vedenie, vedenie dláta, skúšobné implantáty activC
- (TR)** Kullanım Kılavuzu/Teknik açıklama
- activC el tutanağı, freze kılavuzu, keski kılavuzu, deneme implantları

FW647R

FW645R

B | BRAUN

SHARING EXPERTISE


Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 011996 03/14 V6 Änd.-Nr. 46271

Technical alterations reserved

Aesculap® activC

activC handle, reamer guide, chisel guide, trial implants

Legend

- 1 Trial implant
- 2 Handle
- 3 Reamer guide
- 4 Mounting block (for reamer guide)
- 5 Reamer
- 6 Reamer handpiece
- 7 HiLAN pneumatic motor
- 8 activC Chisel insert
- 9 Chisel holder
- 10 Chisel guide
- 11 Locking pipes
- 12 Adjusting wheel
- 13 Locking clip

Symbols on product and packages

	Not for reuse in intended applications as defined by the manufacturer
	Caution, general warning symbol Caution, see documentation supplied with the product
	"Closed" position
	"Open" position
	Ready

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

The trial implants are used in conjunction with the handpiece to verify the size, height and shape of the activC intervertebral disk prosthesis to be implanted, and in conjunction with the reamer guide or the chisel guide to prepare the fin bed for the inferior implant plate of the activC intervertebral disk prosthesis.

Available sizes

Art. no.	Designation
FW789R ¹	Trial implant size XS, H5
FW790R ¹	Trial implant size XS, H6
FW791R ¹	Trial implant size S, H5
FW792R ¹	Trial implant size S, H6
FW793R ¹	Trial implant size M, H5
FW794R ¹	Trial implant size M, H6
FW790R ¹	Trial implant size L, H5
FW790R ¹	Trial implant size L, H6
FW797R ¹	Trial implant size XL, H5
FW798R ¹	Trial implant size XL, H6
FW799R ¹	Trial implant size XXL, H5
FW800R ¹	Trial implant size XXL, H6
FW8701 ¹	Handle
FW871R ²	Reamer guide (with mounting block)
GB771R ²	Hi-Line XS handpiece angled XL-I
GE700SU ⁴	Hi-Line XS activC drill XL-I D1, 5
FW784R ³	activC chisel guide
FW786R ³	activC chisel holder
FW787R ⁴	activC chisels

¹ Trays FW645R and FW647R

² only tray FW647R

³ only tray FW645R

⁴ single-use products

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information about product training, please contact your national B. Braun/ Aesculap agency.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

Art. no.	Designation
GE700SU	Hi-Line XS activC drill XL-I D1, 5
FW787R	activC Chisels

The product is gamma-sterilized and supplied in sterile packaging.
The product must not be reused.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Do not use products from open or damaged sterile packaging.
- Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Do not use the product after its use-by date.

Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

Note
Follow OP-Manual no. 031302 and instructions for use "activC intervertebral disk prosthesis", see TA011995.

Mounting the trial implant of the selected size on the handle

- Turn locking pipes 11 of handle 2 to "open" position see Fig. 1.
- Mount the trial implant with the selected size, height and shape onto the handle 2. Note the "CRANIAL" and "CAUDAL" markings on the handle and on the trial implant.
- Turn locking pipes 11 of handle 2 to "closed" position .
- Check the trial implant for firm seating.

Compression of the spinal canal and of other posterior elements caused by trial implants inserted too deeply!

- Prior to insertion of trial implants, adjust the depth stop to minimum insertion depth.
- Insert trial implants in the intervertebral space under X-ray control.

- Use adjusting wheel 12 to set the depth stop of the trial implant to minimum insertion depth.
- Insert the trial implant in the intervertebral space under X-ray control. Ensure correct cranial and caudal alignment of the trial implant.
- Align the markings of the trial implant center to the midline marking.
- If necessary, correct the position of the depth stop with adjusting wheel 12 until the intended position is reached. Ensure that depth stop abuts to the anterior edge of the vertebral body.

Increased risk of migration if the selected implant size is too small!

- Select the trial implant or implant size, respectively, to achieve maximum coverage of the vertebral body endplates.
- On the X-ray, in AP and lateral view, check to ensure that the size, position and alignment of the trial implant are correct.

Irritation of the facet joints and strain on spinal cord and nerve roots due to excessive or insufficient distraction!

- Select the correct trial implant/ implant height.
- Check for correct height of the trial implant in a lateral X-ray.

- On the intraoperative X-ray image, in AP and lateral view, check the size, height and position of the trial implant.
- With the trial implant properly positioned, turn locking pipes 11 to "open" position and pull off handle 2 carefully in axial direction from the trial implant.
The trial implant remains in the intervertebral space.

Note
The fin bed can either by prepared by reaming or chiseling. This can be carried out either with the reamer guide with handpiece or the chisel guide with chisel holder, depending on which set is used.

Fin bed reaming

Assembling the reamer guide

- Note*
The handle and the reamer guide are fitted with identical locking mechanisms.
- ▶ Turn locking pipes 11 of reamer guide 3 to "open" position 
 - ▶ Mount the reamer guide 3 onto the trial implant in the intervertebral space. Note the "CRANIAL" and "CAUDAL" markings on the reamer guide and on the trial implant.
 - ▶ Turn locking pipes 11 of reamer guide 3 to "closed" position 
 - ▶ Reamer guide 3 is fixed to the trial implant.
 - ▶ Check again, by X-ray imaging, the position of the trial implant.

Assembling the reamer handpiece

- 

WARNING

Damage to nerves, spinal cord and other tissue because of fin bed reamed too deeply due to using too long a reamer handpiece or reamer.

 - ▶ Only use reamer handpiece Hi-Line XS XL-1 (GB771R).
 - ▶ Only use activC reamer (GE700SU).
- ▶ Prior to mounting, verify that the safety stop of the pneumatic motor is in position "0".
 - ▶ Assemble the compressed-air tube, the HiLAN pneumatic motor and reamer handpiece 6 as described and illustrated in the respective instructions for use, see Fig. 2.
 - ▶ Hold locking clip 13 of mounting block 4 pressed down and push the mounting block over the reamer handpiece 6 see Fig. 3.
 - ▶ Slide on mounting block 4 as far as it will go and release locking clip 13, see Fig. 4.
 - ▶ Pull at the mounting block 4 to check whether it is engaged properly. The arrow marking on the locking clip 13 must point to the "closed"  pictogram.
 - ▶ Mount reamer 5 as described in the instructions for use of reamer handpiece 6.

- 

WARNING

Risk of injury to blood vessels, nerves, spinal cord or other tissue if wrong reamer is used!

 - ▶ Only use the activC reamer GE700SU.
 - ▶ Check for technical faults prior to use.
 - ▶ Do not use bent reamers.
 - ▶ Do not reuse reamers after first use.
- 

DANGER

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

 - ▶ Do not reprocess the product.

Note
The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk. The fin bed must only be prepared for the inferior implant plate.

- 

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

 - ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
 - ▶ Make certain that the position of the trial implant does not change during reaming.
- 

WARNING

Risk of injury caused by sharp chisel!

 - ▶ To avoid injuries to critical blood vessels, nerves and other tissue, as well as injuries to the operating surgeon, always work with appropriate care when a reamer is clamped in the machine.
 - ▶ Do not activate the motor before the handpiece with the coupled reamer is properly positioned in the guide hole and the guiding slot of the reamer guide.
- ▶ A: Insert the reamer handpiece 6 with the reamer clamped in into the front guide hole of the reamer guide 3, see Fig. 5.
 - ▶ B: Insert the laterally protruding pins of the mounting block 4 into the rear guide slots of the reamer guide 3, see Fig. 5.
 - ▶ With the hand control not actuated, push the safety stop to . The system is ready for operation.

Note
The reaming of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk.

- 

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

 - ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
 - ▶ Make certain that the position of the trial implant does not change during reaming.
- 

WARNING

Risk of injury to nerves and spinal cord caused by too deep insertion of the reamer!

 - ▶ Use activC instruments only.
 - ▶ Use the instruments according to their intended use.
 - ▶ For reaming the holes, apply the reamer only in combination with a correctly selected and positioned trial implant and with the reamer guide coupled.
 - ▶ Take care that the depth stop of the trial implant abuts the anterior edge of the vertebral body.
 - ▶ Always maintain lateral X-ray control while reaming the holes.

- ▶ To ream the first hole under lateral X-ray control: With the motor running, push reamer handpiece 6 at an even speed and with mild forward pressure into the first guiding slot down to the stop.
- ▶ Pull back reamer handpiece 6 and, after moving the pins of mounting block 4 in the second guiding slot, repeat the reaming procedure under lateral X-ray control.
The reaming depth will be limited by the stop on the reamer guide 3.
- ▶ With the hand control not actuated, shift the safety stop on the pneumatic motor to position "0".
- ▶ Cautiously pull out reamer handpiece 6 first from the rear and then from the front reamer guide 3.
- ▶ Remove reamer guide 3 and the trial implant from the vertebral disk space. Clean the fin bed with nerve hooks from the discectomy set.

Chiseling the fin bed

Assembling the chisel guide

- Note*
The handle and the chisel guide are fitted with identical locking mechanisms.
- ▶ Turn locking pipes 11 of chisel guide 10 to "open" position 
 - ▶ Mount the chisel guide 10 onto the trial implant in the intervertebral space. Note the "CRANIAL" and "CAUDAL" markings on the handle and on the trial implant.
 - ▶ Turn locking pipes 11 of chisel guide 10 to "closed" position 
 - ▶ Chisel guide 10 is fixed to the trial implant.
 - ▶ Check again, by X-ray imaging, the position of the trial implant.

- 

WARNING

Damage to nerves, spinal cord and other tissue because of fin bed reamed too deeply due to using too long a chisel.

 - ▶ Only use chisel holder FW786R.
 - ▶ Only use chisel FW787R.
- ▶ Mount the chisel onto the chisel holder 9, see Fig. 6.
 - ▶ To do this mount the chisel sideways onto the chisel holder and push it backwards until it engages. Pull at the chisel, to check that it is engaged properly.

- 

WARNING

Risk of injury to blood vessels, nerves, spinal cord or other tissue if wrong chisel is used!

 - ▶ Only use chisel FW787R.
 - ▶ Do not use bent chisels. Do not reuse chisels after first use.

Note
The chiseling of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk. The fin bed must only be prepared for the inferior implant plate.

- 

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

 - ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
 - ▶ Make certain that the position of the trial implant does not change during reaming.
- 

WARNING

Risk of injury caused by sharp chisel!

 - ▶ To avoid injuries to critical blood vessels, nerves and other tissue, as well as injuries to the operating surgeon, always work with appropriate care when a chisel is clamped in the chisel holder.

- ▶ A: Carefully push the chisel holder 9 with the chisel through the rear chisel guide, see Fig. 7.
- ▶ B: Carefully push the chisel holder 9 with the chisel through the front chisel guide, see Fig. 7.

Note
The chiseling of the fin bed determines the position and orientation of the disk prosthesis in the plane of the intervertebral disk.

- 

WARNING

Malpositioning of the intervertebral disk prosthesis due to failure to observe the midline!

 - ▶ Check the central position of the trial implant by means of the midline marking and on the AP-view X-ray image.
 - ▶ Make certain that the position of the trial implant does not change during chiseling.
- 

WARNING

Risk of injury to nerves and spinal cord caused by too deep insertion of the reamer!

 - ▶ Use activC instruments only.
 - ▶ Use the instruments according to their intended use.
 - ▶ Only use the fin in combination with a correctly selected and positioned trial implant and with the chisel guide coupled.
 - ▶ Take care that the depth stop of the trial implant abuts the anterior edge of the vertebral body.
 - ▶ Always maintain lateral X-ray control while chiseling the fin.

- ▶ Carefully mount the chisel guide with chisel on the vertebral body.
- ▶ With slight pressure and at an even speed, chisel under lateral X-ray control. The chiseling depth is limited by the stop on the chisel guide 10.
- ▶ Carefully pull the chisel backwards out of the chisel guide.
- ▶ Remove chisel guide 10 and the trial implant from the vertebral disk space.
- ▶ Clean the fin bed with nerve hooks from the discectomy set.

Disassembling

Uncoupling the reamer guide/handle from the trial implant

- ▶ Turn locking pipes 11 of reamer guide 3 or of handle 2 to "open" position 
- ▶ Pull off the trial implant in axial direction.
- ▶ Decouple reamer 5 from reamer handpiece 6 as described in instructions for use TA011267.

Detaching the mounting block from the reamer handpiece

- ▶ Hold locking clip 13 of mounting block 4 pressed down (arrow marking points to "open" position ) and pull off the mounting block from reamer handpiece 6.

Uncoupling the chisel guide/handle from the trial implant

- ▶ Turn locking pipes 11 of chisel guide 10 or of handle 2 to "open" position 
- ▶ Pull off the trial implant in axial direction.
- ▶ Raise the chisel by the tip on the instrument and push out forwards.

Validated reprocessing procedure

General safety instructions

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this. The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Single-use products

WARNING

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

► **Do not reprocess the product!**

Art. no.	Designation
GE700SU	Hi-Line XS activC drill XL-I D1, 5
FW787R	activC chisels

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary. Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for use, for example, on aluminum, plastic materials, and high-grade steel,
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum permitted cleaning temperature of 60 °C.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection FW789R–FW800R	- Cleaning brush: TA011944 20 ml disposable syringe - Drying phase: Use a lint-free cloth or medical compressed air - Clean reamer guide 3, chisel guide 10 and handle 2 with locking pipes 11 in closed position (H) and do not move them during cleaning.	Chapter Manual cleaning/ disinfection and sub-chapter: Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection - FW870 - FW871R - FW784R	- Cleaning brush: TA011944 20 ml disposable syringe - When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. - Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/ disinfection and sub-chapter: Chapter Manual cleaning with ultrasound and immersion disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection FW789R–FW800R	- Cleaning brush: TA011944 20 ml disposable syringe - Mount jaws protection on the product. - Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). - Connect components with lumens and channels directly to the rinsing port of the injector carriage. - Keep working ends open for cleaning. - Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: - Chapter Manual pre-cleaning with a brush - Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection - FW870 - FW871R - FW784R	- Cleaning brush: TA011944 20 ml disposable syringe - Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). - Connect components with lumens and channels directly to the rinsing port of the injector carriage. - Keep working ends open for cleaning. - Place products in the tray with their hinges open. - Clean reamer guide 3, chisel guide 10 and handle 2 with locking pipes 11 in closed position (H) and do not move them during cleaning.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: - Chapter Manual pre-cleaning with ultrasound and brush - Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D–W	-
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D–W:	Drinking water
FD–W:	Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT:	Room temperature
*Recommended:	BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

- Phase III**
- Fully immerse the product in the disinfectant solution.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

- Phase IV**
- Rinse/flush the product thoroughly (all accessible surfaces).
 - Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
 - Rinse lumens with an appropriate disposable syringe at least five times.
 - Drain any remaining water fully.

- Phase V**
- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I**
- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - Drain any remaining water fully.

- Phase III**
- Fully immerse the product in the disinfectant solution.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
 - Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

- Phase IV**
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
 - Rinse lumens with an appropriate disposable syringe at least five times.
 - Drain any remaining water fully.

- Phase V**
- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I**
- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concen- trate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I**
- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II**
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disin- fecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- This product should be assembled prior to sterilization
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note

The product may only be sterilized when dismantled.

Note

To avoid breakage due to stress crack corrosion, sterilize the product with the lock fully open or locked no further than on the first ratchet tooth.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- ▶ Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Technical Service

WARNING

Risk of injury and/or malfunction!

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of warranty rights and applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Implant Systems LLC
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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