

Aesculap Spine activ L

-
- | | |
|------------|--|
| GB | Instructions for use/Technical description |
| USA | activ L insertion instrument |
| D | Gebrauchsanweisung/Technische Beschreibung |
| | activ L-Einsetzinstrument |
| F | Mode d'emploi/Description technique |
| | Instrument d'insertion activ L |
| E | Instrucciones de manejo/Descripción técnica |
| | Instrumento de inserción activ L |
| I | Istruzioni per l'uso/Descrizione tecnica |
| | Impattatore activ L |
| P | Instruções de utilização/Descrição técnica |
| | Instrumento introdutor activ L |
| NL | Gebruiksaanwijzing/Technische beschrijving |
| | activ L-inbrenginstrument |
| CZ | Návod k použití/Technický popis |
| | Nasazovací instrument activ L |

activ L insertion instrument

Legend

- 1 Spacer
- 2 Clamp
- 3 Clamping sleeve
- 4 Insertion instrument handle
- 5 Slotted-hole pin
- 6 Coupling pin
- 7 Superior prosthesis plate
- 8 Polyethylene inlay
- 9 Inferior prosthesis plate
- 10 Unlocking button for spacer
- 11 Twist cap

Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The activ L insertion instrument is used for implanting the intervertebral disc prosthesis through anterior (0°) or lateral (approx. 45°) approach. Different color-coded and labeled insertion instruments are used for implants of different height.

Available sizes

Art. no.	Size
FW961R	8.5 mm
FW962R	10 mm
FW963R	12 mm
FW964R	14 mm

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only for its intended purpose, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or using a machine, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

Coupling the implant

- Select the appropriate insertion instrument according to the height of the intervertebral disc prosthesis.
- Completely assemble the intervertebral disc prosthesis.

- Extend spacer **1** as far as possible. If necessary, change the position of the spacer by turning twist cap **11**.


Fig. 1 Coupling for anterior access (0°)


Fig. 2 Coupling for lateral access (approx. 45°)

- Put the intervertebral disc prosthesis onto clamp **2** of the insertion instrument, observing the marking on the insertion instrument. Inferior prosthesis plate **9** must butt against the part of clamp **2** marked "CAUDAL" of the insertion instrument; superior prosthesis plate **7** must butt against the part marked as "CRANIAL".
- Guide prosthesis plates on spacer **1** across to position pins **5** and **6** of the clamp **2**.
- Make certain that the coupling pins of the clamp fully engage in the coupling holes of the prosthesis plates.
- Turn clamping sleeve **3** clockwise until the intervertebral disc prosthesis is securely couples.

Decoupling the implant

- Release clamping sleeve **3**.
- Carefully detach the insertion instrument from the implant.
- If the insertion instrument cannot be detached without applying force: Turn twist cap **11** anticlockwise to move back spacer **1**.

Disassembly

- Turn twist cap **11** counterclockwise until it can be removed.
- Hold unlocking button **10** pressed down.
- Pull out spacer **1** to the front.

Assembly

- Insert spacer **1** from the front through clamp **2** and, with unlocking button **10** pressed down, push in the clamp as far as it will go.
- Release unlocking button **10**.
- Attach twist cap **11** and screw it on as far as it will go.
- To position spacer **1** as required: Turn twist cap **11** accordingly.

activ L insertion instrument

Validated reprocessing procedure

Note

National laws, national and international standards and directives, and product-specific hygiene regulations for reprocessing must be observed.

Note

Machine reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing procedure. The operator/sterile processing technician is responsible for this.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General information

To prevent increased contamination of loaded instrument trays during use, please ensure that contaminated instruments are collected separately and not returned to the instrument tray.

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion of stainless steel. You should therefore leave no more than 6 hours between use and cleaning, not use pre-clean fixing temperatures >45 °C, and ensure that non-fixing disinfectants (aldehyde-free/alcohol-free) are used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chloride – e.g. in surgical residues, drugs, saline solutions and water for cleaning, disinfection and sterilization – may cause corrosion damage to stainless steel (pinholing, stress corrosion), thus rendering the products unusable.

These must be removed by rinsing thoroughly with demineralized water and then drying.

Only use process chemicals which have been tested and approved (e.g. VAH/DGHMöder FDA-certified or CE marking) and recommended by the chemical manufacturer as being compatible with materials. All the chemical manufacturers' application specifications regarding temperature, concentration and contact time should be strictly observed. Failure to do so may result in the following problems:

- Visible changes to materials e.g. fading or changes in the color of titanium or aluminum. As regards aluminum, visible changes to the surface may already occur at a pH level of >8 in the application/working solution.
- Material damage such as corrosion, cracks, fractures, premature deterioration or swelling.
- Do not use process chemicals which cause stress cracks or brittleness in plastics.
- Clean the product immediately after use.

Please see www.a-k-i.org for more detailed information on hygienically safe reprocessing which is protective of materials and retains instruments' value.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open up instruments with hinges.
- Rinse non-visible surfaces such as those in instruments with concealed crevices, lumens or complex geometries, preferably with distilled water, using a disposable syringe for instance.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry instrument in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Disassemble the product prior to cleaning, see Disassembly.

Cleaning/disinfection

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for e.g. aluminum, plastics and high-grade steel, according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 60 °C.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection,
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting,
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting,
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products using a machine if they can be securely placed in the machine or positioning aid to ensure thorough cleaning.

Manual cleaning/disinfection

- Keep working ends open for cleaning.
- When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Where necessary, repeat the cleaning process.

activ L insertion instrument

Manual Cleaning with Ultrasound and Immersion Disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	5	2	T-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	-	T-W	-
III	Disinfection	RT (cold)	15	2	T-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	1	-	VE-W	-
V	Drying	-	-	-	-	-

D-W: Drinking water

DM-W: Demineralized water

RT: Room temperature

Phase I

- Mount jaws protection on the product.
- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Remove jaws protection.
- Clean the instrument with a suitable cleaning brush for as long as it takes to remove all discernable residues.
- Brush through surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using an appropriate cleaning brush (brush head length: 30 mm/Ø: 4,5 mm, e.g. TA011944).
- Maneuver non-rigid components, such as adjustable screws, hinges etc during cleaning.

- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Drain any remaining water fully.

Phase III

- Fully immerse the instrument in the disinfectant solution.
- Rinse the lumens at least five times at the beginning of the contact period using a disposable syringe (20 ml) and a suitable irrigation adapter. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Rinse the lumens at least five times using a disposable syringe (20ml) and a suitable irrigation adapter.
- Drain any remaining water fully.

Phase V

- Thoroughly dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note

For thermal disinfection, always use demineralized water. Ensure that Ao is >3 000 for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

activ L insertion instrument

Manual Pre-Cleaning with Ultrasound and Brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	15	2	T-W	BBraun Stabimed; aldehyde-free, phenol-free and QAV-free
II	Rinsing	RT (cold)	1	-	T-W	-

D-W: Drinking water

RT: Room temperature

Phase I

- Mount jaws protection on the product.
- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Remove jaws protection.
- Clean the instrument with a suitable cleaning brush for as long as it takes to remove all discernable residues.
- Brush through surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using an appropriate cleaning brush (brush head length: 30 mm/Ø: 4,5 mm, e.g. TA011944).
- Maneuver non-rigid components, such as adjustable screws, hinges etc during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.

Mechanical neutral or mild-alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

- Place the instrument in a screen basket suitable for washing (to ensure that all parts are cleaned).
- Keep working ends open for cleaning.

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	T-W	-
II	Cleaning	55/131	10	VE-W	Neutral: BBRAUN HELIMATIC CLEANER neutral, pH neutral, process solution 0.5 % Mildly alkaline: concentrate, pH = 9.5 < 5 % anionic tensides - 0.5 % solution
III	Intermediate rinse	>10/50	1	VE-W	-
IV	Thermal disinfection	90/194	5	VE-W	-
V	Drying	-	-	-	According to disinfectant program

D-W: Drinking water

DM-W: Demineralized water

- Check visible surfaces for residues after mechanical cleaning/disinfecting.
- Repeat the cleaning process if necessary.

activ L insertion instrument

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).
- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable instruments, see Assembly.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product into its holder or put it on a suitable screen basket. Ensure that all cutting edges are protected.
- Package screen baskets appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization

Note

The product may only be sterilized in disassembled condition.

- Ensure that the sterilant can reach all external and internal surfaces (e.g. by opening valves and stopcocks).
- Validated sterilization process
 - Disassemble the instrument
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer specified in DIN EN 285 and validated in accordance with DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer: Ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F— 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems, LLC.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems, LLC.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG