

Aesculap Spine activ L

-
- | | |
|--------------|--|
| (GB) | Instructions for use/Technical description
activ L insertion instrument |
| (USA) | |
| (D) | Gebrauchsanweisung/Technische Beschreibung
activ L-Einsetzinstrument |
| (F) | Mode d'emploi/Description technique
Instrument d'insertion activ L |
| (E) | Instrucciones de manejo/Descripción técnica
Instrumento de inserción activ L |
| (I) | Istruzioni per l'uso/Descrizione tecnica
Impattore activ L |
| (P) | Instruções de utilização/Descrição técnica
Instrumento introdutor activ L |
| (NL) | Gebruiksaanwijzing/Technische Beschrijving
activ L-inzetinstrument |
| (CZ) | Návod k použití/Technický popis
Nasazovací instrument activ L |

activ L insertion instrument

Legend

- 1 Spacer
- 2 Clamp
- 3 Clamping sleeve
- 4 Insertion instrument handle
- 5 Slotted-hole pin
- 6 Coupling pin
- 7 Superior prosthesis plate
- 8 Polyethylene inlay
- 9 Inferior prosthesis plate
- 10 Unlocking button for spacer
- 11 Twist cap

Intended use

The activ L insertion instrument is used for implanting the intervertebral disc prosthesis through anterior (0°) or lateral (approx. 45°) approach. Different color-coded and labeled insertion instruments are used for implants of different height.

Available sizes

Art. no.	Size
FW961R	8.5 mm
FW962R	10 mm
FW963R	12 mm
FW964R	14 mm

Safe handling and preparation**CAUTION**

Federal law restricts this device to sale by or on order of a physician!

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside damaged products immediately.
- Replace any damaged components immediately with original spare parts.

Safe operation

Coupling the implant

- Select the appropriate insertion instrument according to the height of the intervertebral disc prosthesis.
- Completely assemble the intervertebral disc prosthesis.
- Extend spacer **1** as far as possible. If necessary, change the position of the spacer by turning twist cap **11**.


Fig. 1 Coupling for anterior access (0°)


Fig. 2 Coupling for lateral access (approx. 45°)

- Put the intervertebral disc prosthesis onto clamp **2** of the insertion instrument, observing the marking on the insertion instrument. Inferior prosthesis plate **9** must butt against the part of clamp **2** marked "CAUDAL" of the insertion instrument; superior prosthesis plate **7** must butt against the part marked as "CRANIAL".
- Guide the prosthesis plates along spacer **1** onto coupling pins **5** and **6** of clamp **2**.
- Make certain that the coupling pins of the clamp fully engage in the coupling holes of the prosthesis plates.
- Turn clamping sleeve **3** clockwise until the intervertebral disc prosthesis is securely coupled.

Decoupling the implant

- Release clamping sleeve **3**.
- Carefully detach the insertion instrument from the implant.
- If the insertion instrument cannot be detached without applying force: Turn twist cap **11** anti-clockwise to move back spacer **1**.

activ L insertion instrument

Disassembling

- Turn twist cap **11** counterclockwise until it can be removed.
- Hold unlocking button **10** pressed down.
- Pull out spacer **1** to the front.

Assembling

- Insert clamp **1** from the front through clamp **2** and, with unlocking button **10** pressed down, push in the clamp as far as it will go.
- Release unlocking button **10**.
- Attach twist cap **11** and screw it on as far as it will go.
- To position spacer **1** as required: Turn twist cap **11** accordingly.

Processing

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net.

Preparations

- Carry out non-fixating/NaCl-free pre-cleaning immediately after use.
- Disassemble separable products, see Disassembling.
- Reprocess the product immediately after use.
- If possible, put the product away in dry condition.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.
- If necessary, clean the product through ultrasound treatment, see Cleaning/Disinfecting.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- **Only use cleaning/disinfecting agents that are approved for the product, and follow the manufacturers' instructions.**
- **Do not exceed the maximum allowable cleaning temperature of 60 °C.**
- Use cleaning/disinfecting agents that are suitable for the product. Always follow the manufacturer's instructions with regard to concentration, temperature and exposure time.
- Avoid encrustation of residues/proteins (e.g. caused by aldehyde/alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - to prepare products with encrusted debris for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - as an aftertreatment for products that are still dirty after mechanical cleaning/disinfection.
- After chemical disinfection, rinse the product thoroughly under running water. Always adhere to the manufacturers' instructions.
- Products with fine working tips must be cleaned and disinfected manually.

Mechanical cleaning/disinfecting

- Place the product on a tray that is suitable for cleaning (to avoid rinsing blind spots).
- Ensure that water can flow out of openings.
- Process the product in a cleaning/disinfecting machine. Follow the instructions provided by the manufacturer of the machine.
- Processing cycle:
 - Use a suitable neutral/mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
 - Do not exceed the maximum cleaning temperature of 60 °C.
 - Clean for at least 10 min.
 - Neutralize if necessary.
 - Apply intermediate rinse for at least 1 min.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
 - For thermal disinfection: Rinse for 10 min at 93 °C with distilled, demineralized or fully desalinated water.
 - Complete the program with a drying phase of at least 10 min at a temperature of not more than 100 °C.

- After completion of the mechanical cleaning/disinfecting cycle, inspect surfaces, cavities, lumens, and openings for visible debris.
- Clean manually if necessary.

Manual cleaning/disinfecting

- Use a suitable neutral/mildly alkaline cleaning/disinfecting agent according to the manufacturer's instructions.
- Immerse the product in the cleaning/disinfecting agent in such a way that all surfaces, cavities, lumens, and openings are covered.
- After the end of the disinfection period, thoroughly rinse the product under running water, making sure that water flows through every lumen and channel and all blind holes are repeatedly filled and drained.
- Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Clean lumens, channels and blind holes with soft round plastic brushes of fitting diameter.
- Carry out an intensive final rinse with distilled, demineralized or fully desalinated water.
- Inspect surfaces, cavities, lumens and openings for visible debris. If necessary, repeat the cleaning/disinfection process.
- Use a lint-free cloth or a compressed-air gun for drying the product.
- Make certain that lumens, channels and blind holes are dried, too.

activ L insertion instrument

Control, care and inspection

- Allow the product to cool down to room temperature.

Contamination of the patient caused by unsterile product!

- **Do not lubricate the product with sterilizable, steam-permeable maintenance oil under any circumstances.**

- Assemble the separable product, see Assembling.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside damaged products immediately.

Packaging

- Appropriately protect products with fine working tips.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that any cutting edges are protected. Observe the weight limit for each tray/container.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

Note

The product can be sterilized either in disassembled or in assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Apply steam sterilization, observing the following rules:

Sterilization has to be performed in a validated steam sterilization process (e.g. in a sterilizer according to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, validated according to EN 554/ISO 13683). For the fractionated vacuum process, the sterilization has to be carried out with 134 °C/ 2-bar program for a minimum holding time of 5 minutes.

- When sterilizing several products at the same time in one steam sterilizer: Make sure that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270— 275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store reprocessed products under conditions as germ-free as possible, in a dry, dark, cool and dust-protected room.

Technical Service

For service, maintenance or repairs, contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Designation Product	Art. no.	Designation Replacement part	Ordering no.	Ordering no. (USA)
Insertion instrument	FW961R	Spacer, height: 8.5 mm	FW961801	FW946SU
Insertion instrument	FW962R	Spacer, height: 10 mm	FW962801	FW947SU
Insertion instrument	FW963R	Spacer, height: 12 mm	FW963801	FW948SU
Insertion instrument	FW964R	Spacer, height: 14 mm	FW964801	FW949SU

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG