



Aesculap® activ L

Aesculap Spine

- GB** Instructions for use/Technical description
activ L revision instruments and distraction forceps
- USA** Instructions for use/Technical description
activ L revision instruments and distraction forceps
Note for U.S. users
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
- D** Gebrauchsanweisung/Technische Beschreibung
activ L-Revisionsinstrumente und -Sprezzangen
- F** Mode d'emploi/Description technique
Instruments de révision et pinces d'écartement activ L
- E** Instrucciones de manejo/Descripción técnica
Instrumentos y tenazas de revisión activ L
- I** Istruzioni per l'uso/Descrizione tecnica
Strumenti di revisione e pinze divaricatrici activ L
- P** Instruções de utilização/Descrição técnica
Instrumento de revisão e tenaz activ L
- NL** Gebruiksaanwijzing/Technische beschrijving
activ L-revisie-instrumenten en -spreidtangen
- S** Bruksanvisning/Teknisk beskrivning
activ L-revisionsinstrument och -spridartänger
- RUS** Инструкция по применению/Техническое описание
Инструменты для проведения ревизии и дистракционные щипцы activ L
- CZ** Návod k použití/Technický popis
Revizní nástroje a dистракční kleště activ L
- PL** Instrukcja użytkowania/Opis techniczny
Przyrządy rewizyjne i kleszczyki rozszerzające activ L
- SK** Technical alterations reserved
Revizne nástroje a rozťahovacie kliešte activ L
- TR** Kullanım Kılavuzu/Teknik açıklama
activ L revizyon ekipmanları ve dистраkasyon penseleri

B|BRAUN
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 011904 04/14 V6 Änd.-Nr. 49358

CE 0123 - DIR 93/42/EEC

Technical alterations reserved

Aesculap® activ L

activ L revision instruments and distraction forceps

Legend

- A activ L angled distraction forceps FW960R
- B activ L angled revision distraction forceps FW965R
- C activ L revision instrument for endplates FW966R; size S/M
- D activ L revision instrument for endplates FW967R; size L/XL
- E activ L revision instrument for PE inlay FW968R
- F activ L handle for revision instrument FW998R
- 1 Handle parts
- 2 Lock
- 3 Locking lever
- 4 Working tip
- 5 Pins
- 6 Jaw piece
- 7 Endplates
- 8 Adjustment knob
- 9 Anchoring rods
- 10 Hook
- 11 Recesses
- 12 Lock nut
- 13 Pins
- 14 Slots
- 15 Tubular shaft

Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

The activ L revision instruments and distraction forceps are used for the distraction/mobilization of two adjacent lumbar vertebral bodies or prostheses plates, as well as for the repositioning and revision of activ L intervertebral disk prostheses.

activ L angled distraction forceps FW960R

The activ L angled distraction forceps is used for the distraction/mobilization of two adjacent vertebral bodies.

activ L angled revision distraction forceps FW965R

The activ L angled revision distraction forceps is used for distracting the two activ L metal plates of an implanted activ L prosthesis in order to remove the PE inlay in corrective surgery.

activ L revision instrument for endplates FW966R/FW967 and handle FW998R

The activ L revision instrument for endplates and the corresponding handle are used for the repositioning and revision of activ L intervertebral disk prostheses.

activ L revision instrument for PE inlay FW968R

The activ L revision instrument for PE inlay is used for the revision of the activ L PE inlay.

Safe handling and preparation

Use of the activ L revision instruments and distraction forceps requires precise knowledge about spine surgery and the stabilization and biomechanical situation at the spine.

Surgical application of the activ L revision instruments and distraction forceps is described in detail in the associated operating manual.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

Risk of injury and/or malfunction!
► Always carry out a function check prior to using the product.

Distracting/mobilizing vertebral bodies with the activ L angled distraction forceps FW960R

Risk of fracture of the vertebral body endplates and/or injury to surrounding structures due to overdistracting!
► Apply proper care when distracting.
► Avoid excessive distraction.

Damage to the product caused by actuation of the locking lever!
► To release lock 2, gently press the handle parts 1 together and actuate locking lever 3.

Note

Mobilization of the intervertebral disk space is a critical process and should be performed under X-ray control.

Note

The intervertebral space distraction is often wedge-shaped. Sufficient load relief in the posterior intervertebral disk space is required so that parallel distraction can be performed over the entire intervertebral space.

- Insert the working end 4 into the intervertebral space, see Fig. 2. When doing so, ensure that the instrument is pushed as dorsally as possible into the intervertebral disk space.
- Distract/mobilize the distraction forceps FW960R by pressing the handle parts 1 together as required.

Note

Lock 2 engages automatically when handle parts 1 are released.

Removing the activ L angled distraction forceps FW960R

- Gently press the handle parts 1 together and actuate locking lever 3, see Fig. 1.
This releases the lock 2 of the distraction forceps FW960R.

Applying the activ L angled revision distraction forceps FW965R

Injury to surrounding blood vessels and soft tissue!
Damage to the product or to the activ L disk prosthesis!
► To protect the blood vessels, always use a retractor in conjunction with this product.
► Always insert the product into the two holes at the right-hand side of the prosthesis.

Damage to the product caused by actuation of the locking lever!
► To release the lock, gently press the handle parts together and actuate locking lever 3.

- Insert pins 5 of the distraction forceps FW965R into the holes on the right-hand side of the implant, see Fig. 4.
- Distract the distraction forceps FW965R by pressing the handle parts 1 together as required.

Note

Lock 2 engages automatically when handle parts 1 are released.

Removing the activ L angled revision distraction forceps FW965R

- Gently press the handle parts 1 together and actuate locking lever 3, see Fig. 3.
This releases the lock 2 of the distraction forceps FW965R.

Removing the PE inlay with activ L revision instrument for PE inlay FW968R from the intervertebral disk prosthesis

- Insert the jaw part 6 of the revision instrument FW968R between the two endplates 7 of the intervertebral disk prosthesis, see Fig. 5. Ensure the following:
 - The orientation of the revision instrument FW968R in the caudal direction corresponds with the "caudal" marking on the revision instrument.
 - The edges of the jaw part 6 are pointing in the direction of the inferior endplates.
- Push the revision instrument along the lower endplate, until the jaw part 6 engages on the PE inlay.
- Turn the adjustment knob 8 clockwise to fix the PE inlay in place on jaw part 6, see Fig. 6.
- Lift the revision instrument FW968R with fixated PE inlay over the edge of the inferior endplate and remove from the endplate by pulling in the anterior direction.
- To remove the PE inlay from the revision instrument FW968R:
 - Rotate the adjustment knob 8 in a counterclockwise direction.
 - Pull the PE inlay from the jaw part 6.

Reposition/remove endplates with activ L revision instrument for endplates FW966R/FW967R

Risk of injury to nerves and surrounding structures due to uncontrolled posterior pressure of the hooks of the anchoring rods!
► Apply proper care when manipulating the anchoring rods.
► Do not apply excessive pressure.

- Screw handle FW998R onto handpiece of revision instrument FW966R/FW967R, see Fig. 7.

Note

The pins 13 on the screw connection of the anchoring rods of the revision instrument will point in the same direction as the corresponding hooks on the front part.

Note

In case of an incomplete revision of the implant, the endplates are removed one after another. The slotted hammer is also used for the revision procedure.

Note

The endplates can grow into and adhere strongly to the vertebral bodies. This can considerably complicate the revision procedure. If the revision is of a prosthesis implanted a number of years ago, the endplates may need to be removed from the intervertebral body with a chisel or elevator prior to using the revision instrument.

- Turn the anchoring rods 9 so that both hooks 10 are pointing outwards, see Fig. 8.
- Push the revision instrument FW966R/FW967R between the endplates. When doing so, ensure that the hooks reach over the posterior edges of the plates.
- Turn both anchoring rods 9 90° in the direction of the plate to be removed.
- Pull back the anchoring rods 9. When doing so, ensure that the anchoring rods engage securely in the recesses 11 on the posterior side of the plate, see Fig. 9.
- Fix the anchoring rods 9 in place by turning the lock nut 12 counterclockwise. When doing so, hold the anchoring rods so that they cannot turn anymore.
- Carefully reposition or withdraw the plate.
- To remove the endplate completely, hit it out with the slotted hammer FW579R.

- ▶ Undo the lock nut **12** by turning it clockwise.
- ▶ Turn the anchoring rods **9** 90° in the anterior direction and remove the endplate.
- ▶ If applicable, remove other endplates in the same way.

Disassembly

activ L revision instrument for endplates FW966R/FW967R

- ▶ Unscrew handle FW998R of revision instrument FW966R/FW967R.
- ▶ Turn the anchoring rods **9** so that the hooks **10** are positioned in the direction of the slots **14**.
- ▶ Slide the anchoring rods **9** behind the handpiece until the hooks **10** are positioned on the middle section of the revision instrument FW966R/FW967R. When doing so, slide the hooks **10** through the slots **14**, see Fig. 10.
- ▶ Turn the lock nuts **12** counterclockwise until they are fully undone.
- ▶ Slide the lock nuts **12** along the anchoring rods behind the handpiece.

activ L revision instrument for PE inlay FW968R

- ▶ Turn the adjustment knob **8** counterclockwise and gently pull on the jaw part **6** until the jaw part comes out, see Fig. 11.

Assembly

activ L revision instrument for endplates FW966R/FW967R

- ▶ Slide the lock nuts **12** along the anchoring rods **9** as far as the threaded rods, see Fig. 10.
- ▶ Turn the lock nuts **12** clockwise onto the threaded rods.
- ▶ Turn the anchoring rods **9** so that the hooks **10** are positioned in the direction of the slots **14**.
- ▶ Slide the hooks **10** through the slots **14**.

activ L revision instrument for PE inlay FW968R

- ▶ Insert the jaw part **6** into the tubular shaft **15** of the revision instrument FW968R until it engages, see Fig. 11. When doing so, ensure that the "caudal" marking on the tubular shaft is aligned with the edge of the jaw part **6**.
- ▶ Push the jaw part **6** in the direction of the adjustment knob **8** and turn the adjustment knob clockwise until the jaw part is positioned just in front of the tubular shaft. When doing so, ensure that the jaw part is not compressed by the tubular shaft.

Validated reprocessing procedure

General safety instructions

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.
The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bb Braun.com
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- ▶ Open up instruments with hinges.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparations before cleaning

- FW966R, FW967R, FW968R:
- ▶ Disassemble the product prior to cleaning, see Disassembly.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents approved for, e.g., aluminum, plastics, and high-grade steel, according to the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 60 °C.

- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfection.
- ▶ Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FW960R ■ FW965R ■ FW968R ■ FW998R	■ Cleaning brush: TA011944 ■ 20 ml disposable syringe ■ Keep working ends open for cleaning. ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the hinge while cleaning. ■ Drying phase: Use lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subchapter: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ FW966R ■ FW967R	■ Cleaning brush: TA011944 ■ 20 ml disposable syringe ■ Keep working ends open for cleaning. ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW960R ■ FW965R ■ FW968R ■ FW998R	■ Cleaning brush: TA011944 ■ 20 ml disposable syringe ■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Keep working ends open for cleaning. ■ Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW966R ■ FW967R	■ Cleaning brush: TA011944 ■ 20 ml disposable syringe ■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Keep working ends open for cleaning. ■ Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subchapter: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- ▶ Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	–	FD–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concen- trate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disin- fecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).

- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

FW966R, FW967R, FW968R:

Note

The product may only be sterilized when dismantled.

Note

To avoid breakage due to stress crack corrosion, sterilize the instruments with the lock fully open or locked no further than on the first ratchet tooth.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the product (FW966, FW967, FW998)
 - Steam sterilization using fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

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