
Aesculap Spine

activ L

GB	Instructions for use
USA	activ L distraction forceps revision instrument
D	Gebrauchsanweisung
	activ L-Revisionsinstrument-Spreizgabel
F	Mode d'emploi
	Fourche d'écartement de révision activ L
E	Instrucciones de manejo
	Separador de revisión activ L
I	Istruzioni per l'uso
	Forcella divaricatrice-strumento da revisione activ L
P	Instruções de utilização
	Garfo distractor para instrumento de revisão activ L
NL	Gebruiksaanwijzing/Technische Beschrijving
	activ L revisie-instrument-spreidvork
CZ	Návod k použití
	Rozvírací vidlice activ L-revizního instrumentu

activ L distraction forceps revision instrument

Legend

- 1 Handle parts
- 2 Lock
- 3 Locking lever
- 4 Pins
- 5 Jaw parts

Symbols on product

Attention, see instructions for use

Intended use

The activ L distraction forceps revision instrument FW965R is used for distracting the two activ L metal plates of the implanted activ L prosthesis to enable the removal of the PE inlays in corrective surgery.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Use of the activ L distraction forceps revision instrument requires precise knowledge about spine surgery and the stabilization and biomechanical situation at the spine.

Surgical application of the activ L distraction forceps revision instrument is described, in detail, in the associated operating manual.

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

DANGER

Injury to surrounding blood vessels and soft tissue!

Damage to the product or to the activ L disk prosthesis!

- To protect the blood vessels, always also use a retractor when using this product.
- Always insert the product in the two drill holes at the right-hand side of the prosthesis.

CAUTION

Damage to the product caused by actuation of locking lever 3!

- To release lock 2, slightly compress handle parts 1 and actuate locking lever 3.

Applying the activ L distractor fork revision instrument


Fig. 1

Legend

4 Pins

- Insert pins 4 of the activ L distraction forceps revision instrument in the holes at the right-hand side of the implant.
- Spread open the activ L distraction forceps revision instrument as required, by compressing the handle parts 1.

Note

Lock 2 engages automatically when handle parts 1 are released.

activ L distraction forceps revision instrument

Releasing the activ L distraction forceps revision instrument


Fig. 2

Legend

1 Handle parts

2 Lock

3 Locking lever

- Slightly compress handle parts 1 and actuate locking lever 3.

This releases lock 2 of the activ L distraction forceps revision instrument.

Processing

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net.

Preparations

- Carry out non-fixating/NaCl-free pre-cleaning immediately after use.
- Open jointed products.
- Reprocess the product immediately after use.
- Put away the product in dry condition.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.
- If necessary, clean the product through ultrasound treatment, see Cleaning/Disinfecting.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- **Only use cleaning/disinfecting agents that are approved for the product, and follow the manufacturer's instructions.**
- **Do not exceed the maximum allowable cleaning temperature of 60 °C.**
- Use cleaning/disinfecting agents that are suitable for the product. Always follow the manufacturer's instructions regarding concentration, temperature and exposure time.
- Avoid encrustation of residues/proteins (e.g. caused by aldehyde/alcohol).
- Only use bactericidal, fungicidal and virucidal disinfecting agents.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - to prepare products with encrusted debris for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - as an aftertreatment for products that are still dirty after mechanical cleaning/disinfecting.
- After chemical disinfection, rinse the product thoroughly under running water. Always adhere to the manufacturer's instructions.
- Manually clean and disinfect products with fine working tips.

Mechanical cleaning/disinfecting

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Store any products with hinges or joints on the tray in such a way that the joints are open.
- Ensure that water can flow out of openings.
- Process the product in a cleaning/disinfecting machine. Follow the instructions provided by the manufacturer of the machine.
- Process cycle:
 - Use an appropriate neutral or mildly alkaline cleaning/disinfecting agent, following the manufacturer's instructions.
 - Do not exceed the maximum cleaning temperature of 60 °C.
 - Clean for at least 10 minutes.
 - Neutralize if necessary.
 - Apply intermediate rinse for at least 1 min.
 - Carry out intensive final rinse with distilled, demineralized or fully desalinated water.
 - For thermal disinfection: 10 min rinsing at 93 °C with distilled, demineralized or fully desalinated water.
 - Complete the program with a drying phase of at least 40 min at a temperature of not more than 100 °C.
- After completion of the mechanical cleaning/disinfecting cycle, inspect all surfaces, cavities, lumens, and openings for visible debris.
- Clean manually if necessary.

activ L distraction forceps revision instrument

Manual cleaning/disinfecting

- Use an appropriate neutral or mildly alkaline cleaning/disinfecting agent, following the manufacturer's instructions.
- Immerse the product in the cleaning/disinfecting agent in such a way that all surfaces, cavities, lumens and openings are covered.
- After the end of the disinfection period, thoroughly rinse the product under running water, ensuring that water flows through every lumen and channel, and all blind holes are repeatedly filled and drained.
- Clean hinged or jointed products in open and closed positions.
- Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Clean lumens, channels and blind holes with soft round plastic brushes of fitting diameter.
- Carry out an intensive final rinse with distilled, demineralized or fully desalinated water.
- Inspect surfaces, cavities, lumens and openings for visible debris. If necessary, repeat the cleaning/disinfection process.
- Use a lint-free cloth or a compressed-air gun for drying the product.
- Make certain that lumens, channels and blind holes are dried, too.

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts such as hinges and latches with a sterilizable, steam-permeable and tissue-compatible maintenance oil (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, and not damaged (e.g. insulation), and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside the product if it is damaged.

Packaging

- Appropriately protect products with fine working tips.
- Secure product with retractor in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected. Observe the weight limit for each tray/container.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters*Note*

The product may only be sterilized in fully assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Sterilize with steam, observing the following rules: Sterilization has to be performed through a validated steam sterilization process (e.g. in a sterilizer according to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, validated according to EN 554/ISO 13683). For the fractionated vacuum process, sterilization has to be carried out running the 134 °C/ 2-bar program for a minimum holding time of 5 minutes.
- When sterilizing several products at the same time in one steam sterilizer: Make sure that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270–275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store reprocessed products under conditions as germ-free as possible, in a dry, dark, cool and dust-protected room.

Technical Service

For service, maintenance or repairs, contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420
Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG